

custombiotech.roche.com

Please contact your local CustomBiotech representative

Europe, Middle East, Africa, Latin America

Phone +49 621 759 8580
Fax +49 621 759 6385
mannheim.custombiotech@roche.com

United States

Phone +1 800 428 5433, ext. 14649 (toll-free)
Fax +1 317 521 4065
custombiotech.ussales@roche.com

Canada

Phone +1 450 686 7050
Fax +1 450 686 7012
custombiotech.can@roche.com

Japan

Phone +81 3 6634 1046
Fax +81 3 5479 0585
japan.custombiotech@roche.com

Asia Pacific

Phone +65 6371 6638
Fax +65 6371 6601
apac.custombiotech@roche.com

Published by

Roche Diagnostics GmbH
Sandhofer Straße 116
68305 Mannheim
Germany

© 2017 Roche Diagnostics GmbH
All rights reserved.

May 2017



CustomBiotech Catalog 15th Edition

Clinical Chemistry and Immunology



CUSTOMBIOTECH

We are CustomBiotech from Roche

Making your vision a reality

You are driving a paradigm shift. Scientific and technological advances emerging from creative companies like yours are transforming the way we look at disease and breaking down barriers to novel treatments. The path to these innovations, however, is accompanied by many challenges. How do you secure a supply chain in a global landscape? How do you ensure compliance in a demanding regulatory environment? How do you remain competitive in a fast-paced market? Now more than ever, business success hinges on working with suppliers who understand your needs and allow you to focus on what matters most.

What if all your needs were met by a single partner — a market leader with decades of experience in Biotechnology and *in vitro* Diagnostics working side by side with your team?

Introducing Roche CustomBiotech.

Leveraging the unique know-how of Roche Diagnostics and Roche Pharmaceuticals, we work with you to deliver high quality raw materials, instrumentation, products and services for your biopharmaceutical, cell therapy, or *in vitro* diagnostics business, customized to meet your unique quality and regulatory needs. Our support builds on four core benefits to help you succeed, from development to commercialization.



Innovation

Count on great ideas that work. From research and development to manufacturing and logistics, our experts and facilities cover an unparalleled spectrum of skills and technologies to explore any idea.



Service

Count on accessible and experienced support. From business and regulatory issues to production troubleshooting, we help safeguard your operations and market standing with fast answers to problems, anytime and anywhere.



Customization

Invest time and resources into what you do best. For everything else, rely on us. With customized development, manufacturing, labeling, packaging and filling of components, we streamline the path of your product to market.



Security

Set your mind at ease. Our global reach, stringent standards and state-of-the-art manufacturing mean a secured supply of products and services to drive your business forward, when and where needed.



Security and support with CustomBiotech


*In your operations, behind your decisions,
powering your products*

Rely on a smart partner

As a long-standing industry leader, we have confronted many of the challenges you face as a diagnostics manufacturer. That experience and knowledge allows us to anticipate and overcome roadblocks. For you, this means access to unique insights from specialists who know the industry landscape and support from a team of experts who are passionate about solutions.

Secure your operations

At CustomBiotech, we are committed to support the long-term growth and evolution of your manufacturing — from small to large scale. Operating under stringent quality standards, we offer raw materials and assay components with proven lot-to-lot consistency in narrow specifications so that your products meet regulations. Our solutions enable streamlined workflows and traceable quality to strengthen your operational performance.



Our experts know and understand the industry landscape



World-class infrastructure



Visit us online at custombiotech.roche.com

- Explore all CustomBiotech products and services — also for molecular diagnostics and biopharma products and services
- Connect with your local Roche CustomBiotech representative
- Download product literature
- Obtain comprehensive information on products and applications



The Roche CustomBiotech Team

Your gateway to the worldwide Roche network



Center of Excellence - CustomBiotech headquarters in Penzberg

Our facilities in Penzberg are one of the largest and most productive biotechnology centers in Europe. With an unusual setup merging diagnostics and pharmaceutical know-how under one roof, Penzberg is a multidisciplinary incubator of novel technologies for production and automation processes, and a state-of-the-art manufacturer of high-quality raw materials and biotechnology products.





Your strategic partner in diagnostics manufacturing
Clinical Chemistry assays and Immunoassays powered by Roche



Benefit from the pioneer's expertise in Diagnostics

For decades, Roche has been synonymous with excellent performance in the area of *in vitro* diagnostics. At the core of our market-leading products is an extensive portfolio of carefully designed and expertly manufactured assay components. These components — from enzymes, substrates and cofactors, to antibodies, antigens, streptavidin, biotin products, and immunoassay interference blockers — are available through CustomBiotech to bring your diagnostic product to market.

However, delivering raw materials of excellent quality is not enough. At Roche CustomBiotech, we strive to be your single strategic partner in enabling better diagnostics with leading Clinical Chemistry and Immunology solutions.

This catalog is designed to provide easy access to details about CustomBiotech raw materials for Clinical Chemistry and Immunology assays. For the most up-to-date product information, please visit custombiotech.roche.com, where you will find the entire CustomBiotech product portfolio.



Clinical Chemistry

Biochemicals	10	Cofactors/Nucleotides	30
Buffers	10	Cofactors	30
Detergents	17	Nucleotides	40
<i>Ionic Detergents</i>	17	Enzymes	47
<i>Non-ionic Detergents</i>	20	Substrates	148
<i>Zwitterionic Detergents</i>	22	Colorimetric Substrates	148
Protease Inhibitors	24	Non-Colorimetric Substrates	156
Additional Biochemicals	26		

Diethanolamine 85%

solution

Buffer for enzymatic assays of alkaline phosphatase.

Application

Use Diethanolamine 85% as a buffer in applications that test or use the activity of alkaline phosphatase, such as diagnostic tests for alkaline phosphatase.

CAS: 111-42-2

Properties

Formula: C₄H₁₁O₂N

Molecular weight: 105.14 D

Solubility: Miscible with water and ethanol

Suggested pH range: Approximately 10

Specification

Appearance: Colorless, clear liquid

Solubility: Miscible with water, ethanol

Refractive index: n 20/D: 1.4575-1.4595

Density: D 20/20: 1.092-1.095

Coloration of sample (against water):

A₄₀₅: ≤0.043

A₄₀₅ (10 days at +4°C): ≤0.043

A₄₀₅ (10 days at +35°C): ≤0.051

Diethanolamine (HClO₄ titration): 84.0-86.0%

Thin layer chromatography: Corresponds to reference

Mono-ethanolamine (TLC): ≤0.4%

Fe (calculated on diethanolamine 100%): ≤2 ppm

Heavy metals (as Pb, calculated on Diethanolamine 100%): ≤5 ppm

Water (K. Fischer): 14.0-16.0%

Stability: At +2 to +8°C within specification range 18 months.

Store dry. Protect from light. Keep in tightly sealed containers.

Catalog number

10 201 294 103

Pack size

custom fill

Will be supplied as "Diethanolamine 85%". Unit of measure is "kg".

For further processing only.

Glycylglycine

crystalline powder

Zwitterionic buffer for diagnostic tests or as a substrate for γ-glutamyltransferase tests.

Catalog number

10 002 887 103

Pack size

custom fill

Will be supplied as "Glycylglycine". Unit of measure is "kg".

For further processing only.

Application

Use Glycylglycine as a buffer in diagnostic reagents, or as a substrate in γ -glutamyltransferase tests where γ -glutamyltransferase transfers the gamma-glutamyl group of L- γ -glutamyl-3-carboxy-4-nitroanilide (Glupa-C) to glycylglycine.

CAS: 556-50-3

Properties

Formula: $C_4H_8N_2O_3$

Molecular weight: 132.1 D

Solubility: Easily soluble in water

Suggested pH range: 7.5-8.9

Specification

Appearance: White crystalline powder

Solubility: Easily soluble in water

pH value (c=0.35%, w/v): 5.5-6.5

Glycylglycine (HClO₄ titration, based on anhydrous substance): 99.0-100.5%

TLC: Corresponds to reference

A₄₀₅ (c=10%; w/v, against water): $\leq 0.01\%$

Heavy metals (as Pb): ≤ 5 ppm

Water (K. Fischer): $\leq 0.5\%$

Fe: ≤ 10 ppm

Sulfate ash: $\leq 0.1\%$

Glycine (TLC): $\leq 0.2\%$

Contaminating amino acids (TLC): Not detectable

Microbiological test: Corresponds to specification

IR Spectrum: Corresponds to reference

Stability: At +15 to +40°C within specification range for 36 months.

Store dry in tightly sealed containers.

Hepes

crystalline powder

Buffer for diagnostic tests, such as amylase test.

Application

Use Hepes as a buffer in reagents.

CAS: 7365-45-9

Catalog number

10 172 944 103

Pack size

custom fill

Will be supplied as "Hepes". Unit of measure is "kg".

For further processing only.

Properties**Nomenclature:** 4-(2-Hydroxyethyl)-1-piperazineethanesulfonic acid**Formula:** C₈H₁₈N₂O₄S**Molecular weight:** 238.3 D**Suggested pH range:** 6.8-8.2**Specification****Appearance:** White crystalline powder**Solubility:** Clear, colorless solution in water (c=1 mol/L)**pK value:** 7.21-7.41 (+37°C); 7.45-7.65 (+20°C)**Melting range:** +207 to +213°C**Hepes** (alcalimetric): ≥97%**Hepes** (from N): ≥97%**N** (elementary analysis): ≥11.4%**Thin layer chromatography** (TLC): Chromatographically homogeneous**A**₂₆₀ (against water): ≤0.050**A**₄₀₅ (against water): ≤0.030**Cl** (chloride meter): ≤0.04%**Hepes mA/min** (Purity check, α-amylase contamination): ≤0.1**Exclusion of skin contact and contamination with saliva:**

Corresponds to specification.

Stability: At +15 to +25°C within specification range for 24 months.**Imidazole**

crystalline powder

Buffer for diagnostic tests, such as creatine kinase test.

Application

Use Imidazole as a buffer for diagnostic tests and other reagents, especially enzymatic reactions such as creatine kinase test.

CAS: 288-32-4**Properties****Formula:** C₃H₄N₂**Molecular weight:** 68.08 D**Suggested pH range:** 6.2-7.8**Specification****Appearance:** White or slightly yellowish crystallizate**Melting range:** +86 to +91°C**Imidazol** (titrimetric): ≥99.0%**dA**₂₅₀ - **dA**₃₆₀: ≤0.050**dA**₃₃₄: ≤0.050**Catalog number****10 034 428 103****Pack size**

custom fill

Will be supplied as "Imidazole". Unit of measure is "kg".

For further processing only.

dA_{405} : ≤ 0.010

Stability: At +15 to +25°C within specification range for 24 months.

Pipes

free acid

Buffer for diagnostic tests, such as cholesterol test.

Application

Use Pipes as a buffer in a variety of diagnostic tests, especially in tests for cholesterol and triglycerides.

CAS: 5625-37-6

Properties

Nomenclature: Piperazine-1,4-bis-2-ethane sulfonic acid

Formula: $C_8H_{18}N_2O_6S_2$

Molecular weight: 302.4 D

Suggested pH range: 6.1-7.5

Specification

Appearance: Colorless crystals

Solubility: Clear, colorless solution in water (c=10 mg/mL)

pK value: 6.8 ± 0.1

Pipes (alkalimetric): $\geq 98\%$

Pipes (from N): $\geq 98\%$

A_{260} (c=10 mg/ml water) : ≤ 0.05

Thin layer chromatography (TLC): Chromatographically homogeneous

Br (ion chromatography): $\leq 0.5\%$

Identity (IR spectrum): Corresponds to reference

Stability: At +15 to +25°C within specification range for 24 months.

Catalog number

10 239 500 103

Pack size

custom fill

Will be supplied as "Pipes, Free Acid". Unit of measure is "kg".

For further processing only.

Pipes

disodium salt

Buffer for diagnostic tests, such as cholesterol test.

Application

Use Pipes as a buffer in a variety of diagnostic tests, especially in tests for cholesterol and triglycerides.

CAS: 5625-37-6

Catalog number

10 735 361 103

Pack size

custom fill

Will be supplied as "Pipes, Sodium Salt". Unit of measure is "kg".

For further processing only.

Properties**Nomenclature:** Piperazine-1,4-bis-2-ethane sulfonic acid**Formula:** C₈H₁₆N₂O₆S₂Na₂**Molecular weight:** 346.3 D**Suggested pH range:** 6.1-7.5**Specification****Appearance:** White powder**Solubility:** Clear, colorless solution in water (c=10 mg/mL)**Pipes** (from N): ≥82%**Na** (flame photometric): 11-14%**Water** (K. Fischer): ≤5%**N** (elementary analysis) : ≥7.6%**A₃₄₀** (c=10 mg/ml water): ≤0.010**HPTLC:** Chromatographically homogeneous**Heavy metals:** ≤5 ppm**Identity** (IR spectrum): Corresponds to reference**Stability:** At +15 to +25°C within specification range for 24 months.**Tris****crystallizate**

Buffer for diagnostic tests, such as tests for aminotransferases.

Application

Use Tris as a buffer in diagnostic reagents, especially in tests for aminotransferases or γ-glutamyltransferase.

CAS: 77-86-1**Properties****Nomenclature:** Tris(hydroxymethyl)-aminomethane**Formula:** C₄H₁₁NO₃**Molecular weight:** 121.1 D**Suggested pH range:** 7.0-9.0**Specification****Appearance:** Colorless, odorless crystallizate**Solubility:** Clear, colorless solution in water (c=100 mg/mL), free of fuz**Flow properties:** Passes**Melting range:** +168 to +171°C**Conductivity** (water, 1 μS, +25°C): ≤110 μS**pH value** (c=6 mg/ml, in water): 10.0-11.0**Tris** (titrimetric, based on dry weight): 99.5-100.5%**Water** (K. Fischer): ≤0.2%**Sulfate ash** (with concentrated H₂SO₄ at +600°C): ≤0.05%**Catalog number****10 153 265 103****Pack size**

custom fill

Will be supplied as "Tris-(hydroxymethyl)-aminomethane". Unit of measure is "kg".

For further processing only.

Fe (AAS): ≤ 1 ppm

As (AAS): ≤ 1 ppm

Heavy metals (as Pb): ≤ 1 ppm

Reducing substances (KMnO_4 , 0.002 mol/L): ≤ 3 mL/100 mg

Acetone (GC): $\leq 0.05\%$

Methanol (GC): $\leq 0.05\%$

Cl (turbidimetric test with AgNO_3): ≤ 20 ppm

Bioburden: ≤ 100 CFU/g

A₃₀₀ (against water, c=100 mg/mL): ≤ 0.020

A₄₀₅ (against water, c=100 mg/mL): ≤ 0.004

Stability: At +15 to +25°C within specification range for 24 months.

Protect from light.

MES, free acid

crystalline powder

Buffer for a variety of diagnostic tests, such as glucose test.

Application

Use MES as a buffer in reagents that require a pH of approximately 6.

Properties

Nomenclature: 4-Morpholineethane sulphonic acid

Formula: $\text{C}_6\text{H}_{13}\text{NO}_4\text{S}$

Molecular weight: 195.2 D

Suggested pH range: 5.5-6.7

Specification

Appearance: White, crystalline powder

Solubility: Clear, colorless solution in water (20 mg/mL)

Identity (NIR): Corresponds to reference

PVS/Igl screening: ≤ 1 ppm

MES (alkalimetric, based on dry weight): $\geq 98.0\%$

Water (K. Fischer): $\leq 10.0\%$

Na (flame photometric): ≤ 20 ppm

Stability: At +15 to +25°C within specification range for 24 months.

Catalog number

04 808 177 103

Pack size

custom fill

Will be supplied as "4-Morpholineethane Sulfonic Acid (Mes)". Unit of measure is "kg".

For further processing only.

Detergents

Application

Detergents are used in diagnostic kits to:

- enhance the solubility of test or sample compounds
- activate enzymes such as esterases
- reduce interference from serum lipids
- reduce carryover effects on analyzers and facilitate dispensing processes
- reduce non-specific binding to solid phases in immunoassays
- pretreat samples
- make use of possible antimicrobial effects

Properties

The choice of a specific detergent depends on the protein to be solubilized, the need of removal, toxicity data, interference with UV-VIS-absorption, interference with subsequent isoelectric focusing or ultracentrifugation.

	Ability to disperse protein aggregates	Denaturation of protein	Ease of removal
Non-ionic detergents			
n-Octyl- β -D-glycoside	Low	No	Very easy
Polidocanol (Thesit)	Low	No	Difficult
n-Dodecyl- β -D-maltoside	Low	No	Difficult
Tween 20	Low	No	Difficult
Ionic detergents			
Cholate	Low	(no)	Very easy
Deoxycholate	High	(no)	easy
Zwitterionic detergents			
CHAPS	High	No	Very easy

Dilaurylglycosulfate

powder

Detergent for diagnostic tests.

Application

Use Dilaurylglycosulfate as a co-emulsifier in the diagnostic test for the determination of lipase.

CAS: 99387-94-7

Properties

Formula: C₂₇H₅₆O₆S

Molecular weight: 508.8 D

Specification

Appearance: White powder

Dilaurylglycosulfate (from C): ≥90%

C (elementary analysis): ≥57.3%

H (elementary analysis): ≥10.0%

Water (K. Fischer): ≤5%

Stability: At +2 to +8°C within specification range for 12 months.

Catalog number

11 827 294 103

Pack size

custom fill

Will be supplied as "Dilaurylglycosulfat". Unit of measure is "g".

For further processing only.

Cholate

ionic detergent, sodium salt

Anionic detergent for diagnostic tests.

Application

Use Cholate in diagnostic reagents, such as for the determination of cholesterol and triglycerides.

CAS: 81-25-4

Properties

Formula: C₂₄H₃₉NaO₅

Molecular weight: 430.6 D

Detergent type: Anionic detergent

Solubility: Limited solubility in the presence of Ca²⁺

Handling advice: Harmful if exposed to skin and if inhaled. Adequate precautions as for handling of irrigating products must be taken.

Specification

Appearance: White crystalline powder

Solubility:

Clear, colorless solution in water (c=10 mg/mL)

Catalog number

10 261 084 103

Pack size

custom fill

Will be supplied as "Cholic Acid Sodium Salt". Unit of measure is "kg".

For further processing only.

Clear, colorless to yellowish solution in water (c=150 mg/mL, +20°C)

Identity (NIR): Corresponds to reference

Cholic acid, Na-salt (HPLC): ≥93 area%

Water (K. Fischer): ≤6%

C (elementary analysis): 62-67%

Na (flame photometric): 5.0-5.5%

Heavy metals (as Pb): ≤10 ppm

Flame coloration: Positive

A₃₄₀ (against water): ≤0.100

A₅₀₅ (against water): ≤0.005

A₅₄₆ (against water): ≤0.005

A₅₀₅ to A₅₅₀ (against water): ≤0.025

Hydrophilic contaminants (HPLC): ≤15 area%

Lipophilic contaminants (HPLC): ≤4.0 area%

Reducing substances: ≤0.25 ml (KMnO₄, 0.002 mol/L, per 100 mg)

Oxidizing substances: Negative

Bioburden: ≤100 CFU/g, ≤10 moulds/g

Performance: Corresponds to specification

Stability: At +15 to +25°C within specification range for 36 months.

Protect from light.

Deoxycholate

ionic detergent, sodium salt

Anionic detergent for diagnostic tests.

Application

Use Deoxycholate in diagnostic reagents, such as for the determination of lipase.

CAS: 83-44-3

Properties

Formula: C₂₄H₃₉O₄Na

Molecular weight: 414.6 D

Detergent type: Anionic detergent

Handling advice: Harmful if exposed to skin and if inhaled. Adequate precautions as for handling of irrigating products must be taken.

Specification

Appearance: White, crystalline powder

Solubility: Clear, colorless solution in water (c=50 mg/mL), free of fuzz

Deoxy cholic acid (sodium deoxycholate, HPLC): ≥83.0 area%

Hydrophilic contaminants: ≤15.0 area%

Lipophilic contaminants: ≤5.0 area%

Water (K. Fischer): ≤10.0%

Catalog number

11 434 314 103

Pack size

custom fill

Will be supplied as "Desoxycholat, Mono-NA, Crystal". Unit of measure is "kg".

For further processing only.

Acetone (GC): ≤0.5%

Heavy metals (as Pb): ≤5 ppm

Stability: At +15 to +25°C within specification range for 36 months.

Taurodesoxycholat

sodium salt

Anionic detergent for diagnostic tests.

Application

Use Taurodesoxycholat in diagnostic reagents, such as for the determination of lipase.

CAS: 1180-95-6

Properties

Formula: C₂₆H₄₄NO₆SNa

Molecular weight: 521.7 D

Detergent type: Anionic detergent

Handling advice: Harmful if exposed to skin and if inhaled. Adequate precautions for handling hazardous products must be used.

Specification

Appearance: White lyophilizate

Taurodesoxy cholate, Na (from C): ≥90%

Taurodesoxy cholate, Na (HPLC): ≥89 area%

C (elementary analysis): 53.9-61.0%

H (elementary analysis): 7.9-8.9%

N (elementary analysis): 2.4-3.0%

Na (flame photometric): 4.4-6.6%

Water (K. Fischer): ≤5%

Stability: At +2 to +8°C within specification range for 24 months.

Catalog number	Pack size
11 332 686 103	custom fill

Will be supplied as "Taurodesoxycholic acid, Na, pur.". Unit of measure is "g".

For further processing only.

n-Dodecyl-β-D-maltoside

nonionic detergent, powder

Nonionic detergent for diagnostic tests

Application

Use n-Dodecyl-β-D-maltoside as a mild, nondenaturing detergent for the solubilization of proteins, especially antibodies.

CAS: 69227-93-6**Properties****Nomenclature:** 1-O-n-Dodecyl-β-D-glucopyranosyl(1-4)α-D-glucopyranoside**Formula:** C₂₄H₄₆O₁₁**Molecular weight:** 510.62 D**Detergent type:** Nonionic alkyl maltoside type**Specification****Appearance:** White, crystalline powder**Specific rotation [α]_D²⁵** (in MeOH): +46.0±2.0°**n-Dodecylmaltoside** (from C): ≥98%**C** (elementary analysis): ≥55.20%**Dodecanol** (GC): ≤0.10%**Stability:** At +15 to +25°C within specification range for 24 months.**Catalog number****10 808 342 103****Pack size**

custom fill

Will be supplied as "n-Dodecyl-β-D-maltoside". Unit of measure is "g".

For further processing only.

n-Octyl β-D-glucoside

nonionic detergent, powder

Nonionic detergent for diagnostic tests

Application

Use n-Octyl β-D-glucoside as a mild, nondenaturing detergent for the solubilization of proteins, especially antibodies. n-Octyl β-D-glucoside can be easily removed by dialysis.

CAS: 29836-26-8**Properties****Nomenclature:** 1-O-Octyl-β-D-glucopyranoside**Formula:** C₁₄H₂₈O₆**Molecular weight:** 292.4 D**Detergent type:** Nonionic alkyl glucoside type**pH stability:** Stable in solutions above pH 6.5**Catalog number****10 411 469 103****Pack size**

custom fill

Will be supplied as "n-Octylglucoside". Unit of measure is "g".

For further processing only.

Specification**Appearance:** White powder**n-Octylglucoside** (from C): ≥99%**C** (elementary analysis): ≥56.9%**Octanol** (GC): ≤0.1%**Stability:** At +15 to +25°C within specification range for 36 months.

Store dry.

Polidocanol (Thesit)

Nonionic detergent for use in diagnostic reagents.

Application

Use Polidocanol in diagnostic reagents to enhance solubility, minimize interferences and reduce carryover effects.

CAS: 9002-92-0**Properties****Nomenclature:** Dodecylpolyethyleneglycolether**Formula:** C₃₀H₆₂O₁₀ (n≈approximately 9)**Molecular weight:** Approximately 600 D**Detergent type:** Nonionic polyoxyethylene type**Handling advice:** Polidocanol must be moderately heated (+40 to +50°C) and carefully homogenized by gentle stirring before dispensing.**Specification****Appearance:** White, pasty, fatty substance; clear, colorless to slightly yellow liquid at approximately +30°C**Solubility:** Clear, colorless solution in water (c=100 mg/mL)**Peroxide** (as H₂O₂): ≤1 ppm**Stability:** At +2 to +8°C within specification range for 24 months. Keep under argon or nitrogen. Protect from light.**Catalog number****10 831 620 103****Pack size**

custom fill

Will be supplied as "Polidocanol (PEG Monododecyl Ether)". Unit of measure is "kg".

For further processing only.

Tween 20

purified, solution

Non-ionic detergent for diagnostic tests.

Application

Use Tween 20 as non-denaturing detergent for the solubilization of proteins, especially antibodies.

CAS: 9005-64-5**Catalog number****11 334 000 103****Pack size**

custom fill

Will be supplied as "Tween 20, gereinigt, Lsg". Unit of measure is "L".

For further processing only.

Properties

Formula: $C_{58}H_{114}O_{26}$ (for $w+x+y+z=n=20$)

Molecular weight: 1228 g/mol

Specification

Appearance: Clear, yellow solution

Conductivity: $\leq 100 \mu\text{S/cm}$

Peroxide (as H_2O_2): $\leq 2 \text{ ppm}$

Aldehyde: $\leq 0.02 \text{ mg/ml}$

Stability: At +2 to +8°C within specification range for 24 months. Store under nitrogen. Protect from light.

CHAPS

zwitterionic detergent, crystalline powder

Zwitterionic detergent for diagnostic tests.

Application

Use CHAPS to reduce protein-protein interactions. CHAPS can be easily removed by dialysis.

CAS: 75621-03-3

Properties

Nomenclature: 3-[[3-Cholamidopropyl]dimethylammonio]-1-propanesulfonate

Formula: $C_{32}H_{58}N_2O_7S$

Molecular weight: 614.9 D

Detergent type: Zwitterionic detergent, nondenaturing

Specification

Appearance: White, crystalline powder

CHAPS (from N): $\geq 99\%$

N (elementary analysis): $\geq 4.49\%$

A₂₆₀ (against water): ≤ 0.10

A₂₈₀ (against water): ≤ 0.10

Thin layer chromatography (TLC): Chromatographically homogeneous, corresponds to reference

Stability: At +2 to +8°C within specification range for 24 months. Store dry. Protect from light.

Catalog number

10 810 681 103

Pack size

custom fill

Will be supplied as "CHAPS". Unit of measure is "kg".

For further processing only.

CHAPSO

zwitterionic detergent, crystalline powder

Zwitterionic detergent for diagnostic tests.

Application

Use CHAPSO to reduce protein-protein interactions. CHAPSO can be easily removed by dialysis.

CAS: 82473-24-3

Properties

Nomenclature: 3-[(3-Cholamidopropyl)dimethylammonio]-2-hydroxy-1-propansulfonate

Formula: $C_{32}H_{58}N_2O_8S$

Molecular weight: 630.9 D

Detergent type: Zwitterionic detergent, similar to CHAPS but more soluble.

Handling advice: Harmful if exposed to skin and if inhaled. Adequate precautions as for handling of irrigating products must be taken.

Specification

Appearance: White crystalline powder

CHAPSO (from C): $\geq 95\%$

C (elementary analysis): $\geq 57.6\%$

Water (K. Fischer): $\leq 4\%$

HPTLC: Chromatographically homogeneous

Stability: At +2 to +8°C within specification range for 36 months.

Catalog number

11 112 392 103

Pack size

custom fill

Will be supplied as "CHAPSO". Unit of measure is "g".

For further processing only.

cOplete, EDTA free

7500 tablets in glass vial

Mix of protease inhibitors

Application

cOplete, EDTA-free is used for the inhibition of serine and cysteine proteases in bacterial, yeast, plant, and animal cell extracts. cOplete, EDTA-free Tablets, are used for the inhibition of proteolytic activity in large volumes (up to 50 mL) in which EDTA may interfere with protein stability (e.g., metal-containing proteins) or subsequent assays.

Specification**Appearance:** White tablets**Resolving time:** ≤180 seconds**Abrasion:** ≤1.0%**Diameter of tablets** (average): 6 mm**Height of tablets** (average): 1.75-2.00 mm**Weight of tablets** (average): 56.0-66.0 mg**Homogeneity of lot** (based on weight of tablets): S<3.0**Inhibition of pancrease, chymotrypsin and thermolysin** (after 60 minutes at +20 to +25°C):

pancrease: ≥90%

chymotrypsin: ≥90%

thermolysin: ≤15%

Stability: At +2 to +8°C within specification range for 24 months.**Catalog number****04 574 834 001****Pack size**

7500 tablets in glass vial

Will be supplied as "Complete". Unit of measure is "piece".

Aprotinin

from bovine lung, lyophilizate

Protease inhibitor

Application

Use Aprotinin in reagents to inhibit serine proteases, such as kallikrein, plasmin, trypsin and chymotrypsin.

CAS: 9087-70-1**Specification****Appearance:** White lyophilizate**Activity** (Chromozym TRY, +25°C): ≥630 inhibitor U/mg lyophilizate**Activity** (BAEE, +25°C): ≥200 inhibitor U/mg lyophilizate**Protein** (Lowry): 90-100%**Electrophoresis** (SDS Page): Corresponds to reference**Country of origin:** USA, South Africa, New Zealand, Australia, or Uruguay, respective**Catalog number****10 236 632 103****Pack size**

custom fill

Will be supplied as "Aprotinin from Bovine Lung". Unit of measure is "g".

For further processing only.

Stability: At +2 to +8°C within specification range for 12 months.

Store dry.

Remarks:

Official veterinary certificate of health of the donor animals is available.

Official certificate of the deactivation of animal material including the method (acid treatment at up to pH 5 for up to 5 h) is available.

Pefabloc SC (AEBSF)

4-(2-Aminoethyl)-benzenesulfonyl fluoride hydrochloride, powder

Protease inhibitor

Application

Use Pefabloc in reagents to inhibit serine proteases, such as thrombin in serum or plasma.

CAS: 34284-75-8

Properties

Molecular weight: 239.5 D

Specification

Appearance: White powder

Pefabloc HCl (HPLC): ≥90 area%

Pefabloc HCl (from C): ≥95% (theory 100%)

C (elementary analysis): ≥38.1% (theory 40.08%)

H (elementary analysis): ≥4.2% (theory 4.59%)

N (elementary analysis): ≥5.5% (theory 5.84%)

Thin layer chromatography (TLC): Chromatographically homogeneous

Inhibition chymotrypsin: Corresponds to specification

Stability: At +2 to +8°C within specification range for 24 months.

Catalog number

11 427 393 103

Pack size

custom fill

Will be supplied as "Pefabloc SC AEBSF, Hydrochloride". Unit of measure is "g".



For further processing only.

1,4-Dithiothreitol (DTT)

crystallizate

Reducing agent

Application

Use 1,4-Dithiothreitol primarily to protect free SH-groups from oxidation. Use it routinely in all work with enzymes and proteins during enzyme measurement and the characterization of proteins.

CAS: 3483-12-03**Properties****Nomenclature:** Threo-1,4-dimercapto-2,3-butanediol**Formula:** C₄H₁₀O₂S₂**Molecular weight:** 154.3 D**Specification****Appearance:** White to yellowish crystallizate**DTT** (with Ellman's reagent): ≥97%**Thin layer chromatography:** Chromatographically homogeneous, corresponds to reference**Stability:** At +2 to +8°C within specification range for 24 months.

Store dry. Protect from light.

Catalog number**10 197 785 103****Pack size**

custom fill

Will be supplied as "Dithiothreitol (DTT) Cleland's Reagent". Unit of measure is "g".

For further processing only.

3-Hydroxy-1,2,3,4-tetrahydrobenzo[h]quinoline

crystalline powder

Chemical for dry chemistry diagnostic tests

Application

Use 3-Hydroxy-1,2,3,4-tetrahydrobenzo[h]quinoline in dry chemistry application for the determination of urea.

CAS: 5423-67-6**Properties****Nomenclature:** 3-Hydroxy-1,2,3,4-tetrahydrobenzo[h]quinoline**Formula:** C₁₃H₁₃NO**Molecular weight:** 199.25**Toxicity:** Harmful**Catalog number****10 003 174 103****Pack size**

custom fill

Will be supplied as "BM 32.027". Unit of measure is "kg".

For further processing only.

Specification**Appearance:** White to greyish powder**Solubility:** Clear, colorless solution in methanol (c=0.2%; w/v)**Melting range:** +150 to +154°C**Loss on drying:** ≤0.5%**Sulphate ash:** ≤0.2%**3-Hydroxy-1,2,3,4-tetrahydrobenzo[h]quinoline** (HClO₄ titration, based on dried substance): 98.0-102.0%**UV/VIS spektrum:**Maximum I: 250 to 254 nm (specific absorbance (A_{1%/1cm}): 988-1050)Maximum II: 334 to 338 nm (specific absorbance (A_{1%/1cm}): 286-306)**Thin layer chromatography (TLC):** Corresponds to reference**IR Spectrum:** Corresponds to reference**Stability:** At +2 to +8°C within specification range for 24 months.

Protect from light.

D-Mannitol**reduced sodium**

Excipient (inactive substance) for the production of tablets or granulated material

Application

Use D-Mannitol as an excipient (inactive substance) to produce tablets or granulated material that contain reagents or components used in diagnostic applications.

CAS: 69-65-8**Properties****Formula:** C₆H₁₄O₆**Molecular weight:** 182.2 D**Specification****Appearance:** White, silky crystals or white, crystalline powder**Solubility:** Clear, colorless solution in water (c=10%, w/v), slightly soluble in water and hot ethanol, heavy soluble in ethanol**Melting range:** +166 to +168°C**Specific rotation** (c=10%, w/v, calculated on dry substance): +23.0° to +25.0°**Purity (HPLC):** D-Mannit: ≥97.5 area%**Alkaline impurities** (calculated as NaOH): ≤80 ppm**Acid impurities** (calculated as HCl): ≤45 ppm**Cl:** ≤50 ppm**Fe:** ≤0.1 ppm**Ni:** ≤1 ppm**Catalog number****11 371 754 103****Pack size**

custom fill

Will be supplied as "D-Mannit, Na-arm". Unit of measure is "kg".

For further processing only.

Reducing sugar: Corresponds to reference

Heavy metals (as Pb): ≤2 ppm

Sulphate: ≤100 ppm

Sulphate ash: ≤500 ppm

Loss on drying: ≤0.3%

Na: ≤200 ppm

Microbiological analysis: Corresponds to reference

Stability: At +15 to +25°C within specification range for 36 months.

Store dry. Protect from light.

Kryptofix 221

solution

Cryptant that binds cations in aqueous solutions

Application

Use Kryptofix 221 in enzymatic tests for potassium to decrease the sodium concentration relative to potassium.

CAS: 31364-42-8

Properties

Formula: $C_{16}H_{32}N_2O_5$

Molecular weight: 332.44 D

Specification

Appearance: Clear, yellow liquid

Solubility: Clear yellowish solution in water (c=4 mg/mL)

A₃₄₀ (aqueous solution): ≤0.100

A₄₀₅ (aqueous solution): ≤0.060

Identity (NIR): Corresponds to reference

¹³C-NMR spectrum: Corresponds to masterlot

Kryptofix 221 agent (HPLC; based on masterlot): ≥75%

Na (AES): ≤65 ppm

K (AES): ≤5 ppm

NH₃ (evolution in buffer; after 10 days at +55°C): ≤50 µmol/L

Stability: At +2 to +8°C within specification range for 12 months. Keep under nitrogen. Protect from light.

Catalog number

11 183 958 103

Pack size

custom fill

Will be supplied as "Kryptofix 221". Unit of measure is "kg active ingredient".

For further processing only.

Valinomycin

crystallizate

Potassium selective ionophoric cyclodepsipeptide

Application

Use Valinomycin in diagnostic tests for potassium where it acts as an ion carrier in potassium selective electrodes.

CAS: 2001-95-8

Properties

Formula: C₅₄H₉₀N₆O₁₈

Molecular weight: 1111.4 D

Specification

Appearance: White crystallizate

Solubility: Clear, colorless solution in chloroform (c=10 mg/mL)

Melting point: ≥+183°C

Specific rotation (in chloroform):: +30.0±2.0°

Valinomycin (from N): ≥94%

Valinomycin (HPLC): ≥85.0 area%

C (elementary analysis) : ≥54.86%

H (elementary analysis) : ≥7.67%

N (elementary analysis): ≥7.10%

Thin layer chromatography (HPTLC):

a) UV: Homogeneous

b) to spray with H₂SO₄ (1%); to vaporize with Iodine: Corresponds to reference

Stability: At +2 to +8°C within specification range for 36 months. Store in safety zone dedicated to poisonous agents.

Catalog number

10 161 594 103

Pack size

custom fill

Will be supplied as "Valinomycin". Unit of measure is "g".

For further processing only.

Acetyl-Coenzyme A

trilithium salt

Cofactor for carnitine acetyl transferase

Application

Use Acetyl-Coenzyme A for the determination of L-carnitine.

CAS: 72-89-9

Properties

Formula: $C_{23}H_{35}N_7O_{17}P_3SLi_3$

Molecular weight: 827.4 D (Acetyl-CoA: 809.6 D)

Specification

Appearance: White lyophilizate

Solubility: Clear, colorless solution in water (c=10 mg/mL)

Acetyl-CoA (enzymatic): ≥83%

Acetyl-CoA (A_{260} , $\epsilon=16.0$ [L x mmol⁻¹ x cm⁻¹]): ≥85%

Li (flame photometric): 2±0.3%

Stability: At -15 to -25°C within specification range for 12 months.

Store dry.

Catalog number

10 150 932 103

Pack size

custom fill

Will be supplied as "Acetyl Coenzyme A (Acetyl-CoA) Tri-Li". Unit of measure is "g".



For further processing only.

Coenzyme A, Grade I

free acid

Coenzyme A is a cofactor for some enzymes, e.g., citrate lyase.

Application

Use Coenzyme A in diagnostic tests measuring citrate or citrate lyase. It is also used as an enhancer of luciferase light emission or as a precursor for Acetyl-Coenzyme A or in other chemical or enzymatical reactions.

CAS: 85-61-0

Properties

Formula: $C_{21}H_{36}N_7O_{16}P_3S$

Molecular weight: 767.6 D (CoA: 767.6 D)

Specification

Appearance: White to slightly yellow lyophilizate

CoA, reduced (enzymatic, 10 U phosphotransacetylase): ≥85%

CoA (A_{260} , $\epsilon=16.0$ [L x mmol⁻¹ x cm⁻¹]): ≥88%

Water (K. Fischer): ≤6%

Glutathione, reduced (enzymatic): ≤1%

Catalog number

10 151 009 103

Pack size

custom fill

Will be supplied as "Coenzyme A (CoA), Free Acid, Grade I". Unit of measure is "g".



For further processing only.

Stability: At -15 to -25°C within specification range for 12 months.
Store dry.

Coenzyme A, Grade I

trilithium salt

Coenzyme A is a cofactor for some enzymes, *e.g.*, citrate lyase.

Application

Use Coenzyme A in diagnostic tests measuring citrate or citrate lyase. It is also used as an enhancer of luciferase light emission or as a precursor for Acetyl-Coenzyme A or in other chemical or enzymatical reactions.

CAS: 85-61-0

Properties

Formula: C₂₁H₃₃N₇O₁₆P₃SLi₃

Molecular weight: 785.4 D (CoA: 767.6 D)

Specification

Appearance: White to slightly yellow lyophilizate

CoA, reduced (enzymatic, with 10 U phosphotransacetylase): ≥83%

CoA (A₂₆₀, ε=16.0 [L x mmol⁻¹ x cm⁻¹]): ≥84%

Water (K. Fischer): ≤6%

Glutathione, reduced (enzymatic): ≤1%

Stability: At -15 to -25°C within specification range for 12 months.

Catalog number

10 121 541 103

Pack size

custom fill

Will be supplied as "Coenzyme A (CoA), Tri-Li Salt, Grade I". Unit of measure is "g".



For further processing only.

Coenzyme A, Grade II

trilithium salt

Coenzyme A is a cofactor for some enzymes, *e.g.*, citrate lyase.

Application

Use Coenzyme A in diagnostic tests measuring citrate or citrate lyase. It is also used as an enhancer of luciferase light emission or as a precursor for Acetyl-Coenzyme A or in other chemical or enzymatical reactions.

CAS: 85-61-0

Properties

Formula: C₂₁H₃₃N₇O₁₆P₃SLi₃

Molecular weight: 785.4 D (CoA: 767.6 D)

Catalog number

10 155 969 103

Pack size

custom fill

Will be supplied as "Coenzyme A (CoA), Ttri-Li, Grade II". Unit of measure is "g".



For further processing only.

Specification

Appearance: White to slightly yellow lyophilizate

CoA, reduced (enzymatic, with 10 U phosphotransacetylase): ≥73%

CoA (A_{260} , $\epsilon = 16.0$ [L x mmol⁻¹ x cm⁻¹]): ≥81%

Water (K. Fischer): ≤8%

Stability: At -15 to -25°C within specification range for 12 months.

Store dry.

FAD

disodium salt

Cofactor for dehydrogenases and oxidases.

Application

Use FAD as a cofactor in a variety of enzymatic test or assays to activate enzymes, especially diagnostic tests for triglycerides.

CAS: 146-14-5

Properties

Nomenclature: Flavine-adenine dinucleotide

Formula: C₂₇H₃₁N₉O₁₅P₂Na₂

Molecular weight: 829.6 D (FAD: 785.7 D)

Specification

Appearance: Yellow powder

FAD (A_{450} , $\epsilon = 11.3$ [L x mmol⁻¹ x cm⁻¹]): ≥86%

Na (flame photometric): 5±1%

Water (K. Fischer): ≤9%

P_i: ≤0.6%

Stability: At +2 to +8°C within specification range for 24 months.

Protect from light.

Catalog number

10 154 032 103

Pack size

custom fill

Will be supplied as "Flavine-adenine Dinucleotide (FAD), Di-Na".
Unit of measure is "g".

For further processing only.

NAD, Grade I

free acid

Cofactor for dehydrogenases, *e.g.*, lactate dehydrogenase.

Application

Use NAD, Grade I as a cofactor in a variety of diagnostic tests, such as for the determination of ethanol and lactate dehydrogenase.

CAS: 53-84-9

Catalog number

10 004 618 103

Pack size

custom fill

Will be supplied as "b-Nicotinamide-adenine Dinucleotide, I". Unit of measure is "kg".

For further processing only.

Properties**Formula:** C₂₁H₂₇N₇O₁₄P₂**Molecular weight:** 663.4 D**Specification****Appearance:** Colorless to slightly yellowish lyophilizate**Solubility:** Clear, colorless to slightly yellowish solution in water (c=200 mg/mL)**β-NAD** (from value found enzymatically, based on dry weight): ≥99%**β-NAD** (enzymatic, A₃₄₀): ≥96.5%**β-NAD** (A₂₆₀, ε=17.6 [L x mmol⁻¹ x cm⁻¹]): ≥96.5%**NAD** (HPLC): ≥98 area%**Water** (K. Fischer): ≤3.0%**Fe** (AA): ≤50 ppm**AMP** (enzymatic): ≤0.1%**Ethanol** (GC): ≤30 ppm**Aceton, isopropanol, methanol** (GC): each ≤0.05%**Reaction rates (LDH) based on NAD II, acid:** 95-105%**A₂₅₀/A₂₆₀**: 0.81-0.85**A₂₈₀/A₂₆₀**: 0.20-0.24**Stability:** At +2 to +8°C within specification range for 12 months.

Store dry.

NAD, Grade II

free acid

Cofactor for dehydrogenases, *e.g.*, lactate dehydrogenase.**Application**

Use NAD, Grade II as a cofactor in a variety of diagnostic tests, such as for the determination of ethanol and lactate dehydrogenase.

CAS: 53-84-9**Properties****Formula:** C₂₁H₂₇N₇O₁₄P₂**Molecular weight:** 663.4 D**Specification****Appearance:** Colorless to slightly yellowish lyophilizate**Solubility:** Clear, colorless to slightly yellowish solution in water (c=200 mg/mL)**β-NAD** (from value found enzymatically, based on dry weight): ≥97.5%**β-NAD** (enzymatic, A₃₄₀): ≥94.5%**NAD** (A₂₆₀, ε=17.6 [L x mmol⁻¹ x cm⁻¹]): ≥94.5%**NAD** (HPLC): ≥95 area%**Catalog number****10 004 626 103****Pack size**

custom fill

Will be supplied as "b-Nicotinamide-adenine Dinucleotide, II". Unit of measure is "kg".

For further processing only.

Water (K. Fischer): ≤3.5%

Fe (AAS): ≤25 ppm

AMP (enzymatic): ≤0.1%

Ethanol (GC): ≤40 ppm

Aceton, isopropanol, methanol (GC): ≤0.1%, ≤0.15%, ≤0.15%

Reaction rates (LDH) based on NAD II, acid: 95-105%

A_{250}/A_{260} : 0.81-0.85

A_{280}/A_{260} : 0.20-0.24

Stability: At +2 to +8°C within specification range for 12 months.

Store dry.

NADH, Grade I

disodium salt

Cofactor for a variety of dehydrogenases, *e.g.*, malate dehydrogenase and lactate dehydrogenase.

Application

Use NADH, Grade I, as a cofactor in a variety of diagnostic tests, such as for glutamate dehydrogenase, lactate dehydrogenase, α -hydroxybutyrate dehydrogenase, aminotransferases and urea.

CAS: 58-68-4

Properties

Formula: $C_{21}H_{27}N_7O_{14}P_2Na_2$

Molecular weight: 709.4 D (NADH: 665.4 D)

Specification

Appearance: White to slightly yellowish amorphous powder

Solubility: Clear, colorless to slightly yellowish solution in water (c=50 mg/mL)

NADH- Na_2 (calculated from value found enzymatically, based on dry weight): ≥99%

NADH (enzymatic): ≥85%

NADH (A_{340} , $\epsilon=6.3$ [L x mmol⁻¹ x cm⁻¹]): ≥85%

NADH (A_{260} , $\epsilon=14.3$ [L x mmol⁻¹ x cm⁻¹]): ≥85%

Na (flame photometric): 6.5±0.5%

Water (K. Fischer): ≤5%

NAD (enzymatic): ≤0.5%

AMP (enzymatic): ≤0.2%

Ethanol (GC): ≤4%

Reaction rates (LDH) based on standard: 95-105%

A_{260}/A_{340} : ≤2.35

A_{260}/A_{240} : 1.57-2.17

Catalog number

10 004 634 103

Pack size

250 g or custom fill

Will be supplied as "b-NADH, Reduced, Disodium Salt, Grade I".
Unit of measure is "kg".

For further processing only.

Stability: At +2 to +8°C within specification range for 12 months. Keep under nitrogen. Protect from light.

NADH, Grade II

disodium salt

Cofactor for a variety of dehydrogenases, *e.g.*, malate dehydrogenase and lactate dehydrogenase.

Application

Use NADH, Grade II, as a cofactor in a variety of diagnostic tests, such as for glutamate dehydrogenase, lactate dehydrogenase, α -hydroxybutyrate dehydrogenase, aminotransferases and urea.

CAS: 58-68-4

Properties

Formula: C₂₁H₂₇N₇O₁₄P₂Na₂

Molecular weight: 709.4 D (NADH: 665.4 D)

Specification

Appearance: White to slightly yellowish amorphous powder

Solubility: Clear, colorless to slightly yellowish solution in water (c=50 mg/mL)

NADH-Na₂ (calculated from value found enzymatically, based on dry weight): ≥98%

NADH (enzymatic): ≥84%

NADH (A₃₄₀, $\epsilon=6.3$ [l x mmol⁻¹ x cm⁻¹]): ≥84%

NADH (A₂₆₀, $\epsilon=14.3$ [l x mmol⁻¹ x cm⁻¹]): ≥85%

Na (flame photometric): 6.5±0.5%

Water (K. Fischer): ≤6%

NAD (enzymatic): ≤1%

AMP (enzymatic): ≤0.2%

Ethanol (GC): ≤4%

Reaction rates (LDH) based on standard: 95-105%

A₂₆₀/A₃₄₀: ≤2.40

A₂₆₀/A₂₄₀: 1.57-2.17

Stability: At +2 to +8°C within specification range for 12 months. Keep under nitrogen. Protect from light. Store dry.

Catalog number

10 004 642 103

Pack size

250 g or custom fill

Will be supplied as "b-NADH, Reduced, Disodium Salt, Grade II". Unit of measure is "kg".

For further processing only.

NADH, Grade II

for potassium test, disodium salt

NADH quality for enzymatic potassium test.

Application

Use this special NADH as a cofactor in a diagnostic test for potassium, together with glutamate dehydrogenase to remove ammonia from the reaction mixture.

CAS: 58-68-4**Properties****Formula:** C₂₁H₂₇N₇O₁₄P₂Na₂**Molecular weight:** 709.4 D (NADH: 665.4 D)**Specification****Appearance:** White to slightly yellowish amorphous powder**Solubility:** Clear, colorless to slightly yellowish solution in water (c=50 mg/mL)**NADH-Na₂** (from content found enzymatically, based on dry weight): ≥98%**NADH** (enzymatic): ≥82%**NADH** (A₃₄₀, ε=6.3 [L x mmol⁻¹ x cm⁻¹]): ≥84%**NADH** (A₂₆₀, ε=14.3 [L x mmol⁻¹ x cm⁻¹]): ≥85%**Na** (flame photometric): 6.5±0.5%**K** (AAS): ≤250 ppm**Water** (K. Fischer): ≤6%**NAD** (enzymatic): ≤1%**AMP** (enzymatic): ≤0.2%**Ethanol** (GC): ≤4%**Reaction rates (Lactate dehydrogenase), based on standard:**

Freshly dissolved sample: 95-105%

After 8 days incubation at +45°C: ≥95%

HPLC:

Freshly dissolved sample: ≥96.0 area%

After 8 days incubation at +45°C: ≥75 area%

A₂₆₀/A₃₄₀: ≤2.4**A₂₆₀/A₂₄₀**: 1.57-2.17**Stability:** At +2 to +8°C within specification range for 12 months.**Catalog number****11 333 925 103****Pack size**

custom fill

Will be supplied as "NADH Di-Na, Grade II for Potassium Test". Unit of measure is "kg".



For further processing only.

NADP

disodium salt

Cofactor for dehydrogenases, *e.g.*, glucose-6-phosphate dehydrogenase.

Application

Use NADP as a cofactor in a variety of diagnostic tests, such as for the determination of glucose and creatine kinase.

CAS: 53-59-8

Properties

Formula: $C_{21}H_{26}N_7O_{17}P_3Na_2$

Molecular weight: 787.4 D (NADP: 743.4 D)

Specification

Appearance: White to yellowish amorphous powder

Solubility: Clear, colorless to slightly yellowish solution in water (c=50 mg/mL), free of fuzz

NADP, Na-salt (from value found enzymatically, based on dry weight): $\geq 97\%$

NADP (enzymatic, A_{340} , $\epsilon=6.3$ [L x mmol⁻¹ x cm⁻¹]): $\geq 85\%$

NADP (A_{260} , $\epsilon=18$ [L x mmol⁻¹ x cm⁻¹]): $\geq 85\%$

Na (flame photometric): $4.5 \pm 0.5\%$

Water (K. Fischer): $\leq 6\%$

NAD (enzymatic): $\leq 0.5\%$

Methanol (GC): $\leq 2.9\%$

A_{360} (c=10 mg/mL water, against water): ≤ 0.600

A_{340} (c=0.01 mg/mL phosphate buffer, pH 7.0): ≤ 0.005

Stability: At +2 to +8°C within specification range for 24 months.

Store dry.

Catalog number

10 004 669 103

Pack size

custom fill

Will be supplied as "b-NADP, Disodium Salt". Unit of measure is "kg".

For further processing only.

NADP

monopotassium salt

Cofactor for dehydrogenases, *e.g.*, glucose-6-phosphate dehydrogenase.

Application

Use NADP as a cofactor in a variety of diagnostic tests, such as for the determination of glucose and creatine kinase.

CAS: 53-59-8

Catalog number

10 233 536 103

Pack size

custom fill

Will be supplied as "b-NADP, Monopotassium Salt". Unit of measure is "g".



For further processing only.

Properties**Formula:** $C_{21}H_{27}N_7O_{17}P_3K \times 2 H_2O$ **Molecular weight:** 817.4 D (NADP: 743.4 D)**Specification****Appearance:** White crystalline powder**NADP, K-salt** (calculated from value determined enzymatically, based on dry weight): $\geq 97\%$ **NADP** (enzymatic, A_{340} , $\epsilon=6.3$ [L x mmol⁻¹ x cm⁻¹]): $\geq 88\%$ **NADP** (A_{260} , $\epsilon=18$ [L x mmol⁻¹ x cm⁻¹]): $\geq 88\%$ **K** (flame photometric): 4.0-5.0%**Water** (K. Fischer): 4.5±1.0%**NAD** (enzymatic): $\leq 0.2\%$ **Methanol** (GC): $\leq 2\%$ **Mg** (AAS): ≤ 40 ppm**A₃₆₀** (c=10 mg/mL water, against water): ≤ 0.060 **Stability:** At -15 to -25°C within specification range for 12 months.**NADPH**

tetrasodium salt

Cofactor for glutamate dehydrogenase.

Application

Use NADPH as a cofactor in diagnostic tests for ammonia, urea and creatinine.

CAS: 53-57-6**Properties****Formula:** $C_{21}H_{26}N_7O_{17}P_3Na_4$ **Molecular weight:** 833.4 D (NADPH: 745.4 D)**Specification****Appearance:** White to slightly yellowish amorphous powder**Solubility:** Clear, colorless to slightly yellowish solution in water (c=50 mg/mL)**Particle size** (screen analysis according to Ph.Eur., US mesh <30): $\geq 95\%$ **NADPH-Na₄** (calculated from content found enzymatically, based on dry weight): $\geq 97\%$ **NADPH** (enzymatic, A_{340}): $\geq 79\%$ **NADPH** (A_{340} , $\epsilon=6.3$ [L x mmol⁻¹ x cm⁻¹]): $\geq 79\%$ **NADPH** (A_{260} , $\epsilon=15$ [L x mmol⁻¹ x cm⁻¹]): $\geq 80\%$ **NADPH** (HPLC): ≥ 95 area%**Na** (flame photometric): 10-12%**Catalog number****10 041 939 103****Pack size**

custom fill

Will be supplied as "b-NADPH, Reduced, Tetrasodium Salt". Unit of measure is "g".

For further processing only.

Water (K. Fischer): ≤6%
NADH (HPLC): ≤0.5 area%
NADP (HPLC): ≤0.5 area%
Nicotinic acid amide (HPLC): ≤2 area%
NADP (enzymatic): ≤0.5%
Ethanol (GC): ≤3%
 A_{260}/A_{340} : 2.32-2.65
Stability: At +2 to +8°C within specification range for 12 months.

Thio-NAD

free acid

NAD analog with an absorbance maximum of thio-NADH at 405 nm.

Application

Use Thio-NAD instead of NAD in enzymatic reactions to measure the reaction kinetics at 405 nm, *e.g.*, kinetic enzyme cycling methods using Thio-NAD and NADH.

CAS: 4090-29-3

Properties

Formula: C₂₁H₂₇N₇O₁₃SP₂
Molecular weight: 679.5 D

Specification

Appearance: Yellowish powder
Solubility: Clear, colorless to slightly yellowish solution in water (c=37 mg/mL)
Thio-NAD (A_{259} , $\epsilon=19.7$ [L x mmol⁻¹ x cm⁻¹]): ≥95%
Water (K. Fischer): ≤4%
 A_{398}/A_{340} (against water) : ≤0.124
 A_{236}/A_{259} (against water) : 0.640-0.670
 A_{296}/A_{259} (against water) : 0.244-0.264
Revision of absorption in presence of 0.4 mmol/L NADH: ≤0.001 abs/min
Ca (AAS): No limit
Mg (AAS): No limit
Stability: At -15 to -25°C within specification range for 12 months.

Catalog number

04 635 396 103

Pack size

custom fill

Will be supplied as "Thio-NAD free acid". Unit of measure is "g".



For further processing only.

Adenosine-5'-O-(2-thiodiphosphate)

trilithium salt

Non-hydrolyzable ADP analog.

Application

Use Adenosine-5'-O-(2-thiodiphosphate) to inhibit ADP-binding enzymes.

CAS: 73536-95-5**Properties****Formula:** C₁₀H₁₂N₅O₉P₂SLi₃**Molecular weight:** 461.0 D (ATP-β-S: 443.2 D)**Specification****Appearance:** White powder**ATP-β-S, Li₃(A₂₆₀):** ≥81%**ATP-β-S (A₂₆₀, ε= 15.0 [L x mmol⁻¹ x cm⁻¹]):** ≥78%**ATP-β-S (HPLC):** ≥90 area%**Water (K. Fischer):** ≤12%**ATP (HPLC):** ≤5 area%**ADP (HPLC):** ≤1 area%**AMP (HPLC):** ≤4 area%**A₂₅₀/A₂₆₀:** 0.75-0.83**A₂₈₀/A₂₆₀:** 0.14-0.18**A₂₉₀/A₂₆₀:** 0.00-0.01**Stability:** At -15 to -25°C within specification range for 12 months.**Catalog number****10 200 166 103****Pack size**

custom fill

Will be supplied as "Adenosine-5'-O-(2-thiodiphosphate) Tri-Li". Unit of measure is "g".



For further processing only.

Adenosine-5'-O-(3-thiotriphosphate)

tetralithium salt

Non-hydrolyzable ATP analog.

Application

Use Adenosine-5'-O-(3-thiotriphosphate) to inhibit ATP-binding enzymes.

CAS: 35094-45-2**Properties****Formula:** C₁₀H₁₂N₅O₁₂P₃SLi₄**Molecular weight:** 547.0 D (ATP-γ-S: 523.2 g/mol)**Catalog number****10 122 734 103****Pack size**

custom fill

Will be supplied as "Adenosine-5'-O-(3-thiotriphosphate), Li4". Unit of measure is "g".



For further processing only.

Specification**Appearance:** White powder**ATP- γ -S, Li₄** (A_{260}): $\geq 78\%$ **ATP- γ -S** (A_{260} , $\epsilon=15.0$ [L x mmol⁻¹ x cm⁻¹]): $\geq 74\%$ **ATP- γ -S** (HPLC): ≥ 85 area%**Li** (flame photometric): 3-5%**Water** (K. Fischer): $\leq 12\%$ **ADP** (HPLC): ≤ 12 area%**AMP** (HPLC): ≤ 3 area% **A_{250}/A_{260}** : 0.79 ± 0.02 **A_{280}/A_{260}** : 0.16 ± 0.01 **A_{290}/A_{260}** : ≤ 0.05 **Stability:** At -15 to -25°C within specification range for 6 months.

Protect from light. Store dry.

ADP**potassium salt**

Cofactor for diagnostic tests.

Application

Use ADP in variety of diagnostic tests, such as for the determination of creatine kinase and pyruvate kinase. Use it also for the activation of glutamate dehydrogenase in the determination of for example urea or ammonia.

CAS: 58-64-0**Properties****Formula:** C₁₀H₁₄N₅O₁₀P₂K x 2 H₂O**Molecular weight:** 501.3 D (ADP: 427.2 D)**Remark:** Crystalline ADP-K x 2 H₂O is the purest and most stable form of ADP available.**Specification****Appearance:** Colorless crystals**Solubility:** Clear, colorless solution in water (c=50 mg/mL)**ADP-K x 2 H₂O** (based on value found enzymatically): $\geq 98\%$ **ADP** (enzymatic): $\geq 84\%$ **ADP** (A_{260} , $\epsilon=15$ [L x mmol⁻¹ x cm⁻¹]): $\geq 84\%$ **K** (flame photometric): $7.8 \pm 0.5\%$ **Water** (K. Fischer): $7.2 \pm 1\%$ **P_i** (Fiske and Subbarow): $\leq 0.3\%$ **AMP** (enzymatic): $\leq 1\%$ **ATP** (enzymatic): $\leq 0.2\%$ **NH₄** (enzymatic): $\leq 0.005\%$ **Catalog number****10 233 528 103****Pack size**

custom fill

Will be supplied as "Adenosine-5'-diphosphate (ADP), K-Salt". Unit of measure is "kg".

For further processing only.

A_{250}/A_{260} : 0.78±0.02

A_{280}/A_{260} : 0.16±0.01

A_{290}/A_{260} : ≤0.01

Stability: At +2 to +8°C within specification range for 24 months.

Store dry.

ADP

disodium salt

Cofactor for diagnostic tests.

Application

Use ADP in variety of diagnostic tests, such as for the determination of creatine kinase and pyruvate kinase. Use it also for the activation of glutamate dehydrogenase in the determination of for example urea or ammonia.

CAS: 58-64-0

Properties

Formula: C₁₀H₁₃N₅O₁₀P₂Na₂

Molecular weight: 471.2 D (ADP: 427.2 D)

Remark: ATP and AMP may form during storage.

Specification

Appearance: White lyophilizate

Solubility: Clear, colorless solution in water (c=50 mg/mL)

ADP-Na-salt (calculated on value found enzymatically): ≥90%

ADP (enzymatic): ≥82%

ADP (A_{260} , $\epsilon = 15$ [L x mmol⁻¹ x cm⁻¹]): ≥82%

Na (flame photometric): 9±1%

Water (K. Fischer): ≤7%

P_i (Fiske and Subbarow): ≤0.6%

AMP (enzymatic): ≤3%

ATP (enzymatic): ≤1%

NH₄ (enzymatic): ≤0.01%

A_{250}/A_{260} : 0.76-0.80

A_{280}/A_{260} : 0.15-0.17

A_{290}/A_{260} : ≤0.01

Stability: At -15 to -25°C within specification range for 6 months.

Catalog number

10 129 062 103

Pack size

custom fill

Will be supplied as "ADP, Di-Na". Unit of measure is "kg".

For further processing only.

ADP for potassium test

free acid

Application

Use ADP for potassium testing as an activator for glutamate dehydrogenase in enzymatic potassium tests.

CAS: 58-64-0

PropertiesFormula: $C_{10}H_{15}N_5O_{10}P_2$

Molecular weight: 427.2 D

Remark: ATP and AMP may form during storage.

Specification

Appearance: White crystallate

Solubility: Clear, colorless solution in NaOH, 0.1 mol/L (c=50 mg/mL)

IR-spectrum (KBr-pellet): Corresponds to reference

ADP (enzymatic): $\geq 97\%$ ADP (A_{260} , $\epsilon = 15$ [L x mmol⁻¹ x cm⁻¹]): $\geq 97\%$ Water (K. Fischer): $\leq 2\%$ P_i (Fiske and Subbarow): $\leq 0.6\%$ AMP (enzymatic): $\leq 3\%$ ATP (enzymatic, HK/G6P-DH): $\leq 0.3\%$ Na (AAS): ≤ 750 ppmK (AAS): ≤ 20 ppmNH₄ (enzymatic): $\leq 0.01\%$ A_{250}/A_{260} : 0.78 ± 0.02 A_{280}/A_{260} : 0.16 ± 0.01 A_{290}/A_{260} : ≤ 0.01

Stability: At -15 to -25°C within specification range for 12 months.

Catalog number

11 333 879 103

Pack size

custom fill

Will be supplied as "ADP Free Acid for Potassium Test". Unit of measure is "kg".



For further processing only.

AMP

free acid

Nucleotide for diagnostic tests.

Application

Use AMP for the determination of 5'-nucleotidase and in diagnostic tests for the determination of creatine kinase.

CAS: 61-19-8

PropertiesFormula: $C_{10}H_{14}N_5O_7P \times H_2O$

Molecular weight: 365.2 D (AMP: 347.2 D)

Catalog number

10 000 086 103

Pack size

custom fill

Will be supplied as "Adenosine-5'-monophosphoric Acid (AMP)". Unit of measure is "kg".

For further processing only.

Specification**Appearance:** White, crystalline powder**Solubility:** Clear, colorless solution in NaOH, 1 mol/L (c=50 mg/mL)**AMP x H₂O** (based on value found enzymatically): ≥98%**AMP** (enzymatic): ≥93%**AMP** (A_{260} , $\epsilon=15.0$ [L x mmol⁻¹ x cm⁻¹]): ≥93%**Water** (K. Fischer): 5±2%**P_i**: ≤0.3%**Fe** (bathophenanthrolin): ≤10 ppm**Heavy metals** (as Pb): ≤10 ppm**A₂₅₀/A₂₆₀**: 0.78±0.02**A₂₈₀/A₂₆₀**: 0.15±0.01**A₂₉₀/A₂₆₀**: ≤0.01**Stability:** At +15 to +25°C within specification range for 36 months.**AMP**

disodium salt

Nucleotide for diagnostic tests.

Application

Use AMP for the determination of 5'-nucleotidase and in diagnostic tests for the determination of creatine kinase.

CAS: 61-19-8**Properties****Formula:** C₁₀H₁₂N₅O₇PNa₂**Molecular weight:** 391.2 D (AMP: 347.2 D)**Specification****Appearance:** White crystals**Solubility:** Clear, colorless solution in water (c=50 mg/mL)**Water** (K. Fischer): ≤12%**P_i**: ≤0.3%**A₂₅₀/A₂₆₀**: 0.76-0.80**A₂₈₀/A₂₆₀**: 0.15-0.17**A₂₉₀/A₂₆₀**: ≤0.01**Stability:** At +15 to +25°C within specification range for 36 months.**Catalog number****10 000 094 103****Pack size**

custom fill

Will be supplied as "Adenosine-5'-monophosphate (AMP), Di-Na".
Unit of measure is "kg".

For further processing only.

ATP, Grade I

disodium salt

Cofactor for kinases, *e.g.* glycerokinase and hexokinase.**Application**

Use ATP, Grade I, in a variety of diagnostic tests, such as for the determination of triglycerides, creatine kinase and glucose.

CAS: 56-65-5**Properties****Formula:** $C_{10}H_{14}N_5O_{13}P_3Na_2 \times 3 H_2O$ **Molecular weight:** 605.2 D (ATP: 507.2 D)**Specification****Appearance:** White crystals**Solubility:** Clear, colorless solution in water (c=50 mg/mL)**ATP-Na₂H₂ x 3 H₂O** (based on value found enzymatically): ≥99%**ATP** (enzymatic): ≥84%**ATP** (A_{260} , $\epsilon=15.0$ [L x mmol⁻¹ x cm⁻¹]): ≥84%**Na** (flame photometric): 7.5±0.5%**Water** (K. Fischer): ≤8±1%**P_i**: ≤0.15%**ADP, AMP** (enzymatic): ≤0.5%**GTP** (HPLC): ≤0.01 area%**Fe** (AAS): ≤10 ppm**Mg** (AAS): ≤10 ppm**Ca** (AAS): ≤20 ppm**Zn** (AAS): ≤5 ppm**V** (AAS): ≤1 ppm A_{250}/A_{260} : 0.79±0.02 A_{280}/A_{260} : 0.15±0.01 A_{290}/A_{260} : ≤0.01**Stability:** At +2 to +8°C within specification range for 24 months.**Catalog number****10 422 495 103****Pack size**

custom fill

Will be supplied as "ATP, Di-Na, Special Quality". Unit of measure is "kg".

For further processing only.

ATP, Grade II

disodium salt

Cofactor for kinases, *e.g.* glycerokinase and hexokinase.**Application**

Use ATP, Grade II in a variety of diagnostic tests, such as for the determination of triglycerides and glucose.

CAS: 51963-61-2**Catalog number****10 000 116 103****Pack size**

custom fill

Will be supplied as "Adenosine-5'-triphosphate (ATP), Di-Na". Unit of measure is "kg".

For further processing only.

Properties**Formula:** $C_{10}H_{14}N_5O_{13}P_3Na_2 \times 3 H_2O$ **Molecular weight:** 605.2 D (ATP: 507.2 D)**Specification****Appearance:** White crystals**Solubility:** Clear, colorless solution in water (c=50 mg/mL)**ATP- $Na_2H_2 \times 3 H_2O$** (based on value found enzymatically): $\geq 98\%$ **ATP** (enzymatic): $\geq 82\%$ **ATP** (A_{260} , $\epsilon=15.0$ [$L \times mmol^{-1} \times cm^{-1}$]): $\geq 82\%$ **Na** (flame photometric): $7.5 \pm 0.5\%$ **Water** (K. Fischer): $\leq 10\%$ **P_i** : $\leq 0.3\%$ **ADP, AMP** (enzymatic): $\leq 0.5\%$ **GTP** (HPLC): ≤ 0.01 area%**Fe** (AAS): ≤ 15 ppm**Heavy metals** (as Pb): ≤ 30 ppm **A_{250}/A_{260}** : 0.79 ± 0.02 **A_{280}/A_{260}** : 0.15 ± 0.01 **A_{290}/A_{260}** : ≤ 0.01 **Stability:** At +2 to +8°C within specification range for 24 months.**Guanosine-5'-O-(2-thiodiphosphate)**

trilithium salt

Non-hydrolyzable GDP analog.

Application

Use Guanosine-5'-O-(2-thiodiphosphate) to completely inhibit G protein activation by GTP and GTP analogs.

CAS: 71376-97-1**Properties****Formula:** $C_{10}H_{12}N_5O_{10}P_2SLi_3$ **Molecular weight:** 477.0 D (GDP- β -S: 459.3 D)**Specification****Appearance:** White powder**GDP- β -S** (HPLC): ≥ 85 area%**Li** (flame photometric): $4 \pm 1\%$ **Water** (K. Fischer): $\leq 12\%$ **GMP** (HPLC): ≤ 10 area%**Stability:** At -15 to -25°C within specification range for 12 months.**Catalog number****10 526 134 103****Pack size**

custom fill

Will be supplied as "Guanosine-5'-O-(2-thiodiphosphate) Tri-Li".
Unit of measure is "g".


For further processing only.

Enzymes


Colorimetric Tests	48	Glucose Oxidase (GOD), Grade II	84	L-Lactate Dehydrogenase (L-LDH)	118
UV-Tests	50	Glucose Oxidase (GOD), chemically modified	85	L-Lactate Dehydrogenase (L-LDH)	119
Acetate-CoA Ligase (Acetyl-CoA Synthetase)	52	Glucose-6-phosphate Dehydrogenase (G6P-DH)	86	L-Lactate Dehydrogenase (L-LDH)	120
Acyl-CoA Synthetase (Long-Chain-Fatty-Acid-CoA Ligase)	52	Glucose-6-phosphate Dehydrogenase (G6P-DH), chemically modified	88	L-Lactate Dehydrogenase (L-LDH), chemically modified	121
Acid Phosphatase	53	Glucose-6-phosphate Dehydrogenase (G6P-DH)	89	L-Lactate Dehydrogenase (L-LDH)	122
Acyl-CoA Oxidase	53	Glucose-6-phosphate Dehydrogenase (G6P-DH)	91	Lactate 2-Monooxygenase (Lactate oxidase), Grade I	122
Alanine Aminotransferase (ALT) (GPT)	54	Glucose-6-phosphate Dehydrogenase (G6P-DH)	92	Lactate 2-Monooxygenase (Lactate oxidase), Grade II	124
Alanine Aminotransferase (ALT) (GPT)	54	Glucose-6-phosphate Dehydrogenase (G6P-DH)	93	Lactate 2-Monooxygenase (Lactate oxidase)	125
Alcohol Dehydrogenase	55	Glucose-6-phosphate Dehydrogenase (G6P-DH)	94	Lipoprotein Lipase	125
Alcohol Dehydrogenase, chemically modified	56	Glucose-6-phosphate Dehydrogenase (G6P-DH)	95	Lipoprotein Lipase, chemically modified	127
Alcohol Dehydrogenase	58	Glucose-6-phosphate Dehydrogenase (G6P-DH)	96	Lysozyme	128
Aldehyde Dehydrogenase	59	Glucose-6-phosphate Dehydrogenase (G6P-DH)	97	Malate Dehydrogenase	129
Aldose 1-Epimerase (Mutarotase)	60	Glucose-6-phosphate Dehydrogenase (G6P-DH)	98	Malate Dehydrogenase, chemically modified	129
Ascorbate Oxidase	60	Glucose-6-phosphate Dehydrogenase (G6P-DH)	99	Malate Dehydrogenase, IFCC Quality	130
Ascorbate Oxidase	62	Glucose-6-phosphate Isomerase	94	Malate Dehydrogenase, IFCC Quality	131
Ascorbate Oxidase, chemically modified	63	α -Glucosidase	95	N-Methylhydantoinase (ATP-hydrolyzing)	132
Aspartate Aminotransferase (AST) (GOT)	65	β -Glucuronidase	96	NAD(P)H Dehydrogenase (quinone) (Diaphorase)	133
Aspartate Aminotransferase (AST) (GOT)	65	Glutamate Dehydrogenase (NAD(P))	97	Nitrate Reductase	134
N-Carbamoylsarcosine Amidase	66	Glutamate Dehydrogenase (NAD(P))	98	Oxalate Oxidase	134
Cholesterol Esterase, Grade I	67	Glutamate Dehydrogenase (NAD(P))	100	Peroxidase (POD), Grade II	135
Cholesterol Esterase, chemically modified	68	Glutamate Dehydrogenase (NAD(P))	101	Phosphogluconate Dehydrogenase (decarboxylating)	136
Cholesterol Esterase, Grade II	69	γ -Glutamyltransferase	103	6-Phosphogluconolactonase	137
Cholesterol Esterase	70	Glycerol Kinase (GK)	103	Pyruvate Kinase	139
Cholesterol Esterase	71	Glycerol Kinase (GK)	105	Pyruvate Kinase	140
Cholesterol Oxidase	72	Glycerol Kinase (GK), concentrated	106	Pyruvate Oxidase	140
Cholesterol Oxidase	74	Glycerol-3-phosphate Dehydrogenase	107	Sarcosine Oxidase	141
Citrate Lyase	75	Glycerol-3-phosphate Oxidase	107	Triose-phosphate Isomerase	143
Citrate Synthase	75	Glycerol-3-phosphate Oxidase, chemically modified	109	Thrombin	143
Colipase	76	Hexokinase (HK)	110	Urease	144
Creatinase	77	Hexokinase (HK), chemically modified	112	Uricase	145
Creatininase	78	Hexokinase (HK)	113		
Creatinine Deaminase	79	D-Lactate Dehydrogenase (D-LDH)	115		
Formate Dehydrogenase	80	D-Lactate Dehydrogenase (D-LDH), Grade I	115		
Galactose 1-Dehydrogenase	81	D-Lactate Dehydrogenase (D-LDH), Grade II	117		
Galactose 1-Dehydrogenase	81				
Galactose 1-Dehydrogenase	82				
Glucose Oxidase (GOD), Grade I	82				

Colorimetric Tests


*Indicator reaction




Alkaline phosphatase (ALP)







Cholesterol




α -Amylase




Creatinine




Glucose




γ -Glutamyltransferase




Lactate




Lipase



Triglycerides



Uric acid




UV Tests


Alanine aminotransferase (ALT)




Aspartate aminotransferase (AST)




Creatine kinase




Ethanol




Glucose



Lactate dehydrogenase



Urea



Acetate-CoA Ligase (Acetyl-CoA Synthetase)

from yeast, lyophilizate

Ligase that catalyzes the synthesis of acetyl-CoA from acetate and coenzyme A.

Application

Use Acetate-CoA Ligase (Acetyl-CoA Synthetase) in diagnostic tests for the determination of acetic acid in combination with Citrate Synthase, Catalog No. 10 153 605 103 and Malate Dehydrogenase, Catalog No. 10 200 387 103

EC 6.2.1.1

Specification

Appearance: White lyophilizate

pH value (c=20 mg/mL in water): 6.8-7.8

Specific activity (+37°C, acetate): ≥4 U/mg protein

Protein (Biuret): ≤0.25 mg/mg lyophilizate

Acetate (enzymatic): ≤0.1%

Stability: At +2 to +8°C within specification range for 18 months.

Store dry.

Catalog number

10 128 180 103

Pack size

custom fill

Will be supplied as "Acetyl-CoA Synthetase from Yeast". Unit of measure is "kU".

For further processing only.

Acyl-CoA Synthetase (Long-Chain-Fatty-Acid-CoA Ligase)

from microorganism, lyophilizate

Synthetase that catalyzes the synthesis of acyl-CoA from a long-chain carboxylic acid and coenzyme A.

Application

Use Acyl-CoA Synthetase in diagnostic tests for the determination of free fatty acids in combination with Acyl-CoA Oxidase, Catalog No. 10 885 550 103.

EC 6.2.1.3

Specification

Appearance: White lyophilizate

Activity (+37°C, enzymatic): ≥1.5 U/mg lyophilizate

Specific Activity (+37°C, enzymatic): ≥1.5 U/mg protein

Absorbance of the solution (2 U/mL):

$A_{400} \leq 0.05$

$A_{500} \leq 0.025$

$A_{650} \leq 0.012$

Catalog number

10 885 568 103

Pack size

custom fill

Will be supplied as "Acyl-CoA-Synthetase, Lyo.". Unit of measure is "kU".



For further processing only.

Contaminants:

Lipase (indirect): -15% bis +15%

Stability: At -15 to -25 °C within specification range for 12 months.

Store dry in tightly sealed containers.

Acid Phosphatase

from potato, lyophilizate

Application

Use Acid Phosphatase in your controls or calibrators.

EC 3.1.3.2

Specification

Appearance: Pale brown lyophilizate

Activity (+25°C, 4-nitrophenyl phosphate): ≥2 U/mg lyophilizate

Stability: At +2 to +8°C within specification range for 12 months.

Store dry.

Catalog number

10 154 393 103

Pack size

custom fill

Will be supplied as "Phosphatase, Acid, Grade II from Potato". Unit of measure is "kU".

For further processing only.

Acyl-CoA Oxidase

from microorganisms, lyophilizate

Oxidoreductase that catalyzes the interconversion of acyl-CoA to trans-2,3-dehydroacyl-CoA.

Application

Use Acyl-CoA Oxidase in diagnostic tests for the determination of free fatty acids in combination with Acetate-CoA Ligase (Acetyl-CoA Synthetase), Catalog Nos. 10 885 568 103 or 10 128 180 103.

EC 1.3.3.6

Specification

Appearance: Yellowish lyophilizate powder

Activity (+37°C, enzymatic): ≥20 U/mg lyophilizate

Specific activity (enzymatic): ≥20 U/mg protein

Absorbance of the solution (20 U/mL):

A_{400} : ≤0.08

A_{500} : ≤0.04

A_{650} : ≤0.02

Contaminants:

Catalase: ≤12 U/U Acyl-CoA oxidase

Stability: At -15 to -25 °C within specification range for 12 months.

Catalog number

10 885 550 103

Pack size

custom fill

Will be supplied as "Acyl-CoA-Oxidase, Lyo.". Unit of measure is "kU".



For further processing only.

Alanine Aminotransferase (ALT) (GPT)

from pig heart, lyophilizate

Application

Use Alanine Aminotransferase for designing your calibrator/control reagent and for the synthesis of unnatural L-amino acids from α -keto acids.

EC 2.6.1.2

Specification**Appearance:** White to yellowish lyophilizate**pH value** (c=10 mg/mL in water): 7.0-8.0**Activity** (+25°C, L-alanine, α -oxoglutarate): ≥ 3 U/mg lyophilizate**Activity** (+37°C, ALT (ALAT/GPT)-kit): ≥ 4.8 U/mg lyophilizate**Contaminants** (expressed as percentage of Alanine aminotransferase activity):Contaminating oxidases (FOX): ≤ 0.7 Glutamate dehydrogenase: ≤ 0.01 Aspartate aminotransferase (AST/GOT): ≤ 0.135 Lactate dehydrogenase: ≤ 0.01 Malate dehydrogenase: ≤ 0.01 **SVD free:** Corresponds to specification**pH 5.5 treatment** (30 minutes): Corresponds to specification**Stability:** At +2 to +8°C within specification range for 12 months.

Store dry.

Catalog number**10 170 674 103****Pack size**

custom fill

Will be supplied as "GPT from Pig Heart". Unit of measure is "kU".

For further processing only.

Alanine Aminotransferase (ALT) (GPT)

from pig heart, suspension

Application

Use Alanine Aminotransferase for designing your calibrator/control reagent and for the synthesis of unnatural L-amino acids from α -keto acids.

EC 2.6.1.2

Specification**Appearance:** Slightly yellow suspension in ammonium sulfate, 3.2 mol/L**pH value:** 5.5-6.5**Specific activity** (+25 °C; L-alanine, α -ketoglutarate): ≥ 80 U/mg protein**Protein** (Biuret): ≥ 10 mg/mL (standardized to 10 ± 1 mg/mL)**Ammonium sulfate:** 3.2 ± 0.2 mol/L**Contaminants** (expressed as percentage of Alanine Aminotransferase**Catalog number****10 153 443 103****Pack size**

custom fill

Will be supplied as "GPT from Pig Heart". Unit of measure is "kU".

For further processing only.

activity):

Glutamate dehydrogenase: ≤ 0.01

Aspartate aminotransferase (AST/GOT) : ≤ 0.03

Lactate dehydrogenase: ≤ 0.01

Malate dehydrogenase: ≤ 0.01

SVD free: Corresponds to specification

pH 5.5 treatment (30 minutes): Corresponds to specification

Stability: At +2 to +8°C within specification range for 12 months.

Alcohol Dehydrogenase

from yeast, lyophilizate

Dehydrogenase that catalyzes the interconversion of alcohols to the corresponding aldehydes.

Application

Use Alcohol Dehydrogenase in diagnostic tests for the determination of alcohol or aldehyde (formate).

EC 1.1.1.1

Properties

Nomenclature: Alcohol:NAD⁺ oxidoreductase

Molecular weight: 141 kD (pH 7.0)

Isoelectric point: 5.4-5.8

Michaelis constants (Phosphate buffer, pH 7.15, +25°C):

Ethanol: 1.3×10^{-2} mol/L

NAD: 7.4×10^{-5} mol/L

Acetaldehyde: 7.8×10^{-4} mol/L

NADH: 1.1×10^{-5} mol/L

Inhibitor constants (Phosphate buffer, pH 7.15, +25°C):

Ethanol: 4.3×10^{-2} mol/L

NAD: 6.1×10^{-4} mol/L

Acetaldehyde: 6.7×10^{-4} mol/L

NADH: 1.8×10^{-5} mol/L

Inhibitors:

- SH-reagents and heavy metals, such as derivatives, 4-chloromercuribenzoate, iodoacetic acid, N-substituted maleinimides, Hg²⁺, Ag⁺ and Cu²⁺.
- Complexing agents, *e.g.*, o- phenanthroline, EDTA, oxalate.
- NAD analogs and NAD partial structures, *e.g.*, NADP, NADH, ADP, ADP-ribose.
- Substances, which react with enzyme bound NAD, *e.g.*, sulfite, hydroxylamine, cyanide.
- Substrate analogs, *e.g.*, fluoroethanol.
- Oxidizers, *e.g.*, H₂O₂ and aerial oxygen inactivate by oxidation of

Catalog number

11 452 541 103


Pack size

custom fill


Will be supplied as "Alcohol Dehydrogenase, Yeast". Unit of measure is "MU".



For further processing only.



Temperature dependence



Incubation:
10 min
(NH₄)₂SO₄, 3.2 mol/l;
pH 6.0
18 090 U ADH/ml

Thermal stability

essential groups.

pH optimum: 9.0 (see figure)

Temperature dependence: See figure

pH stability: 6.0-8.0 (see figure)

Thermal stability: Up to +50°C (see figure)

Specificity: Alcohol dehydrogenase oxidizes primary alcohols. Isopropanol and secondary butanol are slowly oxidized, while higher secondary and tertiary alcohols do not react. Numerous aldehydes are reduced in the reverse reaction. The enzyme does not react with NADP.

Remarks: Alcohol dehydrogenase tends to show turbidity in solution at +37°C storage. Modified Alcohol dehydrogenase shows no turbidity for at least 4 weeks in solution at +37°C.

Specification

Appearance: White lyophilizate (50 mg lyophilizate contain approximately 30 mg enzyme protein, 15 mg sucrose, 5 mg phosphate)

pH value (c=50 mg/mL in water): 7.0-8.0

Activity (+25°C, ethanol): No limit

Specific activity: ≥400 U/mg (protein)

Protein (Biuret): ≥0.5 mg/mg lyophilizate


Contaminants (expressed as percentage of Alcohol Dehydrogenase activity):

Lactate dehydrogenase: ≤0.01


Malate dehydrogenase ≤0.01

Stability: At -15 to -25°C within specification range for 12 months.

Store dry.



pH optimum



pH stability

Incubation:
25°C, 120 min
pH 3.0 – 5.0:
citrate buffer, 0.2 mol/l
pH 6.0 – 8.0:
phosphate buffer,
0.2 mol/l
pH 9.0 – 11.0:
glycine buffer, 0.2 mol/l
180 U ADH/ml

Alcohol Dehydrogenase, chemically modified

from yeast, lyophilizate

Dehydrogenase that catalyzes the interconversion of alcohols to the corresponding aldehydes.

Application

Use Alcohol Dehydrogenase in diagnostic tests for the determination of alcohol or aldehyde (formate).

EC 1.1.1.1

Properties

Nomenclature: Alcohol:NAD⁺ oxidoreductase

Molecular weight: 141 kD (pH 7.0)

Isoelectric point: 5.4-5.8

Michaelis constants (Phosphate buffer, pH 7.15, +25°C):

Ethanol: 1.3×10^{-2} mol/L

Catalog number

11 644 980 103

Pack size

custom fill

Will be supplied as "Alcohol Dehydrogenase, Yeast, Modified". Unit of measure is "MU".

For further processing only.

NAD: 7.4×10^{-5} mol/L

Acetaldehyde: 7.8×10^{-4} mol/L

NADH: 1.1×10^{-5} mol/L

Inhibitor constants (Phosphate buffer, pH 7.15, +25°C):

Ethanol: 4.3×10^{-2} mol/L

NAD: 6.1×10^{-4} mol/L

Acetaldehyde: 6.7×10^{-4} mol/L

NADH: 1.8×10^{-5} mol/L

Inhibitors:

- SH-reagents and heavy metals, such as derivatives, 4-chloromercuribenzoate, iodoacetic acid, N-substituted maleinimides, Hg^{2+} , Ag^{+} and Cu^{2+} .
- Complexing agents, *e.g.*, o-phenanthroline, EDTA, oxalate.
- NAD analogs and NAD partial structures, *e.g.*, NADP, NADH, ADP, ADP-ribose.
- Substances, which react with enzyme bound NAD, *e.g.*, sulfite, hydroxylamine, cyanide.
- Substrate analogs, *e.g.*, fluoroethanol.
- Oxidizers, *e.g.*, H_2O_2 and aerial oxygen inactivate by oxidation of essential groups.

pH optimum: 9.0 (see figure for product Cat. No. 11 452 541 103)

Temperature dependence: See figure for product Cat. No. 11 452 541 103

pH stability: 6.0-8.0 (see figure for product Cat. No. 11 452 541 103)

Thermal stability: Up to +50°C (see figure for product Cat. No. 11 452 541 103)

Specificity: Alcohol dehydrogenase oxidizes primary alcohols. Isopropanol and secondary butanol are slowly oxidized, while higher secondary and tertiary alcohols do not react. Numerous aldehydes are reduced in the reverse reaction. The enzyme does not react with NADP.

Specification

Appearance: White lyophilizate

Solubility: Clear, colorless solution in water ($c = 50$ mg/mL)

pH value: 6.5-8.0

Activity (+25°C, ethanol): ≥ 25 U/mg lyophilizate

Contaminants (expressed as percentage of Alcohol Dehydrogenase activity):

Lactate dehydrogenase: ≤ 0.01

Malate dehydrogenase: ≤ 0.01

Stability: At +2 to +8°C within specification range for 12 months.

Alcohol Dehydrogenase

from yeast, suspension

Dehydrogenase that catalyzes the interconversion of alcohols to the corresponding aldehydes.

Application

Use Alcohol Dehydrogenase in diagnostic tests for the determination of alcohol or aldehyde (formate).

EC 1.1.1.1

Properties

Nomenclature: Alcohol:NAD⁺ oxidoreductase

Molecular weight: 141 kD (pH 7.0)

Isoelectric point: 5.4-5.8

Michaelis constants (Phosphate buffer, pH 7.15, +25°C):

Ethanol: 1.3×10^{-2} mol/L

NAD: 7.4×10^{-5} mol/L

Acetaldehyde: 7.8×10^{-4} mol/L

NADH: 1.1×10^{-5} mol/L

Inhibitor constants (Phosphate buffer, pH 7.15, +25°C):

Ethanol: 4.3×10^{-2} mol/L

NAD: 6.1×10^{-4} mol/L

Acetaldehyde: 6.7×10^{-4} mol/L

NADH: 1.8×10^{-5} mol/L

Inhibitors:

- SH-reagents and heavy metals, such as derivatives, 4-chloromercuribenzoate, iodoacetic acid, N-substituted maleinimides, Hg²⁺, Ag⁺ and Cu²⁺.
- Complexing agents, *e.g.*, o-phenanthroline, EDTA, oxalate.
- NAD analogs and NAD partial structures, *e.g.*, NADP, NADH, ADP, ADP-ribose.
- Substances, which react with enzyme bound NAD, *e.g.*, sulfite, hydroxylamine, cyanide.
- Substrate analogs, *e.g.*, fluoroethanol.
- Oxidizers, *e.g.*, H₂O₂ and aerial oxygen inactivate by oxidation of essential groups.

pH optimum: 9.0 (see figure for product Cat. No. 11 452 541 103)

Temperature dependence: See figure for product Cat. No. 11 452 541 103

pH stability: 6.0-8.0 (see figure for product Cat. No. 11 452 541 103)

Thermal stability: Up to +50°C (see figure for product Cat. No. 11 452 541 103)

Specificity: Alcohol dehydrogenase oxidizes primary alcohols.

Isopropanol and secondary butanol are slowly oxidized, while higher secondary and tertiary alcohols do not react. Numerous aldehydes are reduced in the reverse reaction. The enzyme does not react with NADP.

Catalog number

11 531 034 103

Pack size

custom fill

Will be supplied as "ADH-Y, As, new". Unit of measure is "g".

For further processing only.

Remark: Alcohol dehydrogenase tends to show turbidity in solution at +37°C storage.

Modified Alcohol dehydrogenase shows no turbidity for at least 4 weeks in solution at +37°C.

Specification

Appearance: White to yellow-brown crystalline suspension in ammonium sulfate solution, 3.2 mol/L, pH approximately 6.5

pH value: 6.0-7.0

Activity: 9000 to 18000 U/mL

Specific activity (+25°C, ethanol): ≥300 U/mg (protein)

Protein (Biuret): 30±3 mg/mL lyophilizate

Contaminants (expressed as percentage of Alcohol Dehydrogenase activity):

Lactate dehydrogenase: ≤0.01

Malate dehydrogenase: ≤0.01

Stability: At +2 to +8°C within specification range for 9 months. Store under nitrogen.

Aldehyde Dehydrogenase

from yeast, lyophilizate

Dehydrogenase that catalyzes the oxidation of aldehydes using NAD(P)⁺ as acceptor.

Application

Use Aldehyde Dehydrogenase in diagnostic tests that use an NADH/ NADPH recycling system.

EC 1.2.1.5

Specification

Appearance: White lyophilizate

Solubility: Clear, colorless solution in water (c=20 mg/mL)

Activity (+25°C, acetaldehyde): ≥2.0 U/mg lyophilizate

Specific activity: ≥20 U/mg protein

Protein (Biuret): No limit (approximately 10%)

Contaminants (expressed as percentage of Aldehyde Dehydrogenase activity):

Alcohol dehydrogenase: ≤0.01

Lactate dehydrogenase: ≤0.01

"NADH oxidase": ≤0.01

"NADPH oxidase": ≤0.01

Stability: At +2 to +8°C within specification range for 12 months.

Store dry. Store under nitrogen.

Catalog number

10 145 947 103

Pack size

custom fill

Will be supplied as "Aldehyde Dehydrogenase from Yeast". Unit of measure is "kU".

For further processing only.

Aldose 1-Epimerase (Mutarotase)

from pig kidney, suspension

Enzyme for mutarotation of sugars.

Application

Use Aldose 1-Epimerase (Mutarotase) in diagnostic tests for the determination of glucose anomers.

EC 5.1.3.3

Specification

Appearance: White suspension in ammonium sulfate, 3.2 mol/L, pH approximately 6

pH value: 5.5-6.5

Specific activity (+25°C, α-D-glucose): ≥5,000 U/mg protein

Protein (Biuret): 5±0.5 mg/mL

SVD free: Corresponds to specification

pH 5.5 treatment (30 minutes): Corresponds to specification

Stability: At +2 to +8°C within specification range for 12 months.

Catalog number

10 152 331 103

Pack size

custom fill

Will be supplied as "Mutarotase from Hog Kidney". Unit of measure is "MU".

For further processing only.

Ascorbate Oxidase

from *Cucurbita* species, lyophilizate

Oxidoreductase that oxidizes ascorbic acid to dehydroascorbate.

Application

Use Ascorbate Oxidase in a variety of diagnostic tests to eliminate the interference of ascorbic acid, since ascorbic acid interferes with the Trinder reaction that is widely used for the colorimetric determination of analytes. It is useful in liquid as well as dry chemistry test, *e.g.*, for the determination of uric acid, lactate or creatinine.

EC 1.10.3.3

Properties

Nomenclature: L-ascorbate:oxygen oxidoreductase

Molecular weight: Approximately 140 kD

Isoelectric point: 5.0-6.0

Michaelis constant (Phosphate buffer, pH 5.6, +25°C):

L-ascorbate: 3×10^{-4} mol/L

Inhibitors: 4-chloromercuribenzoate, CN⁻, Na₂S, diethyl-dithiocarbamate, 8-hydroxyquinoline, K-ethylxanthate

pH optimum: 5.6-7.0 (see figure)

Temperature dependence: See figure

pH stability: 6.5-9.0 (see figure)

Catalog number

10 199 605 103


Pack size

custom fill

Will be supplied as "Ascorbate Oxidase from *Cucurbita* species". Unit of measure is "kU".



For further processing only.



pH optimum

Thermal stability: Up to +70°C (see figure)

Stability of the lyophilizate: Stable at +35°C for at least 3 weeks (see figure).

Stability in solution: See figure

Specificity: Several analogs of ascorbate react.

Remark: A decrease in activity of approximately 10% may occur.

Specification

Appearance: Turquoise lyophilizate

Solubility: Clear, slightly turquoise solution in water (c=50 mg/mL)

pH value (c=50 mg/mL in water): 7.0-8.0

Activity (+25°C, L-ascorbate): ≥170 U/mg lyophilizate

Specific activity (+25°C): ≥1,700 U/mg protein

Protein (Biuret): 0.07-0.14 mg/mg lyophilizate

Contaminants (+25°C; expressed as percentage of Ascorbate Oxidase activity):


Catalase: ≤0.2

Aspartate aminotransferase (AST/GOT): No limit


Alanine aminotransferase (ALT/GPT): No limit

Peroxidase: ≤0.005


Stability: At -15 to -25°C within specification range for 12 months. Keep tightly sealed.




pH stability




Thermal stability



pH stability



pH optimum



Thermal stability

Ascorbate Oxidase

from *Cucurbita* species, poor of AST/ALT, lyophilizate

Oxidoreductase that oxidizes ascorbic acid to dehydroascorbate.

Application

Use Ascorbate Oxidase, poor of AST/ALT, in a variety of diagnostic tests to eliminate the interference of ascorbic acid, since ascorbic acid interferes with the Trinder reaction that is widely used for the colorimetric determination of analytes. It is useful in liquid as well as dry chemistry test, *e.g.*, for the determination of uric acid, lactate, creatinine or transaminases.

EC 1.10.3.3

Properties

Nomenclature: L-ascorbate:oxygen oxidoreductase

Molecular weight: Approximately 140 kD

Isoelectric point: 5.0-6.0

Michaelis constant (Phosphate buffer, pH 5.6, +25°C):

L-ascorbate: 3×10^{-4} mol/L

Inhibitors: 4-chloromercuribenzoate, CN^- , Na_2S , diethyl-dithiocarbamate, 8-hydroxyquinoline, K-ethylxanthate

pH optimum: 5.6-7.0 (see figure for product Cat. No. 10 199 605 103)

Temperature dependence: See figure for product Cat. No. 10 199 605 103

pH stability: 6.5-9.0 (see figure for product Cat. No. 10 199 605 103)

Thermal stability: Up to +70°C (see figure for product Cat. No. 10 199 605 103)

Stability of the lyophilizate: Stable at +35°C for at least 3 weeks (see figure).

Stability in solution: See figure for product Cat. No. 10 199 605 103

Specificity: Several analogs of ascorbate react.

Remark: A decrease in activity of approximately 10% may occur.

Catalog number

11 136 364 103

Pack size

custom fill

Will be supplied as "Ascorbate Oxidase GOT-deficient". Unit of measure is "MU".



For further processing only.

Specification

Appearance: Turquoise lyophilizate

Solubility: Clear, slightly turquoise solution in water (c=50 mg/mL)

pH value (c=50 mg/mL in water): 7.0-8.0

Activity (+25°C, L-ascorbate): ≥170 U/mg lyophilizate

Specific activity (+25°C): ≥1,700 U/mg protein

Protein (Biuret): 8-14 mg/100mg lyophilizate

Contaminants (+25°C, expressed as percentage of Ascorbate Oxidase activity):

Catalase: ≤0.2

Aspartate aminotransferase (AST/GOT): ≤0.0003

Alanine aminotransferase (ALT/GPT): ≤0.0005

Contaminating oxidases (FOX): ≤0.0002

Stability: At -15 to -25°C within specification range for 12 months. Keep tightly sealed.

Ascorbate Oxidase, chemically modified

from *Cucurbita* species, lyophilizate

Oxidoreductase that oxidizes ascorbic acid to dehydroascorbate.

Application

Use Ascorbate Oxidase, chemically modified, in a variety of diagnostic tests to eliminate the interference of ascorbic acid, since ascorbic acid interferes with the Trinder reaction that is widely used for the colorimetric determination of analytes. It is useful in liquid as well as dry chemistry test, e.g., for the determination of uric acid, lactate or creatinine.

EC 1.10.3.3

Properties

Nomenclature: L-ascorbate:oxygen oxidoreductase

Molecular weight: Approximately 140 kD

Isoelectric point: 5.0-6.0

Michaelis constant (Phosphate buffer, pH 5.6, +25°C):

L-ascorbate: 3×10^{-4} mol/L

Inhibitors: 4-chloromercuribenzoate, CN⁻, Na₂S, diethyl-dithiocarbamate, 8-hydroxyquinoline, K-ethylxanthate

pH optimum: 5.6-7.0 (see figure for product Cat. No. 10 199 605 103)

Temperature dependence: See figure

pH stability: 6.5-9.0 (see figure for product Cat. No. 10 199 605 103)

Thermal stability: Up to +70°C (see figure for product Cat. No. 10 199 605 103)

Catalog number

11 558 668 103

Pack size

custom fill

Will be supplied as "AOD, modified". Unit of measure is "MU".

For further processing only.

Stability of the lyophilizate: Stable at +35°C for at least 3 weeks (see figure for product Cat. No. 10 199 605 103).

Stability in solution: See figure for product Cat. No. 10 199 605 103

Specificity: Several analogs of ascorbate react.

Remark: The modified enzyme is especially suited for liquid stable applications with extended shelf life requirements.

Specification

Appearance: Turquoise lyophilizate

Solubility: Clear, slightly turquoise solution in water (c=50 mg/mL)

pH value (c=50 mg/mL in water): 7.2-8.2

Activity (+25°C, L-ascorbate): ≥120 U/mg lyophilizate

Specific activity (+25°C): ≥1,200 U/mg protein

Activity (+37°C, L-ascorbate): ≥180 U/mg lyophilizate

Specific activity (+37°C): ≥1,800 U/mg protein

Protein (BCA): 0.04-0.10 mg/mg lyophilizate

Contaminants (+25°C, expressed as percentage of Ascorbate Oxidase activity):

Catalase: ≤0.2

Glutamate oxalacetate transaminase (AST): ≤0.0003

Glutamate pyruvate transaminase (ALT): ≤0.0005

Contaminating oxidases (FOX): ≤0.0002

Stability: At -15 to -25°C within specification range for 12 months. Keep tightly sealed.

Aspartate Aminotransferase (AST) (GOT)

from pig heart, lyophilizate

Application

Use Aspartate Aminotransferase for designing your calibrator/control reagent and for the synthesis of unnatural L-amino acids from α -keto acids.

EC 2.6.1.1

Specification

Appearance: Slightly yellow lyophilizate

pH value (c=10 mg/mL in water): 7.0-8.0

Activity (+37°C, AST (ASAT/GOT)-kit): ≥ 45 U/mg lyophilizate

Contaminants (expressed as percentage of Aspartate Aminotransferase activity):

Contaminating oxidases (FOX): ≤ 0.7

Glutamate dehydrogenase: ≤ 0.01

Alanine Aminotransferase (ALT/GPT) : ≤ 0.01

Lactate dehydrogenase: ≤ 0.01

Malate dehydrogenase: ≤ 0.01

Oxaloacetate decarboxylase: ≤ 0.01

SDV free: Corresponds to specification

pH 5.5 treatment (30 minutes): Corresponds to specification

Stability: At +2 to +8°C within specification range for 12 months.

Store dry.

Catalog number

10 170 666 103

Pack size

custom fill

Will be supplied as "GOT from Pig Heart". Unit of measure is "kU".

For further processing only.

Aspartate Aminotransferase (AST) (GOT)

from pig heart, suspension

Application

Use Aspartate Aminotransferase for designing your calibrator/control reagent and for the synthesis of unnatural L-amino acids from α -keto acids.

EC 2.6.1.1

Specification

Appearance: Yellow suspension in ammonium sulfate

pH value: 5.5-6.5

Specific activity (+25°C; L-aspartate, α -ketoglutarate): ≥ 200 U/mg protein

Protein: ≥ 10 mg/ml (standardized to 10 ± 1 mg/mL)

Ammonium sulfate: 3.2 ± 0.2 mol/L

Contaminants (expressed as percentage of Aspartate Aminotransferase activity):

Catalog number

10 153 354 103

Pack size

custom fill

Will be supplied as "GOT from Pig Heart". Unit of measure is "MU".

For further processing only.

Glutamate dehydrogenase: ≤ 0.01

Alanine Aminotransferase (ALT/GPT): ≤ 0.01

Lactate dehydrogenase: ≤ 0.01

Malate dehydrogenase: ≤ 0.01

Oxaloacetate decarboxylase: ≤ 0.01

SDV free: Corresponds

pH 5.5 treatment (30 minutes): Corresponds to specification

Stability: At +2 to +8°C within specification range for 24 months.

Store dry.

N-Carbamoylsarcosine Amidase

from *E.coli* overproducer, lyophilizate

Hydrolase that catalyzes the interconversion of N-carbamoylsarcosine to sarcosine.

Application

Use N-Carbamoylsarcosine Amidase in diagnostic tests for the determination of creatinine in combination with Creatinine Deaminase Catalog No. 11 330 764 103, N-Methylhydantoinase (ATP-hydrolysing), Catalog No. 11 288 555 103, and Sarcosine Oxidase, Catalog No. 11 378 856 103.

EC 3.5.1.59

Specification

Appearance: White lyophilizate

Solubility: Clear, colorless solution in water (c=10 mg/mL)

pH value (c=10 mg/mL in water): 7.3-8.3

Activity (+25°C, carbamoylsarcosine): 0.80-1.30 U/mg lyophilizate

Protein (Biuret): 30-50 mg/100 mg lyophilizate

Contaminants (expressed as percentage of Carbamoylsarcosine Amidase activity):

Creatinase: ≤ 0.013

Creatininase: ≤ 0.01

Catalase: ≤ 30

Uricase: ≤ 0.01

Stability: At -15 to -25°C within specification range for 12 months.

Store dry. Protect from light.

Catalog number

11 248 847 103

Pack size

custom fill

Will be supplied as "N-Carbamoyl-sarcosine Hydrolase". Unit of measure is "kU".



For further processing only.

Cholesterol Esterase, Grade I

from *Pseudomonas species*, lyophilizate

Hydrolase that splits fatty acids from sterols.

Application

Use Cholesterol Esterase in diagnostic tests for the determination of cholesterol in combination with Cholesterol Oxidase, Catalog Nos. 10 634 522 103, or 11 479 709 103.

EC 3.1.1.13

Properties

Nomenclature: Sterol-ester acylhydrolase

Molecular weight: ~129 kD

Isoelectric point: 4.5

Michaelis constant (Phosphate buffer, pH 7.5):

Cholesterol oleate: 7×10^{-5} mol/L

Inhibitors: Heavy metals such as Cu^{2+} , Ag^+ , Zn^{2+}

Activators: Detergents

pH optimum: 6.0–8.0; (maximum at pH 7.6) (see figure)

Temperature dependence: Not possible to determine under assay conditions due to turbidity of Thesit at temperatures above +27°C.

pH stability: 6.0–6.5 (see figure)

Thermal stability: Below +20°C (see figure)

Specificity: Cholesterol esterase is an enzyme of lipid metabolism and gives complete cleavage of all serum cholesterol esters.

Remark: This Cholesterol esterase is especially suited for liquid stable applications with extended shelf life requirements.

Specification

Appearance: Brownish lyophilizate

Solubility: Clear, brown solution in water (c=50 mg/mL)

pH value (c=50 mg/mL in water): 7.0–8.0

Activity (+25°C, cholesterol oleate): ≥ 100 U/mg lyophilizate

Specific Activity: ≥ 100 U/mg protein

Protein (Biuret): No limit

Contaminants (expressed as percentage of Cholesterol Esterase activity):

ATPase: ≤ 0.005

Catalase: ≤ 1.00

Glycerokinase: ≤ 0.001

Glucose oxidase: ≤ 0.001

Hexokinase: ≤ 0.005

“NADH oxidase”: ≤ 0.001

Uricase: ≤ 0.005

Stability: At +2 to +8°C within specification range for 12 months.

Store dry.

Catalog number


11 520 857 103

Pack size


custom fill

Will be supplied as “CE, Ps.species, Lyo., SQ”. Unit of measure is “MU”.

For further processing only.




pH optimum



pH stability

Incubation:
25°C, 25 h
K-phosphate buffer,
0.7 mol/l
42.9 U CE/ml



Thermal stability

Incubation:
10 min
K-phosphate buffer,
0.05 mol/l;
pH 6.5
18 U CE/ml

Cholesterol Esterase, chemically modifiedfrom *Pseudomonas* species, lyophilizate

Hydrolase that splits fatty acids from sterols.

Application

Use Cholesterol Esterase, chemically modified in diagnostic tests for the determination of cholesterol in combination with Cholesterol Oxidase, Catalog Nos. 11 479 709 103, or 10 634 522 103.

EC 3.1.1.13

Properties**Nomenclature:** Sterol-ester acylhydrolase**Molecular weight:** ~129 kD**Isoelectric point:** 4.5**Michaelis constant** (Phosphate buffer, pH 7.5):Cholesterol oleate: 7×10^{-5} mol/L**Inhibitors:** Heavy metals such as Cu^{2+} , Ag^{+} , Zn^{2+} **Activators:** Detergents**pH optimum:** 6.0-8.0; (maximum at pH 7.6) (see figure for product Cat. No. 11 520 857 103)**Temperature dependence:** Not possible to determine under assay conditions due to turbidity of Thesit at temperatures above +27°C.**pH stability:** 6.0-6.5 (see figure for product Cat. No. 11 520 857 103)**Thermal stability:** Below +20°C (see figure for product Cat. No. 11 520 857 103)**Specificity:** Cholesterol esterase is an enzyme of lipid metabolism and gives complete cleavage of all serum cholesterol esters.**Specification****Appearance:** Brownish lyophilizate**Solubility:** Clear, brown solution in water (c=50 mg/mL)**Activity** (+25°C, cholesterol oleate): ≥ 10 U/mg lyophilizate**Contaminants** (expressed as percentage of Cholesterol Esterase activity):ATPase: ≤ 0.005 Catalase: ≤ 1.00 Glycerokinase: ≤ 0.001 Hexokinase: ≤ 0.005 "NADH oxidase": ≤ 0.001 Uricase: ≤ 0.005 **NaCl:** 3 ± 0.2 mol/L**Stability:** At +2 to +8°C within specification range for 12 months. Store dry.**Catalog number****11 641 735 103****Pack size**

custom fill

Will be supplied as "Cholesterol Esterase Modified". Unit of measure is "MU".

For further processing only.

Cholesterol Esterase, Grade II

from *Pseudomonas* species, lyophilizate

Hydrolase that splits fatty acids from sterols.

Application

Use Cholesterol Esterase in diagnostic tests for the determination of cholesterol in combination with Cholesterol Oxidase, Catalog Nos. 10 634 522 103 or 11 479 709 103.

EC 3.1.1.13

Properties

Nomenclature: Sterol-ester acylhydrolase

Molecular weight: ~129 kD

Isoelectric point: 4.5

Michaelis constant (Phosphate buffer, pH 7.5):

Cholesterol oleate: 7×10^{-5} mol/L

Inhibitors: Heavy metals such as Cu^{2+} , Ag^{+} , Zn^{2+}

Activators: Detergents

pH optimum: 6.0-8.0; (maximum at pH 7.6) (see figure for product Cat. No. 11 520 857 103)

Temperature dependence: Not possible to determine under assay conditions due to turbidity of Thesit at temperatures above +27°C.

pH stability: 6.0-6.5 (see figure for product Cat. No. 11 520 857 103)

Thermal stability: Below +20°C (see figure for product Cat. No. 11 520 857 103)

Specificity: Cholesterol esterase is an enzyme of lipid metabolism and gives complete cleavage of all serum cholesterol esters.

Specification

Appearance: Slightly yellowish lyophilizate

Solubility: Clear, colorless or yellowish to brownish solution in water (c=50 mg/mL)

pH value (c=50 mg/mL in water): 7.0-8.0

Ficoll (colorimetric): 40-80 mg/100 mg lyophilizate

Activity (+25°C, cholesterol oleate): ≥ 30 U/mg lyophilizate

Specific Activity: ≥ 100 U/mg protein

Protein (Biuret): 0.1-0.4 mg/mg lyophilizate

Contaminants (expressed as percentage of Cholesterol Esterase activity):

ATPase: ≤ 0.005

Catalase: ≤ 200 U/mg lyophilizate

Glycerokinase: ≤ 0.001

GOD: ≤ 0.001

Hexokinase: ≤ 0.005

"NADH oxidase": ≤ 0.005

Uricase: ≤ 0.005

Stability: At +2 to +8°C within specification range for 12 months.

Catalog number

11 015 923 103

Pack size

custom fill

Will be supplied as "Cholesterol Esterase from Pseudom.species".
Unit of measure is "MU".

For further processing only.

Cholesterol Esterase

from *Candida cylindracea*, lyophilizate

Hydrolase that splits fatty acids from sterols.

Application

Use Cholesterol Esterase in diagnostic tests for the determination of cholesterol in combination with Cholesterol Oxidase, Catalog Nos. 10 634 522 103 or 11 479 709 103.

EC 3.1.1.13

Specification

Appearance: Almost white lyophilizate

Solubility: Clear, colorless solution in phosphate buffer, 0.05 mol/L, pH 6.0 (c=10 mg/mL)

pH value: 5.5-6.5

Activity (+25°C; cholesterol oleate): ≥10.5 U/mg lyophilizate

Protein (Lowry): 0.20-0.30 mg/mg lyophilizate

Contaminants (expressed as percentage of Cholesterol Esterase activity):

ATPase: ≤0.005

Glucose oxidase: ≤0.001

Glycerokinase: ≤0.001

Hexokinase: ≤0.005

Catalase: ≤1 U/mg lyophilizate

“NADH oxidase”: ≤0.005

Proteases: No limit

Uricase: ≤0.005

Stability: At +2 to +8°C within specification range for 12 months. Store dry.

Catalog number


10 129 046 103

Pack size


custom fill

Will be supplied as “Cholesterol Esterase, *Candida cylindracea*”. Unit of measure is “MU”.

For further processing only.




pH optimum



Incubation:
25°C, 20 h
pH 3.0 – 5.0:
Na-acetate buffer,
0.7 mol/l
pH 4.0 – 8.2:
K-phosphate buffer,
0.7 mol/l
48.7 U CE/ml

pH stability



Incubation:
10 min
K-phosphate buffer,
0.05 mol/l;
pH 6.5
48.7 U CE/ml

Thermal stability

Cholesterol Esterase

from *Candida cylindracea*, solution

Hydrolase that splits fatty acids from sterols.

Application

Use Cholesterol Esterase in diagnostic tests for the determination of cholesterol in combination with Cholesterol Oxidase, Catalog Nos. 10 634 522 103 or 11 479 709 103. Apply this ready-to-use enzyme directly in your diagnostic test.

EC 3.1.1.13

Specification

Appearance: Clear to turbid, brownish-yellow solution in NaCl

pH value: 5.7-6.3

Specific activity (+25°C; cholesterol oleate): ≥26 U/mg

Volume activity: For information only [U/mL]

Protein (Lowry): ≥2 mg/mL

NaCl (chloride meter): 3±0.2 mol/L

Contaminants (expressed as percentage of Cholesterol Esterase activity):

ATPase: ≤0.005

Glucose oxidase: ≤0.001

Glycerokinase: ≤0.001

Hexokinase: ≤0.005

Catalase: ≤200 U/mg

"NADH oxidase": ≤0.005

Uricase: ≤0.005

Stability: At +2 to +8°C within specification range for 12 months.

Catalog number


10 262 609 103

Pack size


custom fill

Will be supplied as "Cholesterol Esterase, *Candida cylindracea*". Unit of measure is "MU".

For further processing only.




pH optimum



Incubation:
25°C, 20 h
pH 3.0 – 5.0:
Na-acetate buffer,
0.7 mol/l
pH 4.0 – 8.2:
K-phosphate buffer,
0.7 mol/l
48.7 U CE/ml

pH stability



Incubation:
10 min
K-phosphate buffer,
0.05 mol/l;
pH 6.5
48.7 U CE/ml

Thermal stability

Cholesterol Oxidase

from *Brevibacterium* species, expressed in *E.coli*, lyophilizate

Oxidoreductase that catalyzes the interconversion of cholesterol to cholest-4-en-3-one.

Application

Use Cholesterol Oxidase in diagnostic tests for the determination of cholesterol in combination with Cholesterol Esterase, Catalog Nos. 10 129 046 103, 10 262 609 103, 11 015 923 103, 11 520 857 103 or 11 641 735 103.

EC 1.1.3.6

Properties

Nomenclature: Cholesterol:oxygen oxidoreductase

Molecular weight: 60 kD (native and SDS)

Isoelectric point: ~5.0

Michaelis constant (Phosphate buffer, 0.5 mol/l, pH 7.5; +25°C):

Cholesterol: 1×10^{-4} mol/L

Inhibitors: Hg²⁺, ZnCl₂, SDS

Activators: Non ionic detergents

pH optimum: 5.5-8.0 (see figure)

Temperature dependence: See figure

pH stability: 5.0-10.0 (see figure)

Thermal stability: Up to +55°C (see figure)

Storage and Stability: No decrease in activity over 6 weeks at +35°C (see figure)

Specificity:

cholesterol 100%

pregnenolon 52%

stigmasterol 17%

dehydroandrosterone 0.5%

androsterone 0%

estradiol 0%

cholate 0%

Catalog number

11 479 709 103


Pack size

custom fill


Will be supplied as "ChOD, Brevibacterium rec.". Unit of measure is "MU".




For further processing only.



pH optimum



pH stability



Temperature dependence

Specification

Appearance: Yellow lyophilizate

Solubility: Clear, yellowish solution in water (c=10 mg/mL)

pH value: 6.0-7.0

Activity (+25°C, cholesterol): 10-20 U/mg lyophilizate

Protein (Biuret): 10-30 mg/100 mg lyophilizate

Contaminants (expressed as percentage of Cholesterol Oxidase activity):

Catalase: ≤6.0


Glucose oxidase: ≤0.01

"NADH oxidase": ≤0.01


Uricase: ≤0.01

Stability: At -15 to -25°C within specification range for 12 months.

Store dry.



Thermal stability



Stability of the lyophilizate

Cholesterol Oxidase

from *Streptomyces* species, lyophilizate

Oxidoreductase that catalyzes the interconversion of cholesterol to cholest-4-en-3-one.

Application

Use Cholesterol Oxidase in diagnostic tests for the determination of cholesterol in combination with Cholesterol Esterase, Catalog Nos. 10 129 046 103, 10 262 609 103, 11 015 923 103, 11 520 857 103 or 11 641 735 103.

EC 1.1.3.6

Specification

Appearance: Yellow lyophilizate

Solubility: Clear, yellow solution in water (c=20 mg/mL)

pH value (c=20 mg/mL): 7.0-8.0

Activity (+25°C, cholesterol): ≥3.0 to 4.6 U/mg lyophilizate

Specific activity: ≥40.0 U/mg protein

Protein (Biuret): No limit

Contaminants (expressed as percentage of Cholesterol Oxidase activity):

Glucose oxidase: ≤0.01

Catalase: ≤1.00

Uricase: ≤0.01

Stability: At -15 to -25°C within specification range for 12 months. Store dry.

Catalog number

10 634 522 103

Pack size

custom fill

Will be supplied as "Cholesterol Oxidase, *Streptomyces* species". Unit of measure is "kU".



For further processing only.

Citrate Lyase

from *Klebsiella pneumoniae*, lyophilizate

Enzyme that catalyzes the interconversion of oxalacetate and acetate to citrate.

Application

Use Citrate lyase in tests for citric acid in combination with Malate Dehydrogenase, Catalog Nos. 11 866 109 103 and 10 200 387 103, and Lactate Dehydrogenase, Catalog Nos. 11 291 416 103, 12 235 650 103 or 10 003 557 103.

EC 4.1.3.6

Specification

Appearance: Slightly beige lyophilizate

pH value (hydrous solution): 6.5-7.5

Activity (+25 °C, citrate): ≥0.25 U/mg lyophilizate

Contaminants (expressed as percentage of Citrate Lyase activity):

Isocitrate dehydrogenase (NAD specific): ≤0.05

"NADH-oxidase": ≤0.05

Stability: At +2 to +8°C within specification range for 12 months.

Store dry.

Catalog number

10 354 066 103

Pack size

custom fill

Will be supplied as "Citrate Lyase (CL), *Aerobacter aerogenes*". Unit of measure is "kU".

For further processing only.

Citrate Synthase

from pig heart, suspension

Enzyme that catalyzes the formation of citrate from acetyl-CoA and oxalacetate.

Application

Use Citrate Synthase in reagents for acetic acid testing in combination with Acetate-CoA Ligase (Acetyl-CoA Synthetase), Catalog Nos. 10 128 180 103 and 10 885 568 103, and Malate Dehydrogenase, Catalog Nos. 11 866 109 103 and 10 200 387 103.

EC 2.3.3.1

Catalog number

10 153 605 103

Pack size

custom fill

Will be supplied as "Citrate Synthase (CS) from Pig Heart". Unit of measure is "g".

For further processing only.

Specification

Appearance: Slightly grey-brown suspension in ammonium sulfate, 3.2 mol/L; potassium phosphate, 0.02 mol/L; pH approximately 7

Specific activity (+25°C, oxaloacetic acid): ≥110 U/mg protein

Protein (Biuret): ≥10 mg/mL

Contaminants (expressed as percentage of Citrate Synthase activity):
Oxaloacetate decarboxylase: ≤0.1

SVD free: Corresponds to specification

pH 5.5 treatment (30 minutes): Corresponds to specification

Stability: At +2 to +8°C within specification range for 24 months.

Colipase

from porcine pancreas, lyophilizate

Application

Use Colipase as a co-emulsifier in diagnostic tests for the determination of lipase activity in combination with chromogenic Lipase Substrate, Catalog No. 11 034 618 103.

Properties

Molecular weight: Approximately 10 kD

Isoelectric point: 5.0

pH optimum: 8.8 (see figure)

Temperature dependence: See figure

pH stability: 3.5-11.5 (see figure)

Thermal stability: +25 to +80°C (see figure)

Specificity: Pancreatic colipase consists of 3 forms, colipase101 (procolipase), colipase96 and colipase85 (numbers stand for amino acid residues). colipase96 and colipase85 are trypsin digestion products of colipase101.

Specification

Appearance: White lyophilizate, TEA buffer, pH 6.0

Solubility: Clear, colorless solution in water (c=1 mg/mL)

pH value (c=1 mg/mL in water): 5.0-7.0

Activity (+25°C, tributyrin): ≥70,000 U/mg lyophilizate

Colipase: ≥0.6 mg/mg lyophilizate

Protein (Lowry): 0.8±0.2 mg/mg lyophilizate

Contaminants:

Lipase: ≤0.0005 U/mg lyophilizate

Proteases: ≤180 U/mg lyophilizate

SVD free: Corresponds to specification

pH 5.5 treatment (30 minutes): Corresponds to specification

Stability: At +2 to +8°C within specification range for 24 months.

Store dry.

Catalog number


10 204 307 103

Pack size


custom fill

Will be supplied as "Colipase from Porcine Pancreas". Unit of measure is "g active ingredient".


For further processing only.




Temperature dependence



pH optimum



pH stability



Thermal stability

Creatinase

from microorganism, lyophilizate

Hydrolase for creatinine determination that catalyzes the conversion of creatine to sarcosine and urea.

Application

Use Creatinase in your diagnostic tests for the determination of creatinine in combination with Creatininase, Catalog No. 11 865 471 103 and Sarcosine Oxidase, Catalog No. 11 378 856 103.

EC 3.5.3.3

Specification

- Appearance:** White to slightly yellowish lyophilizate
- Solubility:** Clear, colorless solution in water (c=10 mg/mL)
- pH value** (c=10 mg/mL in water): 5.5-6.5
- Activity** (+25°C, creatine, POD/PAP method): ≥4 U/mg lyophilizate
- Specific activity:** ≥9 U/mg protein
- Protein** (Biuret): 0.3-0.5 mg/mg lyophilizate
- Contaminants** (expressed as percentage of Creatinase activity):
 Creatininase: ≤0.01
 Catalase: ≤2
 Creatinine deaminase: ≤0.01
 Proteases (casein/resorufin, 2 hours stress duration): ≤0.001
 Contaminating oxidases (FOX): ≤0.001
- Stability:** At -15 to -25°C within specification range for 12 months.

Catalog number	Pack size
11 799 142 103	custom fill

Will be supplied as "Creatinase, Microbial Lyophil. Substance". Unit of measure is "MU".



For further processing only.

Creatininase

from *Pseudomonas* species, expressed in *E.coli*, lyophilizate

Hydrolase for creatinine determination that catalyzes the conversion of creatinine to creatine.

Application

Use Creatininase in diagnostic tests for the determination of creatinine in combination with Creatinase, Catalog No. 11 799 142 103 and Sarcosine Oxidase, Catalog No. 11 378 856 103.

EC 3.5.2.10

Properties

Nomenclature: Creatinine amidohydrolase

Molecular weight (gel filtration): 175 kD

Structure (SDS PAGE): 8 equal subunits (23 kD + zinc)

Isoelectric point (IEF): 4.7

Michaelis constants (Glycylglycine buffer, pH 8.0, +25°C):

Creatinine: 3×10^{-2} mol/L

Creatine: 6×10^{-2} mol/L

Inhibitors: Hg^{2+} , Fe^{3+} , Cu^{2+} (1 mmol/L), N-bromosuccinimide, o-phenanthroline, 4-chloromercuribenzoate. The enzyme is sensitive against photooxidation. A stable, inactive apoenzyme free of zinc can be obtained after EDTA incubation which can be reactivated completely with Zn^{2+} , Mn^{2+} , Mg^{2+} , Co^{2+} , Fe^{2+} or Ni^{2+} (1 mmol/L).

Activators: Mn^{2+} , Mg^{2+} . The enzyme requires metal ions.

Phenylmethylsulfonyl fluoride and iodoacetamide do not react.

pH optimum: 7.8 (see figure)

Temperature dependence: See figure

pH stability: 7.5-9.0 (see figure)

Thermal stability: Up to +65°C (see figure)

Specificity: Creatininase is specific for creatinine. It also reacts with glycoyamidine and glycoyamine. It does not react with hydantoin and its derivatives.

Catalog number


11 865 471 103

Pack size


custom fill

Will be supplied as "Creatininase, Recombinant Lyo". Unit of measure is "MU".

For further processing only.



pH stability



Temperature dependence

Specification

Appearance: White to slightly yellowish lyophilizate

Solubility: Clear, colorless to slightly yellowish solution in water (c=10 mg/mL)

pH value (c=10 mg/mL in water): 7.0-8.0

Activity (+25°C, creatinine): ≥55 U/mg lyophilizate

Specific activity: ≥250 U/mg protein

Protein (Biuret): 0.10-0.35 mg/mg lyophilizate

Contaminants (expressed as percentage of Creatininase activity):

ATPase: ≤0.01

Catalase: ≤2.0

Contaminating oxidases (FOX): ≤0.001


Creatinine deaminase: ≤0.0015

Kinase test: ≤0.01


Proteases (casein/resorufin, 2 hours stress duration): ≤0.005

Stability: At +2 to +8°C within specification range for 12 months.

Store dry. Protect from light.



Thermal stability



pH optimum

Creatinine Deaminase

from *Corynebacterium lilium*, lyophilizate

Hydrolase for creatinine determination that catalyzes the conversion of creatinine to N-methylhydantoin and ammonia.

Application

Use Creatinine Deaminase in diagnostic tests for the determination of creatinine in combination with N-Carbamoylsarcosine Amidase, Catalog No. 11 248 847 103, N-Methylhydantoinase (ATP-hydrolysing), Catalog No. 11 288 555 103 and Sarcosine Oxidase, Catalog No. 11 378 856 103.

EC 3.5.4.21

Specification

Appearance: Beige lyophilizate

Solubility: Clear, yellowish solution in water (c=10 mg/mL)

pH value (c=10 mg/mL in water): 8.0-9.0

Activity (+25°C, creatinine, via N-methylhydantoin, UV): 45.0-90.0 U/mg lyophilizate

Catalog number

11 330 764 103

Pack size

custom fill

Will be supplied as "Creatinine Deiminase". Unit of measure is "MU".

For further processing only.

Activity (+25°C, creatinine, via NH₃, UV): 35.0-70.0 U/mg lyophilizate

Protein (BCA): 10-30 mg/100 mg lyophilizate

Contaminants (expressed as percentage of Creatinine Desaminase activity)

ATPase: ≤0.1

Creatinase: ≤0.013

Creatininase: ≤0.01

Catalase: ≤10.0

Urease: ≤0.007

Uricase: ≤0.01

NH₃: ≤0.01μg/U

Stability: At +2 to +8°C within specification range for 12 months.

Store dry. Protect from light.

Formate Dehydrogenase

from yeast, lyophilizate

Dehydrogenase that catalyzes the interconversion of formate to carbon dioxide.

Application

Use Formate Dehydrogenase in diagnostic tests for the determination of oxalate in combination with Oxalate Decarboxylase or for the determination of formic acid. Also used in cofactor recycling systems for NADH.

EC 1.2.1.2

Specification

Appearance: White lyophilizate

pH value (c=10 mg/mL in water): Approximately 7.5

Activity (+25°C, formate): ≥0.40 U/mg lyophilizate

Specific activity: ≥3.0 U/mg protein

Contaminants (expressed as percentage of Formate Dehydrogenase activity):

Alcohol dehydrogenase: ≤0.05

Lactate dehydrogenase: ≤0.05

Malate dehydrogenase: ≤0.1

Stability: At +2 to +8°C within specification range for 12 months.

Store dry.

Catalog number

10 204 226 103

Pack size

custom fill

Will be supplied as "Formate Dehydrogenase from Yeast". Unit of measure is "kU".

For further processing only.

Galactose 1-Dehydrogenase

from *E.coli* overproducer, lyophilizate

Dehydrogenase that catalyzes the oxidation of galactose to D-galactono-1,4-lactone.

Application

Use Galactose 1-Dehydrogenase in diagnostic tests for the determination of total galactose.

EC 1.1.1.48

Specification

Appearance: White lyophilizate

Specific activity (+25°C, galactose): ≥50 U/mg protein

Protein (Biuret): ≥0.3-0.7 mg/mg lyophilizate

Contaminants (expressed as percentage of Galactose 1-Dehydrogenase activity):

Alcohol dehydrogenase: ≤0.01

β-Galactosidase: ≤0.01

Glutamate dehydrogenase (standard): ≤0.5

Lactate dehydrogenase: ≤0.1

"NADH-oxidase": ≤0.05

Stability: At +2 to +8°C within specification range for 12 months.

Catalog number

11 290 983 103

Pack size

custom fill

Will be supplied as "β-Galactose Dehydrogenase S". Unit of measure is "kU".

For further processing only.

Galactose 1-Dehydrogenase

from *E.coli* overproducer, suspension

Dehydrogenase that catalyzes the oxidation of galactose to D-galactono-1,4-lactone.

Application

Use Galactose 1-Dehydrogenase in diagnostic tests for the determination of total galactose.

EC 1.1.1.48

Specification

Appearance: White suspension in ammonium sulfate solution, 3.2 mol/L, pH approximately 6

Specific activity (+25°C, D-galactose): ≥100 U/mg protein

Protein (Biuret): ≥1 mg/mL

Contaminants (expressed as percentage of Galactose 1-Dehydrogenase activity):

Alcohol dehydrogenase: ≤0.01

β-Galactosidase: ≤0.01

Catalog number

10 633 313 103

Pack size

custom fill

Will be supplied as "β-Galactose Dehydrogenase S, *E. coli*". Unit of measure is "kU".

For further processing only.

Lactate dehydrogenase: ≤ 0.1

Malate dehydrogenase: ≤ 1.0

"NADH-oxidase": ≤ 0.05

Stability: At +2 to +8°C within specification range for 12 months.

Galactose 1-Dehydrogenase

from *Pseudomonas fluorescens*, suspension

Dehydrogenase that catalyzes the oxidation of galactose to D-galactono-1,4-lactone.

Application

Use Galactose 1-Dehydrogenase in diagnostic tests for the determination of total galactose.

EC 1.1.1.48

Specification

Appearance: White suspension in ammonium sulfate solution, 3.2 mol/L; EDTA, 1 mmol/L; pH approximately 6

Specific activity (+25°C, D-galactose): ≥ 5 U/mg protein

Protein (Biuret): 5 ± 0.5 mg/mL

Contaminants (expressed as percentage of Galactose 1-Dehydrogenase activity):

Alcohol dehydrogenase: ≤ 0.01

β -Galactosidase: ≤ 0.01

Lactate dehydrogenase: ≤ 0.5

"NADH-oxidase": ≤ 0.1

Stability: At +2 to +8°C within specification range for 12 months.

Catalog number

10 150 959 103

Pack size

custom fill

Will be supplied as "β-Gal-DH from *Pseudomonas fluorescens*". Unit of measure is "kU".

For further processing only.

Glucose Oxidase (GOD), Grade I

from *Aspergillus niger* overproducer, lyophilizate

Oxidoreductase that catalyzes the conversion of D-glucose to D-glucono-1,5-lactone which hydrolyzes spontaneously to gluconate.

Application

Use Glucose Oxidase (GOD), Grade I for the determination of α-amylase and D-glucose or O₂.

EC 1.1.3.4

Catalog number

12 158 566 103

Pack size

custom fill

Will be supplied as "GOD, RG I, rec., Lyo.". Unit of measure is "MU".

For further processing only.

Properties

Nomenclature: β -D-glucose: oxygen 1-oxidoreductase

Molecular weight: 79 kD

Isoelectric point: 4.3

Michaelis constants (Glucose):

Acetate buffer, pH 5.0, +25°C: 3.6×10^{-2} mol/L

Potassium phosphate buffer, 0.2 mol/l, pH 7.5, +25°C: 4.8×10^{-2} mol/L

Inhibitors: Ag^+ , Hg^{2+} , Cu^{2+} , 4-chloromercuribenzoate, D-arabinose (50%). FAD binding is inhibited by several nucleotides.

pH optimum: 7.0 (see figure)

Temperature dependence: See figure

pH stability: See figure

Thermal stability: See figure

Specificity: Glucose oxidase is specific for β -D-glucose. O_2 can be replaced by hydrogen acceptors such as 2,6-dichlorophenol indophenol.

Specification

Appearance: Yellowish lyophilizate

Conductivity (1%, w/v): ≤ 250 $\mu\text{S}/\text{cm}$

Activity (+25°C, glucose): ≥ 300 U/mg lyophilizate

Contaminants (expressed as percentage of Glucose Oxidase activity):


Amylase: ≤ 0.01

Catalase: ≤ 0.5


Saccharase: ≤ 0.01

Stability: At +2 to +8°C within specification range for 24 months.


Store dry.




pH optimum



pH stability



Temperature dependence



Thermal stability

Glucose Oxidase (GOD), Grade IIfrom *Aspergillus niger* overproducer, lyophilizate

Oxidoreductase that catalyzes the conversion of D-glucose to D-glucono-1,5-lactone which hydrolyzes spontaneously to gluconate.

Application

Use Glucose Oxidase (GOD), Grade II for the determination of α -amylase and D-glucose or O₂.

EC 1.1.3.4

Properties**Nomenclature:** β -D-glucose:oxygen 1-oxidoreductase**Molecular weight:** 79 kD**Isoelectric point:** 4.3**Michaelis constants** (Glucose):Acetate buffer, pH 5.0, +25°C: 3.6×10^{-2} mol/LPotassium phosphate buffer, 0.2 mol/L, pH 7.5, +25°C: 4.8×10^{-2} mol/L**Inhibitors:** Ag⁺, Hg²⁺, Cu²⁺, 4-chloromercuribenzoate, D-arabinose (50%). FAD binding is inhibited by several nucleotides.**pH optimum:** 7.0 (see figure for product Cat. No. 12 158 566 103)**Temperature dependence:** See figure for product Cat. No. 12 158 566 103**pH stability:** See figure for product Cat. No. 12 158 566 103**Thermal stability:** See figure for product Cat. No. 12 158 566 103**Specificity:** Glucose oxidase is specific for β -D-glucose. O₂ can be replaced by hydrogen acceptors such as 2,6-dichlorophenol indophenol.**Specification****Appearance:** Yellow brown lyophilizate**Solubility:** Clear, yellow solution in phosphate buffer, 0.1 mol/L, pH 7.0 (c=5 mg/mL)**pH value** (c=10 mg/mL in water): 6.8-7.8**Activity** (+25°C, glucose): ≥ 250 U/mg lyophilizate**Protein** (Pierce): No limit**Contaminants** (expressed as percentage of Glucose Oxidase activity):Amylase: ≤ 0.1 Catalase: ≤ 5 U/mg lyophilizateSaccharase: ≤ 0.1 **Stability:** At +2 to +8°C within specification range for 24 months.

Store dry.

Catalog number**11 939 998 103****Pack size**

custom fill

Will be supplied as "GOD, rec., Lyo.". Unit of measure is "MU".

For further processing only.

Glucose Oxidase (GOD), chemically modified

from *Aspergillus niger* overproducer, lyophilizate

Oxidoreductase that catalyzes the conversion of D-glucose to D-glucono-1,5-lactone which hydrolyzes spontaneously to gluconate.

Application

Use Glucose Oxidase (GOD), chemically modified for the determination of α -amylase and D-glucose or O₂. The chemical modification generates an increased liquid stability of the enzyme.

EC 1.1.3.4

Properties

Nomenclature: β -D-glucose: oxygen 1-oxidoreductase

Molecular weight: 79 kD

Isoelectric point: 4.3

Michaelis constants (Glucose):

Acetate buffer, pH 5.0, +25°C: 3.6×10^{-2} mol/L

Potassium phosphate buffer, 0.2 mol/L, pH 7.5, +25°C: 4.8×10^{-2} mol/L

Inhibitors: Ag⁺, Hg²⁺, Cu²⁺, 4-chloromercuribenzoate, D-arabinose (50%). FAD binding is inhibited by several nucleotides.

pH optimum: 7.0 (see figure for product Cat. No. 12 158 566 103)

Temperature dependence: See figure for product Cat. No. 12 158 566 103

pH stability: See figure for product Cat. No. 12 158 566 103

Thermal stability: See figure for product Cat. No. 12 158 566 103

Specificity: Glucose oxidase is specific for β -D-glucose. O₂ can be replaced by hydrogen acceptors such as 2,6-dichlorophenol indophenol.

Remark: The modified enzyme is especially suited for liquid stable applications with extended shelf life requirements.

Specification

Appearance: Yellowish white lyophilizate

pH value (c=40 mg/mL in water): 6.5-7.5

Activity (+25°C, glucose): ≥ 20 U/mg lyophilizate

Contaminants (expressed as percentage of Glucose Oxidase activity):

Catalase: ≤ 20 U/mg lyophilizate

Stability: At +2 to +8°C within specification range for 12 months.

Store dry.

Catalog number

11 485 938 103

Pack size

custom fill

Will be supplied as "GOD, Asp.niger, Bound to Dextran". Unit of measure is "MU".

For further processing only.

Glucose-6-phosphate Dehydrogenase (G6P-DH)

from *Leuconostoc mesenteroides*, expressed in *E. coli*, lyophilizate

Dehydrogenase, that catalyzes the oxidation of Glucose-6-phosphate.

Application

Use Glucose-6-phosphate Dehydrogenase for the determination of blood glucose or creatine kinase.

EC 1.1.1.49

Properties

Nomenclature: D-glucose-6-phosphate:NAD(P)⁺ 1-oxidoreductase

Molecular weight: 110 kD (1) (2 identical subunits 55,000 D)

Isoelectric point: pH 4.6

Michaelis constants (Tris: 0.1 mol/L; pH 7.8, +25°C):

NAD: 1.4×10^{-4} mmol/L

NADP: 3.7×10^{-5} mmol/L

Glucose-6-P: 3.7×10^{-4} mmol/L (NAD as coenzyme)

Glucose-6-P: 2.0×10^{-4} mmol/L (NADP as coenzyme)

Activators/inhibitors:

Phosphate, 5 mmol/L: 100% (NAD), 80% (NADP)

Phosphate, 50 mmol/L: 100% (NAD), 80% (NADP)

Without Mg²⁺: 90% (NAD), 80% (NADP)

Mg²⁺, 3 mmol/L: 100% (NAD), 100% (NADP)

Mg²⁺, 30 mmol/L: 100% (NAD), 100% (NADP)

HCO³⁻, 3 mmol/L: 100% (NAD), 100% (NADP)

Inhibitors: NADPH is a competitive inhibitor in the NAD-dependent reaction. Unlike the yeast enzyme, myristic acid, dehydroepiandrosterone and palmitoyl CoA do not inhibit.

pH optimum: 7.8 (see figure)

Temperature dependence: See figure

pH stability: 5.0-10.0 (see figure)

Thermal stability: Up to +40°C for native G6P-DH, up to +50°C for modified G6P-DH (see figure)

Buffer stability: Temperature stability can be significantly improved by the increase of ionic strength (see figure).

Stability of the lyophilizate: 100% residual activity after 3 weeks at +35°C

Specificity: G6P-DH is highly specific for glucose-6-phosphate and does not react with fructose-6-P, fructose-1,6-P₂ or glucose-1-P. 2-Deoxyglucose-6-P is slowly oxidized with NAD (5%) and with NADP (4%).

Catalog number


11 293 206 103

Pack size


custom fill

Will be supplied as "G6P-DH, rec., Lyo". Unit of measure is "MU".

For further processing only.




pH optimum



Incubation:
25°C, 180 min
● pH 3.0 – 5.0:
citrate buffer, 0.1 mol/l
■ pH 6.0 – 8.0:
phosphate buffer,
0.1 mol/l
▲ pH 9.0 – 11.0:
glycine buffer, 0.1 mol/l
500 U G6P-DH/ml

pH stability



Temperature dependence

Specification

Appearance: White or slightly yellowish lyophilizate

Solubility: Clear, colorless solution in water (c=10 mg/mL)

pH value (c=10 mg/mL, water): 6.5-7.5

Activity (+25°C, glucose-6-P, NAD): ≥600 U/mg lyophilizate

Specific activity (+25°C): ≥800 U/mg protein

Activity (+30°C): ≥750 U/mg lyophilizate

Activity (+37°C): ≥1,000 U/mg lyophilizate

Activity (+25°C, glucose-6-P, NADP): No limit

Protein (Biuret): 0.7-0.9 mg/mg lyophilizate

Contaminants (expressed as percentage of Glucose-6-phosphate dehydrogenase activity):

ATPase: ≤0.02

Creatine kinase: ≤0.001

Glutamate dehydrogenase: ≤0.01

Glutathione reductase: ≤0.001

Hexokinase and Glucose dehydrogenase: ≤0.05

Myokinase: ≤0.001

"NADH oxidase": ≤0.02

"NADPH oxidase": ≤0.0005

6-Phosphogluconate dehydrogenase (NAD as coenzyme): ≤0.001


6-Phosphogluconate dehydrogenase (NADP as coenzyme): ≤0.001

Phosphogluco isomerase : ≤0.01


Phosphoglucomutase: ≤0.001

Stability: At +2 to +8°C within specification range for 12 months.

Store dry.



Thermal stability



Buffer stability

Glucose-6-phosphate Dehydrogenase (G6P-DH), chemically modifiedfrom *Leuconostoc mesenteroides*, expressed in *E. coli*, lyophilizate

Recombinant dehydrogenase, that catalyzes the oxidation of Glucose-6-phosphate.

Application

Use Glucose-6-phosphate Dehydrogenase for the determination of blood glucose or creatine kinase. The chemical modification generates an increased liquid stability of the enzyme.

EC 1.1.1.49

Properties**Nomenclature:** D-glucose-6-phosphate:NAD(P)⁺ 1-oxidoreductase**Molecular weight:** 110 kD (2 identical subunits 55,000 D)**Isoelectric point:** pH 4.6**Michaelis constants** (Tris: 0.1 mol/L; pH 7.8, +25°C):NAD: 1.4×10^{-4} mmol/LNADP: 3.7×10^{-5} mmol/LGlucose-6-P: 3.7×10^{-4} mmol/L (NAD as coenzyme)Glucose-6-P: 2.0×10^{-4} mmol/L (NAD as coenzyme)**Activators/inhibitors:**

Phosphate, 5 mmol/L: 100% (NAD), 80% (NADP)

Phosphate, 50 mmol/L: 100% (NAD), 80% (NADP)

Without Mg²⁺: 90% (NAD), 80% (NADP)Mg²⁺, 3 mmol/L: 100% (NAD), 100% (NADP)Mg²⁺, 30 mmol/L: 100% (NAD), 100% (NADP)HCO³⁻, 3 mmol/L: 100% (NAD), 100% (NADP)

Inhibitors: NADPH is a competitive inhibitor in the NAD-dependent reaction. Unlike the yeast enzyme, myristic acid, dehydroepiandrosterone and palmitoyl CoA do not inhibit.

pH optimum: 7.8 (see figure for product Cat. No. 11 293 206 103)**Temperature dependence:** See figure for product Cat. No. 11 293 206 103**pH stability:** 5.0-10.0 (see figure for product Cat. No. 11 293 206 103)**Thermal stability:** Up to +40°C for native G6P-DH, up to +50°C for modified G6P-DH (see figure for product Cat. No. 11 293 206 103)**Buffer stability:** Temperature stability can be significantly improved by the increase of ionic strength (see figure for product Cat. No. 11 293 206 103).**Stability of the lyophilizate:** 100% residual activity after 3 weeks at +35°C

Specificity: G6P-DH is highly specific for glucose-6-phosphate and does not react with fructose-6-P, fructose-1,6-P₂ or glucose-1P. 2-Deoxyglucose-6-P is slowly oxidized with NAD (5%) and with NADP (4%).

Catalog number**11 389 343 103****Pack size**

custom fill

Will be supplied as "G6P-DH, rec., Lyo., mod.". Unit of measure is "MU".

For further processing only.

Remark: The modified enzyme is especially suited for liquid stable applications with extended shelf life requirements.

Specification

Appearance: White lyophilizate

Solubility: Clear, colorless solution in water (c=40 mg/mL)

pH value (c=40 mg/mL in water): 6.5-7.5

Activity (+25°C, glucose-6-P, NAD): ≥30 U/mg lyophilizate

Activity (+30°C): ≥39 U/mg lyophilizate

Activity (+37°C): ≥54 U/mg lyophilizate

Contaminants (expressed as percentage of Glucose-6-phosphate Dehydrogenase activity):

ATPase: ≤0.02

Creatine kinase: ≤0.001

Glutamate dehydrogenase: ≤0.01

Glutathione reductase: ≤0.001

Hexokinase and Glucose dehydrogenase: ≤0.05

Myokinase: ≤0.05

"NADH oxidase": ≤0.02

"NADPH oxidase": ≤0.0005

6-Phosphogluconate dehydrogenase (NAD as coenzyme): ≤0.001

6-Phosphogluconate dehydrogenase (NADP as coenzyme): ≤0.001

Phosphoglucose isomerase : ≤0.01

Phosphoglucomutase: ≤0.001

Glucose: ≤0.3 µg/mg lyophilizate

Stability: At +2 to +8°C within specification range for 18 months.

Store dry.

Glucose-6-phosphate Dehydrogenase (G6P-DH)

from *Leuconostoc mesenteroides*, expressed in *E. coli*, solution

Recombinant dehydrogenase, that catalyzes the oxidation of Glucose-6-phosphate.

Application

Use Glucose-6-phosphate Dehydrogenase for the determination of blood glucose or creatine kinase.

EC 1.1.1.49

Properties

Nomenclature: D-glucose-6-phosphate:NAD(P)⁺ 1-oxidoreductase

Molecular weight: 110 kD (1) (2 identical subunits 55,000 D)

Isoelectric point: pH 4.6

Michaelis constants (Tris: 0.1 mol/l; pH 7.8, +25°C):

Catalog number

11 650 742 103

Pack size

custom fill

Will be supplied as "G6P-DH, Recombinant (*E. coli*)". Unit of measure is "MU".

For further processing only.

NAD: 1.4×10^{-4} mmol/L

NADP: 3.7×10^{-5} mmol/L

Glucose-6-P: 3.7×10^{-4} mmol/L (NAD as coenzyme)

Glucose-6-P: 2.0×10^{-4} mmol/L (NAD as coenzyme)

Activators/inhibitors:

Phosphate, 5 mmol/L: 100% (NAD), 80% (NADP)

Phosphate, 50 mmol/L: 100% (NAD), 80% (NADP)

Without Mg^{2+} : 90% (NAD), 80% (NADP)

Mg^{2+} , 3 mmol/L: 100% (NAD), 100% (NADP)

Mg^{2+} , 30 mmol/L: 100% (NAD), 100% (NADP)

HCO_3^- , 3 mmol/L: 100% (NAD), 100% (NADP)

Inhibitors: NADPH is a competitive inhibitor in the NAD-dependent reaction. Unlike the yeast enzyme, myristic acid, dehydroepiandrosterone and palmitoyl CoA do not inhibit.

pH optimum: 7.8 (see figure for product Cat. No. 11 293 206 103)

Temperature dependence: See figure for product Cat. No. 11 293 206 103

pH stability: 5.0-10.0 (see figure for product Cat. No. 11 293 206 103)

Thermal stability: Up to +40°C for native G6P-DH, up to +50°C for modified G6P-DH (see figure for product Cat. No. 11 293 206 103)

Buffer stability: Temperature stability can be significantly improved by the increase of ionic strength (see figure for product Cat. No. 11 293 206 103).

Stability of the lyophilizate: 100% residual activity after 3 weeks at +35°C

Specificity: G6P-DH is highly specific for glucose-6-phosphate and does not react with fructose-6-P, fructose-1,6-P2 or glucose-1P. 2-Deoxyglucose-6-P is slowly oxidized with NAD (5%) and with NADP (4%).

Specification

Appearance: Clear, yellowish solution in glycerol

pH value: 6.0-7.0

Activity (+25°C, glucose-6-P): $\geq 2,500$ U/mL

Activity (+30°C): $\geq 3,000$ U/mL

Contaminants (expressed as percentage of Glucose-6-phosphate Dehydrogenase activity):

ATPase: ≤ 0.0200

Creatine kinase: ≤ 0.001

Glutamate dehydrogenase: ≤ 0.01

Glutathione reductase: ≤ 0.001

Hexokinase and Glucose dehydrogenase: ≤ 0.05

Myokinase: ≤ 0.01

"NADH oxidase": ≤ 0.02

"NADPH oxidase": ≤ 0.0005

Phosphoglucomutase: ≤ 0.001

6-Phosphogluconate dehydrogenase (NAD): ≤ 0.001

6-Phosphogluconate dehydrogenase (NADP): ≤ 0.001

Phosphoglucose isomerase : ≤ 0.01

Glycerol (enzymatic): 45-55% (v/v)

Stability: At +2 to +8°C within specification range for 12 months.

Glucose-6-phosphate Dehydrogenase (G6P-DH)

from *Leuconostoc mesenteroides*, lyophilizate

Dehydrogenase, that catalyzes the oxidation of Glucose-6-phosphate.

Application

Use Glucose-6-phosphate Dehydrogenase for the determination of blood glucose or creatine kinase.

EC 1.1.1.49

Specification

Appearance: White or slightly yellowish lyophilizate

Solubility: Clear, colorless solution in water (c=10mg/mL)

pH value (c=10mg/mL in water): 6.5-7.5

Activity (+25°C, glucose-6-P): ≥ 400 U/mg lyophilizate

Activity (+30°C): ≥ 520 U/mg lyophilizate

Activity (+37°C): ≥ 720 U/mg lyophilizate

Contaminants (expressed as percentage of Glucose-6-phosphate Dehydrogenase activity):

ATPase: ≤ 0.02

Creatine kinase: ≤ 0.001

Glutamate dehydrogenase: ≤ 0.01

Glutathione reductase: ≤ 0.001

Hexokinase and Glucose dehydrogenase: ≤ 0.05

Myokinase: ≤ 0.05

"NADH oxidase": ≤ 0.02

"NADPH oxidase": ≤ 0.0005

Phosphoglucomutase: ≤ 0.001

6-Phosphogluconate dehydrogenase (NAD): ≤ 0.001

6-Phosphogluconate dehydrogenase (NADP): ≤ 0.001

Phosphoglucose isomerase: ≤ 0.01

Stability: At +2 to +8°C within specification range for 12 months.

Store dry.

Catalog number

10 186 783 103

Pack size

custom fill

Will be supplied as "G6P-DH from *Leuconostoc mesenteroides*".
Unit of measure is "kU".

For further processing only.

Glucose-6-phosphate Dehydrogenase (G6P-DH)

from *Leuconostoc mesenteroides*, expressed in *E. coli*, reduced phosphate, lyophilizate

Recombinant dehydrogenase, that catalyzes the oxidation of Glucose-6-phosphate.

Application

Use Glucose-6-phosphate Dehydrogenase, reduced phosphate for the determination of inorganic phosphate in a colorimetric enzymatic reaction.

EC 1.1.1.49

Specification

Appearance: White to slightly yellowish lyophilizate

Solubility: Clear, colorless solution in water (c=10 mg/mL)

pH value (c=10mg/mL in water): 6.5-7.5

Activity (+25°C, glucose-6-P): ≥600 U/mg lyophilizate

Activity (+30°C): ≥750 U/mg lyophilizate

Activity (+37°C): ≥1,000 U/mg lyophilizate

Specific activity (+25°C, glucose-6-P): ≥800 U/mg

Protein (Biuret): 0.8-1.0 mg/mg lyophilizate

Contaminants (expressed as percentage of Glucose-6-phosphate Dehydrogenase activity):

ATPase: ≤0.00005

Creatine kinase: ≤0.001

Glutamate dehydrogenase: ≤0.01

Glutathione reductase: ≤0.001

α-Glucosidase: ≤0.00010

Hexokinase and glucose dehydrogenase: ≤0.05

Myokinase: ≤0.001

"NADH oxidase": ≤0.02

"NADPH oxidase": ≤0.0005

Phosphoglucomutase: ≤0.001

6-Phosphogluconate dehydrogenase (NAD): ≤0.001

6-Phosphogluconate dehydrogenase (NADP): ≤0.0001

Phosphoglucose isomerase: ≤0.01

Phosphate (as P_i): ≤10 µg/mg lyophilizate

Stability: At +2 to +8°C within specification range for 12 months.

Store dry.

Catalog number

11 650 734 103

Pack size

custom fill

Unit of measure is "MU".

For further processing only.

Glucose-6-phosphate Dehydrogenase (G6P-DH)

from *Leuconostoc mesenteroides*, suspension

Dehydrogenase, that catalyzes the oxidation of Glucose-6-phosphate.

Application

Use Glucose-6-phosphate Dehydrogenase for the determination of blood glucose or creatine kinase.

EC 1.1.1.49

Specification

Appearance: Yellowish suspension in ammonium sulfate

pH value: 5.5-6.5

Specific activity (+25°C, glucose-6-P): ≥550 U/mg

Specific activity (+30°C): ≥650 U/mg

Specific activity (+37°C): ≥800 U/mg

Protein (Biuret): ≥5 mg/mL

Contaminants (expressed as percentage of Glucose-6-phosphate Dehydrogenase activity):

ATPase: ≤0.02

Creatine kinase: ≤0.001

Glutamate dehydrogenase: ≤0.01

Glutathione reductase (NADH): ≤0.001

Hexokinase and Glucose dehydrogenase: ≤0.05

Myokinase: ≤0.05

"NADH oxidase": ≤0.02

"NADPH oxidase": ≤0.0005

Phosphoglucomutase: ≤0.001

6-Phosphogluconate dehydrogenase (NAD): ≤0.001

6-Phosphogluconate dehydrogenase (NADP): ≤0.001

Phosphoglucose isomerase: ≤0.01

Stability: At +2 to +8°C within specification range for 18 months.

Catalog number

10 128 171 103

Pack size

custom fill

Will be supplied as "G6P-DH from *Leuconostoc mesenteroides*".
Unit of measure is "MU".

For further processing only.

Glucose-6-phosphate Dehydrogenase (G6P-DH)

from yeast, lyophilizate

Dehydrogenase, that catalyzes the oxidation of Glucose-6-phosphate.

Application

Use Glucose-6-phosphate Dehydrogenase for the determination of blood glucose or creatine kinase.

EC 1.1.1.49

Catalog number

10 190 454 103

Pack size

custom fill

Will be supplied as "Glucose-6-phosphate Dehydrogenase, Yeast".
Unit of measure is "MU".

For further processing only.

Specification

Appearance: White to slightly yellowish lyophilizate
Solubility: Clear, colorless solution in water (c=10 mg/mL)
pH value (c=10 mg/mL): 6.0-7.0
Activity (+25°C, glucose-6-P): ≥15.0 U/mg lyophilizate
Contaminants (expressed as percentage of Glucose-6-phosphate Dehydrogenase activity):
 Creatine kinase: ≤0.001
 Glutathione reductase: ≤0.05
 Hexokinase: ≤0.02
 Phosphoglucomutase: ≤0.01
 6-Phosphogluconate dehydrogenase: ≤0.01
 Phosphoglucose isomerase: 0.002
Bioburden: ≤10,000 CFU/g
Stability: At +2 to +8°C within specification range for 12 months.
 Store dry.

Glucose-6-phosphate Isomerase

from yeast, suspension

Isomerase, that catalyzes the interconversion of glucose-6-phosphate into fructose 6-phosphate.

Application

Use Glucose-6-phosphate Isomerase for the isomerization of ketoses to aldoses and can be used for the determination of fructose.

EC 5.3.1.9

Specification

Appearance: White suspension in ammonium sulfate, 3.2 mol/L
Specific activity (+25°C, fructose-6-P): ≥350 U/mg protein
Protein (Biuret): 10±1 mg/mL
Contaminants (expressed as percentage of Glucose-6-phosphate Isomerase activity):
 Fructose-6-phosphate kinase: ≤0.01
 β-Fructosidase: ≤0.2
 Glutathione reductase: ≤0.01
 Phosphoglucomutase: ≤0.01
 6-Phosphogluconate dehydrogenase: ≤0.01
Stability: At +2 to +8°C within specification range for 24 months.

Catalog number

10 154 334 103

Pack size

custom fill

Will be supplied as "Phosphoglucose Isomerase (PGI) from Yeast".
 Unit of measure is "MU".

For further processing only.

α -Glucosidase

from yeast overproducer, multifunctional, lyophilizate

Recombinant glucosidase, that hydrolyzes 1,4-linked α -D-glucose residues with release of α -D-glucose.

Application

Use α -Glucosidase in diagnostic tests for the determination of α -amylase and pancreatic α -amylase activity according to the recommendations of the International Federation of Clinical Chemistry and Laboratory Medicine (IFCC). It catalyzes the 100% liberation of p-nitrophenol residues from the amylase substrate EPS (Catalog No. 10 880 078 103) once it has been cleaved by α -amylase.

EC 3.2.1.20

Specification

Appearance: White lyophilizate

Solubility: Clear, colorless solution in water (c=10 mg/mL)

pH value (c=10 mg/mL): 6.8-7.4

Activity (+37°C, 4-NP- α -D-glucoside): ≥ 60 U/mg lyophilizate

Specific activity: ≥ 130 U/mg protein

Protein (Biuret): 25-45 mg/100 mg lyophilizate

Contaminants (expressed as percentage of α -Glucosidase activity):

α -Amylase: ≤ 0.00001

Stability: At -15 to -25°C within specification range for 12 months.

Store dry.

Catalog number


11 626 329 103

Pack size

custom fill


Will be supplied as "α-Glucosidase Multifunctional". Unit of measure is "MU".

For further processing only.




37°C
 pH 4.0 – 6.0:
 acetate buffer,
 50 mmol/l
 pH 6.0 – 8.0:
 phosphate buffer,
 50 mmol/l
 pH 8.0 – 9.0:
 Tris-borate buffer,
 50 mmol/l

pH optimum




Incubation:
 25°C, 20 h
 pH 4.0 – 6.0:
 acetate buffer,
 50 mmol/l
 pH 6.0 – 8.0:
 phosphate buffer,
 50 mmol/l
 pH 8.0 – 9.0:
 Tris-borate buffer,
 50 mmol/l

pH stability




Incubation:
 15 min
 phosphate buffer,
 50 mmol/l; pH 7.0

Thermal stability




Incubation:
 15 min
 phosphate buffer,
 50 mmol/l
 pH 7.0

Temperature dependence



Stability in solution



Substrate specificity

β -Glucuronidase

from *E.coli*, solution

Hydrolase that cleaves β -linked terminal glucuronic acid.

Application

Use β -Glucuronidase in reagents for drug monitoring and doping analysis where it catalyzes the hydrolysis of steroid conjugates to detect various steroids in urine.

EC 3.2.1.31

Specification

Appearance: Clear, colorless solution, partially with small particles

Specific activity (+25°C, 4-NP-glucuronide): ≥ 80 U/mg protein

Specific activity (+37°C, 4-NP-glucuronide): ≥ 140 U/mg protein

Protein (Biuret): ≥ 0.5 mg/mL solution

Stability: At +2 to +8°C within specification range for 18 months.

Catalog number

03 708 446 103

Pack size

custom fill

Will be supplied as " β -Glucuronidase *E.coli* K12 Glycerol". Unit of measure is "MU".

For further processing only.

Glutamate Dehydrogenase (NAD(P))from *E.coli* overproducer, lyophilizate

Recombinant glutamate dehydrogenase.

Application

Use recombinant Glutamate Dehydrogenase in diagnostic tests for the determination of ammonia, urea, L-glutamate, glutamate pyruvate transaminase and leucine aminopeptidase.

EC 1.4.1.3

Properties**Nomenclature:** L-glutamate:NAD(P)⁺ oxidoreductase (deaminating)**Molecular weight:** ~2 200 kD for the associated enzyme with 8 subunits; 280 kD for one subunit.**Michaelis constants** (Tris buffer, pH 8.0, +23°C):L-glutamate: 1.8×10^{-3} mol/LNADP: 4.7×10^{-5} mol/L α -ketoglutarate: 7.0×10^{-4} mol/LNH₄⁺: 3.2×10^{-3} mol/LNADPH: 2.6×10^{-5} mol/L

K_m values for NAD or NADH are difficult to obtain due to their inhibitory action.

Inhibitors: 4-chloromercuribenzoate, Na₂S, diethyldithiocarbamate, 1,10-phenanthroline,8-hydroxyquinoline, NaN₃, thyroxine, heparin, sulfonylcarbamides, Cu²⁺, Hg²⁺, Ag²⁺, Fe³⁺, Zn²⁺, K⁺, PO₄²⁻, NO₃⁻**Activators:** Thioglycolic acid, b-mercaptoethylamine, EDTA, α , α' -dipyridyl**pH optimum:** 8.0 (see figure for product Cat. No. 10 190 462 103)**Temperature dependence:** See figure for product Cat. No. 10 190 462 103**pH stability:** 5.5-6.5 (see figure for product Cat. No. 10 190 462 103)**Thermal stability:** Up to +60°C (see figure for product Cat. No. 10 190 462 103)**Specificity:** The oxidation of L-glutamate is stimulated by ADP and inhibited by GTP. In contrast, the oxidation of alanine, leucine, isoleucine, methionine, valine, norleucine, norvaline and 2-aminobutyrate is stimulated by GTP and inhibited by ADP.**Catalog number****11 745 727 103****Pack size**

custom fill

Will be supplied as "GIDH, Lyo., rec.". Unit of measure is "MU".

For further processing only.

Specification

Appearance: White lyophilizate

Solubility: Clear, colorless solution in water (c=20 mg/mL)

pH value (c=20 mg/L in water): 6.5-7.5

Activity (+25°C, α-ketoglutarat): ≥80 U/mg lyophilizate

Contaminants (expressed as percentage of Glutamate Dehydrogenase activity):

Alcohol dehydrogenase: ≤0.005

Lactate dehydrogenase: ≤0.005

Malate dehydrogenase: ≤0.005

"NADH-Oxidase": ≤0.005

NH₄: ≤0.05 µg/mg lyophilizate

Stability: At +2 to +8°C within specification range for 12 months.

Store dry.

Glutamate Dehydrogenase (NAD(P))

from beef liver, lyophilizate

Dehydrogenase that catalyzes the conversion of glutamate to α-ketoglutarate.

Application

Use Glutamate Dehydrogenase in diagnostic tests for the determination of ammonia, urea, L-glutamate, glutamate pyruvate transaminase and leucine aminopeptidase.

EC 1.4.1.3

Properties

Nomenclature: L-glutamate:NAD(P)⁺ oxidoreductase (deaminating)

Molecular weight: ~2 200 kD for the associated enzyme with 8 subunits; 280 kD for one subunit.

Michaelis constants (Tris buffer, pH 8.0, +23°C):

L-glutamate: 1.8×10^{-3} mol/L

NADP: 4.7×10^{-5} mol/L

α-ketoglutarate: 7.0×10^{-4} mol/L

NH₄⁺: 3.2×10^{-3} mol/L

NADPH: 2.6×10^{-5} mol/L

K_m values for NAD or NADH are difficult to obtain due to their inhibitory action.

Inhibitors: 4-chloromercuribenzoate, Na₂S, diethyldithiocarbamate, 1,10-phenanthroline, 8-hydroxyquinoline, NaN₃, thyroxine, heparin, sulfonylcarbamides, Cu²⁺, Hg²⁺, Ag²⁺, Fe³⁺, Zn²⁺, K⁺, PO₄²⁻, NO₃⁻

Activators: Thioglycolic acid, b-mercaptoethylamine, EDTA, α, α'-dipyridyl

pH optimum: 8.0 (see figure)

Catalog number


10 190 462 103

Pack size


custom fill

Will be supplied as "Glutamate Dehydrogenase from Beef Liver". Unit of measure is "MU".

For further processing only.



pH optimum



pH stability

Incubation:
25°C, 180 min
pH 3.0 – 5.0:
citrate buffer, 0.1 mol/l
pH 6.0 – 8.0:
phosphate buffer,
0.1 mol/l
pH 9.0 – 11.0:
glycine buffer, 0.1
mol/l
120 U GIDH/ml

Temperature dependence: See figure

pH stability: 5.5-6.5 (see figure)

Thermal stability: Up to +60°C (see figure)

Specificity: The oxidation of L-glutamate is stimulated by ADP and inhibited by GTP. In contrast, the oxidation of alanine, leucine, isoleucine, methionine, valine, norleucine, norvaline and 2-aminobutyrate is stimulated by GTP and inhibited by ADP.

Remarks:

- Glutamate dehydrogenase suspension or solution can be dialyzed against phosphate buffer, 10 mmol/L. Glutamate dehydrogenase molecules have the tendency to associate in some test formulations, modified Glutamate dehydrogenase minimizes this effect.

Specification

Appearance: White lyophilizate

Solubility: Clear, colorless to slightly opalescent solution in water (c=20 mg/mL)

pH value (c=20 mg/mL in water): 6.5-7.5

Activity (+25°C, α-oxoglutarat): ≥10 U/mg lyophilizate

Contaminants (expressed as percentage of Glutamate Dehydrogenase activity):

Alcohol dehydrogenase: ≤0.005

Lactate dehydrogenase: ≤0.005


Malate dehydrogenase: ≤0.005

NH₄⁺: ≤0.1 µg/mg lyophilizate


pH 5.5 treatment (30 minutes): Corresponds to specification

Stability: At +2 to +8°C within specification range for 18 months.

Store dry.



Temperature dependence



Incubation:
10 min
(NH₄)₂ SO₄, 2.0 mol/l,
pH 7.0
2400 U GIDH/ml

Thermal stability

Glutamate Dehydrogenase (NAD(P))

from beef liver, chemically modified, lyophilizate

Dehydrogenase that catalyzes the conversion of glutamate to α -ketoglutarate.

Application

Use modified Glutamate Dehydrogenase in diagnostic tests for the determination of ammonia, urea, L-glutamate, glutamate pyruvate transaminase and leucine aminopeptidase.

EC 1.4.1.3

Properties

Nomenclature: L-glutamate:NAD(P)⁺ oxidoreductase (deaminating)

Molecular weight: ~2 200 kD for the associated enzyme with 8 subunits; 280 kD for one subunit.

Michaelis constants (Tris buffer, pH 8.0, +23°C):

L-glutamate: 1.8×10^{-3} mol/L

NADP: 4.7×10^{-5} mol/L

α -ketoglutarate: 7.0×10^{-4} mol/L

NH₄⁺: 3.2×10^{-3} mol/L

NADPH: 2.6×10^{-5} mol/L

K_m values for NAD or NADH are difficult to obtain due to their inhibitory action.

Inhibitors: 4-chloromercuribenzoate, Na₂S, diethyldithiocarbamate, 1,10-phenanthroline, 8-hydroxyquinoline, NaN₃, thyroxine, heparin, sulfonylcarbamides, Cu²⁺, Hg²⁺, Ag²⁺, Fe³⁺, Zn²⁺, K⁺, PO₄²⁻, NO₃⁻

Activators: Thioglycolic acid, b-mercaptoethylamine, EDTA, α , α' -dipyridyl

pH optimum: 8.0 (see figure for product Cat. No. 10 190 462 103)

Temperature dependence: See figure for product Cat. No. 10 190 462 103

pH stability: 5.5-6.5 (see figure for product Cat. No. 10 190 462 103)

Thermal stability: Up to +60°C (see figure for product Cat. No. 10 190 462 103)

Specificity: The oxidation of L-glutamate is stimulated by ADP and inhibited by GTP. In contrast, the oxidation of alanine, leucine, isoleucine, methionine, valine, norleucine, norvaline and 2-aminobutyrate is stimulated by GTP and inhibited by ADP.

Remarks:

- Glutamate dehydrogenase suspension or solution can be dialyzed against phosphate buffer, 10 mmol/L. Glutamate dehydrogenase molecules have the tendency to associate in some test formulations, modified Glutamate dehydrogenase minimizes this effect.
- The modified enzyme is especially suited for liquid stable applications with extended shelf life requirements.

Catalog number**11 434 993 103****Pack size**

custom fill

Will be supplied as "GIDH, Modified from Beef Liver". Unit of measure is "MU".

For further processing only.

Specification

Appearance: White lyophilizate, stabilized with RPLA 4 and ADP

Solubility: Clear, slightly opalescent solution in water (c=40 mg/mL)

Activity (+25°C, α-oxoglutarat): ≥7 U/mg lyophilizate

Contaminants (expressed as percentage of Glutamate Dehydrogenase activity):

Alcohol dehydrogenase: ≤0.005

Lactate dehydrogenase: ≤0.005

Malate dehydrogenase: ≤0.005

NH₄: ≤0.16 μmol/KU Glutamate Dehydrogenase

K (flame photometric): ≤0.1 μmol/KU Glutamate Dehydrogenase

Na (flame photometric): ≤2.0 μmol/KU Glutamate Dehydrogenase

pH 5.5 treatment (30 minutes): Corresponds to specification

Stability: At +2 to +8°C within specification range for 18 months.

Store dry.

Glutamate Dehydrogenase (NAD(P))

from beef liver, solution

Dehydrogenase that catalyzes the conversion of glutamate to α-ketoglutarate.

Application

Use Glutamate Dehydrogenase in diagnostic tests for the determination of ammonia, urea, L-glutamate, glutamate pyruvate transaminase and leucine aminopeptidase.

EC 1.4.1.3

Properties

Nomenclature: L-glutamate:NAD(P)⁺ oxidoreductase (deaminating)

Molecular weight: ~2 200 kD for the associated enzyme with 8 subunits; 280 kD for one subunit.

Michaelis constants (Tris buffer, pH 8.0, +23°C):

L-glutamate: 1.8×10^{-3} mol/L

NADP: 4.7×10^{-5} mol/L

α-ketoglutarate: 7.0×10^{-4} mol/L

NH₄⁺: 3.2×10^{-3} mol/L

NADPH: 2.6×10^{-5} mol/L

K_m values for NAD or NADH are difficult to obtain due to their inhibitory action.

Inhibitors: 4-chloromercuribenzoate, Na₂S, diethyldithiocarbamate, 1,10-phenanthroline,

8-hydroxyquinoline, NaN₃, thyroxine, heparin, sulfonylcarbamides, Cu²⁺, Hg²⁺, Ag²⁺, Fe³⁺, Zn²⁺, K⁺, PO₄²⁻, NO₃⁻

Catalog number

10 120 847 103

Pack size

custom fill

Will be supplied as "Glutamate Dehydrogenase, Beef Liver". Unit of measure is "L".

For further processing only.

Activators: Thioglycolic acid, b-mercaptoethylamine, EDTA, α , α' -dipyridyl

pH optimum: 8.0 (see figure for product Cat. No. 10 190 462 103)

Temperature dependence: See figure for product Cat. No. 10 190 462 103

pH stability: 5.5-6.5 (see figure for product Cat. No. 10 190 462 103)

Thermal stability: Up to +60°C (see figure for product Cat. No. 10 190 462 103)

Specificity: The oxidation of L-glutamate is stimulated by ADP and inhibited by GTP. In contrast, the oxidation of alanine, leucine, isoleucine, methionine, valine, norleucine, norvaline and 2-aminobutyrate is stimulated by GTP and inhibited by ADP.

Remarks: GIDH suspension or solution can be dialyzed against phosphate buffer, 10 mmol/L. GIDH molecules have the tendency to associate in some test formulations, modified GIDH minimizes this effect.

Specification

Appearance: Clear, colourless solution in glycerol

pH value: 7.0-7.8

Specific Activity: ≥ 120 U/mg

Protein (Biuret): 30 ± 3 mg/mL

Contaminants (expressed as percentage of Glutamate Dehydrogenase activity):

Alcohol dehydrogenase: ≤ 0.01

Lactate dehydrogenase: ≤ 0.01

Malate dehydrogenase: ≤ 0.01

NH₄: ≤ 0.16 $\mu\text{g}/\text{mg}$ protein

Glycerol: 560-680 mg/mL (45-55% (v/v))

EDTA (complexometric): 12.2-13.4 mmol/L

pH 5.5 treatment (30 minutes): Corresponds to specification

Stability: At +2 to +8°C within specification range for 18 months.

γ-Glutamyltransferase

from hog kidney, lyophilizate

Application

Use γ-Glutamyltransferase for designing your calibrator or control reagent.

EC 2.3.2.2

Specification

Appearance: White/off white to buff lyophilizate

Solubility: Clear, colorless to brownish solution in water (c=10 mg/mL)

Activity (+37°C, with γ-GT kit): >23 U/mg lyophilizate

Contaminants (expressed as percentage of γ-Glutamyltransferase activity):

Leucine aminopeptidase: <0.10

Phosphatase, alkaline: <2

"NADH oxidase": ≤0.01

SVD free: Corresponds to specification

pH 5.5 treatment (30 minutes): Corresponds to specification

Stability: At +2 to +8°C within specification range for 12 months.

Store dry.

Catalog number

10 445 363 103

Pack size

custom fill

Will be supplied as "γ-Glutamyltransferase from Hog Kidney". Unit of measure is "kU".

For further processing only.

Glycerol Kinase (GK)

from *Bacillus stearothermophilus*, lyophilizate

Enzyme that catalyzes the phosphorylation of glycerol to glycerol-3-phosphate.

Application

Use Glycerol Kinase for diagnostic tests for the determination of triglycerides together with Glycerol-3-phosphate Oxidase, Catalog Nos. 11 582 003 103 and 11 654 730 103, and Lipoprotein Lipase, Catalog No. 11 145 991 103.

EC 2.7.1.30

Properties

Nomenclature: ATP:glycerol 3-phosphotransferase

Molecular weight: 230 000 D (Sephadex G 200), 4 x 58 000 D (SDS-gel electrophoresis)

Michaelis constants (Glycine buffer, pH 9.8; +30°C):

Glycerol: 4.4×10^{-5} mol/L

Inhibitors: Unknown; Inhibitors of glycerokinase from *Candida mycoderma* do not inhibit the glycerokinase from *Bacillus stearothermophilus*.

Catalog number


11 499 530 103

Pack size

custom fill

Will be supplied as "GK, B.stearot., Lyo., w. Lactose". Unit of measure is "MU".

For further processing only.




pH optimum

pH optimum: 10.0-10.5 (see figure)
Temperature dependence: See figure
pH stability: 5.0-11.0 (see figure)
Thermal stability: Up to +60° C (see figure)


Specification

Appearance: White to slightly yellowish lyophilizate
Solubility: Clear, colorless solution in water (c=10 mg/mL)
pH value (c=10 mg/mL in water): 5.0-7.0
Activity (+25°C, glycerol): 18-25 U/mg lyophilizate
Specific activity: ≥80 U/mg protein
Protein (Biuret): 0.18-0.26 mg/mg lyophilizate
Contaminants (expressed as percentage of Glycerol Kinase activity):
 ATPase: ≤0.005
 Hexokinase: ≤0.01
 "NADH oxidase": ≤0.005
Glycerol (enzymatic): ≤40 µg/kU
Stability: At +2 to +8°C within specification range for 12 months.
 Store dry.




Incubation:
 25°C, 180 min
 pH 3.0 - 5.0:
 citrate buffer, 0.1 mol/l
 pH 6.0 - 8.0:
 phosphate buffer,
 0.1 mol/l
 pH 9.0 - 11.0:
 glycine buffer, 0.1
 mol/l
 85 U GK/ml

pH stability



Temperature dependence



Incubation:
 10 min
 Tris buffer, 0.1 mol/l;
 pH 9.0
 510 U GK/ml

Thermal stability

Glycerol Kinase (GK)

from *Bacillus stearothermophilus*, solution

Enzyme that catalyzes the phosphorylation of glycerol to glycerol-3-phosphate.

Application

Use Glycerol Kinase for diagnostic tests for the determination of triglycerides together with Glycerol-3-phosphate Oxidase, Catalog Nos. 11 582 003 103 and 11 654 730 103, and Lipoprotein Lipase, Catalog No. 11 145 991 103.

EC 2.7.1.30

Properties

Nomenclature: ATP:glycerol 3-phosphotransferase

Molecular weight: 230 000 D (Sephadex G 200), 4 x 58 000 D (SDS-gel electrophoresis)

Michaelis constants (Glycine buffer, pH 9.8; +30°C):

Glycerol: 4.4×10^{-5} mol/L

Inhibitors: Unknown; Inhibitors of glycerokinase from *Candida mycoderma* do not inhibit the glycerokinase from *Bacillus stearothermophilus*.

pH optimum: 10.0-10.5 (see figure for product Cat. No. 11 499 530 103)

Temperature dependence: See figure for product Cat. No. 11 499 530 103

pH stability: 5.0-11.0 (see figure for product Cat. No. 11 499 530 103)

Thermal stability: Up to +60° C (see figure for product Cat. No. 11 499 530 103)

Specification

Appearance: Clear, colorless to slightly yellowish solution in Tris buffer; pH approximately 7.3, stabilized, potential particles as result of recrystallized salts

Activity (+25°C, glycerol): ≥ 500 U/mL

Specific activity: ≥ 85 U/mg protein

Protein (Biuret): No limit

Contaminants (expressed as percentage of Glycerol Kinase activity):

Hexokinase: ≤ 0.01

"NADH oxidase": ≤ 0.005

Stability: At +2 to +8°C within specification range for 12 months.

Catalog number

10 691 666 103

Pack size

custom fill

Will be supplied as "Glycerokinase from Bac.stearothermophil.". Unit of measure is "kU".

For further processing only.

Glycerol Kinase (GK), concentratedfrom *Bacillus stearothermophilus*, solution

Enzyme that catalyzes the phosphorylation of glycerol to glycerol-3-phosphate.

Application

Use Glycerol Kinase for diagnostic tests for the determination of triglycerides together with Glycerol-3-phosphate Oxidase, Catalog Nos. 11 582 003 103 and 11 654 730 103, and Lipoprotein Lipase, Catalog No. 11 145 991 103.

EC 2.7.1.30

Properties**Nomenclature:** ATP:glycerol 3-phosphotransferase**Molecular weight:** 230 000 D (Sephadex G 200), 4 x 58 000 D (SDS-gel electrophoresis)**Michaelis constants** (Glycine buffer, pH 9.8; +30°C):Glycerol: 4.4×10^{-5} mol/L**Inhibitors:** Unknown; Inhibitors of glycerokinase from *Candida mycoderma* do not inhibit the glycerokinase from *Bacillus stearothermophilus*.**pH optimum:** 10.0-10.5 (see figure for product Cat. No. 11 499 530 103)**Temperature dependence:** See figure for product Cat. No. 11 499 530 103**pH stability:** 5.0-11.0 (see figure for product Cat. No. 11 499 530 103)**Thermal stability:** Up to +60° C (see figure for product Cat. No. 11 499 530 103)**Specification****Appearance:** Clear, colorless to slightly yellowish solution in Tris buffer; pH approximately 7.3, stabilized, potential particles as result of recrystallized salts**pH value** (c=10 mg/mL, in water): 7.1-7.5**Activity** (+25°C, glycerol): $\geq 2,200$ U/mL**Specific activity:** ≥ 85 U/mg protein**Protein** (Biuret): No limit**Contaminants** (expressed as percentage of Glycerol Kinase activity):Hexokinase: ≤ 0.01 "NADH oxidase": ≤ 0.005 **Stability:** At +2 to +8°C within specification range for 12 months.**Catalog number****10 539 937 103****Pack size**

custom fill

Will be supplied as "Glycerokinase, Bac. stearothermophilus". Unit of measure is "MU".

For further processing only.

Glycerol-3-phosphate Dehydrogenase

from rabbit muscle, suspension

Dehydrogenase that catalyzes the interconversion of dihydroxyacetone phosphate to glycerol 3-phosphate.

Application

Use Glycerol-3-phosphate Dehydrogenase in diagnostic reagents for the determination of aldolase in combination with Triose-phosphate Isomerase, Catalog No. 10 153 338 103.

EC 1.1.1.8

Specification

Appearance: White suspension in ammonium sulfate, 3.2 mol/L, pH approximately 6

pH value: 5.5-6.5

Activity: $\geq 2,000$ U/mL solution

Specific activity (+25°C, glyceraldehyde-3-phosphate): ≥ 170 U/mg protein

Protein (Biuret): 10 mg/mL

Ammonium sulfate: 3.2 ± 0.2 mol/L

Contaminants (expressed as percentage of Glycerol-3-phosphate Dehydrogenase activity):

Aldolase: ≤ 0.001

Glyceraldehyde-3-phosphate dehydrogenase: ≤ 0.001

Lactate dehydrogenase: ≤ 0.01

Triose-phosphate isomerase: ≤ 0.01

Stability: At +2 to +8°C within specification range for 18 months.

Catalog number

10 151 351 103

Pack size

custom fill

Will be supplied as "GDH from Rabbit Muscle". Unit of measure is "MU".

For further processing only.

Glycerol-3-phosphate Oxidase

from *E.coli* overproducer, lyophilizate

Recombinant oxidoreductase that catalyzes the interconversion of glycerol 3-phosphate to dihydroxyacetone phosphate.

Application

Use Glycerol-3-phosphate Oxidase in diagnostic tests for the determination of triglycerides together with Glycerol Kinase, Catalog Nos. 10 539 937 103 or 11 499 530 103 and Lipoprotein Lipase, Catalog No. 11 145 991 103.

EC 1.1.3.21

Catalog number

11 654 730 103

Pack size

custom fill

Will be supplied as "L-a-Glycerol-phosphate Oxidase, rec., Lyo.". Unit of measure is "MU".

Additional formulation: Lyophilizate, Catalog No. 11 582 003 103

For further processing only.

Properties

Nomenclature: Glycerol-3-phosphate : oxygen oxidoreductase

Molecular weight: 75 kD (SDS-PAGE); 74 kD (gel filtration, Sephadex G 150)

Isoelectric point: ~4.2 (230 000 D (Sephadex G 200), 4 x 58 000 D (SDS-gel electrophoresis))

Michaelis constants (L-glycerol phosphate):

K-phosphate buffer, 0.1 mol/L; pH 7.5: 1.36×10^{-2} mol/L (o-dianisidine assay)

Tris buffer, 0.1 mol/L; pH 7.6: 2.90×10^{-3} mol/L (o-dianisidine assay)

Tris buffer, 0.1 mol/L; pH 8.1: 1.40×10^{-3} mol/L (PAP assay)

Structure: Monomeric protein with FAD as cofactor

Inhibitors: Ag, Hg-salts and SDS

pH optimum: 8.0-8.5 (see figure)

Temperature dependence: See figure

pH stability: 6.5-8.5 (see figure)

Thermal stability: See figure

Specificity: Glycerol phosphate oxidase reacts highly specific with L- α -glycerol phosphate.

Specification

Appearance: Greenish yellow lyophilizate

Solubility: Clear yellow solution in Tris/HCl, 150 mmol/l, pH 7.6 (c=10 mg/ml)

Activity (+25°C, L- α -glycerol-3-phosphate): ≥ 50 U/mg lyophilizate

Activity (+37°C): ≥ 90 U/mg lyophilizate


Contaminants (expressed as percentage of Glycerol-3-phosphate Oxidase activity):

Cholesterol oxidase: ≤ 0.001


Lactate oxidase: ≤ 0.002

Uricase: ≤ 0.001


Stability: At +2 to +8°C within specification range for 12 months. Store dry.




pH optimum



pH stability



Temperature dependence



Thermal stability

Glycerol-3-phosphate Oxidase, chemically modified

from *E.coli* overproducer, lyophilizate

Recombinant oxidoreductase that catalyzes the interconversion of glycerol 3-phosphate to dihydroxyacetone phosphate.

Application

Use Glycerol-3-phosphate Oxidase in diagnostic tests for the determination of triglycerides together with Glycerol Kinase, Catalog Nos. 10 539 937 103 or 11 499 530 103 and Lipoprotein Lipase, Catalog No. 11 145 991 103.

EC 1.1.3.21

Properties

Nomenclature: Glycerol-3-phosphate : oxygen oxidoreductase

Molecular weight: 75 kD (SDS-PAGE); 74 kD (gel filtration, Sephadex G 150)

Isoelectric point: ~4.2 (230 000 D (Sephadex G 200), 4 x 58 000 D (SDS-gel electrophoresis))

Michaelis constants (L-glycerol phosphate):

K-phosphate buffer, 0.1 mol/L; pH 7.5: 1.36×10^{-2} mol/L (o-dianisidine assay)

Tris buffer, 0.1 mol/L; pH 7.6: 2.90×10^{-3} mol/L (o-dianisidine assay)

Tris buffer, 0.1 mol/L; pH 8.1: 1.40×10^{-3} mol/L (PAP assay)

Structure: Monomeric protein with FAD as cofactor

Inhibitors: Ag, Hg-salts and SDS

pH optimum: 8.0-8.5 (see figure for product Cat. No. 11 654 730 103)

Temperature dependence: See figure for product Cat. No. 11 654 730 103

pH stability: 6.5-8.5 (see figure for product Cat. No. 11 654 730 103)

Thermal stability: See figure for product Cat. No. 11 654 730 103

Catalog number

11 582 003 103

Pack size

custom fill

Will be supplied as "L-Glycerol-3-phosphate Oxidase rec. mod.". Unit of measure is "MU".

Additional formulation: Lyophilizate, Catalog No. 11 654 730 103

For further processing only.

Specificity: Glycerol phosphate oxidase reacts highly specific with L- α -glycerol phosphate.

Specification

Appearance: Green-yellow amorphous lyophilizate

Solubility: Clear yellow solution in water (c=10 mg/mL)

pH value (c=10 mg/mL in water): 6.8-7.8

Activity (+25°C, L- α -glycerol phosphate): ≥ 5 U/mg lyophilizate

Activity (+37°C, L- α -glycerol phosphate): ≥ 10 U/mg lyophilizate

Specific activity (+25°C): ≥ 40 U/mg protein

Protein (BCA): ≥ 0.1 mg/mg lyophilizate

Contaminants (expressed as percentage of Glycerol-3-phosphate Oxidase activity):

Cholesterol oxidase: ≤ 0.001

Lactate oxidase: ≤ 0.002

Uricase: ≤ 0.05

Stability: At +2 to +8°C within specification range for 12 months.

Store dry.

Hexokinase (HK)

from yeast overproducer, lyophilizate

Recombinant enzyme that converts hexose to hexose-6-phosphate.

Application

Use Hexokinase in diagnostic tests for blood glucose using the hexokinase method and for the determination of creatine kinase.

EC 2.7.1.1

Properties

Nomenclature: ATP:D-hexose 6-phosphotransferase

Molecular weight: 57 kD (SDS-PAGE)

Isoelectric point: 4.5-5.0

Michaelis constants (D-glucose):

Phosphate buffer, 0.1 mol/L, pH 7.0; +25°C: 3.05×10^{-4} mol/L

Phosphate buffer, 0.1 mol/L, pH 7.4; +30°C: 1.90×10^{-4} mol/L

Tea buffer, 0.1 mol/L, pH 7.6; +25°C: 2.30×10^{-4} mol/L

Michaelis constants (ATP):

Tris buffer, 0.1 mol/L, pH 7.6; +28°C: 1.60×10^{-4} mol/L

Tea buffer, 0.1 mol/L, pH 7.6; +25°C: 1.90×10^{-4} mol/L

Inhibitors: EDTA, SH-blocking compounds, sorbose-1-phosphate, polyphosphates, 6-deoxy-6-fluoroglucose, 2-C-hydroxymethylglucose, lyxose.

Activators: Mg^{2+} , catecholamines

pH optimum: 7.0-10.0 (see figure)

Catalog number


11 119 796 103

Pack size

custom fill

Will be supplied as "Hexokinase (HK) from Yeast Overproducer". Unit of measure is "MU".

For further processing only.



pH optimum

Temperature dependence: See figure

pH stability: 5.0-9.0 (see figure)

Thermal stability: 100% for 20 minutes at +37°C, 50% for 20 minutes at +45°C, 75% for 5 hours at +37°C (see figures)

Specificity: Hexokinase phosphorylates D-glucose, D-fructose, D-mannose, D-glucosamin, 2-deoxyglucose. L-Arabinose, D-xylose, L-rhamnose, D-galactose, D-lactose, sucrose, maltose, trehalose, raffinose, N-acetyl glucosamine do not react. ATP can be partially replaced by other nucleotides.

Specification

Appearance: Yellowish lyophilizate

Solubility: Clear, colorless solution in water (c=10 mg/mL)

pH value (c=10 mg/mL in water): 6.5-7.5

Activity (+25°C, glucose): ≥70 U/mg lyophilizate

Activity (+30°C): ≥98 U/mg lyophilizate

Activity (+37°C): ≥115 U/mg

Protein (Biuret): 0.15±0.05 mg/mg lyophilizate

Contaminants (expressed as percentage of Hexokinase activity):

Alcohol dehydrogenase: ≤0.001

ATPase: ≤0.05

Creatine kinase: ≤0.001

G6P-DH: ≤0.005

Glutamate dehydrogenase: ≤0.05

Glutathione reductase: ≤0.005

Myokinase: ≤0.001

"NADH oxidase": ≤0.001

"NADPH oxidase": ≤0.001


6-Phosphogluconate dehydrogenase: ≤0.001

Phosphoglucose isomerase: ≤0.002


Phosphoglucomutase: ≤0.02

Stability: At +2 to +8°C within specification range for 18 months.


Store dry.




pH stability



Temperature dependence



Thermal stability



Thermal stability

Hexokinase (HK), chemically modified

from yeast overproducer, lyophilizate

Recombinant enzyme that converts hexose to hexose-6-phosphate.

Application

Use Hexokinase in diagnostic tests for blood glucose using the hexokinase method and for the determination of creatine kinase. The chemical modification generates an increased liquid stability of the enzyme.

EC 2.7.1.1

Properties**Nomenclature:** ATP:D-hexose 6-phosphotransferase**Molecular weight:** 57 kD (SDS-PAGE)**Isoelectric point:** 4.5-5.0**Michaelis constants (D-glucose):**Phosphate buffer, 0.1 mol/L, pH 7.0; +25°C: 3.05×10^{-4} mol/LPhosphate buffer, 0.1 mol/L, pH 7.4; +30°C: 1.90×10^{-4} mol/LTea buffer, 0.1 mol/L, pH 7.6; +25°C: 2.30×10^{-4} mol/L**Michaelis constants (ATP):**Tris buffer, 0.1 mol/L, pH 7.6; +28°C: 1.60×10^{-4} mol/LTea buffer, 0.1 mol/L, pH 7.6; +25°C: 1.90×10^{-4} mol/L**Inhibitors:** EDTA, SH-blocking compounds, sorbose-1-phosphate, polyphosphates, 6-deoxy-6-fluoroglucose, 2-C-hydroxymethylglucose, lyxose.**Activators:** Mg²⁺, catecholamines**pH optimum:** 7.0-10.0 (see figure for product Cat. No. 11 119 796 103)**Temperature dependence:** See figure for product Cat. No. 11 119 796 103**pH stability:** 5.0-9.0 (see figure for product Cat. No. 11 119 796 103)**Thermal stability:** See figure for product Cat. No. 11 119 796 103**Specificity:** Hexokinase phosphorylates D-glucose, D-fructose, D-mannose, D-glucosamin, 2-deoxyglucose. L-Arabinose, D-xylose, L-rhamnose, D-galactose, D-lactose, sucrose, maltose, trehalose, raffinose, N-acetyl glucosamine do not react. ATP can be partially replaced by other nucleotides.**Catalog number****11 370 600 103****Pack size**

custom fill

Will be supplied as "Hexokinase (HK) from Rec.Yeast, Modif.". Unit of measure is "MU".

For further processing only.

Specification

Appearance: White lyophilizate

Solubility: Clear, colorless solution in water (c=40 mg/mL)

pH value (c=40 mg/mL in water): 6.5-7.5

Activity (+25°C, glucose): ≥40 U/mg lyophilizate

Contaminants (expressed as percentage of hexokinase activity):

Alcohol dehydrogenase: ≤0.001

ATPase: ≤0.05

Creatine kinase: ≤0.001

G6P-DH: ≤0.005

Glutamate dehydrogenase: ≤0.05

Glutathione reductase: ≤0.005

Myokinase: ≤0.001

"NADH oxidase": ≤0.001

6-Phosphogluconate dehydrogenase: ≤0.001

Phosphoglucose isomerase: ≤0.002

Phosphoglucomutase: ≤0.02

Glucose: ≤0.3 µg/mg lyophilizate

Stability: At +2 to +8°C within specification range for 18 months.

Remark: This enzyme is especially suited for liquid stable applications with extended shelf life requirements.

Hexokinase (HK)

from yeast overproducer, solution

Recombinant enzyme that converts hexose to hexose-6-phosphate.

Application

Use Hexokinase in diagnostic tests for blood glucose using the hexokinase method and for the determination of creatine kinase.

EC 2.7.1.1

Properties

Nomenclature: ATP:D-hexose 6-phosphotransferase

Molecular weight: 57 kD (SDS-PAGE)

Isoelectric point: 4.5-5.0

Michaelis constants (D-glucose):

Phosphate buffer, 0.1 mol/L, pH 7.0; +25°C: 3.05×10^{-4} mol/L

Phosphate buffer, 0.1 mol/L, pH 7.4; +30°C: 1.90×10^{-4} mol/L

Tea buffer, 0.1 mol/L, pH 7.6; +25°C: 2.30×10^{-4} mol/L

Michaelis constants (ATP):

Tris buffer, 0.1 mol/L, pH 7.6; +28°C: 1.60×10^{-4} mol/L

Tea buffer, 0.1 mol/L, pH 7.6; +25°C: 1.90×10^{-4} mol/L

Catalog number	Pack size
11 149 130 103	custom fill

Will be supplied as "Hexokinase (HK) from Recombinant Yeast".
Unit of measure is "MU".

For further processing only.

Inhibitors: EDTA, SH-blocking compounds, sorbose-1-phosphate, polyphosphates, 6-deoxy-6-fluoroglucose, 2-C-hydroxymethylglucose, lyxose.

Activators: Mg²⁺, catecholamines

pH optimum: 7.0-10.0 (see figure for product Cat. No. 11 119 796 103)

Temperature dependence: See figure for product Cat. No. 11 119 796 103

pH stability: 5.0-9.0 (see figure for product Cat. No. 11 119 796 103)

Thermal stability: 100% for 20 minutes at +37°C, 50% for 20 minutes at +45°C, 75% for 5 hours at +37°C (see figures)

Specificity: Hexokinase phosphorylates D-glucose, D-fructose, D-mannose, D-glucosamin, 2-deoxyglucose. L-Arabinose, D-xylose, L-rhamnose, D-galactose, D-lactose, sucrose, maltose, trehalose, raffinose, N-acetyl glucosamine do not react. ATP can be partially replaced by other nucleotides.

Specification

Appearance: Clear, yellowish solution, in 50% glycerol (v/v)

pH value: 6.0-7.0

Activity (+25°C, glucose): ≥1200 U/mL

Activity (+30°C): ≥1680 U/mL

Protein (Biuret): ≥75 mg/mL

Contaminants (expressed as percentage of Hexokinase activity):

Alcohol dehydrogenase: ≤0.001

ATPase: ≤0.05

Creatine kinase: ≤0.001

G6P-DH: ≤0.005

Glutamate dehydrogenase: ≤0.01

Glutathione reductase: ≤0.005

Myokinase: ≤0.001

6-Phosphogluconate dehydrogenase: ≤0.001

Phosphoglucose isomerase: ≤0.002

Phosphoglucomutase: ≤0.02

Glucose: ≤0.125 µg/mg lyophilizate

Glycerol (enzymatic): 560-680 mg/mL ≅ 45-55%

Stability: At +2 to +8°C within specification range for 18 months.

D-Lactate Dehydrogenase (D-LDH)

from microorganism, lyophilizate

Recombinant dehydrogenase that catalyzes the interconversion of D(-)-lactate to pyruvate.

Application

Use D-Lactate Dehydrogenase in a variety of diagnostic tests, *e.g.*, in the determination of alanine aminotransferases, lactate or pyruvate. Used for the removal of pyruvate in determinations working with NADH (*i.e.*, triglycerides, lipase, aldolase, aspartate aminotransferases, glutamate dehydrogenase).

EC 1.1.1.28

Specification

Appearance: White to slightly yellow powder or lyophilizate

Solubility: Soluble in water

Activity (+25°C, lyophilizate): ≥340 U/mg

Stability: At -15 to -25°C within specification range for 18 months.

Store dry in tightly sealed containers.

Catalog number

12 235 650 103

Pack size

custom fill

Will be supplied as "LDH". Unit of measure is "MU".



For further processing only.

D-Lactate Dehydrogenase (D-LDH), Grade I

from *Lactobacillus delbrückii*, lyophilizate

Dehydrogenase that catalyzes the interconversion of D(-)-lactate to pyruvate.

Application

Use D-Lactate Dehydrogenase (D-LDH), Grade I, in a variety of diagnostic tests, *e.g.*, in the determination of alanine aminotransferases, lactate or pyruvate. Used for the removal of pyruvate in determinations working with NADH (*i.e.*, triglycerides, lipase, aldolase, aspartate aminotransferases, glutamate dehydrogenase).

EC 1.1.1.28

Properties

Nomenclature: D-lactate:NAD⁺ oxidoreductase

Michaelis constants (Tris maleate buffer, pH 8.0, +25°C):

D-lactate: 0.7×10^{-1} mol/L (NAD, 2 mmol/L)

Pyruvate: 1.2×10^{-3} mol/L (NADH, 0.1 mmol/L)

NADH: 7.1×10^{-5} mol/L (pyruvate, 20 mmol/L)

pH optimum: 7.0 (see figure)

Temperature dependence: See figure

pH stability: 4.0-10.0 (see figure)

Catalog number


11 291 416 103

Pack size

custom fill

Will be supplied as "D(-)-LDH, Special Quality". Unit of measure is "MU".

For further processing only.



pH optimum

pH stability: 4.0-10.0 (see figure)

Thermal stability: Up to +50°C (see figure)

Specificity: Lactate dehydrogenase is specific for D(-)-lactate, L(+)-lactate does not react.

Remark: Lactate dehydrogenase, Grade I is especially suited for liquid stable applications with extended shelf life requirements.

Specification

Appearance: White to yellowish lyophilizate

Solubility: Clear to yellowish solution in water (c=10 mg/mL)

pH value (c=10 mg/mL in water): 6.0-7.0

Activity (+25°C, pyruvate): ≥180 U/mg lyophilizate

Specific activity: ≥450 U/mg protein

Protein (Biuret) : No limit, 0.3-0.8 mg/mg lyophilizate

Contaminants (expressed as percentage of D-Lactate Dehydrogenase activity):

Alcohol dehydrogenase: ≤0.01

Malate dehydrogenase: ≤0.1

"NADH oxidase": ≤0.0005

Succinate dehydrogenase: ≤0.01


NH₄: ≤0.01 μmol/KU

Na (flame photometric): ≤0.5 μmol/KU

K (flame photometric): ≤0.007 μmol/KU


Stability: At +2 to +8°C within specification range for 12 months.

Store dry.




Incubation:
25°C, 60 min
pH 3.0 - 5.0:
citrate buffer, 0.2 mol/l
pH 6.0 - 8.0:
phosphate buffer,
0.2 mol/l
pH 9.0 -11.0:
glycine buffer, 0.2 mol/l
2400 U D-LDH/ml

pH stability



Temperature dependence



Incubation:
10 min
water; pH 5.2
2400 U D-LDH/ml

Thermal stability

D-Lactate Dehydrogenase (D-LDH), Grade IIfrom *Lactobacillus delbrückii*, lyophilizate

Dehydrogenase that catalyzes the interconversion of D(-)-lactate to pyruvate.

Application

Use D-Lactate Dehydrogenase (D-LDH), Grade II, in a variety of diagnostic tests, *e.g.*, in the determination of alanine aminotransferases, lactate or pyruvate. Used for the removal of pyruvate in determinations working with NADH (*i.e.*, triglycerides, lipase, aldolase, aspartate aminotransferases, glutamate dehydrogenase).

EC 1.1.1.28

Properties**Nomenclature:** D-lactate:NAD⁺ oxidoreductase**Michaelis constants** (Tris maleate buffer, pH 8.0, +25°C):D-lactate: 0.7×10^{-1} mol/L (NAD, 2 mmol/L)Pyruvate: 1.2×10^{-3} mol/L (NADH, 0.1 mmol/L)NADH: 7.1×10^{-5} mol/L (pyruvate, 20 mmol/L)**pH optimum:** 7.0**pH stability:** 4.0-10.0**Thermal stability:** Up to +50°C**Specificity:** Lactate dehydrogenase is specific for D(-)-lactate, L(+)-lactate does not react.**Specification****Appearance:** White to yellowish lyophilizate**Solubility:** Clear, colorless to slightly yellowish solution in water (c=10 mg/mL)**pH value** (c=10 mg/mL in water): 5.7-6.7**Activity** (+25°C, pyruvate): ≥150 U/mg lyophilizate**Specific activity:** ≥300 U/mg protein**Protein** (Biuret): No limit, approximately 0.4-0.7 mg/mg lyophilizate**Contaminants** (expressed as percentage of D-Lactate Dehydrogenase activity):

Alcohol dehydrogenase: ≤0.01

Glucose dehydrogenase: ≤0.01

Malate dehydrogenase: ≤0.1

Succinate dehydrogenase: ≤0.01

Stability: At +2 to +8°C within specification range for 12 months.

Store dry.

Catalog number**10 679 666 103****Pack size**

custom fill

Will be supplied as "D(-)-Lactate Dehydrogenase (D-LDH)". Unit of measure is "MU".

For further processing only.

L-Lactate Dehydrogenase (L-LDH)

from pig muscle, for use of AST/GOT-Determination according to IFCC recommendations, lyophilizate

Dehydrogenase that catalyzes the interconversion of specific for L(+)-lactate to pyruvate.

Application

Use L-Lactate Dehydrogenase in a variety of diagnostic tests for the removal of pyruvate in determinations working with NADH (*i.e.*, triglycerides, lipase, aldolase, aminotransferases, glutamate dehydrogenase).

EC 1.1.1.27

Properties

Nomenclature: L-lactate:NAD⁺ oxidoreductase

Molecular weight: 140 kD

Isoelectric point: 4.6

Michaelis constants (Phosphate buffer, pH 7.5; +25°C):

Pyruvate: 1.5×10^{-4} mol/L (NADH: 0.18 mmol/L)

L-lactate: 3.3×10^{-3} mol/L (NAD: 0.5 mmol/L)

NADH: 1.1×10^{-5} mol/L (Pyruvate: 0.6 mmol/L)

NAD: 6.7×10^{-5} mol/L (L-lactate: 34 mmol/L)

Inhibitors: Oxamate, pyruvate (excess), oxalate, Ag⁺, Hg²⁺, Cu²⁺

pH optimum: 3.0-7.0 (see figure)

Temperature dependence: See figure

pH stability: 5.5-8.5 (see figure)

Thermal stability: Up to +40°C (see figure)

Specificity: Lactate dehydrogenase is specific for L(+)-lactate, D(-)-lactate does not react. Glyoxylate is also a lactate dehydrogenase substrate. Apart from pyruvate some 2-oxoacids are reduced. NAD analogs (*e.g.*, APAD) react at similar rates.

Specification

Appearance: White lyophilizate

Solubility: Clear colorless solution in water (c=10 mg/mL)

pH value (c=10 mg/mL in water): 6.0-7.0

Activity (+30°C, pyruvate; according to IFCC recommendations):

for Aspartate aminotransferase (AST/GOT) determination: ≥ 50 U/mg lyophilizate

for Alanine aminotransferase (ALT/GPT) determination: ≥ 50 U/mg lyophilizate

Activity (mean value of both determinations): ≥ 50 U/mg lyophilizate

Contaminants (expressed as percentage of Lactate Dehydrogenase activity, assayed according to the IFCC recommendations):

Aspartate aminotransferase (AST/GOT): ≤ 0.001

unspecificity of Lactate dehydrogenase: ≤ 0.005

Catalog number

10 254 754 103


Pack size

custom fill


Will be supplied as "LDH IFCC-quality from Hog Muscle". Unit of measure is "MU".

Additional formulation: Suspension in glycerol solution, Catalog No. 10 417 718 103

For further processing only.




pH optimum



Incubation:
25°C, 6 h resp. 24 h
● pH 3.0 – 5.5:
citrate buffer, 0.1 mol/l
■ pH 5.5 – 8.5:
phosphate buffer,
0.1 mol/l
▲ pH 8.5 – 11.0:
glycine buffer, 0.1 mol/l
10 U LDH/ml

pH stability



Temperature dependence


Alanine aminotransferase (ALT/GPT): ≤ 0.001

unspecificity of Lactate dehydrogenase: ≤ 0.005

SVD free: Corresponds to specification

pH 5.5 treatment (30 minutes): Corresponds to specification

Stability: At +2 to +8°C within specification range for 12 months.



Thermal stability

L-Lactate Dehydrogenase (L-LDH)

from pig muscle, for use of AST/GOT-Determination according to IFCC recommendations, solution

Dehydrogenase that catalyzes the interconversion of specific for L(+)-lactate to pyruvate.

Application

Use L-Lactate Dehydrogenase in a variety of diagnostic tests for the removal of pyruvate in determinations working with NADH (*i.e.*, triglycerides, lipase, aldolase, aminotransferases, glutamate dehydrogenase).

EC 1.1.1.27

Properties

Nomenclature: L-lactate:NAD⁺ oxidoreductase

Molecular weight: 140 kD

Isoelectric point: 4.6

Michaelis constants (Phosphate buffer, pH 7.5; +25°C):

Pyruvate: 1.5×10^{-4} mol/L (NADH: 0.18 mmol/L)

L-lactate: 3.3×10^{-3} mol/L (NAD: 0.5 mmol/L)

NADH: 1.1×10^{-5} mol/L (Pyruvate: 0.6 mmol/L)

NAD: 6.7×10^{-5} mol/L (L-lactate: 34 mmol/L)

Inhibitors: Oxamate, pyruvate (excess), oxalate, Ag⁺, Hg²⁺, Cu²⁺

pH optimum: 3.0-7.0 (see figure for product Cat. No. 10 254 754 103)

Temperature dependence: See figure for product Cat. No. 10 254 754 103

pH stability: 5.5-8.5 (see figure for product Cat. No. 10 254 754 103)

Thermal stability: Up to +40°C (see figure for product Cat. No. 10 254 754 103)

Specificity: Lactate dehydrogenase is specific for L(+)-lactate, D(-)-lactate does not react. Glyoxylate is also a lactate dehydrogenase

Catalog number

10 417 718 103

Pack size

custom fill

Will be supplied as "LDH, IFCC-quality from Hog Muscle". Unit of measure is "MU".

For further processing only.

substrate. Apart from pyruvate some 2-oxoacids are reduced. NAD analogs (*e.g.*, APAD) react at similar rates.

Specification

Appearance: Clear, colorless solution in glycerol, 50% (v/v), pH approximately 7

pH value: 6.5-7.5

Activity (+30°C, pyruvate, according to the IFCC recommendations): ≥9,600 U/mL solution

Specific activity: ≥480 U/mg protein

Protein (Biuret): ≥20 mg/mL solution

Glycerol: 560-680 mg/L $\hat{=}$ 45-55 % (v/v)

Contaminants (expressed as percentage of Lactate Dehydrogenase specific activity):

Glutamate dehydrogenase: ≤0.003

Aspartate aminotransferase (AST/GOT): ≤0.005

Alanine aminotransferase (ALT/GPT): ≤0.005

Reagent blank for determination of aspartate aminotransferase (AST/GOT): ≤0.9 mA/min

SVD free: Corresponds to specification

pH 5.5 treatment (30 minutes): Corresponds to specification

Stability: At +2 to +8°C within specification range for 12 months.

L-Lactate Dehydrogenase (L-LDH)

from pig muscle, suspension

Dehydrogenase that catalyzes the interconversion of specific for L(+)-lactate to pyruvate.

Application

Use L-Lactate Dehydrogenase in a variety of diagnostic tests for the removal of pyruvate in determinations working with NADH (*i.e.*, triglycerides, lipase, aldolase, aminotransferases, glutamate dehydrogenase).

EC 1.1.1.27

Specification

Appearance: White suspension in ammonium sulfate, 3.2 mol/l; Tris, 10 mmol/L, pH approximately 6.5

pH value: 6.0-7.0

Specific activity (+25°C, pyruvate): ≥550 U/mg protein

Protein (Biuret): ≥10 mg/mL

Ammonium sulfate: 3.2±0.2 mol/L

Contaminants (expressed as percentage of Lactate Dehydrogenase activity):

Catalog number

10 021 415 103

Pack size

custom fill

Will be supplied as "Lactate Dehydrogenase (LDH), Hog Muscle". Unit of measure is "MU".

For further processing only.

Aldolase: ≤ 0.001

Glutamate dehydrogenase: ≤ 0.01

Aspartate aminotransferase (AST/GOT): ≤ 0.005

Alanine aminotransferase (ALT/GPT): ≤ 0.005

Malate dehydrogenase: ≤ 0.01

Myokinase: ≤ 0.01

Pyruvate kinase: ≤ 0.001

SVD free: Corresponds to specification

pH 5.5 treatment (30 minutes): Corresponds to specification

Stability: At +2 to +8°C within specification range for 12 months.

L-Lactate Dehydrogenase (L-LDH), chemically modified

from pig heart, lyophilizate

Dehydrogenase that catalyzes the interconversion of L(+)-lactate to pyruvate.

Application

Use L-Lactate Dehydrogenase (L-LDH), chemically modified, in a variety of diagnostic tests for the removal of pyruvate in determinations working with NADH (*i.e.*, triglycerides, lipase, aldolase, aminotransferases, glutamate dehydrogenase). The chemical modification generates an increased liquid stability of the enzyme.

EC 1.1.1.27

Specification

Appearance: White lyophilizate

Solubility: Clear, colorless solution in water (c=10 mg/mL)

pH value: 7.1-8.1

Activity (+25°C, pyruvate): ≥ 25 U/mg lyophilizate

Specific activity: ≥ 150 U/mg protein

Protein (BCA): 0.15-0.25 mg/mg lyophilizate

Contaminants (expressed as percentage of Lactate Dehydrogenase activity):

Aspartate aminotransferase (AST/GOT): ≤ 0.005

Unspecificity of Lactate dehydrogenase: ≤ 0.05

Alanine aminotransferase (ALT/GPT): ≤ 0.01

Unspecificity of Lactate dehydrogenase: ≤ 0.05

"NADH-Oxidase": ≤ 0.001

SVD free: Corresponds to specification

pH 5.5 treatment (30 minutes): Corresponds to specification

Stability in CAPSO (pH 9.4, at +60°C for 1 hour): $\geq 50\%$

Stability: At +2 to +8°C within specification range for 12 months.

Catalog number

11 866 117 103

Pack size

custom fill

Will be supplied as "Lactate Dehydrogenase (LDH)". Unit of measure is "MU".

For further processing only.

L-Lactate Dehydrogenase (L-LDH)

from rabbit muscle, suspension

Dehydrogenase that catalyzes the interconversion of L(+)-lactate to pyruvate.

Application

Use L-Lactate Dehydrogenase in a variety of diagnostic tests for the removal of pyruvate in determinations working with NADH (*i.e.*, triglycerides, lipase, aldolase, aminotransferases, glutamate dehydrogenase).

EC 1.1.1.27

Specification

Appearance: Almost white suspension in ammonium sulfate, 3.2 mol/L, pH approximately 6.5

pH value: 6.0-7.0

Specific activity (+25°C, pyruvate): ≥550 U/mg protein

Protein (Biuret): ≥10 mg/mL

Ammonium sulfate: 3.2±0.2 mol/L

Contaminants (expressed as percentage of Lactate Dehydrogenase activity):

Aldolase: ≤0.001

Aspartate aminotransferase (AST/GOT): ≤0.01

Alanine aminotransferase (ALT/GPT): ≤0.01

Malate dehydrogenase: ≤0.01

Myokinase: ≤0.01

Pyruvate kinase: ≤0.001

Stability: At +2 to +8°C within specification range for 12 months.

Catalog number

10 003 557 103

Pack size

custom fill

Will be supplied as "Lactate Dehydrogenase, Rabbit Muscle". Unit of measure is "MU".

For further processing only.

Lactate 2-Monooxygenase (Lactate oxidase), Grade I

from *Aerococcus viridans*, expressed in *E. coli*, lyophilizate

Recombinant oxidoreductase that catalyzes the conversion of lactate to pyruvate.

Application

Use Lactate 2-Monooxygenase, Grade I in diagnostic tests for the determination of L-lactate.

EC 1.13.12.4

Catalog number

04 822 277 103

Pack size

custom fill

Will be supplied as "Lactat-OD, SQ, rec., lyo". Unit of measure is "MU".



For further processing only.

Properties

Nomenclature: L-lactate:oxigen oxidoreductase

Michaelis constant: L-lactate: 5×10^{-4} mol/L

V_{maximum} : L-lactate: 0.2 mol/L

pH optimum: 6.5-7.5 (see figure)

Temperature dependence: See figure

pH stability: 6.0-9.0 (see figure)

Thermal stability: Up to +65°C (see figure)

Specification

Appearance: Yellow lyophilizate

Activity: (+25°C, L-lactate): ≥ 40 U/mg lyophilizate

Specific activity: ≥ 55 U/mg protein

Contaminants (expressed as percentage of Lactate 2-Monooxygenase activity):

Catalase: ≤ 0.2


Glucose oxidase: ≤ 0.001

Pyruvate oxidase: ≤ 0.001

Uricase: ≤ 0.001


Stability: At -15 to -25°C within specification range for 12 months.

Store dry.




pH 4.0 – 5.5:
citrate, 0.1 mol/l
pH 5.5 – 7.5:
K-phosphate, 0.1 mol/l
pH 7.5 – 9.0:
Tris, 0.1 mol/l

pH optimum




Incubation:
25°C, 3 h
pH 4.0 – 5.5:
citrate, 0.1 mol/l
pH 5.5 – 7.5:
K-phosphate, 0.1 mol/l
pH 7.5 – 9.0:
Tris, 0.1 mol/l

pH stability



Temperature dependence



Thermal stability

Incubation:
15 min
K-phosphate buffer;
pH 7.0

Lactate 2-Monooxygenase (Lactate oxidase), Grade IIfrom *Aerococcus viridans*, expressed in *E. coli*, lyophilizate

Recombinant oxidoreductase that catalyzes the conversion of lactate to pyruvate.

Application

Use Lactate 2-Monooxygenase, Grade II in diagnostic tests for the determination of L-lactate.

EC 1.13.12.4

Properties**Nomenclature:** L-lactate:oxigen oxidoreductase**Michaelis constant:** L-lactate: 5×10^{-4} mol/L **V_{maximum} :** L-lactate: 0.2 mol/L**pH optimum:** 6.5-7.5 (see figure for product Cat. No. 04 822 277 103)**Temperature dependence:** See figure for product Cat. No. 04 822 277 103**pH stability:** 6.0-9.0 (see figure for product Cat. No. 04 822 277 103)**Thermal stability:** Up to +65°C (see figure for product Cat. No. 04 822 277 103)**Specification****Appearance:** Yellow lyophilizate**Activity** (25°C; with L-lactate) = ≥ 70 U/mg lyophilized material**Specific activity:** ≥ 55 U/mg protein**Protein** (BCA): 0.3-0.7 mg/mg lyophilizate**Contaminants** (expressed as percentage of Lactate 2-Monooxygenase activity):Catalase: ≤ 0.2 Glucose oxidase: ≤ 0.001 Pyruvate oxidase: ≤ 0.001 Uricase: ≤ 0.001 **Stability:** At -15 to -25°C within specification range for 12 months.

Store dry.

Catalog number**11 798 197 103****Pack size**

custom fill

Will be supplied as "Lactat-OD, rec., Lyo.". Unit of measure is "MU".



For further processing only.

Lactate 2-Monooxygenase (Lactate oxidase)

from *Pediococcus* species, lyophilizate

Oxidoreductase that catalyzes the conversion of lactate to pyruvate.

Application

Use Lactate 2-Monooxygenase in diagnostic tests for the determination of lactate.

EC 1.13.12.4

Specification

Appearance: Yellow lyophilizate

Activity (25°C, L-lactate): ≥20 U/mg lyophilizate

Specific activity: ≥55 U/mg protein

Protein (Lowry): 0.2-0.4 mg/mg lyophilizate

Stability: At -15 to -25°C within specification range for 18 months.

Store dry.

Catalog number

10 980 927 103

Pack size

custom fill

Will be supplied as "Lactate Oxidase from *Pediococcus* species".
Unit of measure is "kU".



For further processing only.

Lipoprotein Lipase

from *Pseudomonas* species, lyophilizate

Enzyme that hydrolyzes triglycerides into three free fatty acids and glycerol.

Application

Use Lipoprotein Lipase in diagnostic tests for the determination of triglycerides together with Glycerol Kinase, Catalog No. 10 539 937 103 or 11 499 530 103 and Glycerol-3-phosphate Dehydrogenase, Catalog No. 11 654 730 103 or 11 582 003 103.

EC 3.1.1.34

Properties

Nomenclature: Triacylglycerol-protein acylhydrolase

Molecular weight: 47 kD

Effectors: Hg²⁺, Ag⁺, Cr²⁺, Sn²⁺, Cu²⁺ and ionic detergents inhibit. Mg²⁺, sodium cholate and BSA stabilize the enzyme. 4-Chloromercuribenzoate (2 mmol/L), monoiodoacetate (2 mmol/L), NaF (20 mmol/L), NaN₃ (20 mmol/L), EDTA (5 mmol/L) and 2-phenanthroline (2 mmol/L) do not affect the enzyme activity while SDS (0.1% (w/v)) is inactivating.

pH optimum: 7.5 (see figure)

pH stability: 6.0-10.0 (see figure)

Thermal stability: Up to +50°C (see figure)

Specificity: Lipoprotein Lipase has both lipolytic and sterol ester hydrolytic activities. It hydrolyzes triacylglycerols in chylomicrons,

Catalog number

10 734 284 103

Pack size

custom fill

Will be supplied as "Lipoprotein Lipase from *Pseudomonas* spec.". Unit of measure is "MU".

For further processing only.

lipoproteins and diacylglycerols. With human plasma as substrate triglycerides are hydrolyzed more rapidly than cholesterol esters. The effects of pH and ionic strength on the enzymatic activity are somewhat different between the hydrolysis of triglyceride and of cholesterol ester depending on the different states of these substrates in the plasma or the transfer of the reaction products at the interface of substrates.

Lipolytic activity (Substrate, Number of C-atoms to number of double bonds, Relative rate):

olive oil:	94%
triolein (18:1):	100%
tripalmitin (16:0):	2%
trimyristin (14:0):	7%
trilaurin (12:0):	4%
tricaprin (10:0):	17%
tricaprylin (8:0):	64%
tricaproin (6:0):	2%
tributylin (4:0):	2%
tripropionin (3:0):	2%
triacetin (2:0):	1%

Specification

Appearance: Brownish lyophilizate

Solubility: Clear, brown solution in water (c=50 mg/mL)

Activity (+25°C, cholesterol oleate): ≥100 U/mg lyophilizate

Contaminants (expressed as percentage of Lipoprotein Lipase activity):

ATPase: ≤0.005

Catalase: ≤1.0

Glycerokinase: ≤0.001

Glucose oxidase: ≤0.001


Hexokinase: ≤0.005

"NADH oxidase": ≤0.001


Uricase: ≤0.005

Stability: At +2 to +8°C within specification range for 12 months.

Store dry.




pH optimum



Incubation:
25°C, 60 min
pH 6.0 – 8.0:
phosphate buffer,
0.7 mol/l
pH 9.0 – 11.0:
glycine buffer, 2.1 mol/l
50 U LPL/ml

pH stability



Incubation:
30 min
Tris buffer, 0.1 mol/l;
pH 7.7
50 U LPL/ml

Thermal stability

Lipoprotein Lipase, chemically modified

from *Pseudomonas* species, lyophilizate

Enzyme that hydrolyzes triglycerides into three free fatty acids and glycerol.

Application

Use Lipoprotein lipase in diagnostic tests for the determination of triglycerides together with Glycerol Kinase, Catalog Nos. 10 539 937 103 or 11 499 530 103 and Glycerol-3-phosphate Dehydrogenase, Catalog Nos. 11 654 730 103 or 11 582 003 103.

The chemical modification generates an increased liquid stability of the enzyme.

EC 3.1.1.34

Properties

Nomenclature: Triacylglycerol-protein acylhydrolase

Molecular weight: 47 kD

Effectors: Hg²⁺, Ag⁺, Cr²⁺, Sn²⁺, Cu²⁺ and ionic detergents inhibit. Mg²⁺, sodium cholate and BSA stabilize the enzyme. 4-Chloromercuribenzoate (2 mmol/L), monoiodoacetate (2 mmol/L), NaF (20 mmol/L), NaN₃ (20 mmol/L), EDTA (5 mmol/L) and 2-phenanthroline (2 mmol/L) do not affect the enzyme activity while SDS (0.1% (w/v)) is inactivating.

pH optimum: 7.5 (see figure)

pH stability: 6.0-10.0 (see figure)

Thermal stability: Up to +50°C (see figure)

Specificity: Lipoprotein Lipase has both lipolytic and sterol ester hydrolytic activities. It hydrolyzes triacylglycerols in chylomicrons, lipoproteins and diacylglycerols. With human plasma as substrate triglycerides are hydrolyzed more rapidly than cholesterol esters. The effects of pH and ionic strength on the enzymatic activity are somewhat different between the hydrolysis of triglyceride and of cholesterol ester depending on the different states of these substrates in the plasma or the transfer of the reaction products at the interface of substrates.

Lipolytic activity (Substrate, Number of C-atoms to number of double bonds, Relative rate):

- olive oil: 94%
- triolein (18:1): 100%
- tripalmitin (16:0): 2%
- trimyristin (14:0): 7%
- trilaurin (12:0): 4%
- tricaprin (10:0): 17%
- tricaprylin (8:0): 64%
- tricaproin (6:0): 2%
- tributylin (4:0): 2%
- tripropionin (3:0): 2%
- triacetin (2:0): 1%

Catalog number


11 145 991 103

Pack size


custom fill

Will be supplied as "Lipoprotein Lipase Modified". Unit of measure is "MU".

For further processing only.




pH optimum



pH stability

Incubation:
25°C, 60 min
pH 6.0 – 8.0:
phosphate buffer,
0.7 mol/l
pH 9.0 – 11.0:
glycine buffer, 2.1 mol/l
50 U LPL/ml



Thermal stability

Incubation:
30 min
Tris buffer, 0.1 mol/l;
pH 7.7
50 U LPL/ml

Remark: Chemically modified Lipoprotein Lipase (LPL) is more hydrophilic than native LPL. Carryover effect is therefore reduced.

Specification

Appearance: Brownish lyophilizate

Solubility: Clear, brown solution in water (c=50 mg/mL)

Activity (+25°C, cholesterol oleate): ≥10 U/mg lyophilizate

Contaminants (expressed as percentage of Lipoprotein Lipase activity):

ATPase: ≤0.005

Catalase: ≤1.0

Glycerokinase: ≤0.001

Hexokinase: ≤0.005

"NADH oxidase": ≤0.001

Uricase: ≤0.005

Stability: At +2 to +8°C within specification range for 12 months.

Store dry.

Lysozyme

from hen egg white, crystalline powder

Glucosidic bond hydrolyzing enzyme

Application

Use Lysozyme for bacteriolysis, preparation of protoplasts and sample preparation prior to isolation of nucleic acids.

Product description

Lysozyme from chicken egg. During purification and processing of this enzyme, steps included pH treatment at pH 3.5 or less for at minimum 30 minutes.

EC 3.2.1.17

Specification

Appearance: White, crystalline powder

Activity (+25°C, with *Micrococcus luteus*): ≥12,200 U/mg substance

Activity (+25°C, with *M. luteus*, previous Roche-substrate; calculated): ≥50,000 U/mg substance

Proteases (azocoll): ≤0.5 U/mg substance

Stability: At +2 to +8°C within specification range for 36 months.

Store dry.

Catalog number

10 153 516 001

10 153 516 103

Pack size

custom fill

custom fill

Will be supplied as "Lysozyme (Muramidase) from Hen Egg White". Unit of measure is "g".

For further processing only.

Malate Dehydrogenase

from pig heart, lyophilizate

Dehydrogenase that catalyzes the interconversion of malate to oxaloacetate.

Application

Use Malate Dehydrogenase in diagnostic tests for the determination of aspartate aminotransferase or in applications for citric and acetic acid testing.

EC 1.1.1.37

Specification

Appearance: White lyophilizate, stabilized

Solubility: Clear, colorless solution in water (c=10 mg/mL)

pH value: 7.0-8.0

Activity (+25°C, oxaloacetate): ≥70 U/mg lyophilizate

Contaminants (expressed as percentage of Malate Dehydrogenase activity):

Fumarase: ≤0.01

Aspartate aminotransferase (AST/GOT): ≤0.002

Lactate dehydrogenase: ≤0.01

SVD free: Corresponds to specification

pH 5.5 treatment (30 minutes): Corresponds to specification

Stability: At +2 to +8°C within specification range for 12 months.

Store dry.

Catalog number

10 200 387 103

Pack size

custom fill

Will be supplied as "Malate Dehydrogenase, Pig Heart (Mitochon.)". Unit of measure is "MU".

For further processing only.

Malate Dehydrogenase, chemically modified

from pig heart, lyophilizate

Dehydrogenase that catalyzes the interconversion of malate to oxaloacetate.

Application

Use Malate Dehydrogenase in diagnostic tests for the determination of aspartate aminotransferase or in applications for citric and acetic acid testing. The chemical modification generates an increased liquid stability of the enzyme.

EC 1.1.1.37

Catalog number

11 866 109 103

Pack size

custom fill

Will be supplied as "MDH, Lyo., mod.". Unit of measure is "MU".



For further processing only.

Specification

Appearance: White lyophilizate
pH value: 7.5-8.5
Activity (+25°C, oxaloacetate): ≥20 U/mg lyophilizate
Specific activity: ≥400 U/mg protein
Protein (BCA): ≥0.02 mg/mg lyophilizate
Contaminants (expressed as percentage of Malate Dehydrogenase activity):
 Aspartate aminotransferase (AST/GOT): ≤0.01
 "NADH-Oxidase": ≤0.005
SVD free: Corresponds to specification
pH 5.5 treatment (30 minutes): Corresponds to specification
Stability: At -15 to -25°C within specification range for 18 months.

Malate Dehydrogenase, IFCC Quality

from pig heart, lyophilizate

Dehydrogenase that catalyzes the interconversion of malate to oxaloacetate.

Application

Use Malate Dehydrogenase in diagnostic tests for the determination of aspartate aminotransferase or in applications for citric and acetic acid testing.

EC 1.1.1.37

Properties

Nomenclature: L-malate:NAD⁺ oxidoreductase
Molecular weight: 70 kD
Isoelectric point: 6.1-6.4
Michaelis constants (Phosphate buffer, 95 mmol/L, pH 8.3, +25°C):
 L-malate: 4.0 x 10⁻⁴ mol/L
 L-tartrate: 9.0 x 10⁻³ mol/L
 meso-tartrate: 1.2 x 10⁻³ mol/L
 oxaloacetate: 3.3 x 10⁻⁵ mol/L
Inhibitors: Iodinated compounds such as iodine cyanide, thyroxine and molecular iodine, phenols, 1,10-phenanthroline, 8-hydroxyquinoline, sulfide, nicotinic acidamide, adenine, AMP, ATP; oxaloacetate (excess).
Activators: Phosphate, arsenate, Zn²⁺
pH optimum: 7.5 (see figure)
Temperature dependence: See figure
pH stability: 7.0-9.0 (see figure)
Thermal stability: Up to +40°C (see figure)
Specificity: L-configuration of malate and tartrate. NAD can be replaced by its analogs, but not by NADP.

Catalog number


10 267 155 103

Pack size


custom fill

Will be supplied as "MDH IFCC-quality, Pig Heart (Mitochon)". Unit of measure is "MU".

For further processing only.



Temperature dependence



Thermal stability

Incubation:
 10 min
 (NH₄)₂ SO₄, 3.2 mol/l;
 pH 6.0
 12 000 U MDH/ml

Specification

Appearance: White lyophilizate

Solubility: Clear, colorless solution in water (c=10 mg/mL)

Activity (+30°C, oxaloacetate; according to the IFCC recommendations): ≥70 U/mg lyophilizate

Contaminants (expressed as percentage of Malat Dehydrogenase activity; assayed according to the IFCC recommendations):

Aspartate aminotransferase (AST/GOT): ≤0.001

Alanine aminotransferase (ALT/GPT): ≤0.001


Glutamate dehydrogenase: ≤0.005

SVD free: Corresponds to specification


pH 5.5 treatment (30 minutes): Corresponds to specification

Stability: At +2 to +8°C within specification range for 12 months.

Store dry.



pH optimum



Incubation:
25°C, 180 min
pH 3.0 – 5.0:
citrate buffer, 0.1 mol/l
pH 6.0 – 8.0:
phosphate buffer,
0.1 mol/l
pH 9.0 – 11.0:
glycine buffer, 0.1 mol/l
1200 U MDH/ml

pH stability

Malate Dehydrogenase, IFCC Quality

from pig heart, solution

Dehydrogenase that catalyzes the interconversion of malate to oxaloacetate.

Application

Use Malate Dehydrogenase in diagnostic tests for the determination of aspartate aminotransferase or in applications for citric and acetic acid testing.

EC 1.1.1.37

Properties

Nomenclature: L-malate:NAD⁺ oxidoreductase

Molecular weight: 70 kD

Isoelectric point: 6.1-6.4

Michaelis constants (Phosphate buffer, 95 mmol/L, pH 8.3, +25°C):

L-malate: 4.0 x 10⁻⁴ mol/L

L-tartrate: 9.0 x 10⁻³ mol/L

Catalog number

10 417 726 103

Pack size

custom fill

Will be supplied as "MDH, IFCC-quality, Pig Heart (Mitochon.)".
Unit of measure is "MU".

For further processing only.

meso-tartrate: 1.2×10^{-3} mol/L

oxaloacetate: 3.3×10^{-5} mol/L

Inhibitors: Iodinated compounds such as iodine cyanide, thyroxine and molecular iodine, phenols, 1,10-phenanthroline, 8-hydroxyquinoline, sulfide, nicotinic acidamide, adenine, AMP, ATP; oxaloacetate (excess).

Activators: Phosphate, arsenate, Zn^{2+}

pH optimum: 7.5 (see figure for product Cat. No. 10 267 155 103)

Temperature dependence: See figure for product Cat. No. 10 267 155 103

pH stability: 7.0-9.0 (see figure for product Cat. No. 10 267 155 103)

Thermal stability: Up to $+40^{\circ}\text{C}$ (see figure for product Cat. No. 10 267 155 103)

Specificity: L-configuration of malate and tartrate. NAD can be replaced by its analogs, but not by NADP.

Specification

Appearance: Clear, colorless solution in glycerol (50% (v/v))

pH value: 6.5-7.5

Specific activity ($+30^{\circ}\text{C}$, oxaloacetate): ≥ 900 U/mg protein

Protein (Biuret): ≥ 10 mg/mL

Contaminants (expressed as percentage of Malate Dehydrogenase activity):

Aspartate aminotransferase (AST/GOT): ≤ 0.005

Alanine aminotransferase (ALT/GPT): ≤ 0.005

Glutamate dehydrogenase: ≤ 0.003

Reagent blank for determination of aspartate aminotransferase

(AST/GOT): ≤ 0.009 ($\Delta A_{334}/10$ minutes)

SVD free: Corresponds to specification

pH 5.5 treatment (for at minimum 30 minutes): Corresponds to specification

Stability: At $+2$ to $+8^{\circ}\text{C}$ within specification range for 12 months.

N-Methylhydantoinase (ATP-hydrolyzing)

from *Arthrobacter* species, expressed in *E. coli*, lyophilizate

Hydrolase for creatinine determination that uses ATP to catalyze the conversion of N-methylhydantoin to N-carbamoylsarcosine.

Application

Use N-Methylhydantoinase (ATP-hydrolyzing) in diagnostic tests for the determination of creatinine in combination with Creatinine Deaminase Catalog No. 11 330 764 103, N-Carbamoylsarcosine Amidase, Catalog No. 11 248 847 103 and Sarcosine Oxidase, Catalog No. 11 378 856 103.

EC 3.5.2.14

Catalog number

11 288 555 103

Pack size

custom fill

Will be supplied as "N-Methylhydantoin Hydrolase". Unit of measure is "kU".



DRY ICE

For further processing only.

Specification

Appearance: White lyophilizate

Solubility: Clear, colorless solution in water (c=10 mg/mL, +25 °C)

pH value (c=100 mg/mL in water): 7.8-8.8

Activity (+25°C, N-methylhydantoin): 0.6-1.0 U/mg lyophilizate

Protein (Biuret): 20-43 mg/100 mg lyophilizate

Contaminants (expressed as percentage of N-Methylhydantoinase activity):

Creatinase: ≤0.013

Creatininase: ≤0.01

Catalase: ≤100

Uricase: ≤0.01

Stability: At -15 to -25°C within specification range for 12 months.

Store dry. Protect from light.

NAD(P)H Dehydrogenase (quinone) (Diaphorase)

from pig heart, suspension

Dehydrogenase that catalyzes the oxidation of dihydrolipoyl groups and has diaphorase activity.

Application

Use the diaphorase activity of NAD(P)H Dehydrogenase (quinone) for the determination of NAD(P)H and many dehydrogenases when coupled with various dyes which act as hydrogen acceptors from NAD(P)H, *e.g.* tetrazolium salts.

EC 1.6.5.2

Specification

Appearance: Yellow suspension in ammonium sulfate, 3.2 mol/L

pH value: 5.5-6.5

Specific activity (+25°C, lipoate): ≥25 U/mg protein

Protein (Biuret): 10±1 mg/mL

Ammonium sulfate: 3.2±0.2 mol/L

SVD free: Corresponds to specification

pH 5.5 treatment (30 minutes): Corresponds to specification

Stability: At +2 to +8°C within specification range for 12 months.

Catalog number

10 153 427 103

Pack size

custom fill

Will be supplied as "Diaphorase, Grade I from Pig Heart". Unit of measure is "g".

For further processing only.

Nitrate Reductase

from *Aspergillus* species, lyophilizate

Oxidoreductase that catalyzes the reduction of nitrate to nitrite.

Application

Use Nitrate Reductase in diagnostic tests for the determination of nitrate.

EC 1.7.1.2

Specification

Appearance: Yellow lyophilizate

Activity (+25°C, nitrate): ≥ 0.4 U/mg lyophilizate

Specific activity: ≥ 10 U/mg protein

Protein (Biuret): No limit

Contaminants (expressed as percentage of Nitrate Reductase activity):

Alcohol dehydrogenase (NADPH dependent): ≤ 0.8

"NADPH oxidase": ≤ 0.5

Nitrite reductase: ≤ 0.15

Stability: At -15 to -25°C within specification range for 12 months.

Store dry. Protect from light.

Catalog number

10 918 202 103

Pack size

custom fill

Will be supplied as "Nitrate Reductase (*Aspergillus* species)". Unit of measure is "kU".



For further processing only.

Oxalate Oxidase

from barley seedings, lyophilizate

Oxidoreductase that catalyzes the interconversion of oxalate to carbon dioxide and hydrogen peroxide.

Application

Use Oxalate Oxidase in diagnostic tests for the determination of oxalate.

EC 1.2.3.4

Specification

Appearance: White lyophilizate

Solubility: Clear, colorless to brownish solution in water (c=40 mg/mL)

Activity (+37°C, oxalate): ≥ 0.25 U/mg lyophilizate

Specific activity: ≥ 5 U/mg protein

Protein (Biuret): No limit

Stability: At -15 to -25°C within specification range for 12 months.

Store dry.

Catalog number

10 570 524 103

Pack size

custom fill

Will be supplied as "Oxalate Oxidase from Barley Seedlings". Unit of measure is "U".



For further processing only.

Peroxidase (POD), Grade II

from horse radish, lyophilizate

Application

Use Peroxidase (POD), Grade II, for the oxidation of reduced dyes in the indicator reaction of many diagnostic tests, *e.g.*, for the determination of blood glucose, triglycerides or lactate. It may also be used as a marker enzyme for enzyme immunoassays (EIA).

EC 1.11.1.7

Properties

Nomenclature: Donor:hydrogen-peroxide oxidoreductase

Molecular weight: 40 kD

Structure: Glycoprotein with one mole of protoheme IX

Isoelectric point: 7.2

Rate constants:

a) Hydrogen acceptors:

H₂O₂ 9 x 10⁸ [L x mol⁻¹ x s⁻¹]

methyl peroxide 1.5 x 10⁶ [L x mol⁻¹ x s⁻¹]

ethyl peroxide 3.6 x 10⁶ [L x mol⁻¹ x s⁻¹]

b) Hydrogen donors: Many

Inhibitors: Cyanide, sulfide, fluoride, azide, hydroxylamine, hydroxyl ions

Activators: Peroxidation of o-dianisidine is accelerated by ammonia, pyridine, imidazole at pH values >7.0

pH optimum: 6.0-6.5 (see figure)

Temperature dependence: See figure

pH stability: 4.0-10.0 (see figure)

Thermal stability: Up to +40°C (see figure)

Specificity: Peroxidase is specific for the hydrogen acceptor; only H₂O₂, methyl- and ethylperoxides are active. In contrast the enzyme is not specific for the hydrogen donor. A large number of phenols, aminophenols, diamines, indophenols, leucocyte dyes, ascorbate and several amino acids react.

Specification

Appearance: Red-brown lyophilizate

Solubility: Clear, red-brown solution in water (c=10 mg/mL)

pH value (c=10 mg/mL in water): 6.6-7.6

Activity: ≥200 U/mg lyophilizate

Purity number (A₄₀₃/A₂₇₅): 2.0-3.5

A₅₀₀ (100 U/mL): ≤0.120

Contaminants (expressed as percentage of Peroxidase activity):

ATPase: ≤0.001

Catalase: ≤0.7

Contaminating oxidases: ≤0.00005

Phosphatase, acidic: ≤0.001

Catalog number


11 378 783 103

Pack size


custom fill

Will be supplied as "Peroxidase (POD), Grade II, Horse-radish". Unit of measure is "MU".

For further processing only.




Temperature dependence



Incubation:
10 min
phosphate buffer,
0.1 mol/l; pH 7.0
2 000 U POD/ml


Thermal stability



pH optimum

Glucose: ≤0.25 µg/mg lyophilizate

Stability: At +2 to +8°C within specification range for 24 months.
Store dry.



pH stability

Phosphogluconate Dehydrogenase (decarboxylating)

from yeast, lyophilizate

Dehydrogenase that catalyzes the formation of ribulose 5-phosphate from 6-phosphogluconate.

Application

Use Phosphogluconate Dehydrogenase in diagnostic tests for the determination of creatine kinase or glucose in the combination with Hexokinase, Catalog Nos. 11 119 796 103, Glucose-6-phosphate Dehydrogenase, Catalog Nos. 10 186 783 103, 11 389 343 103, 11 293 206 103 or 10 190 454 103, and 6-Phosphogluconolactonase, Catalog No. 11 373 129 103.

EC 1.1.1.44

Properties

Nomenclature: 6-phospho D-gluconate:NADP⁺ 2-oxidoreductase (decarboxylating)

Molecular weight: 150 kD (native), 47 kD (SDS-PAGE)

Michaelis constants (TEA, 0.1 mol/L, pH 7.6, +25°C):

6-Phosphogluconate: 7.1×10^{-5} mol/L

NADP: 1.3×10^{-4} mol/L

Inhibitor constant (Phosphate buffer, pH 7.5):

Pyridoxal-5-P: 4.3×10^{-5} mol/L competitive

Inhibitors: Pyridoxal-5-P, iodoacetate and 4-hydroxymercuribenzoate

Activators: Chelators (EDTA, cysteine) plus metal ions (Mg²⁺); NaCl (0.2 mol/L), KCl (0.2 mol/L).

pH optimum: 7.8 (see figure)

Temperature dependence: See figure

pH stability: 5.0-8.0 (see figure)

Thermal stability: Up to +45°C (see figure)

Catalog number


11 126 482 103

Pack size


custom fill

Will be supplied as "6-PGDH from Yeast, Lyophilizate". Unit of measure is "kU".

For further processing only.



pH optimum



pH stability

Specificity: Phosphogluconate dehydrogenase is specific for NADP; NAD does not react.

Specification

Appearance: White lyophilizate

Solubility: Clear, colorless solution in water (c=10 mg/mL)

pH value (c=10 mg/mL in water): 7.0-8.0

Protein (Biuret): 0.08-0.16 mg/mg lyophilizate

Activity (+25°C, gluconate-6-P): ≥2 U/mg lyophilizate

Specific activity: ≥12 U/mg protein

Contaminants (expressed as percentage of 6-Phosphogluconate Dehydrogenase activity):

Creatine kinase: ≤0.006


G6P-DH: ≤0.01

Glutathione reductase: ≤0.01


Hexokinase: ≤0.01

Phosphoglucose isomerase: ≤0.03

Stability: At +2 to +8°C within specification range for 12 months.



Temperature dependence



Incubation:
10 min
Tea buffer, 0.1 mol/l;
pH 7.6
25 U 6-PGDH/ml

Thermal stability

6-Phosphogluconolactonase

from *Leuconostoc mesenteroides*, lyophilizate

Hydrolase that catalyzes the conversion of 6-phosphogluconolactone to 6-phosphogluconate.

Application

Use 6-Phosphogluconolactonase in diagnostic tests for the determination of creatine kinase or glucose in the combination with Hexokinase, Catalog Nos. 11 119 796 103, Glucose-6-phosphate Dehydrogenase, Catalog Nos. 10 186 783 103, 11 389 343 103, 11 293 206 103 or 10 190 454 103, and Phosphogluconate Dehydrogenase, Catalog No. 11 126 482 103.

EC 1.1.1.31

Properties

Nomenclature: 6-phosphogluconolactonase

Molecular weight: 38 kD (SDS)

Isoelectric point: 6.0

Catalog number

11 373 129 103

Pack size

custom fill

Will be supplied as "6-Phosphogluconolactonase". Unit of measure is "kU".



For further processing only.

Michaelis constants (MES buffer, pH 6.5; +25°C):

6-Phosphogluconalactone: $<1 \times 10^{-7}$ mol/L

Inhibitors: $(\text{NH}_4)_2\text{SO}_4$ (> 20 mmol/L), Mg^{2+} (>10 mmol/l), NaCl (>10 mmol/L). The enzyme is not inhibited by Cu^{2+} , Zn^{2+} , EDTA, 5,5'-dithiobis-2-nitrobenzoic acid, octanol (0.01%), Triton X-100 (1%) and Thesit (1%).

pH optimum: 6.0-7.5 (see figure)

Temperature dependence: See figure

pH stability: 7.0-9.0 (at +4°C, see figure)

Thermal stability: Up to +50°C (see figure)

Stability at different ionic strength: See figure

Specificity: 6-Phosphogluconolactone 100%, gluconolactone 0.5%

Specification

Appearance: White lyophilizate

Activity (+25°C, 6-phosphogluconolactone): ≥ 50 U/mg lyophilizate

Contaminants (expressed as percentage of

6-Phosphogluconolactonase activity):

Creatine kinase: ≤ 0.001

G6P-DH (with NAD): ≤ 0.02

Myokinase: ≤ 0.001


"NADPH oxidase": ≤ 0.001

6- Phosphogluconate dehydrogenase (with NAD): ≤ 0.01


Function testing (with G6P-DH, reaction time up to 5 minutes): $\geq 98\%$

Signal stability (δA_{340} , 30 minutes): ≤ 0.02

Stability: At -15 to -25°C within specification range for 12 months.




pH optimum




Incubation:
 ● 4°C, 24 h
 ■ 25°C, 24 h
 pH 3.0 – 5.0: citrate buffer, 0.05 mol/l
 pH 6.0 – 8.0: phosphate buffer, 0.05 mol/l
 pH 9.0 – 11.0: glycine buffer, 0.05mol/l
 1 U 6-PGL/ml


pH stability



Temperature dependence



Thermal stability



Stability at different ionic strength

Pyruvate Kinase

from *Bacillus stearotherophilus*, lyophilizate

Application

Use Pyruvate Kinase to catalyze the transfer of a phosphate group from phosphoenolpyruvate (PEP) to ADP, e.g., for the enzymatic determination of potassium or triglycerides.

EC 2.7.1.40

Specification

Appearance: White lyophilizate

Solubility: Clear, colorless solution in water (c=10 mg/mL)

pH value (+25°C, c=10 mg/mL): 7.6-8.7

Activity (+37°C, PEP): ≥120 U/mg lyophilized material

Specific activity: ≥150 U/mg protein

Protein (Biuret): 0.55-0.95 mg/mg lyophilizate

Contaminants (expressed as percentage of Pyruvate Kinase activity):

“NADH oxidase” (dA₃₆₅, 48 hours): ≤0.060

Na (flame photometric): ≤30 µmol/KU

K (flame photometric): ≤0.6 µmol/KU

NH₄ (enzymatic): ≤0.13 mg/KU

Mg (AAS): ≤50 µmol/KU

Mn (AAS): ≤2.4 µmol/KU

Stability: At +2 to +8°C within specification range for 12 months.

Store dry.

Catalog number

11 462 652 103

Pack size

custom fill

Will be supplied as “PK, Bacillus stearothermoph., Lyo.”. Unit of measure is “MU”.

For further processing only.

Contents

100 mg lyophilized material contains 60 mg protein and 40 mg Tris

Pyruvate Kinase

from rabbit muscle, suspension

Application

Use Pyruvate Kinase to catalyze the transfer of a phosphate group from phosphoenolpyruvate (PEP) to ADP, *e.g.*, for the enzymatic determination of potassium or triglycerides.

EC 2.7.1.40

Specification

Appearance: White suspension in ammonium sulfate, 3.2 mol/L, pH approximately 6

pH value: 5.5 - 6.5

Activity (+25°C, PEP): $\geq 2,000$ U/mL solution

Specific activity: ≥ 200 U/mg protein

Protein (Biuret): ≥ 10 mg/mL

Ammonium sulfate: 3.2 ± 0.2 mol/L

Contaminants (expressed as percentage of Pyruvate Kinase activity):

ATPase: ≤ 0.002

Enolase: ≤ 0.01

Glycerokinase: ≤ 0.001

Hexokinase: ≤ 0.002

Lactate dehydrogenase: ≤ 0.01

Myokinase: ≤ 0.01

"NADH oxidase": ≤ 0.002

Glycerol (enzymatic): ≤ 10 $\mu\text{g}/10$ mg

Stability: At +2 to +8°C within specification range for 24 months.

Catalog number

10 005 533 103

Pack size

custom fill

Will be supplied as "Pyruvate Kinase (PK) from Rabbit Muscle". Unit of measure is "MU".

For further processing only.

Pyruvate Oxidase

from *E.coli* overproducer, lyophilizate

Recombinant oxidoreductase that catalyzes the interconversion of pyruvate to acetyl phosphate.

Application

Use Pyruvate Oxidase in a variety of diagnostic tests, such as for the determination of pyruvate, lactate or aminotransferases.

EC 1.2.3.3

Catalog number

11 418 912 103

Pack size

custom fill

Will be supplied as "Pyruvate Oxidase Recombinant (*E. coli*)". Unit of measure is "kU".



For further processing only.

Specification

Appearance: Yellow lyophilizate

Solubility: Clear, yellowish solution in potassium phosphate buffer, 0.1 M, pH 6.5 (c=10 mg/mL)

Activity (+25°C, pyruvate, O₂, P_i): ≥1.5 U/mg lyophilizate

Specific activity: ≥3 U/mg protein

Protein (Biuret): ≥0.4 mg/mg lyophilizate

Contaminants (expressed as percentage of Pyruvate Oxidase activity):

ATPase: No limit

Glucose oxidase: ≤0.001

Aspartate aminotransferase (AST/GOT): ≤0.01

Alanine aminotransferase (ALT/GPT): ≤0.01

apo-Alanine aminotransferase (apo-ALT/apo-GPT) : ≤0.005

Impurities, total: ≤0.02

Lactate oxidase: ≤0.002

α-Ketoglutarate oxidase: ≤0.02

"NADH oxidase": ≤0.02

Stability: At -15 to -25°C within specification range for 12 months.

Store dry.

Sarcosine Oxidase

from *E.coli* overproducer, lyophilizate

Oxidoreductase that catalyzes the demethylation of sarcosine to glycine.

Application

Use Sarcosine Oxidase in diagnostic tests for the determination of creatinine. This can be done using one of two methods:

(1) In combination with Creatinase, Catalog No. 11 799 142 103 and Creatininase, Catalog No. 11 865 471 103.

(2) In combination with Creatinine Deaminase, Catalog No. 11 330 764 103, N-Carbamoylsarcosine Amidase, Catalog No. 11 248 847 103 and N-Methylhydantoinase (ATP-hydrolyzing), Catalog No. 11 288 555 103.

EC 1.5.3.1

Catalog number


11 378 856 103

Pack size

custom fill

Will be supplied as "Sarcosine Oxidase, Recombinant (E. coli)". Unit of measure is "kU".

For further processing only.



Temperature dependence

Properties

Nomenclature: Sarcosine:oxygen oxidoreductase (demethylating)

Molecular weight: 40 kD (PAGE, native Phast®-System)

Isoelectric point: 5.3 (Phast®-System)

Michaelis constants (Tris buffer, 0.1 mol/L, pH 8.0; Sarcosine):

at +25°C: 3.7×10^{-3} mol/L

at +37°C: 6.3×10^{-3} mol/L

Inhibitors: Completely inhibited by $ZnCl_2$ (7 mmol/L), $CdCl_2$ (7 mol/L), heavy metals and NaN_3 . Chloroacetic amine (0.2%) does not inhibit.

pH optimum: 8.0 (see figure)

Temperature dependence: See figure

pH stability: 7.0-10.0 (see figure)

Thermal stability: Up to +50°C (see figure)

Specificity: Sarcosine Oxidase reacts with sarcosine (100%), N-ethylglycine, 2 mmol/L (4%), L(-)-proline (0.28%), carbamoylsarcosine (0%), and glycine (0%).

Specification

Appearance: Yellow lyophilizate

Solubility: Clear, yellow solution in water (c=10 mg/mL)

pH value (c=10 mg/mL in water): 7.5-8.5

Activity (+25°C, sarcosine): 22-40 U/mg lyophilizate

Specific activity: ≥ 45 U/mg protein

Protein (Biuret): 0.4-0.6 mg/mg lyophilizate

Contaminants (expressed as percentage of Sarcosine Oxidase activity):

ATPase: ≤ 0.01

Catalase: ≤ 10.0

Contaminating oxidases (FOX): ≤ 0.005

Creatinase: ≤ 0.001

Creatininase: ≤ 0.01


Creatinine deaminase: ≤ 0.001

N-Carbamoylsarcosine amidohydrolase: ≤ 0.001


N-Methylhydantoin hydrolase: ≤ 0.001

Stability: At -15 to -25°C within specification range for 12 months.

Store dry. Protect from light.




pH optimum



Incubation:
25°C, 360 min
● pH 3.0 - 5.0: citrate buffer, 50 mmol/l
■ pH 6.0 - 8.0: phosphate buffer, 50 mmol/l
▲ pH 9.0 - 11.0: glycine buffer, 50 mmol/l
10 U sarcosine OD/ml

pH stability



Incubation:
10 min
phosphate buffer, 0.1 mol/l; pH 8.0
10 U sarcosine OD/ml

Thermal stability

Triose-phosphate Isomerase

from rabbit muscle, suspension

Isomerase that interconverts dihydroxyacetone phosphate and D-glyceraldehyde 3-phosphate.

Application

Use Triose-phosphate Isomerase in diagnostic reagents for the determination of aldolase in combination with Glycerol-3-phosphate Dehydrogenase, Catalog No. 10 151 351 103.

EC 5.3.1.1

Specification

Appearance: White suspension in ammonium sulfate

pH value: 5.5-6.5

Specific activity (+25°C, glyceraldehyde-3-phosphate): ≥5,000 U/mg protein

Protein (Biuret): 10±1 mg/mL

Contaminants (expressed as percentage of Triose-phosphate Isomerase activity):

Aldolase: ≤0.01

Glyceraldehyde-3-phosphate dehydrogenase: ≤0.001

Glycerol-phosphate dehydrogenase: ≤0.01

Stability: At +2 to +8°C within specification range for 24 months.

Catalog number

10 153 338 103

Pack size

custom fill

Will be supplied as "Triosephosphate Isomerase, Rabbit Muscle". Unit of measure is "MU".

For further processing only.

Thrombin

from human plasma, lyophilizate

Plasma derived coagulation factor II a that selectively cleaves the Arg-Gly bonds of fibrinogen to form fibrin.

Application

Use Thrombin to generate reference antigens for anti D-dimer antibodies.

EC 3.4.21.5

Specification

Appearance: White lyophilizate

Specific activity (+25°C, Chromozym TH): ≥120 U/mg protein

Protein (Lowry): Approximately 0.004 mg/mg lyophilizate

Factor Xa: ≤3%

Anti HIV: Negative

HBsAg: Negative

Stability: At +2 to +8°C within specification range for 24 months.

Catalog number

10 582 514 103

Pack size

custom fill

Will be supplied as "Thrombin (Coagulation Factor II a)". Unit of measure is "U".

For further processing only.

Urease

from jack bean, lyophilizate

Hydrolase that catalyzes the breakdown of urea in carbon dioxide and ammonia.

Application

Use Urease in diagnostic tests for the determination of urea in combination with Glutamate Dehydrogenase, Catalog No. 11 434 993 103.

EC 3.5.1.5

Properties

Nomenclature: Urea amidohydrolase

Molecular weight: 480 kD

Isoelectric point: 5.0-5.1

Michaelis constant (Phosphate buffer, pH 7.0; +25°C):

Urea: 1.05×10^{-2} mol/L

Inhibitors: Na^+ , K^+ , NH_4^+ ; suramin and thiourea are competitive inhibitors.

Activators: P_i

pH optimum: 7.5 (see figure)

Temperature dependence: See figure

pH stability: 6.0-9.5 (see figure)

Thermal stability: Up to +70° C (see figure)

Specificity: Urease is specific for urea.

Specification

Appearance: Almost white lyophilizate

Solubility: Clear, colorless solution in water (c=20 mg/mL)

pH value (c=20 mg/mL in water): 6.0-7.0

Activity (+25° C, urea): ≥ 45 U/mg lyophilizate

Specific activity: ≥ 600 U/mg protein

Protein (Biuret): ≤ 0.15 mg/mg lyophilizate

Contaminants (expressed as percentage of Urease activity):

L-Arginase: ≤ 0.002

NH_4 : ≤ 1.5 $\mu\text{g}/\text{KU}$

Stability: At +2 to +8°C within specification range for 12 months.

Store dry.

Catalog number


11 759 132 103

Pack size


custom fill

Will be supplied as "Urease, Lyo., SQ". Unit of measure is "MU".

For further processing only.




pH optimum




Incubation:
25°C, 24 h
K-phosphate buffer,
0.1 mol/l
urease,
approx. 45 U/ml

pH stability



Temperature dependence



Thermal stability

Uricase

from *Arthrobacter protophormiae*, lyophilizate

Oxidase that catalyzes the oxidation of uric acid to 5-hydroxyisourate which decomposes to allantoin under *in vitro* conditions.

Application

Use Uricase in diagnostic tests for the determination of uric acid and for the elimination of uric acid interferences.

EC 1.7.3.3

Properties

Nomenclature: Urate:oxygen oxidoreductase

Molecular weight: ~170 kD, with four subunits of ~40 kD

Michaelis constant (Phosphate buffer, 0.1 mol/L, pH 8.0; +25°C):

Urate: 6.6×10^{-5} mol/L

Stabilizer/activators: EDTA is good for stabilization. DTT or DTE may show a stabilizing effect depending on reagent composition. Triton X-100 (1-2 mL/L) may show an activating effect.

Inhibitors: Zn^{2+} , Cl^- (Tris-HCl buffer is not suitable) and borate inhibit strongly. NaN_3 , 0.1% does not inhibit.

pH optimum: 9.0 (see figure). Roche uric acid reagent contains phosphate buffer, pH 7.8. For these conditions, high activity and higher stability of the Uricase are achieved.

Temperature dependence: See figure, above +50°C there is a decrease in activity due to lower O_2 concentration.

pH stability: 6.5-10.0 (see figure)

Thermal stability: +20 to +60°C (see figure)

Specificity: Uricase is specific for urea.

Catalog number


10 828 475 103

Pack size


custom fill

Will be supplied as "Uricase from *Arthrobact. protophormiae*". Unit of measure is "MU".

For further processing only.



pH optimum



pH stability

Specification

Appearance: White lyophilizate

pH value (c=10 mg/mL in water): 6.7-7.5

Activity (+25°C, urate): ≥20 U/mg lyophilizate

Specific activity: ≥50 U/mg protein


Protein (Biuret): No limit

Contaminant (expressed as percentage of Uricase activity):


Catalase: ≤0.5

Stability: At +2 to +8°C within specification range for 12 months.

Store dry.



Temperature dependence



Incubation:
15 min
phosphate buffer,
0.1 mol/l;
pH 8.0
40 U uricase/ml

Thermal stability

3,5-Dichlorophenolsulfonic Acid

disodium salt

Color reagent for diagnostic tests

Application

Use 3,5-Dichlorophenolsulfonic Acid instead of phenol as a component in the trinder reaction.

CAS: 95041-38-6

Properties

Nomenclature: 3,5-Dichloro-2-hydroxy-benzosulphonic acid disodium salt

Formula: $C_6H_2O_4Cl_2Na_2$

Molecular weight: 287.0 D

Specification

Appearance: White powder

Solubility: Clear, colorless solution in water (c=20 mg/mL)

Dichlorophenolsulphonic acid-Na₂ (from C): ≥98.0%

C (elementary analysis): 24.6-25.6%

H (elementary analysis): 0.69-0.80%

Thin layer chromatography (TLC, silica gel 60 F₂₅₄, 1-butanol/acetic acid/water=50/15/25): Chromatographically homogeneous; corresponds to reference

Stability: At +15 to +25°C within specification range for 36 months.

Catalog number

10 667 536 103

Pack size

custom fill

Will be supplied as "3,5-Dichlorophenol Sulfonic Acid, Di-Na". Unit of measure is "g".

For further processing only.

4-Aminoantipyrine (4-APP)

crystalline powder

Substrate for peroxidase

Application

Use 4-Aminoantipyrine in a variety of diagnostic tests that use the Trinder reaction for the colorimetric determination of analytes, such as for the determination of cholesterol, glucose, creatinine or uric acid.

CAS: 83-07-8

Properties

Formula: $C_{11}H_{13}N_3O$

Molecular weight: 203.25 D

Catalog number

10 073 474 103

Pack size

custom fill

Will be supplied as "4-Aminoantipyrine". Unit of measure is "kg".

For further processing only.

Specification**Appearance:** Yellow to redish brown, crystalline powder**Solubility:** Clear, colorless solution in water (c=0.1%, w/v)**Melting range:** +106 to +109°C**A 370** (c=0.1%, w/v): ≤0.04**UV-spectrum** (c=0.002%, w/v): Corresponds to reference

Maximum: 242-246 nm

Minimum: 217-219 nm

Shoulder: 274 nm

Heavy metals (as Pb): ≤5 ppm $\hat{=}$ 0.0005%**IR-spectrum:** Corresponds to reference**4-Aminoantipyrine** (HClO₄-titration, based on undried substance):

≥98.0%

Purity (HPLC): ≥99.0 area%**Stability:** At +15 to +40°C within specification range for 36 months.

Store dry in tightly closed containers.

4-Nitrophenyl- α -D-maltohexaoside

powder

Nitrophenyl substrate

ApplicationUse 4-Nitrophenyl- α -D-maltohexaoside in diagnostic tests for the determination of α -amylase.**CAS:** 74173-30-1**Properties****Formula:** C₄₂H₆₅NO₃₃**Molecular weight:** 1112.1 D**Specification****Appearance:** White to slightly yellowish, amorphous powder**Solubility:** Clear, slightly yellowish solution in water (c=70 mg/mL)**4-Nitrophenyl-maltohexaoside** (enzymatically): ≥90%**4-Nitrophenyl-maltohexaoside** (HPLC): ≥96.0 area%**Water** (K. Fischer): ≤3.0%**4-Nitrophenyl-maltopentaoside** (HPLC): ≤1.0 area%**4-Nitrophenyl-maltoheptaoside** (HPLC): ≤2.0 area%**4-Nitrophenol, free:** ≤0.05%**2-Propanol** (GC): ≤6%**Stability:** At +2 to +8°C within specification range for 18 months.

Store dry.

Catalog number**10 691 682 103****Pack size**

custom fill

Will be supplied as "4-Nitrophenyl- α -D-malto-hexaoside". Unit of measure is "g".

For further processing only.

4-Nitrophenyl Phosphate (4-NPP)

di-Tris salt

Substrate for alkaline phosphatase

Application

Use 4-Nitrophenyl Phosphate in diagnostic test for the determination of alkaline phosphatase according to the recommendations of the International Federation of Clinical Chemistry and Laboratory Medicine (IFCC).

CAS: 91-81-6

Properties

Formula: $C_6H_4NO_6P(C_4H_{12}NO_3)_2$

Molecular weight: 461.3 D (4-NPP: 219.1 D)

Specification

Appearance: White to slightly yellow, crystalline powder

Solubility: Clear, colorless to slightly yellow solution in water

4-NPP, di-Tris (from content found enzymatically): $\geq 88\%$

4-NPP (enzymatic): $\geq 42\%$

Tris (titrimetric): $\geq 46\%$

Water (K. Fischer): $\leq 6\%$

4-NP free: $\leq 0.07\%$

Reaction rates (alkaline phosphatase): $100 \pm 5\%$

Stability: At +2 to +8°C within specification range for 24 months.

Store dry. Protect from light.

Catalog number

10 270 857 103

Pack size

custom fill

Will be supplied as "4-Nitrophenyl Phosphate, Di-Tris Salt". Unit of measure is "kg".

For further processing only.

Benzylidene-4-NP-G7

4,6-Benzylidene-4-nitrophenyl- α -D-maltoheptaoside, lyophilizate

Nitrophenyl substrate

Application

Use Benzylidene-4-NP-G7 in diagnostic tests for the determination of α -amylase.

CAS: 109055-07-4

Properties

Formula: $C_{55}H_{79}NO_{38}$

Molecular weight: 1362.1 D

Catalog number

11 378 872 103

Pack size

custom fill

Will be supplied as "Benzylidene-4-NP-G7". Unit of measure is "kg".



For further processing only.

Specification**Appearance:** White to slightly yellowish lyophilizate**Benzylidene-4-NP-G7** (enzymatic): $\geq 90\%$ **Water** (K. Fischer): $\leq 3\%$ **4-NP-maltoheptaoside** (HPLC): ≤ 1.0 area%**4-Nitrophenol, free:** $\leq 0.01\%$ **Reaction rates** (α -amylase):In Precinorm® U: $100 \pm 5\%$ In Precipath® U: $100 \pm 5\%$ **Stability:** At -15 to -25°C within specification range for 18 months.

Store dry.

Chromogenic Substrate for Lipase

Substrate for lipase

Application

Use Chromogenic Substrate for Lipase in diagnostic tests for the determination of lipase activity.

CAS: 195833-46-6**Properties****Nomenclature:** 1,2-O-Dilauryl-rac-glycero-3-glutaric acid-(6'-methylresorufin)ester**Formula:** $\text{C}_{45}\text{H}_{69}\text{NO}_8$ **Molecular weight:** 752.05 D $\lambda_{\text{max,substrate}}$: 470 nm (Tris-HCl, pH 8.4) ϵ_{470} : 57.94 [L x mmol⁻¹ x cm⁻¹] $\lambda_{\text{max,methylresorufin}}$: 581 nm (Tris-HCl, pH 8.4)**Melting range:** +29 to +31°C**pH optimum:** 7.0-9.5**Solubility:** Soluble in polar organic solvents, e.g., n-propanol, ethyl acetate, dioxane, methanol, dimethyl sulfoxide. The limit of solubility in n-propanol is 42.9 mg/mL.**Specification****Appearance:** Red, smear substance**Chromogenic Lipase Substrate** (from C): $\geq 95\%$ **Chromogenic Lipase Substrate** (HPLC): ≥ 95 area%**C** (elementary analysis): 68.2-72.5%**H** (elementary analysis): 8.7-9.7%**N** (elementary analysis): 1.4-2.4%**Methylresorufin, free** (HPLC): ≤ 0.5 area%**Isomer** (HPLC): ≤ 2 area%**Stability:** At $+2$ to $+8^\circ\text{C}$ within specification range for 36 months.**Catalog number****11 034 618 103****Pack size**

custom fill

Will be supplied as "Chromogenic Substrate for Lipase". Unit of measure is "g".

For further processing only.

Background information**Reagent proposal for lipase test**

The sensitivity of this lipase test is especially influenced by the extinction of the substrate solution, the concentration of taurodesoxycholate, the pH value, and the molarity of the Tris buffer.

Final test concentration of the substrate solution:

chrom. Lipase Substrate: 0.24 mmol/L (=180 mg/L)

colipase: 0.98 mg/L

taurodesoxycholate *(see buffer solution)

CaCl₂: 0.1 mol/L

tartrate buffer, pH 4.0: 1.6 mmol/L

stabilizers

The lipase substrate has to be dissolved in a small quantity of an organic solvent (*e.g.* n-propanol) first. Under vigorous stirring this organic solution should be injected into the tartare buffered aqueous solution with a thin beam. (The lipase substrate starts to hydrolyze at alkaline pH values.) The lipase substrate containing solution should be a micro-emulsion with an extinction of about 0.5 E. (Lower extinction of the reagent results in measurement of nonspecific serum esterases.) Stabilizers like mannitol, polywax 4000 and co-emulsifiers like lecithin, phosphoryl choline or dilauryl-glycerol-sulfate improve the stability of the micro-emulsion.

Final test concentration of the buffer solution:

Tris-HCl, pH 8.4: 41 mmol/L

taurodesoxycholate: *7.2 mmol/L (total)

desoxycholate: 1.77 mmol/L

Wavelength: 578 nm or 580 nm

Temperature: +25°C, +30°C or +37°C, respective

Buffer solution: 1 mL

Substrate solution: 0.2 mL

Sample volume: 0.02 mL

Assay time: 2 to 10 minutes

Start of reaction: with substrate

Ethylidene-4-NP-G7

Ethyliden-4-nitrophenyl- α -D-maltoheptaosid (EPS), powder

Nitrophenyl substrate

Application

Use Ethylidene-4-NP-G7 in diagnostic tests for the determination of α -amylase and pancreatic α -amylase according to the recommendations of the International Federation of Clinical Chemistry and Laboratory Medicine (IFCC).

Catalog number

10 880 078 103

Pack size

custom fill

Will be supplied as "Ethylidene-4-NP-G7". Unit of measure is "kg". Additional products: OEM reagents for the determination of α -amylase and pancreatic amylase, as well as specific inhibitory antibodies. Catalog Nos. 11 543 598 103 and 11 543 601 103 for the detection of pancreatic amylase.

For further processing only.

CAS: 96597-16-9

Properties

Formula: C₅₀H₇₇NO₃₈

Molecular weight: 1300.1 D

Specification

Appearance: White to slightly yellowish, amorphous powder

Solubility: Clear, slightly yellowish solution in water (c=70 mg/mL)

EPS (enzymatic): ≥90%

Water (K. Fischer): ≤3%

pNP-G7 (enzymatic): ≤0.1%

pNP, free: ≤0.01%

Reaction rates (α-amylase):

In Precinorm® U: 95-105%

In Precipath® U: 95-105%

Stability: At +2 to +8°C within specification range for 24 months.

Store dry.

Glupa-carboxylate

monoammonium salt

Substrate for γ-glutamyltransferase

Application

Use Glupa-carboxylate in diagnostic tests for the determination of γ-glutamyltransferase, according to the recommendations of the International Federation of Clinical Chemistry and Laboratory Medicine (IFCC).

CAS: 63699-78-5

Properties

Formula: C₁₂H₁₂N₃O₇NH₄

Molecular weight: 328.3 D

Specification

Appearance: White to yellowish crystalline powder

Solubility: Clear, yellow solution in water (c=100 mg/mL), free of fuzzi

pH value: 4.0-6.0

Molar rotation [α] 25/D: +32.0±2.0°

Melting range (Kofler): Approximately +170 to +180°C

Glupa-carboxylate, free acid (enzymatic): ≥87%

Glupa-carboxylate (HPLC): ≥99 area%

Water (K. Fischer): ≤6.2%

NH₄ (Nessler's reagent): 5.2±1%

Catalog number

10 413 151 103

Pack size

custom fill

Will be supplied as "Glupa-carboxylate, Monoammonium Salt". Unit of measure is "kg".

For further processing only.

5-Amino-2-nitrobenzoate (HPLC): ≤ 0.1 area%
 α -Glupa-carboxylate (HPLC): ≤ 0.4 area%
Thin layer chromatography (silica gel F; n-butanol/glacial acetic acid/
 $H_2O = 50/15/25$; UV, with Nihydrin): Chromatographically homogeneous
 A_{405} (Glupa-carboxylate, 6 mmol/L): 0.65-0.80
Stability: At +2 to +8°C within specification range for 24 months.
 Protect from light.

TOOS

(N-Ethyl-N-(2-hydroxy-3-sulfopropyl)-m-toluidine, monosodium salt, dihydrate

Substrate for peroxidase

Application

Use TOOS together with 4-Aminoantipyrine in an indicator reaction using peroxidase to form a quinoneimine dye.

CAS: 679787-10-1

Properties

Nomenclature: Dihydrate (N-Ethyl-N-(2-hydroxy-3-sulfopropyl)-m-toluidine

Formula: $C_{12}H_{18}NO_4SNa \times 2 H_2O$

Molecular weight: 331.37 D

Specification

Appearance: White to slightly bluish crystallizate

TOOS, mono-Na x 2 H₂O (from C): $\geq 98.0\%$

C (elementary analysis): 42.6-46.5%

H (elementary analysis): 6.3-6.8%

N (elementary analysis): 3.8-4.7%

Water (K. Fischer): 7.0-12.0%

Heavy metals (as Pb): ≤ 20 ppm

Stability: At +15 to +25°C within specification range for 24 months.

Catalog number

11 650 670 103

Pack size

custom fill

Will be supplied as "TOOS". Unit of measure is "kg".

For further processing only.

Tribromo-hydroxybenzoic acid

crystallizate

Color reagent for diagnostic tests

Application

Use Tribromo-hydroxybenzoic acid instead of phenol as a component in the trinder reaction.

CAS: 14348-40-4**Properties****Formula:** $C_7H_3O_3Br_3$ **Molecular weight:** 374.8 D**Specification****Appearance:** White crystallizate**Solubility:** Clear, colorless solution in NaOH, 0.1 mol/L (c=20.9 mg/mL)**A₄₀₅** (c=20.9 mg/mL, in NaOH, 0.1 mol/L): ≤0.020**Melting range:** +143 to +148°C**Thin-layer chromatography** (HPTLC): Corresponds to reference**2,4,6-Tribromo-hydroxybenzoic acid** (alkalimetric): ≥98.0%**Stability:** At +15 to +25°C within specification range for 36 months.**Catalog number****10 755 745 103****Pack size**

custom fill

Will be supplied as "Tribrom-Hydroxybenzoic acid". Unit of measure is "g".

For further processing only.

L(+)-Alanine

crystalline powder

Substrate for alanine aminotransferase

Application

Use L(+)-Alanine in diagnostic tests for the determination of alanine aminotransferase (ALT).

CAS: 56-41-7**Properties****Formula:** C₃H₇NO₂**Molecular weight:** 89.09 D**Specification****Appearance:** White, crystalline powder or crystals**Solubility:** Easily soluble in water and mineral acids, insoluble in organic solvents**Appearance in buffer solution:** Clear, colorless solution in phosphate buffer (c=0.9%, w/v, pH 7.4)**Microbiological test:** Corresponds**Heavy metals** (as Pb): ≤20 ppm ≙ 0.002%**Sulfate ash:** ≤0.1%**Thin layer chromatography:** Corresponds to reference**Water** (K.Fischer): ≤1.0%**L-alanine** (HClO₄ titration, based on anhydrous substance): 98.5-100.5%**L-alanine** (enzymatic, based on anhydrous substance): 97.0-105.0%**Stability:** At +15 to +40°C within specification range for 36 months. Store dry in tightly closed containers.**Catalog number****10 136 921 103****Pack size**

custom fill

Will be supplied as "L(+)-Alanin". Unit of measure is "kg".

For further processing only.

α-Ketoglutarate (2-Oxoglutarate)

free acid

Substrate for transaminases and glutamate dehydrogenase

Application

Use α-Ketoglutarate in a variety of diagnostic tests, such as for the determination of alanine aminotransferase, aspartate aminotransferase, ammonia, urea and glutamate dehydrogenase.

CAS: 328-50-7**Catalog number****10 156 736 103****Pack size**

custom fill

Will be supplied as "α-Ketoglutaric Acid, Free Acid". Unit of measure is "kg".

For further processing only.

Properties**Formula:** C₅H₆O₅**Molecular weight:** 146.1 D**Specification****Appearance:** White crystallizate**Solubility:** Clear, colorless solution in water (c=50 mg/mL)**A₄₀₅** (c=50 mg/ml in water, against water): ≤0.020**Melting range:** +113 to +117°C**α-Ketoglutaric acid** (enzymatic): ≥98%**Water** (K. Fischer): ≤1%**NH₄** (enzymatic): ≤0.1%**Heavy metals** (as Pb): ≤10 ppm**Bioburden:** ≤100 CFU/g**Reaction rates** (Glutamate pyruvate transaminase (ALT)): ≥95%**Reaction rates** (Glutamate oxalacetate transaminase (AST)): ≥95%**Stability:** At +15 to +25°C within specification range for 36 months.**α-Ketoglutarate (2-Oxoglutarate)**

disodium salt, dihydrate

Substrate in enzymatic reactions with glutamate dehydrogenase or transaminases

Application

Use α-Ketoglutarate in a variety of diagnostic tests, such as for the determination of glutamate dehydrogenase, ammonia, alanine- and aspartate aminotransferases and urea. The dihydrate formulation is well suited for dry chemistry tests.

CAS: 305-72-6**Properties****Formula:** C₅H₈O₇Na₂**Molecular weight:** 226.1 D (α-KG: 146.1 D)**Specification****Appearance:** White, crystalline powder**Solubility:** Clear, colorless solution in water, pH 7.3 (c=200 mg/mL)**A₄₀₅** (against water): ≤0.020**α-Ketoglutarate, salt** (based on value found enzymatically): ≥97%**α-Ketoglutarate, free acid** (enzymatic): ≥63%**Na** (flame photometric): 20.5±1%**Water** (K. Fischer): 15±2%**Heavy metals** (as Pb): ≤20 ppm**Stability:** At +15 to +25°C within specification range for 24 months.**Catalog number****Pack size****10 040 584 103**

custom fill

Will be supplied as "α-Ketoglutarate (α-Oxoglutarate), Di-Na". Unit of measure is "kg".

Additional formulation: Crystallized free acid, Catalog No. 10 156 736 103

For further processing only.

α -Ketoglutarate (2-Oxoglutarate)

disodium salt

Substrate in enzymatic reactions with glutamate dehydrogenase or transaminases

Application

Use α -Ketoglutarate in a variety of diagnostic tests, such as for the determination of glutamate dehydrogenase, ammonia, alanine- and aspartate aminotransferases and urea. The disodium formulation is well suited for liquid tests.

CAS: 305-72-6**Properties****Formula:** $C_5H_4O_5Na_2$ **Molecular weight:** 190.1 D (α -KG: 146.1 D)**Specification****Appearance:** White, crystalline powder**Solubility:** Clear, colorless solution in water, pH 7.3 (c=200 mg/mL)**A₄₀₅** (c=10 mg/mL in water, against water): ≤ 0.020 **α -Ketoglutarate, salt** (based on value found enzymatically): $\geq 97.5\%$ **α -Ketoglutarate, free acid** (enzymatic): $\geq 74\%$ **Na** (flame photometric): $24 \pm 2\%$ **Water** (K. Fischer): $\leq 2\%$ **Heavy metals** (as Pb): ≤ 20 ppm**Stability:** At +15 to +25°C within specification range for 24 months.**Catalog number****10 266 400 103****Pack size**

custom fill

Will be supplied as " α -Ketoglutarate, Di-Na, (M 190.1 g/mol)". Unit of measure is "kg". Additional formulation: Crystallized free acid, Catalog No. 10 156 736 103

For further processing only.

 α -Ketoglutarate (2-Oxoglutarate) for potassium test

free acid

Substrate for transaminases and glutamate dehydrogenase

Application

Use α -Ketoglutarate for enzymatic potassium tests especially to remove ammonia from the reaction.

CAS: 328-50-7**Properties****Formula:** $C_5H_6O_5$ **Molecular weight:** 146.1 D**Catalog number****11 332 775 103****Pack size**

custom fill

Will be supplied as " α -Ketoglutaric Acid for Potassium Test". Unit of measure is "kg".

For further processing only.

Specification**Appearance:** White crystallize**Solubility:** Clear, colorless solution in water (c=50 mg/mL)**A₄₀₅** (c=50 mg/mL in water, against water): ≤0.020**Melting range:** +113 to +117°C**α-Ketoglutaric acid** (enzymatic): ≥98%**Water** (K. Fischer): ≤1%**NH₄** (enzymatic): ≤0.1%**Na** (AES): ≤500 ppm**K** (AES): ≤10 ppm**Heavy metals** (as Pb): ≤10 ppm**Bioburden:** ≤100 CFU/g**Reaction rates** (Glutamate pyruvate transaminase (ALT)): ≥95%**Reaction rates** (Glutamate oxalacetate transaminase (AST)): ≥95%**Stability:** At +15 to +25°C within specification range for 36 months.**Creatine Phosphate**

disodium salt

Substrate for creatine kinase (reverse reaction)

Application

Use Creatine Phosphate in diagnostic tests for the determination of creatine kinase, according to the recommendations of the International Federation of Clinical Chemistry and Laboratory Medicine (IFCC).

CAS: 922-32-7**Properties****Molecular weight:** 327.2 D (CP: 211.1 D)**Specification****Appearance:** White crystals**Solubility:** Clear, colorless solution in water (c=150 mg/mL), free from fuzz**pH value** (c=10 mg/mL in water): 7.7-8.7**Creatine-P-Na₂ x 4 H₂O** (based on value found enzymatically): ≥97%**Creatine-P** (enzymatic): ≥63%**Creatine-P** (from P_{organic}): ≥63%**Na** (flame photometric): 14±1%**Water** (K. Fischer): 22±2%**P_{organic}** (P_{total} - P_i): ≥9.25%**P_{total}**: ≥9.25%**P_i** (acid labile): ≤0.5%**P_i** (Fiske and Subbarow): ≤1.5%**PP_i** (enzymatic): ≤0.02%**Catalog number****10 003 506 103****Pack size**

custom fill

Will be supplied as "Creatine Phosphate, Disodium Salt". Unit of measure is "kg".

For further processing only.

ATP (enzymatic, with hexokinase/G6P-DH): $\leq 0.002\%$
Sulfate (qualitative): Negative
Creatine, free: $\leq 0.5\%$
Glucose-6-P (enzymatic): $\leq 0.006\%$
PEP (enzymatic): $\leq 0.05\%$
Pyruvate (enzymatic): $\leq 0.02\%$
Kinetic of creatine kinase reaction: Corresponds to standard
Reaction rates (creatine kinase): 95-105%
A₃₃₄ (c=9 mL/mL water): ≤ 0.005
A₃₃₄ (against reaction mixture CK NAC active): ≤ 0.040
A₃₄₀ (hydrous solution): ≤ 0.120
Stability: At +2 to +8°C within specification range for 36 months.

D(-)-Lactate

monolithium salt

Substrate for D-lactate dehydrogenase

Application

Use D(-)-Lactate as a standard in tests for lactic acid.

CAS: 27848-80-2

Properties

Formula: $C_3H_5O_3Li$

Molecular weight: 96.0 D (Lactate: 89.1 D, Lactic acid: 90.1 D)

Specification

Appearance: White, crystalline powder

D(-)-Lactate (enzymatic, as anion): $\geq 91\%$

Li (flame photometric): $7.0 \pm 1.0\%$

L(+)-Lactate (enzymatic, as anion): $\leq 0.2\%$

Stability: At +15 to +25°C within specification range for 36 months.

Catalog number

10 151 874 103

Pack size

custom fill

Will be supplied as "D(-)-Lactate, Monolithium Salt". Unit of measure is "g".

For further processing only.

Di(adenosine-5'-)penta-phosphate

trilithium salt

Inhibitor of adenylate kinase

Application

Use Di(adenosine-5'-)penta-phosphate in diagnostic reagents for the determination of creatine kinase to inhibit adenylate kinase in the reaction.

Catalog number

10 161 624 103

Pack size

custom fill

Will be supplied as "Di(adenosine-5'-)penta-phosphate, Tri-Li". Unit of measure is "g".

For further processing only.

CAS: 75522-97-3

Properties

Formula: $C_{20}H_{26}N_{10}O_{22}P_5Li_3$

Molecular weight: 934.2 D (Ap5A: 916.4 D)

Specification

Appearance: White to slightly yellowish lyophilizate

Solubility: Clear, colorless to slightly yellowish solution in water (c=10 mg/mL)

Ap5A-Li₃ (from P_{organic}): ≥91%

Ap5A (A_{260} , $\epsilon=26.4$ [L x mmol⁻¹ x cm⁻¹]): ≥90%

Ap5A (from P_{organic}): ≥90%

Ap5A (HPLC): ≥95 area%

Li (flame photometric): 2.1±0.3%

Water (K. Fischer): ≤5%

P_{organic} (P_{total}-P_i): ≥15.2%

P_i: ≤1.5%

Thin layer chromatography (PEI-cellulose, KH₂PO₄, 0.75 mol/L):

Chromatographically homogeneous

A₂₅₀/A₂₆₀: 0.79±0.04

A₂₈₀/A₂₆₀: 0.21±0.03

A₂₉₀/A₂₆₀: 0.02±0.02

Stability: At +2 to +8°C within specification range for 24 months.

Glucose-1,6-diphosphate

tetra(cyclohexylammonium) salt

Substrate in diagnostic tests

Application

Use Glucose-1,6-diphosphate in diagnostic tests for the determination of inorganic phosphate and sucrose.

CAS: 10139-18-1

Properties

Formula: $C_6H_{14}O_{12}P_2 \times (C_6H_{14}N)_4 \times 4 H_2O$

Molecular weight: 808.9 D (Glucose-1,6-P₂: 340.1 D)

Specification

Appearance: Yellowish crystallizate

Glucose-1,6-P₂ (CHA)₄ x 4 H₂O: 93.0-105.0%

Glucose-1,6-P₂ (from P_{organic}): 39.0-44.0%

CHA (titrimetric): 46.-50.0%

Water (K. Fischer): 6.0-10.0%

Catalog number

10 150 827 103

Pack size

custom fill

Will be supplied as "Glucose-1,6-diphosphate, Tetra-CHA Salt".
Unit of measure is "g".

For further processing only.

P_{organic} ($P_{\text{total}} - P_i$): 7.10-8.00%

P_i: ≤0.30%

Stability: At +15 to +25°C within specification range for 36 months.

Glucose-6-phosphate

disodium salt

Substrate for glucose-6-phosphate dehydrogenase

Application

Use Glucose-6-phosphate in diagnostic tests for the determination of glucose-6-phosphate dehydrogenase.

CAS: 3671-99-6

Properties

Formula: C₆H₁₁O₉PNa₂

Molecular weight: 304.2 D (Glucose-6-P: 260.2 D)

Specification

Appearance: White lyophilizate

Solubility: Clear solution in water (c= 50 mg/mL)

Glucose-6-P (enzymatic): ≥77%

Glucose-6-P (from P_{organic}): ≥77%

Na (flame photometric): 12.5±1%

Water (K. Fischer): 8±2%

P_{organic} ($P_{\text{total}} - P_i - P_{\text{fructose-6-P}}$): ≥8.9%

P_i: ≤0.6%

Fructose-6-P (enzymatic): ≤2%

Glucose (enzymatic): ≤0.2%

Stability: At +15 to +25°C within specification range for 24 months.

Store dry.

Catalog number

10 153 079 103

Pack size

custom fill

Will be supplied as "Glucose-6-phosphate, Disodium Salt". Unit of measure is "g".

For further processing only.

N-Acetyl-L-Cysteine

crystallizate

Activator of creatine kinase

Application

Use N-Acetyl-L-Cysteine in diagnostic tests for the determination of creatine kinase, where it is used to reactivate creatine kinase as recommended by the International Federation of Clinical Chemistry and Laboratory Medicine (IFCC).

Catalog number

10 068 365 103

Pack size

custom fill

Will be supplied as "N-Acetyl-L-Cystein". Unit of measure is "kg".

For further processing only.

CAS: 616-91-1

Properties

Formula: $C_5H_9NO_3S$

Molecular weight: 163.19 D

Specification

Appearance: White crystals

Solubility: Clear, colorless solution in water (c=50 mg/mL)

Melting range: +107 to +113°C

Specific rotation [α] 25/D: +3.0° to +5.0°

IR-spectrum: Corresponds to reference

Purity (HPLC): ≥99.0 area%

N-Acetyl-L-Cysteine (content from SH-groups): ≥99.0%

N-Acetyl-L-Cysteine (alkalimetric): ≥99.0%

N-Acetyl-L-Cysteine (content from N): ≥99.0%

Content of nitrogen (elementary analysis): ≥8.5%

Screening analysis:

Particle size ≥60 μm: ≤15%

Particle size ≥140 μm: ≤50%

Heavy metals (as Pb): ≤5 ppm

Fe (AAS): ≤2.0 ppm

Cu (AAS): ≤1.0 ppm

Mn (AAS): ≤1.0 ppm

Bioburden: 100 CFU/g

Country of origin: China, USA, Canada

Heat treatment (only chinese goods): Corresponds to reference

Stability: At +15 to +40°C within specification range for 24 months.

Phosphoenolpyruvate (PEP)

tri(cyclohexylammonium) salt

Substrate for phosphoenolpyruvate carboxylase

Application

Use Phosphoenolpyruvate in diagnostic tests for the determination of carbon dioxide, creatinine or pyruvate kinase.

CAS: 138-08-9

Properties

Formula: $C_3H_2O_6P(C_6H_{14}N)_3 \times H_2O$

Molecular weight: 483.3 D (PEP: 168.0 D)

Catalog number

10 005 185 103

Pack size

custom fill

Will be supplied as "Phosphoenolpyruvate (PEP), CHA-Salt". Unit of measure is "kg".

For further processing only.

Specification**Appearance:** Colorless, crystallize**PEP salt** (based on value found enzymatically): ≥96%**PEP** (enzymatic): ≥34.5%**CHA** (titrimetric with perchloric acid): 57-67%**Water** (K. Fischer): ≤4.5%**P_i** (Fiske and Subbarow): ≤0.6%**Pyruvate** (enzymatic): ≤0.1%**Stability:** At +2 to +8°C within specification range for 24 months.

Additional formulation crystallized monosodium salt, Catalog No. 10 152

960 103

Phosphoenolpyruvate (PEP)

monosodium salt

Substrate for phosphoenolpyruvate carboxylase

Application

Use Phosphoenolpyruvate in diagnostic tests for the determination of carbon dioxide, creatinine or pyruvate kinase.

CAS: 138-08-9**Properties****Formula:** C₃H₄O₆PNa x H₂O**Molecular weight:** 208.0 D (PEP: 168.0 D)**Specification****Appearance:** White, crystalline powder**PEP-Na x H₂O** (based on value found enzymatically): ≥94%**PEP** (enzymatic): ≥76.0%**Na** (flame photometric): 9-13%**Water** (K. Fischer): 8-10%**P_i**: ≤0.6%**Pyruvate** (enzymatic): ≤0.1%**Stability:** At +2 to +8°C within specification range for 24 months.**Catalog number****10 152 960 103****Pack size**

custom fill

Will be supplied as "Phosphoenolpyruvate (PEP), Mono-Na Salt".
Unit of measure is "g".

For further processing only.

Phosphoenolpyruvate (PEP), for potassium test

tri(cyclohexylammonium) salt

Substrate for phosphoenolpyruvate carboxylase

Application

Use Phosphoenolpyruvate as a substrate for pyruvate kinase, stimulated by potassium, for the enzymatic determination of potassium.

CAS: 138-08-9

Properties

Formula: $C_3H_2O_6P(C_6H_{14}N)_3 \times H_2O$

Molecular weight: 483.3 D (PEP: 168.0 D)

Specification

Appearance: Colorless crystallate

PEP-(CHA)₃ (from content found enzymatically): ≥96%

PEP (enzymatic): ≥34.5%

CHA (titrimetric): 57-67%

Water (K. Fischer): ≤4.5%

P_i: ≤0.6%

Pyruvate (enzymatic): ≤0.1%

Na (AES): ≤100 ppm

K: ≤10 ppm

Stability: At +2 to +8°C within specification range for 24 months.

Catalog number

11 333 968 103

Pack size

custom fill

Will be supplied as "PEP, tri-CHA for Potassium Test". Unit of measure is "kg".

Additional formulation: Crystallized monosodium salt, Catalog No. 10 152 960 103

For further processing only.

Pyruvate

monosodium salt

Substrate for many enzymes, such as lactate dehydrogenase and pyruvate kinase.

Application

Use Pyruvate in diagnostic tests for the determination of lactate dehydrogenase.

CAS: 57-60-3

Properties

Formula: $C_3H_3O_3Na$

Molecular weight: 110.0 D

Catalog number

10 005 525 103

Pack size

custom fill

Will be supplied as "Pyruvate Monosodium Salt". Unit of measure is "g".

For further processing only.

Specification**Appearance:** White, crystalline powder**Pyruvate-Na** (from content found enzymatically): 96-103%**Pyruvate** (enzymatic, based on anion): 77-81%**Na** (flame photometric): 20.5-21.5%**Bioburden:** ≤100 CFU/g**Heavy metals** (as Pb): ≤10 ppm**Stability:** At +15 to +25°C within specification range for 18 months.**S-Butyrylthiocholine Iodide**

crystallizate

Substrate for cholinesterase

Application

Use S-Butyrylthiocholine Iodide in diagnostic tests for the determination of cholinesterase.

CAS: 1866-16-6**Properties****Formula:** C₉H₂₀NOSJ**Molecular weight:** 317.2 D**Specification****Appearance:** Colorless crystallizate**Solubility:** Clear, colorless solution in water (c=70 mg/mL)**Melting range:** +172 to +174°C**Butyrylthiocholine iodide** (titrimetric) : ≥98.0%**Thiocholine iodide, free:** ≤0.15%**Reaction rates** (inhibitors of choline esterase): 100±5%**Stability:** At +2 to +8°C within specification range for 24 months.**Catalog number****10 034 614 103****Pack size**

custom fill

Will be supplied as "S-Butyrylthiocholine Iodide". Unit of measure is "kg".

For further processing only.



Immunology

Antibodies/Antigens	170	Interference Blocker	217
Monoclonal Antibodies	170	Antibody Interference Blocker	217
Polyclonal Antibodies	192	Antibody Interference Control	224
Antigens	195	Surface Interference Blocker	226
Biotin/Streptavidin System	204	Detection Interference Blocker	233
Streptavidin	204	Marker Enzymes and Substrates	235
Biotin Labels	206	Enzymes	235
Fluorescent Labels	208	Substrates	242
Solid Phases	210	Serum	246
Dyes	215		

MAB<CK-MB>M-7.4.5 IgG

lyophilizate

Qualified for **cobas**[®] platforms**Application**

MAB<CK-MB>M-7.4.5 IgG is qualified for heterogeneous immunoassays (HetIAs).

Product description**Antibody class:** IgG1, kappa**Properties**

MAB<CK-MB>M-7.4.5 IgG is a monoclonal antibody directed to creatine kinase MB. It is lyophilized from a solution containing protein, potassium phosphate buffer and sodium chloride. No preservatives are added. Sandwich partner to MAB<CK-MB>M-6.12.47 IgG.

Specification**Appearance:** White lyophilizate**Solubility:** Clear to slightly opalescent solution in NaCl, 0.9% (c=10 mg/mL)**Protein** (Biuret): ≥0.6 mg/mg lyophilizate**Purity** (HPLC / Mono Q): ≥90 area%**pH 5.5 treatment** (30 minutes): Corresponds to specification**Stability:** At -15 to -25°C within specification range for 60 months.

Avoid repeated freezing and thawing.

Catalog number**11 719 815 103****Pack size**

5 mg (samples), ≥50 mg (custom fill)

Will be supplied as "MAK<CK-MB>M-7.4.5 IGG". Unit of measure is "mg".

For further processing only.

MAB<CK-MB>M-6.12.47 IgG

lyophilizate

Qualified for **cobas**[®] platforms**Application**

MAB<CK-MB>M-6.12.47 IgG is qualified for heterogeneous immunoassays (HetIAs).

Product description**Antibody class:** IgG1, kappa**Properties**

MAB<CK-MB>M-6.12.47 IgG is a monoclonal antibody directed to creatine kinase MB. It is lyophilized from a solution containing protein, potassium phosphate buffer and sodium chloride. No preservatives are added. Sandwich partner for MAB<CK-MB>M-7.4.5 IgG.

Catalog number**11 719 823 103****Pack size**

5 mg (samples), ≥50 mg (custom fill)

Will be supplied as "MAK<CK-MB>M-6.12.47 IGG". Unit of measure is "mg".

For further processing only.

Specification**Appearance:** White lyophilizate**Solubility:** Clear to slightly opalescent solution in NaCl, 0.9% (c=10 mg/mL)**Protein** (Biuret): ≥0.6 mg/mg lyophilizate**Purity** (HPLC / Mono Q): ≥90 area%**pH 5.5 treatment** (30 minutes): Corresponds to specification**Stability:** At -15 to -25°C within specification range for 60 months.

Avoid repeated freezing and thawing.

MAB<CK-MM>Mix

frozen solution

For measurement of human creatine kinase isoenzyme (CK-MB), which is a well established tool in confirming the diagnosis of acute myocardial infarction.

Application

MAB<CK-MM>Mix is a main compound of the CK-MB assay.

Product description**Immunogen:** Human creatinine kinase isoenzyme MM (h-CK-MM)**Inhibitor capacity** (for information only):

4500 U/L CK-MM: ≥99.6%

300 U/L CK-BB: ± 5%

Properties

The MAB<CK-MM>Mix consists of four highly specific monoclonal mouse antibodies directed to human creatinin kinase isoenzyme MM (h-CK-MM). The frozen solution contains protein, potassium phosphate buffer and sodium chloride. No preservative are added.

Specification**Appearance:** Clear to slightly opalescent solution after thaw**pH value** (+25°C): 7.4-7.6**HPLC** (Mono Q, basic material): ≥80%**Protein** (BCA): 10.8-13.2 mg/mL**Aggregates** (TSK 3000): ≤10%**Inhibitor capacity:**

CK-MM 4,500 U/l (+37°C): ≥96%

CK-BB 300 U/l (+37°C): ±5%

pH 5.5 treatment (30 minutes): Corresponds to specification**Stability:** At -15 to -25°C within specification range for 36 months.

Avoid repeated freezing and thawing.

Catalog number**04 688 457 103****Pack size**

1 mL, 10 mL, 50 mL, 100 mL, 1000 mL

Will be supplied as "Mab<CK-MM>Mix". Unit of measure is "L".



For further processing only.

Background information

Human creatine kinase isoenzyme CK-MB consist of two subunits: CK-M and CK-B. By inhibiting the CK-M subunit, the creatine kinase reaction is triggered exclusively by the β subunit of creatine kinase, which accounts for one-half of the activity of CK-MB.

MAB<AFP>M-LJ738 IgG

lyophilizate

Qualified for **cobas**[®] platforms

Application

MAB<AFP>M-LJ738 IgG is qualified for heterogeneous immunoassays (HetIAs).

Product description

Antibody class: IgG1, kappa

Properties

MAB<AFP>M-LJ738 IgG is a monoclonal antibody directed to alpha fetoprotein. It is lyophilized from a solution containing protein, potassium phosphate buffer and sodium chloride. No preservatives are added. Sandwich partner for MAB<AFP>M-TU11 IgG.

Specification

Appearance: White lyophilizate

Solubility: Clear to slightly opalescent solution in NaCl, 0.9% (c=10 mg/mL)

Protein (Biuret): ≥ 0.6 mg/mg lyophilizate

Purity (HPLC / Mono Q): ≥ 90 area%

pH 5.5 treatment (30 minutes): Corresponds to specification

Stability: At -15 to -25°C within specification range for 60 months.

Avoid repeated freezing and thawing.

Catalog number

11 492 101 103

Pack size

5 mg (samples), ≥ 50 mg (custom fill)

Will be supplied as "MAK<AFP>M-LJ738-IGG(DE)". Unit of measure is "g active ingredient".

For further processing only.

MAB<AFP>M-TU11 IgG

lyophilizate

Qualified for **cobas**[®] platforms

Application

MAB<AFP>M-TU11 IgG is qualified for heterogeneous immunoassays (HetIAs).

Catalog number

11 492 080 103

Pack size

5 mg (samples), ≥ 50 mg (custom fill)

Will be supplied as "MAK<AFP>M-TU11-IGG(DE)". Unit of measure is "g active ingredient".

For further processing only.

Product description

Antibody class: IgG2a, kappa

Properties

MAB<AFP>M-TU11 IgG is a monoclonal antibody directed to alpha fetoprotein. It is lyophilized from a solution containing protein, potassium phosphate buffer and sodium chloride. No preservatives are added. Sandwich partner for MAB<AFP>M-LJ738 IgG.

Specification

Appearance: White lyophilizate

Solubility: Clear to slightly opalescent solution in NaCl, 0.9% (c=10 mg/mL)

Protein (Biuret): ≥0.6 mg/mg lyophilizate

Purity (HPLC / Mono Q): ≥90 area%

pH 5.5 treatment (30 minutes): Corresponds to specification

Stability: At -15 to -25°C within specification range for 60 months.

Avoid repeated freezing and thawing.

MAB<CEA>M-TU2 IgG

lyophilizate

Application

MAB<CEA>M-TU2 IgG is qualified for heterogeneous immunoassays (HetIAs).

Product description

Immunogen: Carcinoembryonic antigen (CEA)

Spleen donor: Mouse Balb/c

Antibody class: IgG1, kappa

Properties

MAB<CEA>M-TU2 IgG is a monoclonal antibody directed to carcinoembryonic antigen. It is lyophilized from a solution containing protein, potassium phosphate buffer and sodium chloride; pH 7.5. No preservatives are added. Sandwich partner for MAB<CEA>M-TU3 IgG.

Specification

Appearance: White lyophilizate

Solubility: Clear to slightly opalescent solution in NaCl, 0.9% (c=5mg/mL)

Purity (HPLC / Mono Q): ≥90 area%

pH 5.5 treatment (30 minutes): Corresponds to specification

Stability: At -15 to -25°C within specification range for 24 months.

Avoid repeated freezing and thawing.

Catalog number

11 353 713 103

Pack size

10 mg (samples), 100 mg

Will be supplied as "MAK<CEA>M-TU2-IGG *SQ". Unit of measure is "mg active ingredient".

For further processing only.

MAB<CEA>M-TU3 IgG

lyophilizate

Qualified for **cobas**[®] platforms**Application**

MAB<CEA>M-TU3 IgG is qualified for heterogeneous immunoassays (HetIAs).

Product description

Antibody class: IgG1, kappa

Properties

MAB<CEA>M-TU3 IgG is a monoclonal antibody directed to carcinoembryonic antigen. It is lyophilized from a solution containing protein, potassium phosphate buffer and sodium chloride; pH 7.5. No preservatives are added. Sandwich partner for MAB<CEA>M-TU2 IgG.

Specification

Appearance: White lyophilizate

Solubility: Clear, to slightly opalescent solution in NaCl, 0.9% (c=10 mg/mL)

Protein (Biuret): ≥0.6 mg/mg lyophilizate

Purity (HPLC / Mono Q): ≥90 area%

pH 5.5 treatment (30 minutes): Corresponds to specification

Stability: At -15 to -25°C within specification range for 24 months.

Avoid repeated freezing and thawing.

Catalog number**10 777 498 103****Pack size**

5 mg (samples), ≥50 mg (custom fill)

Will be supplied as "MAK<CEA>M-TU3-IGG(DE)". Unit of measure is "mg active ingredient".

For further processing only.

MAB<DD>M-1.2.57 IgG

lyophilizate

Qualified for **cobas**[®] platforms**Application**

MAB<DD>M-1.2.57 IgG is qualified for heterogeneous immunoassays (HetIAs).

Product description

Immunogen: Human fibrinogen cleavage product D-Dimer

Spleen donor: Mouse Balb/c

Antibody class: IgG1, kappa

Properties

MAB<DD>M-1.2.57 IgG is a monoclonal antibody directed to human fibrinogen cleavage product D-Dimer. It is lyophilized from a solution containing protein, potassium phosphate buffer and sodium chloride.

Catalog number**12 156 903 103****Pack size**

5 mg (samples), ≥50 mg (custom fill)

Will be supplied as "MAK<DD>M-1.2.57-IgG(SP/Q)". Unit of measure is "g active ingredient".

For further processing only.

No preservatives are added. Sandwich partner for MAB<DD>M-2.1.16 IgG.

Specification

Appearance: White lyophilizate

Solubility: Clear to slightly opalescent solution in NaCl, 0.9% (c=10 mg/mL)

Protein (Biuret) : ≥0.6 mg/mg lyophilizate

Aggregated IgG (HPLC / TSK 3000): ≤10 area%

pH 5.5 treatment (30 minutes): Corresponds to specification

Stability: At -15 to -25°C within specification range for 24 months.

Avoid repeated freezing and thawing.

MAB<DD>M-2.1.16 IgG

lyophilizate

Qualified for **cobas**[®] platforms

Application

MAB<DD>M-2.1.16 IgG is qualified for heterogeneous immunoassays (HetIAs).

Product description

Immunogen: Human fibrinogen cleavage product D-Dimer

Spleen donor: Mouse Balb/c

Antibody class: IgG1, kappa

Properties

MAB<DD>M-2.1.16 IgG is a monoclonal antibody directed to human fibrinogen cleavage product D-Dimer. It is lyophilized from a solution containing protein, potassium phosphate buffer and sodium chloride. No preservatives are added. Sandwich partner for MAB<DD>M-1.2.57 IgG.

Specification

Appearance: White lyophilizate

Solubility: Clear to slightly opalescent solution in NaCl, 0.9% (c=10 mg/mL)

Protein (Biuret): ≥0.6 mg/mg lyophilizate

Aggregated IgG (HPLC / TSK 3000): ≤10 area%

pH 5.5 treatment (30 minutes): Corresponds to specification

Stability: At -15 to -25°C within specification range for 24 months.

Avoid repeated freezing and thawing.

Catalog number

12 045 206 103

Pack size

5 mg (samples), ≥50 mg (custom fill)

Will be supplied as "MAK<DD>M-2.1.16-IgG(SP/Q)". Unit of measure is "g active ingredient".

For further processing only.

MAB<Ferr>M-3.170 IgG

lyophilizate

Qualified for **cobas**[®] platforms**Application**

MAB<Ferr>M-3.170 IgG is qualified for heterogeneous immunoassays (HetIAs).

Product description**Immunogen:** Human liver ferritin**Spleen donor:** Mouse Balb/c**Antibody class:** IgG1, kappa**Cross reactivity to:** Spleen ferritin 74%; heart ferritin 11%**Properties**

MAB<Ferr>M-3.170 IgG is a monoclonal antibody directed to human liver ferritin. It is lyophilized from a solution containing protein, potassium phosphate buffer and sodium chloride. No preservatives are added. Sandwich partner for MAB<Ferr>M-4.184 IgG.

Specification**Appearance:** White lyophilizate**Solubility:** Clear to slightly opalescent solution in NaCl, 0.9% (c=10 mg/mL)**Protein** (Biuret): ≥0.6 mg/mg lyophilizate**Purity** (HPLC / Mono Q): ≥90 area%**pH 5.5 treatment** (30 minutes): Corresponds to specification**Stability:** At -15 to -25°C within specification range for 24 months.

Avoid repeated freezing and thawing.

Catalog number**11 547 089 103****Pack size**

5 mg (samples), ≥50 mg (custom fill)

Will be supplied as "MAK<Ferr>M-3.170-IgG". Unit of measure is "g active ingredient".

For further processing only.

MAB<Ferr>M-4.184 IgG

lyophilizate

Qualified for **cobas**[®] platforms**Application**

MAB<Ferr>M-3.170 IgG is qualified for heterogeneous immunoassays (HetIAs).

Product description**Immunogen:** Human liver ferritin**Spleen donor:** Mouse Balb/c**Antibody class:** IgG2a, kappa**Cross reactivity to:** Spleen ferritin 60%; heart ferritin 11%**Catalog number****11 547 119 103****Pack size**

5 mg (samples), ≥50 mg (custom fill)

Will be supplied as "MAK<Ferr>M-4.184-IgG". Unit of measure is "g active ingredient".

For further processing only.

Properties

MAB<Ferr>M-4.184 IgG is a monoclonal antibody directed to human liver ferritin. It is lyophilized from a solution containing protein, potassium phosphate buffer and sodium chloride. No preservatives are added. Sandwich partner for MAB<Ferr>M-3.170 IgG.

Specification

Appearance: White lyophilizate

Solubility: Clear to slightly opalescent solution in NaCl, 0.9% (c=10 mg/mL)

Protein (Biuret): ≥0.6 mg/mg lyophilizate

Purity (HPLC / Mono Q): ≥90 area%

pH 5.5 treatment (30 minutes): Corresponds to specification

Stability: At -15 to -25°C within specification range for 24 months.

Avoid repeated freezing and thawing.

MAB<FSH>M-1.303 IgG

lyophilizate

Qualified for **cobas**[®] platforms

Application

MAB<FSH>M-1.303 IgG is qualified for heterogeneous immunoassays (HetIAs).

Product description

Immunogen: Human follicle stimulating hormone (FSH)

Spleen donor: Mouse Balb/c

Antibody class: IgG1, kappa

Cross reactivity to: Human chorionic gonadotropin (HCG) <0.1%; Luteinizing hormone (LH) <1.5%; Thyroid stimulating hormone (TSH) <2.5%

Properties

MAB<FSH>M-1.303-IgG is a monoclonal antibody directed to human follicle stimulating hormone. It is lyophilized from a solution containing protein, potassium phosphate buffer and sodium chloride. No preservatives are added. Sandwich partner for MAB<FSH>M-W3 IgG.

Specification

Appearance: White lyophilizate

Solubility: Clear to slightly opalescent solution in NaCl, 0.9% (c=10 mg/mL)

Protein (Biuret): ≥0.6 mg/mg lyophilizate

Purity (HPLC / Mono Q): ≥90 area%

pH 5.5 treatment (30 minutes): Corresponds to specification

Catalog number

11 493 540 103

Pack size

5 mg (samples), ≥50 mg (custom fill)

Will be supplied as "MAK<FSH>M-1.303-IGG". Unit of measure is "mg active ingredient".

For further processing only.

Stability: At -15 to -25°C within specification range for 24 months.
Avoid repeated freezing and thawing.

MAB<FSH>M-W3g IgG

lyophilizate

Qualified for **cobas**[®] platforms

Application

MAB<FSH>M-W3 IgG qualified for heterogeneous immunoassays (HetIAs).

Product description

Immunogen: Human follicle stimulating hormone (FSH)

Spleen donor: Mouse Balb/c

Antibody class: IgG2a, kappa

Cross reactivity to: Human chorionic gonadotropin (HCG) not detectable; Luteinizing hormone (LH) <1.0%; Thyroid stimulating hormone (TSH) <3.5%

Properties

MAB<FSH>M-W3 IgG is a monoclonal antibody directed to human follicle stimulating hormone. It is lyophilized from a solution containing protein, potassium phosphate buffer and sodium chloride. No preservatives are added. Sandwich partner for MAB<FSH>M-1.303 IgG.

Specification

Appearance: White lyophilizate

Solubility: Clear to slightly opalescent solution in NaCl, 0.9% (c=10 mg/mL)

Protein (Biuret): ≥0.6 mg/mg lyophilizate

Purity (HPLC / Mono Q): ≥90 area%

pH 5.5 treatment (30 minutes): Corresponds to specification

Stability: At -15 to -25°C within specification range for 24 months.
Avoid repeated freezing and thawing.

Catalog number

11 493 531 103

Pack size

5 mg (samples), ≥50 mg (custom fill)

Will be supplied as "MAK<FSH>M-W3-IGG". Unit of measure is "mg active ingredient".

For further processing only.

MAB<HCG>M-INN2 IgG

lyophilizate

Qualified for **cobas**[®] platforms

Application

MAB<HCG>M-INN2 IgG is qualified for heterogeneous immunoassays (HetIAs).

Catalog number

03 116 263 103

Pack size

5 mg (samples), ≥50 mg (custom fill)

Will be supplied as "MAK<HCG>M-INN2-IgG". Unit of measure is "mg active ingredient".

For further processing only.

Product description**Immunogen:** Human chorionic gonadotropin b-chain (HCG)**Spleen donor:** Mouse Balb/c**Antibody class:** IgG1, kappa**Cross reactivity to:** Luteinizing hormone (LH) <0,3%; Follicle stimulating hormone (FSH) <0.1%; Thyroid stimulating hormone (TSH) <0.1%**Properties**

MAB<HCG>M-INN2 IgG is a monoclonal antibody directed to human chorionic gonadotropin β -chain. It is lyophilized from a solution containing protein, potassium phosphate buffer and sodium chloride. No preservatives are added. Sandwich partner for MAB<HCG>M-INN2.

Specification**Appearance:** White lyophilizate**Solubility:** Clear, to slightly opalescent solution in NaCl, 0.9% (c=10 mg/mL)**Protein** (Biuret): ≥ 0.6 mg/mg lyophilizate**Purity** (HPLC / Mono Q): ≥ 90 area%**pH 5.5 treatment** (30 minutes): Corresponds to specification**Stability:** At -15 to -25°C within specification range for 24 months.

Avoid repeated freezing and thawing.

MAB<HCG>M-INN2 IgG

lyophilizate

Qualified for **cobas**[®] platforms**Application**

MAB<HCG>M-INN2 IgG is qualified for heterogeneous immunoassays (HetIAs).

Product description**Immunogen:** Human chorionic gonadotropin β -chain (HCG)**Spleen donor:** Mouse Balb/c**Antibody class:** IgG1, kappa**Cross reactivity to:** Luteinizing hormone (LH) <5.0%; Follicle stimulating hormone (FSH) <0.2%; Thyroid stimulating hormone (TSH) <0.5%**Catalog number****11 812 564 103****Pack size**5 mg (samples), ≥ 50 mg (custom fill)

Will be supplied as "MAK<HCG>M-INN22-IgG". Unit of measure is "mg active ingredient".



For further processing only.

Properties

MAB<HCG>M-INN22 IgG is a monoclonal antibody directed to human chorionic gonadotropin β -chain. It is lyophilized from a solution containing protein, potassium phosphate buffer and sodium chloride. No preservatives are added. Sandwich partner for MAK<HCG>M-INN2.

Specification

Appearance: White lyophilizate

Solubility: Clear, to slightly opalescent solution in NaCl, 0.9% (c=10 mg/mL)

Protein (Biuret): ≥ 0.6 mg/mg lyophilizate

Purity (HPLC / Mono Q): ≥ 90 area%

pH 5.5 treatment (30 minutes): Corresponds to specification

Stability: At -15 to -25°C within specification range for 24 months. Avoid repeated freezing and thawing.

MAB<H-S-Amy>M-Tu66C7 IgG

lyophilizate

For measurement of acute pancreatitis (pancreatic α -amylase) in human serum and urine selective blocking of salivary α -amylase isoenzyme is achieved in the presence of the pancreatic h- α -amylase. Qualified for the **cobas**[®] platforms.

Application

The combination of MAB<H-S-Amy>Tu88E8 and MAB<H-S-Amy>Tu66C7 inhibits the human salivary α -amylase $\geq 97\%$ while maintaining the activity of the pancreatic h- α -amylase.

Product description

Immunogen: Human salivary amylase

Spleen donor: Mouse Balb/c

Antibody class: IgG1, kappa

Clone: Tu66C7

Cross reactivity to h-pancreas α -amylase: ≤ 1 U/gW

Properties

MAB<H-S-Amy>M-Tu66C7 IgG fraction is purified by chromatography and lyophilized from a solution containing protein (≥ 20 mg/mL), potassium-phosphate buffer and NaCl. No preservatives are added.

Specification

Appearance: White lyophilizate

Solubility: Clear colorless solution in NaCl, 0.9% (c=10 mg/mL)

Protein (Biuret): ≥ 0.7 mg/mg lyophilizate

Purity (HPLC / Mono Q): ≥ 90 area%

Catalog number

11 543 601 103

Pack size

custom fill

Will be supplied as "MAK<H-S-Amy>M-Tu66C7-IgG(BR)SQ". Unit of measure is "g active ingredient".

For further processing only.

Function testing (synergetic effects at +37°C):

h salivary amylase + MAB : ≤3% amylase activity

h pancreas amylase + MAB <S-AMY>: ≥98% amylase activity

Stability: At -15 to -25°C within specification range for 24 months.

Avoid repeated freezing and thawing.

MAB<H-S-Amy>M-Tu88E8 IgG**lyophilizate**

For measurement of acute pancreatitis (pancreatic α-amylase) in human serum and urine selective blocking of salivary α-amylase isoenzyme is achieved in the presence of the pancreatic h-α-amylase.

Qualified for the **cobas**[®] platforms.**Application**

The combination of MAB<H-S-Amy>Tu88E8 and MAB<H-S-Amy>Tu66C7 inhibits the human salivary α-amylase ≥97% while maintaining the activity of the pancreatic h-α-amylase.

Product description**Immunogen:** Human salivary amylase**Spleen donor:** Mouse Balb/c**Antibody class:** IgG2a, kappa**Clone:** Tu88E8**Cross reactivity to h-pancreas α-amylase:** ≤1 U/gW**Properties**

MAB<H-S-Amy>M-Tu88E8 IgG fraction is purified by chromatography and lyophilized from a solution containing protein (≥20 mg/mL), potassium-phosphate buffer and NaCl. No preservatives are added.

Specification**Appearance:** White lyophilizate**Solubility:** Clear, colorless solution in NaCl, 0.9% (c=10 mg/mL)**Protein** (Biuret): ≥0.7 mg/mg lyophilizate**Purity** (HPLC / Mono Q): ≥90 area%**Function testing** (synergetic effects at +37°C):

h salivary amylase + MAB : ≤3% amylase activity

h pancreas amylase + MAB : ≥98% amylase activity

Stability: At -15 to -25°C within specification range for 24 months.

Avoid repeated freezing and thawing.

Catalog number**11 543 598 103****Pack size**

custom fill

Will be supplied as "MAK<H-S-Amy>M-Tu88E8-IgG(BR)SQ". Unit of measure is "g active ingredient".

For further processing only.

MAB<IGE>M-323 IgG

lyophilizate

Qualified for **cobas**[®] platforms**Application**

MAB<IGE>M-323 IgG is qualified for heterogeneous immunoassays (HetIAs).

Product description**Immunogen:** Human IgE**Spleen donor:** Mouse Balb/c**Antibody class:** IgG1, light chain kappa**Cross reactivity to:** Human IgM, IgG and IgA: Not detectable**Properties**

MAB<IGE>M-323 IgG is a monoclonal antibody directed to human IgE. It is lyophilized from a solution containing protein, potassium phosphate buffer and sodium chloride. No preservatives are added. Sandwich partner for MAB<IGE>M-7H8.

Specification**Appearance:** White lyophilizate**Solubility:** Clear to slightly opalescent solution in NaCl, 0.9% (c=10 mg/mL)**Protein** (Biuret): ≥0.6 mg/mg lyophilizate**Purity** (HPLC / Mono Q): ≥90 area%**pH 5.5 treatment** (30 minutes): Corresponds to specification**Stability:** At -15 to -25°C within specification range for 24 months.

Avoid repeated freezing and thawing.

Catalog number**11 543 393 103****Pack size**

5 mg (samples), ≥50 mg (custom fill)

Will be supplied as "MAK<IGE>M-323-IgG". Unit of measure is "g active ingredient".

For further processing only.

MAB<IGE>M-7H8 IgG

lyophilizate

Qualified for **cobas**[®] platforms**Application**

MAB<IGE>M-7H8 IgG is qualified for heterogeneous immunoassays (HetIAs).

Product description**Immunogen:** Human IgE**Spleen donor:** Mouse Balb/c**Antibody class:** IgG1, kappa**Cross reactivity to:** Human IgM, IgG and IgA not detectable**Catalog number****11 988 204 103****Pack size**

5 mg (samples), ≥50 mg (custom fill)

Will be supplied as "MAK<IGE>M-7H8-IGG". Unit of measure is "g active ingredient".

For further processing only.

Properties

MAB<IGE>M-7H8 IgG is a monoclonal antibody directed to human IgE. It is lyophilized from a solution containing protein, potassium phosphate buffer and sodium chloride. No preservatives are added. Sandwich partner for MAB<IGE>M-323 IgG.

Specification

Appearance: White lyophilizate

Solubility: Clear to slightly opalescent solution in NaCl, 0.9% (c=10 mg/mL)

Protein (Biuret): ≥0.6 mg/mg lyophilizate

Purity (HPLC / Mono Q): ≥90 area%

pH 5.5 treatment (30 minutes): Corresponds to specification

Stability: At -15 to -25°C within specification range for 24 months. Avoid repeated freezing and thawing.

MAB<INSULIN>M-BM1 IgG

lyophilizate

Qualified for **cobas**[®] platforms

Application

MAB<INSULIN>M-BM1 IgG is qualified for heterogeneous immunoassays (HetIAs).

Product description

Antibody class: IgG1, kappa

Properties

MAB<INSULIN>M-BM1 IgG is a monoclonal antibody to insulin. It is lyophilized from a solution containing protein, potassium phosphate buffer and sodium chloride. No preservatives are added. Sandwich partner for MAB<INSULIN>M-ST3 IgG.

Specification

Appearance: White lyophilizate

Solubility: Clear to slightly opalescent solution in NaCl, 0.9% (c=10 mg/mL)

Protein (Biuret): ≥0.6 mg/mg lyophilizate

Purity (HPLC / Mono Q): ≥90 area%

pH 5.5 treatment (30 minutes): Corresponds to specification

Stability: At -15 to -25°C within specification range for 60 months. Avoid repeated freezing and thawing.

Catalog number

12 208 725 103

Pack size

5 mg (samples), ≥50 mg (custom fill)

Will be supplied as "MAK<INSULIN>M-BM1-IgG". Unit of measure is "mg active ingredient".

For further processing only.

MAB<INSULIN>M-ST3 IgG

lyophilizate

Qualified for **cobas**[®] platforms**Application**

MAB<INSULIN>M-ST3 IgG is qualified for heterogeneous immunoassays (HetIAs).

Product description**Antibody class:** IgG1, kappa**Cross reactivity to:** Human pro insulin 1.2%; porcine insulin (strong recognition, determined by radio immuno assay (RIA)); bovine insulin (weak recognition, determined by RIA)**Properties**

MAB<INSULIN>M-ST3 IgG is a monoclonal antibody directed to insulin. It is lyophilized from a solution containing protein, potassium phosphate buffer and sodium chloride. No preservatives are added. Sandwich partner for MAB<INSULIN>M-BM1 IgG.

Specification**Appearance:** White lyophilizate**Solubility:** Clear to slightly opalescent solution in NaCl, 0.9% (c=10 mg/mL)**Protein** (Biuret): ≥0.6 mg/mg lyophilizate**Purity** (HPLC / Mono Q): ≥90 area%**pH 5.5 treatment** (30 minutes): Corresponds to specification**Stability:** At -15 to -25°C within specification range for 60 months. Avoid repeated freezing and thawing.**Catalog number****12 208 750 103****Pack size**

5 mg (samples), ≥50 mg (custom fill)

Will be supplied as "MAK<INSULIN>M-ST3-IgG". Unit of measure is "mg active ingredient".

For further processing only.

MAB<LH>M-11412 IgG

lyophilizate

Qualified for **cobas**[®] platforms**Application**

MAB<LH>M-11412 IgG is qualified for heterogeneous immunoassays (HetIAs).

Product description**Immunogen:** Luteinizing hormone (LH)**Antibody class:** IgG1, kappa**Cross reactivity to:** Human chorionic gonadotropin <0.1%**Catalog number****11 547 925 103****Pack size**

5 mg (samples), ≥50 mg (custom fill)

Will be supplied as "MAK<LH>M-11412-IgG". Unit of measure is "g active ingredient".

For further processing only.

Properties

MAB<LH>M-11412 IgG is a monoclonal antibody directed to luteinizing hormone. It is lyophilized from a solution containing protein, potassium phosphate buffer and sodium chloride. No preservatives are added. Sandwich partner for MAB<LH>M-2.406 IgG.

Specification

Appearance: White lyophilizate

Solubility: Clear to slightly opalescent solution in NaCl, 0.9% (c=10 mg/mL)

Protein (Biuret): ≥0.6 mg/mg lyophilizate

Purity (HPLC / Mono Q): ≥90 area%

pH 5.5 treatment (30 minutes): Corresponds to specification

Stability: At -15 to -25°C within specification range for 60 months. Avoid repeated freezing and thawing.

MAB<LH>M-2.406 IgG**lyophilizate**

Qualified for **cobas**[®] platforms

Application

MAB<LH>M-2.406 IgG is qualified for heterogeneous immunoassays (HetIAs).

Product description

Immunogen: Luteinizing hormone (LH)

Antibody class: IgG1, kappa

Cross reactivity to: Human chorionic gonadotropin (HCG) not detectable; Follicle stimulating hormone (FSH) <0.3%; Thyroid stimulating hormone (TSH) <0.3%

Properties

MAB<LH>M-2.406-IgG is a monoclonal antibody directed to luteinizing hormone. It is lyophilized from a solution containing protein, potassium phosphate buffer and sodium chloride. No preservatives are added. Sandwich partner for MAB<LH>M-11412 IgG.

Specification

Appearance: White lyophilizate

Solubility: Clear to slightly opalescent solution in NaCl, 0.9% (c=10mg/mL)

Protein (Biuret): ≥0.6 mg/mg lyophilizate

Purity (HPLC / Mono Q): ≥90 area%

pH 5.5 treatment (30 minutes): Corresponds to specification

Catalog number

11 547 038 103

Pack size

5 mg (samples), ≥50 mg (custom fill)

Will be supplied as "MAK<LH>M-2.406-IgG". Unit of measure is "g active ingredient".

For further processing only.

Stability: At -15 to -25°C within specification range for 60 months.
Avoid repeated freezing and thawing.

MAB<TSH>M-A8 IgG

lyophilizate

Application

MAB<TSH>M-A8 IgG is qualified for heterogeneous immunoassays (HetIAs).

Product description

Immunogen: Human thyroid stimulating hormone

Spleen donor: Mouse Balb/c

Antibody class: IgG1, kappa

Cross reactivity to: Luteinizing hormone (LH) <0.1%; Follicle stimulating hormone (FSH) <0.1%

Properties

MAB<TSH>M-A8 IgG is a monoclonal antibody directed to thyroid stimulating hormone. It is lyophilized from a solution containing protein, potassium phosphate buffer and sodium chloride. No preservatives are added. Sandwich partner for MAB<TSH>M-TU1.20 IgG.

Specification

Appearance: White lyophilizate

Solubility: Clear to slightly opalescent solution in NaCl, 0.9% (c=10 mg/mL)

Protein (Biuret): ≥0.6 mg/mg lyophilizate

Purity (HPLC / Mono Q): ≥90 area%

pH 5.5 treatment (30 minutes): Corresponds to specification

Stability: At -15 to -25°C within specification range for 24 months.

Avoid repeated freezing and thawing.

Catalog number

11 367 978 103

Pack size

5 mg (samples), ≥50 mg (custom fill)

Will be supplied as "MAK<TSH>M-A8-IGG(BR)". Unit of measure is "mg active ingredient".

For further processing only.

MAB<TSH>M-TU1.20 IgG

lyophilizate

Qualified for **cobas**[®] platforms

Application

MAB<TSH>M-Tu1.20 IgG is qualified for heterogeneous immunoassays (HetIAs).

Catalog number

10 767 778 103

Pack size

5 mg (samples), ≥50 mg (custom fill)

Will be supplied as "MAK<TSH>M-TU1.20-IGG(DE),LYO.". Unit of measure is "mg active ingredient".

For further processing only.

Product description

Immunogen: Human thyroid stimulating hormone
Spleen donor: Mouse Balb/c
Antibody class: IgG1, kappa
Cross reactivity to: Luteinizing hormone (LH) <0.2%; Follicle stimulating hormone (FSH) <3.0%.

Properties

MAB<TSH>M-Tu1.20 IgG is a monoclonal antibody directed to thyroid stimulating hormone. It is lyophilized from a solution containing protein, potassium phosphate buffer and sodium chloride. No preservatives are added. Sandwich partner for MAB<TSH>M-A8 IgG.

Specification

Appearance: White lyophilizate
Solubility: Clear to slightly opalescent solution in NaCl, 0.9% (c=10 mg/mL)
Protein (Biuret): ≥0.6 mg/mg lyophilizate
Purity (HPLC / Mono Q): ≥90 area%
pH 5.5 treatment (30 minutes): Corresponds to specification
Stability: At -15 to -25°C within specification range for 24 months.
 Avoid repeated freezing and thawing.

MAB<Dig>M-19-11 IgG

powder

Qualified for **cobas®** platforms

Application

MAB<Dig>M-19-11 IgG is qualified for heterogeneous immunoassays (HetIAs).

Product description

Immunogen: Digoxin coupled to carrier
Spleen donor: Mouse, Balb/c
Antibody class: IgG1, kappa

Properties

MAB<Dig>M-19-11 IgG is a monoclonal antibody directed to Digoxigenin. It is lyophilized from a solution containing protein, potassium phosphate buffer and sodium chloride. No preservatives are added. Sandwich partner for MAB<TSH>M-A8 IgG.

Catalog number	Pack size
05 910 714 103	5 mg (samples), ≥50 mg (custom fill)

Will be supplied as "MAB<Dig>M-19-11-IgG *SQ". Unit of measure is "mg active ingredient".

For further processing only.

Specification**Appearance:** White powder**Protein** (OD₂₈₀): ±20% (based on manufacturer's data)**Aggregates** (HPLC; TSK 3000): ≤10%**pH 5.5 treatment** (30 minutes): Corresponds to specification**Stability:** At -15 to -25°C within specification range for 36 months.**MAB<C-Peptid>M-RPEP001 IgG**

frozen solution

Qualified for **cobas**[®] platforms**Application**

MAB<C-Peptid>M-RPEP001 IgG is qualified for heterogeneous immunoassays (HetIAs).

Product description**Immunogen:** Human C-Peptide**Spleen donor:** Mouse, Balb/c**Properties**

MAB<C-Peptid>M-RPEP001 IgG is a monoclonal antibody directed to human C-Peptid. It is lyophilized from a solution containing protein, potassium phosphate buffer and sodium chloride. No preservatives are added. Sandwich partner for MAB<C-Peptid>M-RCPT3F11 IgG.

Specification**Appearance:** Clear, colorless solution**Storage buffer:** 0.01mol/L phosphate buffer, 0.15 mol/L NaCl, pH 7.2**Protein** (OD₂₈₀): 3.0-6.0 mg/mL**Purity** (HPLC / TSK 3000): ≥90 area%**pH-5.5 treatment** (30 minutes): Corresponds to specification**Stability:** At -60 to -90°C within specification range for 24 months.

Avoid repeated freezing and thawing.

Catalog number**03 732 584 103****Pack size**5 mg, 10 mg (samples), ≥50 mg
(custom fill)

Will be supplied as "MAB<C-Peptid>M-RPEP001-IgG". Unit of measure is "mg active ingredient".



For further processing only.

MAB<C-Peptid>M-RCPT3F11 IgG

frozen solution

Qualified for **cobas**[®] platforms**Application**

MAB<C-Peptid>M-RCPT3F11 IgG is qualified for heterogeneous immunoassays (HetIAs).

Catalog number**03 732 606 103****Pack size**5 mg, 10 mg (samples), ≥50 mg
(custom fill)

Will be supplied as "MAB<C-Peptid>M-RCPT3F11-IgG". Unit of measure is "mg active ingredient".



For further processing only.

Product description**Immunogen:** Human C-Peptide**Spleen donor:** Mouse, Balb/c**Properties**

MAB<C-Peptid>M-RCPT3F11 IgG is a monoclonal antibody directed to human C-Peptid. It is lyophilized from a solution containing protein, potassium phosphate buffer and sodium chloride. No preservatives are added. Sandwich partner for MAB<C-Peptid>M-RPEP001 IgG.

Specification**Appearance:** Clear, colorless solution**Storage buffer:** 0.01mol/L phosphate buffer, 0.15 mol/L NaCl, pH 7.2**Protein** (OD₂₈₀): 3.0-6.0 mg/mL**Purity** (HPLC / TSK 3000): ≥90 area%**pH-5.5 treatment** (30 minutes): Corresponds to specification**Stability:** At -60 to -90°C within specification range for 24 months.

Avoid repeated freezing and thawing.

MAB<HCT>M-11.3.59 IgG

solution

Qualified for **cobas**[®] platforms**Application**

MAB<HCT>M-11.3.59 IgG is qualified for heterogeneous immunoassays (HetIAs).

Product description**Immunogen:** Peptide derived from human calcitonin, coupled to carrier protein**Spleen donor:** Mouse, Balb/c**Properties**

MAB<HCT>M-11.3.59 IgG is a monoclonal antibody directed to human Calcitonin. It is lyophilized from a solution containing protein, potassium phosphate buffer and sodium chloride. No preservatives are added.

Specification**Appearance:** Clear, colorless solution**Storage buffer:** 50 mM potassium phosphate, 0,15 M NaCl, pH 8,5**Protein** (OD₂₈₀): ≥10 mg/mL**Purity** (HPLC / TSK 3000): ≥90 area%**pH-5.5 treatment** (30 minutes): Corresponds to specification**Stability:** At -60 to -90°C within specification range for 24 months.

Avoid repeated freezing and thawing.

Catalog number**06 548 776 103****Pack size**

5 mg, 10 mg (samples), ≥50 mg (custom fill)

Will be supplied as "MAB<HCT>M-11.3.59-IgG". Unit of measure is "mg active ingredient".



For further processing only.

MAB<Digit>M-145A41 IgG

solution

Qualified for **cobas**[®] platforms

Application

MAB<Digit>M-145A41 IgG is qualified for heterogeneous immunoassays (HetIAs).

Product description

Immunogen: Digitoxin, coupled to carrier protein

Spleen donor: Mouse, Balb/c

Antibody class: IgG1, kappa

Properties

MAB<Digit>M-145A41 IgG is a monoclonal antibody directed to Digitoxin. It is lyophilized from a solution containing protein, potassium phosphate buffer and sodium chloride.

Related product: Alpha-Biotin-Epsilon-Dig-Lys hapten.

Specification

Appearance: Clear, yellowish solution

Storage buffer: potassium phosphate, 10 mM; NaCl, 70 mol/L; Tween, 0.2%; sodium azide, 0.09%; pH 7.5

Protein (OD₂₈₀): ≥10 mg/mL

Purity (HPLC / Mono Q): ≥90% of total protein

pH-5.5 treatment (30 minutes): Corresponds to specification

Stability: At -60 to -90°C within specification range for 24 months.

Avoid repeated freezing and thawing.

Catalog number**12 041 740 103****Pack size**

5 mg, 10 mg (samples), ≥50 mg (custom fill)

Will be supplied as "MAB<Digit>M-145A41-IgG". Unit of measure is "mg active ingredient".



For further processing only.

MAB<Keratin>M-AE3 IgG

frozen solution

Application

Qualified for immunohistochemistry

Product description

Spleen donor: Mouse, Balb/c

Antibody class: IgG1, kappa

Properties

MAB<Keratin>M-AE3 IgG is a monoclonal antibody directed to human keratin AE3. It is lyophilized from a solution containing protein, potassium phosphate buffer and sodium chloride. Related product: MAB<Keratin>M-AE1 IgG.

Catalog number**03 808 513 103****Pack size**

1 mg, 5 mg, 10 mg, 50 mg, 100 mg

Will be supplied as "MAB<Keratin>M-AE3-IgG". Unit of measure is "mg ingredient".



For further processing only.

Specification

Appearance: Clear to opaque solution
pH value: 7.4-7.6
Protein (OD₂₈₀): 9.0-11.0 mg/mL
HPLC, MONOQ: ≥90%
Function test (IC 50%, ≥0 µg/mL): Corresponds to specification
pH 5.5 treatment (30 minutes): Corresponds to specification
Stability: At -15 to -25°C within specification range for 24 months.

MAB<Keratin>M-AE1 IgG

frozen solution

Application

Qualified for immunohistochemistry

Product description

Spleen donor: Mouse, Balb/c
Antibody class: IgG1, kappa

Properties

MAB<Keratin>M-AE1 IgG is a monoclonal antibody directed to human Keratin AE1. It is lyophilized from a solution containing protein, potassium phosphate buffer and sodium chloride. Related product: MAB<Keratin>M-AE3 IgG.

Specification

Appearance: Clear to opaque solution
pH value: 7.4-7.6
Protein (OD₂₈₀): 9.0-11.0 mg/mL
HPLC, MONOQ: ≥90%
Function test (IC 50%, ≥0 µg/mL): Corresponds to specification
pH 5.5 treatment (30 minutes): Corresponds to specification
Stability: At -15 to -25°C within specification range for 24 months.

Catalog number

03 808 521 103

Pack size

1 mg, 5 mg, 10 mg, 50 mg, 100 mg

Will be supplied as "MAB<Keratin>M-AE1-IgG". Unit of measure is "mg ingredient".



For further processing only.

PAB<E2>K IgG

frozen solution

Qualified for **cobas**[®] platforms**Application**

PAB<E2>K IgG is qualified for heterogeneous immunoassays (HetIAs).

Product description**Spleen donor:** Mouse, Balb/c**Antibody class:** IgG1, kappa**Properties**

PAB<E2>K IgG is a polyclonal antibody directed to estradiol produced in rabbits. No preservative added.

Specification**Appearance:** Frozen solution**Appearance of solution:** Clear to slightly opalescent, colorless**Storage buffer:** Potassium phosphate buffer, pH 8**Protein** (OD₂₈₀): ≥30.0 mg/mL**Purity** (HPLC / TSK3000): ≥90% of total protein**Country of origin:** Germany**pH 5.5 treatment** (30 minutes): Corresponds to specification**Stability:** At -60 to -90°C within specification range for 24 months.

Avoid multiple freeze/thaw cycles.

Catalog number**05 344 255 103****Pack size**

10 mg (samples), ≥50 mg (custom fill)

Will be supplied as "PAK<E2>K-IgG(DE)". Unit of measure is "mg active ingredient".



For further processing only.

PAB<CRP>S IgG

frozen solution

For measurement of c-reactive protein (CRP).

Application

Use PAB<CRP>S IgG for turbidimetric / nephelometric assays. It can be coupled to latex surfaces and used in the respective assays.

Product description**Immunogen:** Human C-reactive protein**Related product:** MAB<Digit>M-145A41 IgG**Properties**The polyclonal antibody IgG directed to c-reactive protein is produced in sheep. It is prepared as solution containing protein (≥40 g/L); Tris buffer with NaN₃, 0.09% (w/v).**Recommended working concentration:** 15 mg/mL**Catalog number****11 888 714 103****Pack size**

0.010 L, 0.100 L, 0.5 L, 1 L

Will be supplied as "PAK<CRP>S-IgG *SQ". Unit of measure is "L".



For further processing only.

Remark: When stored over longer periods at +4°C, a slight turbidity may occur which can easily be removed by centrifugation. No alteration of antibody properties occurs thereby.

Specification

Appearance: Clear to slightly opalescent yellowish solution

pH value (+25°C): 7.8-8.2

Protein (Biuret): 50-60 mg/mL

Function: Calibration curve characterization defined by turbidimetric measurement A_{340}/A_{700}

δA (Standard 1mg/dL): ≥ 0.038

δA (Standard 10mg/dL): ≥ 0.228

δA (Standard 25mg/dL): ≥ 0.370

pH 5.5 treatment (30 minutes): Corresponds to specification

Stability: At -15 to -25°C within specification range for 24 months.

Avoid repeated freezing and thawing.

Background information

C-reactive protein (CRP) is part of the β -globulin family found in human plasma. Increased levels of CRP are involved in a variety of inflammatory diseases. Furthermore elevated CRP serum levels indicate tissue injury, transplant rejection, carcinogenesis and acute myocardial infarction.

PAB<T3>S IgG

lyophilizate

Qualified for **cobas**[®] platforms

Application

PAB<T3>S IgG is qualified for heterogeneous immunoassays (HetIAs).

Product description

PAB<T3>S IgG is a polyclonal antibody directed to triiodothyronine produced in sheep.

Immunogen: Triiodothyronine derivative

Specification

Appearance: White lyophilizate

Protein (OD_{280}): ≥ 0.7 mg/mg lyophilizate

HPLC (HPLC / TSK 3000): ≥ 90 area%

pH 5.5 treatment (30 minutes): Corresponds to specification

Stability: At -15 to -25°C within specification range for 24 months.

Catalog number

10 907 332 103

Pack size

10 mg (samples), ≥ 50 mg (custom fill)

Will be supplied as "PAK<T3>S-IGG(DE), (ES 3G)". Unit of measure is "g active ingredient".

For further processing only.

PAB<T4>S IgG

lyophilizate

Qualified for **cobas**[®] platforms**Application**

PAB<T4>S IgG is qualified for heterogeneous immunoassays (HetIAs).

Product description

PAB<T4>S IgG is a polyclonal antibody directed to triiodothyronine produced in sheep.

Immunogen: Triiodothyronine derivative**Specification****Appearance:** White lyophilizate**Protein** (OD₂₈₀): ≥0.7 mg/mg lyophilizate**HPLC** (HPLC /TSK 3000): ≥90 area%**pH 5.5 treatment** (30 minutes): Corresponds to specification**Stability:** At -15 to -25°C within specification range for 24 months.**Catalog number****10 767 794 103****Pack size**

10 mg (samples), ≥50 mg (custom fill)

Will be supplied as "PAK<T4>S-IGG(DE),(ES)". Unit of measure is "mg active ingredient".

For further processing only.

Alpha-Biotin-Epsilon-Dig-Lys

powder

Application

Alpha-Biotin-Epsilon-Dig-Lys is qualified for heterogenous immunoassays (HetIAs).

Properties

Related product: MAB<Dig>M-145A41 IgG

Specification

Appearance: Colorless powder

Purity (HPLC): ≥95 area%

Stability: At -60 to -90°C within specification range for 12 months.

Catalog number**11 831 658 103****Pack size**

0.5 mg, custom fill

Will be supplied as "alpha-Biotin-epsilon-Dig-Lys". Unit of measure is "mg ingredient".



For further processing only.

CA 15-3 Ag

liquid

Application

For development and manufacturing of target-specific assay reagents.

Properties

Highly purified native breast cancer (CA 15-3) antigen derived from human fluids.

Specification

Identification (ELISA Calbiotech): Positive

CA 15-3 content (Roche Elecsys®): >10 kU/mL

Purity Ratio (Roche Elecsys®): Activity/ml/OD280 nm

Co product 19-9: <25%

Co product 125: <25%

Co product AFP: <25%

Co product 72-4: <25%

Co Product Ferritin: <25%

Co product CEA: <25%

Bioburden: <10 CFU/mL

Viral Marker Test (ELISA) for HBsAg, Anti-HCV and Anti-HIV I + II (US-FDA approved BioRad & Ortho Clinical Diagnostics): Negative

Absence of viral DNA/RNA (HBV, HIV I and HCV; PCR): Negative

Stability: At +2 to +8°C within specification range for 10 months.

Catalog number**07 190 999 001****Pack size**

custom fill

Will be supplied as "CA 15-3 Ag". Unit of measure is "kU".

For further processing only.

Quality

Extensively tested for absence of HBsAg, Anti-HCV and Anti-HIV I + II (ELISA), and for viral DNA/RNA (PCR: HBV, HIV I and HCV) by applying FDA-approved methods.

Background information

Liquid (aseptically filtered) in 50 mM PBS (pH 7.2 -7.6) and 0.09% w/v sodium azide

CA 19-9 Ag

liquid

Application

For development and manufacturing of target-specific assay reagents.

Properties

Highly purified native ovarian cancer (CA 125) antigen derived from human fluids.

Specification

Identification (ELISA Calbiotech): Positive

CA 19-9 content (Roche Elecsys®): >10 kU/mL

Purity Ratio (Roche Elecsys®): Activity/ml/OD280 nm

Co product 15-3: <25%

Co product 125: <25%

Co product AFP: <25%

Co product 72-4: <25%

Co Product Ferritin: <25%

Co product CEA: <25%

Bioburden: <10 CFU/mL

Viral Marker Test (ELISA) for HBsAg, Anti-HCV and Anti-HIV I + II (US-FDA approved BioRad & Ortho Clinical Diagnostics): Negative

Absence of viral DNA/RNA for (HBV, HIV I and HCV) (PCR): Negative

Stability: At +2 to +8°C within specification range for 10 months.

Quality

Extensively tested for absence of HBsAg, Anti-HCV and Anti-HIV I + II (ELISA), and for viral DNA/RNA (PCR: HBV, HIV I and HCV) by applying FDA-approved methods.

Background information

Liquid (aseptically filtered) in 50 mM PBS (pH 7.2 -7.6) and 0.09% w/v sodium azide

Catalog number

07 191 006 001

Pack size

custom fill

Will be supplied as "CA 19-9 Ag". Unit of measure is "kU".

For further processing only.

CA 72-4 Ag

liquid

Application

For development and manufacturing of target-specific assay reagents.

Properties

Highly purified native gastric cancer (CA 72-4) antigen derived from human fluids.

Specification

Identification (ELISA): Positive

CA 72-4 content (Roche Elecsys®): >10 kU/mL

Purity Ratio (Roche Elecsys®): Activity/ml/OD280 nm

Co product 15-3: <25%

Co product 125: <25%

Co product AFP: <25%

Co product 19-9: <25%

Co Product Ferritin: <25%

Co product CEA: <25%

Bioburden: <10 CFU/mL

Viral Marker Test (ELISA) for HBsAg, Anti-HCV and Anti-HIV I + II (US-FDA approved BioRad & Ortho Clinical Diagnostics): Negative

Absence of viral DNA/RNA for (HBV, HIV I and HCV; PCR): Negative

Stability: At +2 to +8°C within specification range for 10 months.

Quality

Extensively tested for absence of HBsAg, Anti-HCV and Anti-HIV I + II (ELISA), and for viral DNA/RNA (PCR: HBV, HIV I and HCV) by applying FDA-approved methods.

Background information

Liquid in salt buffer (50 mM phosphate, 0.15 M NaCl; pH 7.2-7.6), preservative 0.09% w/v sodium azide and 0.2% Proclin 300.

Catalog number

07 191 014 001

Pack size

custom fill

Will be supplied as "CA 72-4 Ag". Unit of measure is "kU".

For further processing only.

CA 125 Ag

liquid

Application

For development and manufacturing of target-specific assay reagents.

Properties

Highly purified native carcinoembryonic (CEA) antigen derived from cancer fluids.

Catalog number

07 190 930 001

Pack size

custom fill

Will be supplied as "CA 125 Ag". Unit of measure is "kU".

For further processing only.

Specification**Identification** (ELISA Calbiotech): Positive**CA 125 content** (Roche Elecsys®): >10 kU/mL**Purity Ratio** (Roche Elecsys®): Activity/mL/OD280 nm

Co product 19-9: <25%

Co product 15-3: <25%

Co product AFP: <25%

Co product 72-4: <25%

Co Product Ferritin: <25%

Co product CEA: <25%

Bioburden: <10 CFU/mL**Viral Marker Test** (ELISA) for HBsAg, Anti-HCV and Anti-HIV I + II (US-FDA approved BioRad & Ortho Clinical Diagnostics): Negative**Absence of viral DNA/RNA** (HBV, HIV I and HCV; PCR): Negative**Stability:** At +2 to +8°C within specification range for 10 months.**Quality**

Extensively tested for absence of HBsAg, Anti-HCV and Anti-HIV I + II (ELISA), and for viral DNA/RNA (PCR: HBV, HIV I and HCV) by applying FDA-approved methods.

Background information

Liquid (aseptically filtered) in 50 mM phosphate, 0.15 M NaCl (pH 7.2-7.6) and 0.09% w/v sodium azide

CEA Ag

liquid

Application

For development and manufacturing of target-specific assay reagents.

Product description

Antigen derived from cancer fluids.

Properties

Highly purified native carcinoembryonic (CEA) antigen derived from cancer fluids.

Related products: MAB<CEA>M-TU3 IgG, MAB<CEA>M-TU2 IgG**Specification****Identification** (ELISA): Positive**CEA content** (Roche **cobas**® e411): 1-2 mg/mL**Purity** (SDS PAGE): >98%**Purity profile** (SDS PAGE): Single band at 180 kD**Bioburden:** <10 CFU/mL**Catalog number****07 190 948 001****Pack size**

custom fill

Will be supplied as "CEA Ag". Unit of measure is "mg".

For further processing only.

Viral Marker Test (ELISA) for HBsAg, Anti-HCV and Anti-HIV I + II (US-FDA approved BioRad & Ortho Clinical Diagnostics): Negative
Absence of viral DNA/RNA for (HBV, HIV I and HCV; PCR): Negative
Stability: At +2 to +8°C within specification range for 10 months.

Quality

Extensively tested for absence of HBsAg, Anti-HCV and Anti-HIV I + II (ELISA), and for viral DNA/RNA (PCR: HBV, HIV I and HCV) by applying FDA-approved methods.

Background information

Liquid (aseptically filtered) 50 mM sodium citrate and 150 mM NaCl (pH 7.2-7.6) containing 20% glycerol and 0.09% w/v sodium azide

AFP Ag

liquid

Application

For development and manufacturing of target-specific assay reagents.

Properties

Highly purified native α feto protein (AFP) antigen derived from human cord blood.

Related products: MAB<AFP>M-TU11 IgG, MAB<AFP>M-LJ738 IgG

Specification

Identification (ELISA): Positive

AFP content (Roche Elecsys®): 1-2 mg/mL

Total protein (OD₂₈₀): >AFP content

Purity (Activity/TP): ≥98%

Purity profile (SDS PAGE): Single band at 66 kD

Bioburden: <10 CFU/mL

Viral Marker Test (ELISA) for HBsAg, Anti-HCV and Anti-HIV I + II (US-FDA approved BioRad & Ortho Clinical Diagnostics): Negative

Absence of viral DNA/RNA for (HBV, HIV I and HCV; PCR): Negative

Stability: At +2 to +8°C within specification range for 10 months.

Quality

Extensively tested for absence of HBsAg, Anti-HCV and Anti-HIV I + II (ELISA), and for viral DNA/RNA (PCR: HBV, HIV I and HCV) by applying FDA-approved methods.

Background information

Liquid (aseptically filtered) in 50 mM Tris and 150 mM NaCl, (pH 7.2±0.2) containing 0.09% w/v sodium azide

Catalog number

07 190 905 001

Pack size

custom fill

Will be supplied as "AFP Ag". Unit of measure is "mg".

For further processing only.

PSA Ag

liquid

Application

For development and manufacturing of target-specific assay reagents.

Product description

Antigen derived from human seminal plasma

Properties

Highly purified native prostate specific (PSA) antigen derived from human seminal plasma.

Specification

Identification (ECLIA): Positive

PSA activity (ECLIA, Roche Elecsys®): >0.5 mg/ml

Purity (SDS PAGE): >96%, cluster of bands at 29 kD

Bioburden: <10 CFU/ml

Viral Marker Test (ELISA) for HBsAg, Anti-HCV and Anti-HIV I + II (US-FDA approved BioRad & Ortho Clinical Diagnostics): Negative

Absence of viral DNA/RNA for (HBV, HIV I and HCV; PCR): Negative

Stability: At +2 to +8°C within specification range for 10 months.

Quality

Extensively tested for absence of HBsAg, Anti-HCV and Anti-HIV I + II (ELISA), and for viral DNA/RNA (PCR: HBV, HIV I and HCV) by applying FDA-approved methods.

Background information

Liquid (aseptically filtered) in 0.1 M Tris acetate (pH 6 ±0.2) containing 20% glycerol and 0.09% w/v sodium azide

Catalog number

07 190 956 001

Pack size

custom fill

Will be supplied as "PSA Ag". Unit of measure is "mg".

For further processing only.

CRP Ag

liquid

Application

For development and manufacturing of target-specific assay reagents.

Properties

Highly purified native C reactive protein (CRP) antigen derived from human fluids.

Specification

Identification (Western Blot): Positive

CRP content (Turbilatex): 2-3 mg/mL (Spinreact)

Total protein (OD₂₈₀): >CRP content

Catalog number

07 191 022 001

Pack size

custom fill

Will be supplied as "CRP Ag". Unit of measure is "mg".

For further processing only.

Purity (CRP content/TP): $\geq 99\%$

Purity profile (SDS PAGE): Single band at 21 kD

Bioburden: < 10 CFU/mL

Viral Marker Test (ELISA) for HBsAg, Anti-HCV and Anti-HIV I + II

(US-FDA approved BioRad & Ortho Clinical Diagnostics): Negative

Absence of viral DNA/RNA for (HBV, HIV I and HCV; PCR): Negative

Stability: At +2 to +8°C within specification range for 10 months.

Quality

Extensively tested for absence of HBsAg, Anti-HCV and Anti-HIV I + II (ELISA), and for viral DNA/RNA (PCR: HBV, HIV I and HCV) by applying FDA-approved methods.

Background information

Liquid (aseptically filtered) in 20 mM Tris and 280 mM NaCl, 5 mM CaCl₂ (pH 7.9-8.1) containing 0.09%w/v sodium azide

B2M Ag

liquid

Application

For development and manufacturing of target-specific assay reagents.

Properties

Highly purified native B2M antigen derived from human urine.

Specification

Identification (Western Blot): Positive

B2M content (Turbilatex): 1-2 mg/mL (Spinreact)

Total protein (OD₂₈₀): As observed

Purity (B2M content/TP): $> 98\%$

Purity profile (SDS PAGE): Single band at 12 kD

Bioburden: < 10 CFU/mL

Viral Marker Test (ELISA) for HBsAg, Anti-HCV and Anti-HIV I + II

(US-FDA approved BioRad & Ortho Clinical Diagnostics): Negative

Absence of viral DNA/RNA for (HBV, HIV I and HCV; PCR): Negative

Stability: At +2 to +8°C within specification range for 10 months.

Quality

Extensively tested for absence of HBsAg, Anti-HCV and Anti-HIV I + II (ELISA), and for viral DNA/RNA (PCR: HBV, HIV I and HCV) by applying FDA-approved methods.

Background information

Liquid (aseptically filtered) 40 mM PBS, 0.15 M NaCl, 1 mM EDTA, sodium salt, (pH 7.2 \pm 0.2) containing 0.09% w/v sodium azide

Catalog number

07 191 073 001

Pack size

custom fill

Will be supplied as "B2M Ag". Unit of measure is "mg".

For further processing only.

HBsAg (ad)

liquid

Application

For development and manufacturing of target-specific assay reagents.

Product description

Antigen derived from human blood.

Properties

Highly purified native hepatitis B surface -ad type (HBsAg) (ad) antigen derived from human blood.

Related product: HBsAg (ay)

Specification

Identification (ELISA): Positive

HBsAg content (BioRad, against PEI standard): 2-4 mg/mL

Total protein (OD_{280-320nm}): 2-4 mg/mL

'Ad' subtype purity (ELISA) >99%

Purity profile (SDS PAGE): Corresponds to internal reference standard

Aggregation (gel filtration): <9%

Bioburden: <10 CFU/mL

Viral Marker Test (ELISA) for HBsAg, Anti-HCV and Anti-HIV I + II

(US-FDA approved BioRad & Ortho Clinical Diagnostics): Negative

Absence of viral DNA/RNA for (HBV, HIV I and HCV; PCR): Negative

Stability: At +2 to +8°C within specification range for 10 months.

Quality

Extensively tested for absence of HBsAg, Anti-HCV and Anti-HIV I + II (ELISA), and for viral DNA/RNA (PCR: HBV, HIV I and HCV) by applying FDA-approved methods.

Background information

Liquid (aseptically filtered) in 50 mM phosphate, 150 mM NaCl (pH 7.2-7.6) containing 0.09% w/v sodium azide.

Heat inactivated at +60°C for 15 hours.

Catalog number

07 191 057 001

Pack size

custom fill

Will be supplied as "HBsAg (ad)". Unit of measure is "mg".

For further processing only.

HBsAg (ay)

liquid

Application

For development and manufacturing of target-specific assay reagents.

Properties

Highly purified native hepatitis B surface -ay type (HBsAg) (ay) antigen derived from human blood.

Related product: HBsAg (ad)

Specification

Identification (ELISA): Positive

HBsAg content (BioRad, against PEI standard): 2-4 mg/mL

Total protein (OD_{280-320nm}): 2-4 mg/mL

'Ay' subtype purity (ELISA): >99%

Aggregation (gel filtration): <9%

Purity profile (SDS PAGE): Corresponds to internal reference standard

Bioburden: <10 CFU/mL

Viral Marker Test (ELISA) for HBsAg, Anti-HCV and Anti-HIV I + II

(US-FDA approved BioRad & Ortho Clinical Diagnostics): Negative

Absence of viral DNA/RNA for (HBV, HIV I and HCV) (PCR): Negative

Stability: At +2 to +8°C within specification range for 10 months.

Quality

Extensively tested for absence of HBsAg, Anti-HCV and Anti-HIV I + II (ELISA), and for viral DNA/RNA (PCR: HBV, HIV I and HCV) by applying FDA-approved methods.

Background information

Liquid (aseptically filtered) in 50 mM phosphate, 150 mM NaCl (pH 7.2-7.6) containing 0.09% w/v sodium azide.

Heat inactivated at +60°C for 15 hours.

Catalog number

07 191 065 001

Pack size

custom fill

Will be supplied as "HBsAg (ay)". Unit of measure is "mg".

For further processing only.

Streptavidin, recombinant

from *Streptomyces avidinii*, expressed in *E. coli*, lyophilizate

Application

Use Streptavidin, recombinant as a tool for solid phase technology and universal detection systems in immunology and molecular diagnostics.

Specification

Appearance: White lyophilizate

Protein (A_{282} ; factor 3.1): 0.6-0.8 mg/mg lyophilizate

Specific activity/Biotin binding capacity: ≥ 17 U/mg protein

Proteasen (incubation with Azocoll for up to 24 hours at +25°C): ≤ 0.001 U/mg lyophilizate

Absorption (A_{405} , against repurified water): ≤ 0.1

Water (K. Fischer): $\leq 12\%$

IEF (pH 6-9): Two main bands between 6.8 and 7.5

SDS-PAGE: Chromatographically homogeneous

Stability: At +2 to +8°C within specification range for 24 months.

Background information

Streptavidin consists of four subunits with a molecular weight of 13 kD, each containing a single biotin binding site. Each subunit has six tyrosine residues. The protein is carbohydrate free.

Streptavidin + 4 biotin \rightarrow streptavidin · (biotin)₄

The formation of the complex is measured at 233 nm.

Catalog number

11 520 679 103

Pack size

custom fill

Will be supplied as "Streptavidin Special Quality". Unit of measure is "g active ingredient".

For further processing only.

Streptavidin-Mutein-Sepharose

suspension

Sepharose-Mutein-Sepharose provides a recombinant streptavidin mutein with reduced binding affinity to biotin, immobilized onto cross-linked agarose beads in a highly stable form.

Application

Streptavidin-Mutein-Sepharose allows purification of biotinylated antibodies and other proteins, resulting in excellent purity and recovery.

Properties

Gel suspension (50%), ready to use, bulk format.

Specification

Appearance: White suspension

Binding capacity of Fab-Bi: ≥ 4.5 mg/mL

Shipping category: Cool packs

Stability: At+2 to +8°C within specification range 36 months.

Catalog number

03 314 731 103

Pack size

0.010 L (samples), 0.2 L (bulk)

Will be supplied as "SA-MUTEIN-SEPHAROSE". Unit of measure is "L".

For further processing only.

Streptavidin POD

frozen solution

Application

Universal reagent for the detection of biotinylated compounds.

Specification**Appearance:** Clear, brownish solution**POD activity** (+25°C, ABTS, pH 5.0): 2,000-8,000 U/mL**Performance test** (recovery Streptavidin POD): Corresponds to reference**Stability:** At -60 to -90°C within specification range for 24 months.**Catalog number****11 096 346 103****Pack size**

custom fill

Will be supplied as "Streptavidin-POD Conjugate". Unit of measure is "kU".



DRY ICE

For further processing only.

D-Biotin-N-hydroxysuccinimide ester

crystalline powder

Application

Use D-Biotin-N-hydroxysuccinimide ester as biotinylating reagent for proteins and amino labeled oligonucleotides.

CAS: 35013-72-0**Properties****Nomenclature:** D-Biotinyl-N-hydroxy-succinimide ester**Formula:** C₁₄H₁₉N₃O₅S**Molecular weight:** 341.4 D

Remark: Under mild conditions the activated ester reacts with amino groups. The aminocaproic acid spacer is useful if biotinylated macromolecules are coupled because steric hindrance is minimized. Gel suspension (50%), ready to use, bulk format.

Specification**Appearance:** White crystallize**Biotin ester** (from N): 98-103%**N** (elementary analysis): 12.0-12.73%

Purity (TLC: silica gel, 1-butanol/glacial acetic acid/H₂O= 2/1/1; a) in UV; b) with 4-Dimethylamino cinnamaldehyde): Chromatographically homogeneous

Stability: At +2 to +8°C within specification range for 24 months.**Catalog number****10 734 250 103****Pack size**

custom fill

Will be supplied as "D-Biotin-N-hydroxy-succinimide Ester". Unit of measure is "g".

For further processing only.

D-Biotinoyl-ε-aminocaproic acid-N-hydroxysuccinimide ester

powder

Application

Use D-Biotinoyl-ε-aminocaproic acid-N-hydroxysuccinimide ester as biotinylating reagent for proteins and amino labeled oligonucleotides.

CAS: 72040-63-2**Properties****Nomenclature:** D-Biotinyl-ε-amido caproic acid N-hydroxysuccinimide ester**Formula:** C₂₀H₃₀N₄O₆S**Molecular weight:** 454.5 D

Remark: Under mild conditions the activated ester reacts with amino groups. The aminocaproic acid spacer is useful if biotinylated macromolecules are coupled because steric hindrance is minimized.

Catalog number**11 003 933 103****Pack size**

custom fill

Will be supplied as "D-Biotinyl-ε-aminocap. Acid-N-Hydroxy Succ.". Unit of measure is "g".



For further processing only.

Specification

Appearance: White to beige powder

Biotin ester (from N): ≥ 97 -103%

N (elementary analysis): 11.9-12.8%

C (elementary analysis): 51.0-54.7%

H (elementary analysis): 6.4-6.9%

Purity (TLC: silica gel, 1-butanol/glacial acetic acid/H₂O= 50/15/25, iodide stream/UV): Chromatographically homogeneous

Hydrolysis product (NMR): ≤ 20 %

Stability: At -15 to -25°C within specification range for 24 months.

Streptavidin R-Phycoerythrin LumiGrade Reagent**solution**

Standard for highly sensitive fluorescent detection.

Application

Conjugated reporter dyes such as Streptavidin R-Phycoerythrin (SA-PE) are well established for Luminex's xMAP Assay Kits or array-based applications due to their excellent spectral characteristics. Their dedicated and reproducible design ensure high performance multiplex assays using antibodies, receptors, peptides and oligonucleotides.

Product description

The Streptavidin R-Phycoerythrin LumiGrade conjugate contains Phycoerythrin from red algae (RPE) and Streptavidin (SA), recombinant.

Properties

High molecular weight conjugate size distribution (700-1200 kD).

Specification

Appearance: Reddish solution

Streptavidin (A_{280}): 16-21 weight%

Absorption ratio $A_{566/280}$: >3.3

Protein: 1.00 ± 0.10 mg/ml

Content of color relating to SA: 0.95-1.40 (molar ratio RPE : SA)

Purity (HPLC/ TSK 3000 XL): free SA <1%

Fluorescence emission: Maximum / intensity

488 nm excitation: 576 nm ± 5 / ≥200

545 nm excitation : 576 nm ± 5 / ≥250

Stability: At +2 to +8°C within specification range for 36 months.

Catalog number**05 065 925 103****Pack size**

1 mL, 5 mL, 100 mL, custom fill

Will be supplied as "Streptavidin R-Phycoerythrin Lumi Grade Reagent". Unit of measure is "mg active ingredient".

For further processing only.

Streptavidin R-Phycoerythrin LumiGrade Ultrasensitive Reagent**solution**

Standard for ultrasensitive fluorescent detection.

Application

Conjugated reporter dyes such as Streptavidin R-Phycoerythrin (SA-PE) are well established for Luminex's xMAP® Assay Kits and array-based applications due to their excellent spectral characteristics. Their dedicated and reproducible design ensure high performance multiplex assays using antibodies, receptors, peptides and oligonucleotides.

Catalog number**05 351 693 103****Pack size**

1 mL, 100 mL, custom fill

Will be supplied as "Streptavidin R-Phycoerythrin Lumi Grade ultrasensitive". Unit of measure is "mg active ingredient".

For further processing only.

Product description

The Streptavidin R-Phycoerythrin LumiGrade conjugate contains Phycoerythrin (PE) from red algae and Streptavidin (SA), recombinant.

Properties

High molecular weight conjugate size distribution (1500-50000 kD).

Specification

Appearance: Reddish solution

A₅₆₆: 0.73-0.81

SA-R-PE from R-Phycoerythrin (A₅₆₆/7.7): 0.95-1.05 mg/mL

Purity (HPLC / TSK 6000): ≥ 99.7 area%

Contamination (HPLC / TSK 6000): ≤ 0.3 area%

Stability: At +2 to +8°C within specification range for 24 months.

Multi Analyte Stripe

universal device

Universal test stripe

Application

Use the Multi Analyte Stripes for a range of applications, such as classical sandwich immunoassays to NA/Oligo-detection.

Product description

Test stripes are delivered in boxes, each containing 50 stripes.

Specification

Components of the test stripe (4.6 nm):

MAB<Dig>IgG on gold conjugate: 0.2 µg per test stripe

Poly-Streptavidin (result line): 0.8 µg per test stripe

PAB<MouseFc>IgG (control line): 0.1 µg per test stripe

Sensitivity (analytical):

A biotin/-digoxigenin-peptide in 700 µL buffer solution at a concentration of 50 pg/mL is visually detected as positive on the basis of the result line in the read out zone after chromatography.

Specificity (analytical):

The test stripes don't show a visible result line after chromatography of 700 µL buffer without adding biotin/-digoxigenin-peptide (negative control).

Stability: At +2 to +8°C from date of manufacturing for 18 months.

Remark: The Multi Analyte Stripe is an immunoassay test stripe employing anti-biotin/anti-digoxigenin.

Catalog number

05 354 358 103

Pack size

50 stripes in 1 box

Will be supplied as "Multi Analyte Stripe". Unit of measure is "piece".

Minimum order size: 85 000 single strips.

For further processing only.

Streptavidin Magnetic Particles

suspension

Streptavidin-coated magnetic particles

Application

Use Streptavidin Magnetbeads for the fast and simple separation of a variety of biotin-labeled molecules in solution. This includes single mRNA, oligonucleotides, DNA, DNA fragments, glycoconjugates and protein isolation. Streptavidin Magnetbeads can also be used as a solid phase within liquid ELISA systems.

Properties

The beads are offered in a protease- and DNase/RNase-free format.

Catalog number

11 636 502 103

Pack size

custom fill

Will be supplied as "Streptavidin Magnetic-Particles". Unit of measure is "g active ingredient".

For further processing only.

Specification**Appearance:** Brown suspension**Content:** 9-11 mg/mL (solid binding)**Performance test in mRNA-HS-kit:** Function corresponds**Specific activity/Biotin binding capacity:** ≥ 1700 pmol/mg**Performance test** (in mRNA-HS-Kit): Function corresponds**Proteases** (incubation for up to 30 minutes at +37°C, casein-resorufin-marked): Not detectable**RNases** (incubation for up to 4 hours at +37°C): Not detectable**DNases** (incubation for up to 4 hours at +37°C): Not detectable**Unspecific binding of protein:** Not detectable**Stability:** At +2 to +8°C within specification range for 24 months.**StreptaWell, 384/C1**

transparent, coated with recombinant streptavidin

Application

Use Streptavidin-coated microplates in all kinds of immunoassays as well as in assay systems for detecting protein-nucleic acid interactions and nucleic acid amplification and hybridization. Various types of transparent and white streptavidin-coated plates (96- and 384-well) are available for colorimetric, chemiluminescent and fluorescent applications. The minimum order size is 500 plates/order with a delivery time of 8 –10 weeks.

Specification**Type of coating / number of cups:** C1 / 384**Type of plate** (384 NUNC): Corresponds to specification**Color** (transparent): Corresponds to specification**Biotin binding capacity:** ≥ 1.5 ng/well**Homogeneity [VK] of plate:** $\leq 8\%$ **Homogeneity [VK] of series:** $\leq 15\%$ **Leaching:** < 5 ngSA/well**Loading volume:** ≥ 90 μ L/well**NBL volume:** > 90 μ L/well**Stability:** At +2 to +8°C within specification range for 36 months.**Catalog number****Pack size****11 974 998 103**

1 plate

Will be supplied as "TRSA-SA MTP 384-well, clear". Unit of measure is "piece".

For further processing only.

StreptaWell, N-breakapart C8/C1

transparent, coated with recombinant streptavidin

Application

Use Streptavidin-coated microplates in all kinds of immunoassays as well as in assay systems for detecting protein-nucleic acid interactions and nucleic acid amplification and hybridization. Various types of transparent and white streptavidin-coated plates (96- and 384-well) are available for colorimetric, chemiluminescent and fluorescent applications. The minimum order size is 500 plates/order with a delivery time of 8–10 weeks.

Specification

Kind of coating/ number of cups (C1 / 96): Corresponds to specification

Plate type: C8 / breakapart

Color: Transparent

Biotin binding capacity: >5 ng/well

Homogeneity [VK] of plate: <5%

Homogeneity [VK] of series: <10%

Leaching: <5 ngSA/well

Loading volume: ≥250 µL/well

NBL volume: >250 µL/well

Stability: At +2 to +8°C within specification range for 36 months.

Catalog number**03 246 507 103****Pack size**

1 plate

Will be supplied as "SA-MTP (N-breakap. transp./C1)". Unit of measure is "piece".

For further processing only.

StreptaWell, N-breakapart C8/C2 plus, high binding capacity

transparent, coated with recombinant streptavidin

Application

Use Streptavidin-coated microplates in all kinds of immunoassays as well as in assay systems for detecting protein-nucleic acid interactions and nucleic acid amplification and hybridization. Various types of transparent and white streptavidin-coated plates (96- and 384-well) are available for colorimetric, chemiluminescent and fluorescent applications. The minimum order size is 500 plates/order with a delivery time of 8–10 weeks.

Specification

Kind of coating/ number of cups (C2 / 96): Corresponds to specification

Plate type (C8 / breakapart): Corresponds to specification

Color: Transparent

Biotin binding capacity: ≥20 ng/well

Homogeneity [VK] of plate: ≤5%

Homogeneity [VK] of series: ≤10%

Catalog number**11 986 694 103****Pack size**

1 plate

Will be supplied as "SA-MTP (N-breakap.C8/C2 plus)". Unit of measure is "piece".

For further processing only.

Leaching: <5 ngSA/well

Loading volume: ≥250 µL/well

NBL volume: >250 µL/well

Stability: At +2 to +8°C within specification range for 36 months.

StreptaWell, F8/C1 module

transparent, coated with recombinant streptavidin

Application

Use Streptavidin-coated microplates in all kinds of immunoassays as well as in assay systems for detecting protein-nucleic acid interactions and nucleic acid amplification and hybridization. Various types of transparent and white streptavidin-coated plates (96- and 384-well) are available for colorimetric, chemiluminescent and fluorescent applications. The minimum order size is 500 plates/order with a delivery time of 8–10 weeks.

Specification

Kind of coating/ number of cups (C1 / 96): Corresponds to specification

Plate type (F8): Corresponds to specification

Color: Transparent

Biotin binding capacity: ≥5 ng/well

Homogeneity [VK] of plate: <5%

Homogeneity [VK] of series: <10%

Leaching: <2 ngSA/well

Loading volume: ≥300 µL/well

NBL volume: >300 µL/well

Stability: At +2 to +8°C within specification range for 36 months.

Catalog number

11 940 279 103

Pack size

1 plate

Will be supplied as "SA coated MTP Nunc F8". Unit of measure is "piece".

For further processing only.

StreptaWell, F8/C2 plus, high binding capacity

transparent, coated with recombinant streptavidin

Application

Use Streptavidin-coated microplates in all kinds of immunoassays as well as in assay systems for detecting protein-nucleic acid interactions and nucleic acid amplification and hybridization. Various types of transparent and white streptavidin-coated plates (96- and 384-well) are available for colorimetric, chemiluminescent and fluorescent applications. The minimum order size is 500 plates/order with a delivery time of 8–10 weeks.

Catalog number

11 965 875 103

Pack size

1 plate

Will be supplied as "SA-MTP (Nunc F8 transp./C2+)". Unit of measure is "piece".

For further processing only.

Specification

Kind of coating/ number of cups (C2 / 96): Corresponds to specification

Plate type (F8): Corresponds to specification

Color: Transparent

Biotin binding capacity: ≥ 25 ng/well

Homogeneity [VK] of plate: $\leq 5\%$

Homogeneity [VK] of series: $\leq 10\%$

Leaching: < 5 ngSA/well

Loading volume: ≥ 300 μ L/well

NBL volume: > 300 μ L/well

Stability: At +2 to +8°C within specification range for 36 months.

Colloidal Gold 20 nm

suspension

Colloidal Gold is well known as an established labeling tool for a broad range of blotting and diagnostic applications and for electron- /light microscopy.

Application

Use Colloidal Gold, 20 nm as conjugation partner for all kind of antibodies, proteins and macromolecules. It is recommended for low to medium sensitive assays.

Product description

The 20 nm Goldsol quality is red and spherical. The ready to use “gold suspension” is pH adjusted.

Properties

“Citrate Gold” obtained from reduction of Tetrachloro-auric acid (HAuCl₄) with citric acid.

Specification

Appearance: Clear, light red liquid

Particle size: 19-23 nm

Particle concentration (A₅₂₀ unit): 0.85-1.00

λ_{max}: 516.0-518.5 nm

Stability: At +2 to +8°C within specification range for 12 months.

Background information

Due to its intense red color Colloidal Gold is one of the basic components for test strip development and manufacturing.

Catalog number	Pack size
05 418 291 103	1 L, 5 L, 25 L

Will be supplied as “Colloidal Gold 20 nm”. Unit of measure is “L”.

For further processing only.

Colloidal Gold 40 nm

suspension

Colloidal Gold is well known as an established labeling tool for a broad range of blotting and diagnostic applications and for electron- /light microscopy.

Application

Use Colloidal Gold, 40 nm as conjugation partner for all kind of antibodies, proteins and macromolecules. It is recommended for high sensitive assays due to its unique shape and color.

Product description

The 40 nm Goldsol quality is red/violet (“potato shape”). The ready to use “gold suspension” is pH adjusted.

Catalog number	Pack size
05 416 744 103	1 L, 5 L, 25 L

Will be supplied as “Colloidal Gold 40 nm”. Unit of measure is “L”.

For further processing only.

Properties

"Citrate Gold" obtained from reduction of Tetrachloro-auric acid (HAuCl_4) with citric acid.

Specification

Appearance: Turbid, slightly opalescent raspberry red liquid

Particle size: 38.0-43.0 nm

Particle concentration (A_{520} unit): 1 ± 0.2

λ_{max} : 531 ± 1 nm

Stability: At +2 to +8°C within specification range for 12 months.

Background information

Due to its intense red/violet color Colloidal Gold is one of the basic components for test strip development and manufacturing.

Framework IEP

lyophilizate

Application

The monoclonal Framework IEP eliminates monomeric (“immuno-response”) and highly specific interferences against the framework regions of the antibodies. Its interference eliminating strength is based on identical immuno-reactive conformational epitopes in the framework region of the test antibody.

Product description**Immunogen:** PTH**Spleen donor:** Mouse Balb/c**Antibody class:** IgG1, kappa

Preparation: Framework IEP is lyophilized from a solution containing protein, potassium phosphate and NaCl. No further preservatives are added.

Properties**Molecular structure:** IgG1, monomer**Remark:** Cannot be used in test systems for determination of PTH.

Serum concentrations of CK-MM: ≤3 U/mL have shown neither influence on interference elimination properties nor on recovery of analyte.

Specification**Appearance:** White lyophilizate

Solubility: Clear, colorless to slightly opalescent solution in NaCl, 0.9% (c=10 mg/mL)

Purity (HPLC): ≥90% IgG of total protein**Functional activity (relative titer based on masterlot):** ≥80%

Recommended working concentration: 20-5,000 µg/ml incubation buffer

pH 5.5 treatment (30 minutes): Corresponds to specification

Stability: At -15 to -25°C within specification range for 24 months. Avoid repeated freezing and thawing.

Catalog number**03 369 846 103****Pack size**

custom fill

Will be supplied as “Framework IEP *SQ”. Unit of measure is “mg active ingredient”.

For further processing only.

MAB IgG2b/Fab2a Poly

lyophilizate

Application

MAB IgG2b/Fab2a Poly mainly covers polymeric interference against IgG2a and/or IgG2b antibodies. It also covers Fab neo-epitopes.

Catalog number**11 355 830 103****Pack size**

5 mg, 50 mg, 250 mg

Will be supplied as “MAB-IgG(2b)/Fab(2a) Polymer, PolyMAB2b/2a”. Unit of measure is “mg active ingredient”.

For further processing only.

Product description**Immunogen:** DPH (2b), human S-AMY (2a)**Spleen donor:** Mouse Balb/c**Antibody class:** IgG2b, kappa/IgG2a, kappa**Preparation:** Lyophilized from a solution containing potassium phosphate and NaCl and 6% sucrose. No further preservatives are added.**Properties****Molecular structure:** IgG2a-Fab, polymerized with IgG2b; defined molecular range distribution**Remarks:** Cannot be used in test systems for determination of DPH, human S-AMY.**Specification****Appearance:** White lyophilizate**Solubility:** Clear, colorless to slightly opalescent solution in NaCl, 0.9% (c=5 mg/mL)**Turbidity properties:** Corresponds to specification**Functional activity** (relative titer based on master lot determined by MTP assay): ≥80%**Recommended working concentration:** 0.5-500 µg/ml incubation buffer**pH 5.5 treatment** (30 minutes): Corresponds to specification**Stability:** At -15 to -25°C within specification range for 24 months. Avoid repeated freezing and thawing.**MAB33 IgG1**

lyophilizate

Application

MAB33 IgG1 is used for test formulations employing intact IgG1. MAB33 IgG1 is especially suitable for the elimination of monomeric and specific interference.

Product description**Immunogen:** h CK-MM**Spleen donor:** Mouse Balb/c**Antibody class:** IgG1, kappa**Preparation:** Lyophilized from a solution containing potassium phosphate and NaCl. No further preservatives are added.**Properties****Molecular structure:** IgG1, monomer**Remarks:** Cannot be used in test systems for determination of CK-MM and CK-MB.**Catalog number****11 200 941 103****Pack size**

custom fill

Will be supplied as "MABM-33-IgG(DE),SQ MAB 33". Unit of measure is "g active ingredient".

For further processing only.

Specification**Appearance:** White lyophilizate**Solubility:** Clear, colorless to slightly opalescent solution in NaCl, 0.9% (c=10 mg/mL)**Protein** (Biuret): ≥0.7 mg protein/mg lyophilizate**Purity** (HPLC / Mono Q): ≥90 area% IgG of total protein**Functional activity** (relative titer based on master lot determined by MTP assay): ≥80%**Recommended working concentration:** 50-5,000 µg/ml incubation buffer**pH 5.5 treatment** (30 minutes): Corresponds to specification**Stability:** At -15 to -25°C within specification range for 24 months. Avoid repeated freezing and thawing.**Background information**

MAB33 IgG1 is a monoclonal antibody with defined specificity.

MAB33 IgG1/Fab1 Poly

lyophilizate

Application

MAB33-IgG1/Fab1 Poly is used for assays working with FAB-conjugates. MAB33-IgG1/Fab1 Poly is more efficient for polymeric and less specific types of interferences.

Product description**Immunogen:** h CK-MM**Spleen donor:** Mouse Balb/c**Antibody class:** IgG1**Preparation:** Lyophilized from a solution containing potassium phosphate and NaCl and 6% sucrose. No further preservatives are added.**Properties****Molecular structure:** IgG1-Fab, polymerized with IgG1; defined molecular range distribution**Remarks:** Cannot be used in test systems for determination of CK-MM and CK-MB.**Specification****Appearance:** White lyophilizate**Solubility:** Clear, colorless to slightly opalescent solution in NaCl, 0.9% (c=10 mg/mL)**Functional activity** (relative titer based on master lot determined by MTP assay): ≥80%**Recommended working concentrations:** 0.5-500 µg/ml incubation buffer**Catalog number****11 368 338 103****Pack size**

5 mg, 50 mg, 250 mg

Will be supplied as "MAB-IgG/Fab (Polymer), SQ Poly MAB 33". Unit of measure is "g active ingredient".



For further processing only.

pH 5.5 treatment (30 minutes): Corresponds to specification

Stability: At -15 to -25°C within specification range for 24 months.

Avoid repeated freezing and thawing.

Background information

The polymer MAB33 IgG1/Fab1 Poly incorporates elements comparable with those of MAB33 IgG1/IgG1 Poly, and also covers interferences against Fab neo-epitopes.

MAB33 IgG1/IgG1 Poly

frozen solution

Application

MAB33-IgG1/IgG1 Poly is used for test formulations employing intact IgG1. MAB33-IgG1/IgG1 Poly is more efficient for polymeric and less specific types of interferences.

Product description

Immunogen: h CK-MM

Spleen donor: Mouse Balb/c

Antibody class: IgG1

Preparation: Lyophilized from a solution containing potassium phosphate and NaCl and 6% sucrose. No further preservatives are added.

Properties

Molecular structure: Molecular structure: IgG1, polymerized with IgG1; defined molecular range distribution

Remarks: Cannot be used in test systems for determination of CK-MM and CK-MB.

Specification

Appearance: Frozen liquid

Solubility: Yellowish clear to slightly opalescent solution, containing K-phosphate buffer, NaCl and 4% sucrose, pH 7,5

Protein (Biuret): ≥ 30 mg protein/mg lyophilizate

Purity (HPLC): ≥ 90 area% IgG of total protein

Functional activity (relative titer based on master lot determined by MTP assay): $\geq 80\%$

Recommended working concentration: 0.5-500 $\mu\text{g/mL}$ incubation buffer

Turbidity properties δA_{334} : Corresponds to specification

Bioburden: ≤ 250 CFU/mL

pH 5.5 treatment (30 minutes): Corresponds to specification

Catalog number

11 939 661 103

Pack size

5 mg, 50 mg, 250 mg, 1 g, 5 g

Will be supplied as "MAB-33-IgG-Polymer *SQ". Unit of measure is "g active ingredient".



For further processing only.

Stability: At -60 to -90°C within specification range for 24 months.
Avoid repeated freezing and thawing.

Background information

MAB33 IgG1/IgG1 Poly is the polymerized chemical version of MAB33 IgG1.

PAB<->H-IgG/Fab Poly

frozen solution

PAB<->H-IgG/Fab Poly is designed as a specific blocker to eliminate false positive/negative antibody-directed interactions.

Application

PAB<->H-IgG/Fab Poly eliminates antibody-directed interactions derived from patient sera in immunoassays employing human and/or chimeric antibody conjugates.

Product description

Preparation: Lyophilized from a solution containing acetic acid, sodium chloride, and 6% sucrose.

Properties

Co-polymerization of human IgG-Fab with IgG ensures highly consistent molecular weight distribution.

Specification

Appearance: White lyophilizate

Solubility: Clear, to slightly opalescent colorless solution in NaCl, 0.9% (c=5 mg/mL)

Protein (Biuret): 0.25-0.60 mg/mg lyophilizate

Relative titer (refers to master lot): Corresponds to specification

Stability: At -15 to -25°C within specification range for 24 months.

Catalog number

11 668 544 103

Pack size

10 mg (samples), custom fill

Will be supplied as "PAB<->H-IgG/Fab Polymer". Unit of measure is "g active ingredient".

For further processing only.

Sheep IgG (PAB<->S-IgG)

lyophilizate

Application

Sheep IgG reduces nonspecific antibody interferences in assays employing sheep antibodies.

Catalog number

10 717 606 103

Pack size

custom fill

Will be supplied as "PAB<->S-IgG". Unit of measure is "g".

For further processing only.

Product description

The polyclonal antibody Sheep IgG is lyophilized from a solution containing potassium phosphate and sodium chloride. No preservatives are added. IgG fraction is produced in sheep, purified by anion-exchange chromatography.

Properties

Recommended working concentration: 0.5-1.5 mg/mL

Specification

Appearance: White lyophilizate

Solubility: Clear to slightly opalescent solution in NaCl, 0.9% (c= 5 mg/mL)

Turbidimetric measurement (A_{546} , against water): ≤ 100 mE

Protein (Biuret): ≥ 0.8 mg/mg lyophilizate

Aggregated IgG (TSK 3000): $\leq 10\%$

Stability: At -15 to -25°C within specification range for 24 months.

Rabbit IgG (PAB<->K-IgG)

lyophilizate

Application

Rabbit IgG reduces nonspecific antibody interferences in assays employing rabbit antibodies.

Product description

The polyclonal antibody Rabbit IgG is lyophilized from a solution containing potassium phosphate and sodium chloride. No preservative are added. IgG fraction is produced in rabbit, purified by anion-exchange chromatography.

Specification

Appearance: White lyophilizate

Solubility: Clear to slightly opalescent solution in NaCl 0.9% (c=10 mg/mL)

Protein (Biuret) : $\geq 0,7$ mg/mg lyophilizate

Purity (HPLC / TSK 3000): Corresponds to specification (in comparison to master lot)

Stability: At -15 to -25°C within specification range for 36 months.

Avoid repeated freezing and thawing.

Catalog number

10 912 280 103

Pack size

custom fill

Will be supplied as "PAB<->K-IgG(DE-FF)". Unit of measure is "g active ingredient".

For further processing only.

Bovine IgG (PAB<->R-IgG)

lyophilizate

Application

Bovine IgG reduces nonspecific adsorption of antibodies to the solid phase and other cross-reactive, nonspecific antibody interactions.

Product description

The polyclonal antibody Bovine IgG is lyophilized from a solution containing potassium phosphate and sodium chloride. No preservatives are added. IgG fraction is produced in bovine, purified by anion-exchange chromatography.

Properties

Recommended working concentration: 0.5-2.5 mg/mL

Specification

Appearance: White lyophilizate

Solubility: Clear, colorless solution in NaCl, 0.9% (c=10 mg/mL)

Turbidimetric measurement (A_{546} , against water): ≤ 100 mE

Protein (Biuret): ≥ 0.8 mg/mg lyophilizate

Aggregated IgG (HPLC / TSK 3000): $\leq 5\%$

Country of origin: USA or NZL

pH 5.5 treatment (30 minutes): Corresponds to specification

Stability: At -15 to -25°C within specification range for 24 months.

Catalog number**11 293 621 103****Pack size**

custom fill

Will be supplied as "PAB<->R-IgG(DE) Bovine IgG". Unit of measure is "g active ingredient".



For further processing only.

HAMA Serum, Type I

lyophilizate

Application

HAMA Serum Type I is primarily intended for polyvalent and spontaneous HAMA interference occurring in healthy donors.

Product description

HAMA Serum Type I is a lyophilized human serum with serum ingredients within normal range. No preservative are added. The product must be handled just as carefully as patient specimens.

Properties

pH value (+25°C): 7-8

Protein (Biuret): ≥65 mg/vial

Dilution: Up to 1:20 recommended

Solubility: Reconstitution of lyophilizate/vial in 1 mL water results in serum typical solution

Typical concentrations determined by Roche cobas® assays:

Alpha fetoprotein (AFP): 2.0 IU/mL

Carcinoembryonic antigen (CEA): 1.3 ng/mL

Follicle stimulating hormone (FSH): 10.0 mIU/mL

Human chorionic gonadotropin (HCG): 2.8 mIU/mL

Luteinizing hormone (LH): 8.1 mIU/mL

Prolactin: 160.7 µIU/mL

Prostate specific antigen (PSA): 0.3 ng/mL

PSA free: 0.06 ng/mL

Thyroid stimulating hormone (TSH): 2.5 µIU/mL

Troponin T (TN-T): Not detectable

Specification

Appearance: Yellowish lyophilizate

Interference effect: Corresponds to specification

Infectious parameters (determined by FDA approved methods):

Non reactive HbsAg: Corresponds to specification

Anti HIV 1+2: Negative

Anti HCV: Negative

HIV 1 NAT, non reactive: Corresponds to specification

Stability: At -15 to - 25°C within specification range for 24 months.

Background information

HAMA serum interferences in immunoassays can vary within a broad range depending on the person's immune system.

Catalog number

11 767 275 103

Pack size

1 mL vial

Will be supplied as "Hama-Serum I - Qual. Standard *SQ". Unit of measure is "piece".

For further processing only.

HAMA Serum, Type II

lyophilizate

Application

HAMA Serum Type II primarily represents mono-/bivalent and specific HAMA interference occurring after treatment with monoclonal antibodies.

Product description

HAMA Serum Type II is a lyophilized human serum with serum ingredients within normal range. No preservative are added. The product must be handled just as carefully as patient specimens.

Properties

pH value (+25°C): 7-8

Protein (Biuret): ≥65 mg/vial

Dilution: Up to 1:20 recommended

Solubility: Reconstitution of lyophilizate/vial in 1 mL water results in serum typical solution

Typical concentrations determined by Roche cobas® assays:

Alpha fetoprotein (AFP): 3.5 IU/mL

Carcinoembryonic antigen (CEA): 1.2 ng/mL

Follicle stimulating hormone (FSH): 7.0 mIU/mL

Human chorionic gonadotropin (HCG): 3.0 mIU/mL

Luteinizing hormone (LH): 6.7 mIU/mL

Prolactin: 114 µIU/mL

Prostate specific antigen (PSA): 0.3 ng/mL

PSA free: 0.05 ng/mL

Thyroid stimulating hormone (TSH): 1.6 µIU/mL

Troponin T (TN-T): 0.01 ng/mL

Specification

Appearance: Yellowish lyophilizate

Interference effect: Corresponds to specification

Infectious parameters:

Non reactive HbsAg: Corresponds to specification

Anti HIV 1+2: Negative

Anti HCV: Negative

HIV 1 NAT, non reactive: Corresponds to specification

Stability: At -15 to -25°C within specification range for 24 months.

Background information

HAMA serum interferences in immunoassays can vary within a broad range depending on the person's immune system.

Catalog number

05 167 060 103

Pack size

1 mL vial

Will be supplied as "HAMA Serum 2L *SQ". Unit of measure is "piece".

For further processing only.

Bovine Serum Albumin I

lyophilizate

Application

Bovine Serum Albumin I reduces nonspecific adsorption to the solid phase or saturates unoccupied binding sites. Bovine Serum Albumin I represents a very high purity grade ($\geq 95\%$) and is recommended for assays with higher demands on sensitivity.

Product description

BPLA Type I is lyophilized from a solution containing potassium phosphate and sodium chloride. No preservatives are added.

Specification

Appearance: Yellowish lyophilizate

Solubility: Clear to slightly turbid yellowish solution in water ($c=60$ mg/mL)

A₄₀₅ (against water): ≤ 0.250

pH value: 6.5-7.5

Protein (from N, elementary analysis, factor 6.25): $\geq 95\%$

Purity (HPLC / TSK 3000): $\geq 95\%$ (monomer)

Water (K. Fischer): $\leq 5\%$

Bioburden: ≤ 50 CFU/g

Ca: $\leq 0.1\%$

Fe: $\leq 0.005\%$

Cu: $\leq 0.002\%$

Complexing agent:

Recovery of Fe: $100 \pm 20\%$

Recovery of Cu: $100 \pm 20\%$

Heavy metals (as Pb): $\leq 0.002\%$

P_i: $\leq 0.005\%$

Octanoic acid (GC): $\leq 0.5\%$

Analysis of T3, T4, Estradiol, Testosterone, Progesterone (for information only): Values stated on certificate of analysis.

Country of origin: USA or NZL

pH 4.5 treatment (up to 3 hours): Corresponds to specification

Stability: At $+2$ to $+8^\circ\text{C}$ within specification range for 24 months.

Store dry.

Catalog number**11 726 536 103****Pack size**

custom fill

Will be supplied as "Albumin RPLA 1 Assay Quality". Unit of measure is "kg".

For further processing only.

Bovine Serum Albumin IV

lyophilizate

Application

Bovine Serum Albumin IV reduces nonspecific adsorption to the solid phase or saturates unoccupied binding sites. Bovine Serum Albumin IV represents a very high purity grade ($\geq 95\%$) and is recommended for assays with higher demands on sensitivity.

Product description

BPLA Type IV is lyophilized from a solution containing potassium phosphate and sodium chloride. No preservatives are added.

Specification

Appearance: Yellowish lyophilizate

Solubility: Clear to slightly turbid yellowish solution in water (c=60 mg/mL)

A₄₀₅ (against water): ≤ 0.200

pH value: 6.5-7.5

Protein (from N, elementary analysis, factor 6.25): $\geq 95\%$

Purity (HPLC / TSK 3000): $\geq 95\%$ (monomer)

Water (K. Fischer): $\leq 5\%$

Bioburden: ≤ 50 CFU/g

Ca: $\leq 0.1\%$

Fe: $\leq 0.005\%$

Cu: $\leq 0.002\%$

pH 4.5 treatment (up to 3 hours): Corresponds to specification

Country of origin: USA, NZL

Stability: At +2 to +8°C within specification range for 24 months. Store dry.

Catalog number**11 726 544 103****Pack size**

custom fill

Will be supplied as "Albumin RPLA 4 Assay Quality". Unit of measure is "kg".

For further processing only.

Bovine Serum Albumin (BSA), Fraction V

lyophilizate

Serum albumin protein that has numerous biochemical applications.

Application

Use Bovine Serum Albumin (BSA) as a buffering agent, stabilizer, standard and for blending. Bovine Serum Albumin (BSA) is also a versatile tool against non-specific solid phase interference. As blocking reagent Bovine Serum Albumin (BSA) saturates unoccupied binding sites on the solid phase. Use Bovine Serum Albumin (BSA) typically at a concentration of 0.5 to 3% within the reagent buffer.

CAS: 9048-46-8

Catalog number**10 738 328 103****Pack size**

custom fill

Will be supplied as "Albumin, Fraction V from Bovine Serum". Unit of measure is "kg".

For further processing only.

Properties**Molecular weight:** 68 kD

Bovine Serum Albumin (BSA) contains no detectable IgG.

Bovine Serum Albumin (BSA) is controlled for low molecular weight contaminants.

Bovine Serum Albumin (BSA) consists primarily of monomeric albumin.

Specification**Appearance:** Slightly yellow lyophilizate**Solubility:** Clear, odourless solution in water**A₄₀₅** (against water): ≤0.200**Albumin** (gel electrophoresis): ≥98%**Protein** (from N according to elementary analysis): ≥95%**pH value:** 6.8-7.2**Water** (K. Fischer): ≤5%**Heavy metals** (as Pb): ≤ 0.003%**P_i**: ≤0.003%**Chloride** (chloride meter): ≤0.15%**Glucose** (enzymatic): ≤0.05%**Glycerol** (enzymatic): ≤0.005%**L-Lactate** (enzymatic): ≤0.1%**Na** (flame photometric): ≤0.8%**K** (flame photometric): ≤0.015%**Li** (flame photometric): ≤0.0005%**Ca:** ≤0.05%**Mg:** ≤0.005%**Fe:** ≤0.002%**Bioburden:** ≤100 CFU/g lyophilizate**Country of origin:** New Zealand, USA**Stability:** At +2 to +8°C within specification range for 24 months.

Store dry.

Remarks:

Official veterinary certificate of health of the donor animals is available.

Official certificate of the deactivation of animal material including the method (acid treatment at pH 5 for 3 hours) is available.

Bovine Serum Albumin (BSA), Fraction Vfatty acids ≤ 0.2 mg/g, lyophilizate

Highly purified serum albumin protein that has numerous biochemical applications.

Application

Use Bovine Serum Albumin (BSA) as a buffering agent, stabilizer, standard and for blending. Bovine Serum Albumin (BSA) is also a versatile tool against non-specific solid phase interference. As blocking reagent Bovine Serum Albumin (BSA) saturates unoccupied binding sites on the solid phase. Use Bovine Serum Albumin (BSA) typically at a concentration of 0.5 to 3% within the reagent buffer.

CAS: 9048-46-8**Properties****Molecular weight:** 68 kD

Bovine Serum Albumin (BSA) contains no detectable IgG.

Bovine Serum Albumin (BSA) is controlled for low molecular weight contaminants.

Bovine Serum Albumin (BSA) consists primarily of monomeric albumin.

Specification**Appearance:** Slightly yellow lyophilizate**Protein** (from N, according to elementary analysis): $\geq 97.0\%$ **Water** (K. Fischer): $\leq 5.0\%$ **Na** (flame photometric): $\leq 0.5\%$ **K** (flame photometric): $\leq 0.01\%$ **Fe** (AAS): $\leq 0.001\%$ **Cu** (AAS): $\leq 0.002\%$ **Fatty acids, total** (GC): ≤ 0.2 mg/g**Triglycerides** (enzymatic): Not detectable**Immunoglobulines** (ELISA): Not detectable**Country of origin:** USA**Stability:** At +2 to +8°C within specification range for 24 months.**Remarks:**

Official veterinary certificate of health of the donor animals is available.

Official certificate of the deactivation of animal material including the method (acid treatment at pH 5 for 3 hours) is available.

Catalog number**10 774 111 103****Pack size**

custom fill

Unit of measure is "kg".

For further processing only.

Bovine Serum Albumin (BSA), reduced sodium and potassium

lyophilizate

Serum albumin protein for tests that require a strongly reduced concentration of sodium and potassium.

Application

Use Bovine Serum Albumin (BSA) as a buffering agent, stabilizer, standard and for blending. Bovine Serum Albumin (BSA) is also a versatile tool against non-specific solid phase interference. As blocking reagent Bovine Serum Albumin (BSA) saturates unoccupied binding sites on the solid phase. Use Bovine Serum Albumin (BSA) typically at a concentration of 0.5 to 3% in the reagent buffer.

CAS: 9048-46-8

Properties

Molecular weight: 68 kD

Bovine Serum Albumin (BSA) contains no detectable IgG.

Bovine Serum Albumin (BSA) is controlled for low molecular weight contaminants.

Bovine Serum Albumin (BSA) consists primarily of monomeric albumin.

Specification

Appearance: Yellowish lyophilizate

A₄₀₅ (against water): ≤0.300

pH value: 4.3-5.3

Protein (Biuret): ≥80%

Water (K. Fischer): ≤5%

Na (AAS): ≤35 ppm

K (AAS) : ≤4ppm

Li (AS) : ≤50ppm

Ca (AAS): ≤500ppm

Fe (bathophenanthroline): ≤25 ppm

Cu (bathocuproin): ≤15ppm

Heavy metals (as Pb): ≤50ppm

P_i: ≤150ppm

Bioburden: ≤100 CFU/g lyophilizate

NH₄ (enzymatic): ≤10ppm

Glucose (enzymatic): ≤0.02%

Complex creator:

Recovery of Fe: 80-120%

Recovery of Cu: 80-120%

Electrophoresis: Chromatographically homogeneous

Country of origin: Germany

Stability: At +2 to +8°C within specification range for 36 months.

Catalog number

11 297 368 103

Pack size

custom fill

Unit of measure is "kg".

For further processing only.

Remarks:

Official veterinary certificate of health of the donor animals is available.
 Official certificate of the deactivation of animal material including the method is available.

Poly BSA Type I

frozen solution

Catalog number**11 866 737 103****Pack size**

1 g, 5 g, 20 g

Application

Polymeric BSA Type I is used for the elimination of nonspecific polymeric interferences (e.g. for antibodies, conjugates or antigens integrated into infectious disease assays and tumor marker tests).

Product description

Poly BSA Type I is produced from bovine serum albumine (BSA) by polymerization. Its surface is chemically modified by acetylation, and the negative charge of the Poly BSA Type I eliminates hydrophobic interactions. The Poly BSA Type I solution contains protein, potassium-phosphate buffer and preservatives (chloroacetamide, methylisothiazolone and sucrose).

Properties

Recommended working concentration: 0.1-20 mg/mL incubation buffer.

Specification

Appearance: Yellowish, clear to slightly opalescent solution (frozen)

pH value (+25°C): 6.8-7.2

Protein (Biuret): ≥40 mg/mL

Particle size (Photon correlations spectrometry): 15-45 nm

Country of origin: USA

pH 4.5 treatment (up to 3 hours): Corresponds to specification

Stability: At -15 to -25°C within specification range for 24 months.

Avoid repeated freezing and thawing.

Will be supplied as "Poly BSA Type I *SQ". Unit of measure is "g active ingredient".



For further processing only.

Poly BSA Type II

frozen solution

Application

Polymeric BSA Type II is used for the elimination of nonspecific polymeric interferences (e.g. for antibodies, conjugates or antigens integrated into infectious disease assays and tumor marker tests).

Product description

Poly BSA Type II is produced from bovine serum albumine (BSA) by polymerization. Its surface is chemically modified by succinylation, and the negative charge of the Poly BSA Type II eliminates hydrophobic interactions. The Poly BSA Type II solution contains protein, potassium-phosphate buffer, and preservatives (chloroacetamide, methylisothiazolone and sucrose).

Properties

Recommended working concentration: 0.1–20 mg/ml incubation buffer

Specification

Appearance: Yellowish, clear to slightly opalescent solution (frozen)

pH value (+25°C): 6.8–7.2

Protein (Biuret): ≥40 mg/mL

Particle size (Photon correlations spectrometry): 25–55 nm

Country of origin: USA

pH 4.5 treatment (up to 3 hours): Corresponds to specification

Stability: At -15 to -25°C within specification range for 24 months. Avoid repeated freezing and thawing.

Catalog number

11 816 438 103

Pack size

1 g, 5 g, 20 g

Will be supplied as "poly BSA Type II *SQ". Unit of measure is "g active ingredient".



For further processing only.

Streptavidin rec. inactive, Poly

frozen solution

Streptavidin rec. inactive, Poly is designed as a specific blocker for streptavidin directed interactions derived from patient sera.

Application

Blocking of interferences in immunoassay formats employing streptavidin/biotin-system.

Properties

It is prepaed as frozen solution containing 50 mM KPO₄, 8% saccharose, pH 6.6.

Specification

Appearance: Opalescent suspension

Function test on Elecsys® TnT hs: Corresponds

Unencumbered based on reference (1 week, +35°C, only loaded R2): Corresponds

Stability: At -60 to -90°C within specification range for 12 months.

Catalog number

11 922 122 103

Pack size

0.01, 0.4 g, custom fill

Will be supplied as "Streptavidin rec. inactive, Poly". Unit of measure is "g active ingredient".



For further processing only.

Alkaline Phosphatase Mutein, recombinantfrom calf intestine, expressed in *Pichia pastoris*, lyophilizate**Application**

Use Alkaline Phosphatase Mutein (AP Mutein) to eliminate human serum derived AP directed assay interferences.

Product description

In addition to the AP Mutein, recombinant the lyophilizate contains a mixture of proteins from the *Pichia* system supporting the interference elimination. AP Mutein is lyophilized from a solution containing in NaCl, 0.2 mol/L; ZnCl₂, 0.1 mmol/L; Tea, 30 mmol/L; MgCl₂, 1 mmol/L; raffinose, 50% (w/v); pH approximately 7.6. Production is established according to the procedures of the active enzyme.

EC 3.1.3.1

Specification

Appearance: White to yellowish lyophilizate

Solubility: Clear to light yellowish solution in 50% glycerol solution (c=10 mg/mL)

pH value: 7.0-8.0

Protein (A₂₈₀, 1 mg/mL=1, against water): ≥0.2 mg protein/mg lyophilizate

Specific activity (+37°C, 4-NPP): ≤10 U/mg protein

Catalog number

04 781 007 103

Pack size

custom fill

Will be supplied as "AP-Mutein, rec.". Unit of measure is "g".

For further processing only.

SDS-gel (qualitative comparison of the gel bands in reference to the bands of a standard): Corresponds to specification

Stability: At +2 to +8°C within specification range for 24 months.

Background information

Alkaline Phosphatase Mutein (AP Mutein), recombinant, is the inactive form of recombinant highly active AP, expressed in *Pichia pastoris*. The inactivation of AP Mutein is based on one single point mutation located in the active site of the alkaline phosphatase.

β-Galactosidase Mutein

from *E. coli* overproducer, lyophilizate

Application

Use β-Galactosidase Mutein to eliminate β-galactosidase directed interferences in immunoassays derived from human sera.

EC 3.2.1.23

Properties

β-Galactosidase Mutein is identical to native β-galactosidase with respect to immuno-reactivity, conjugation, properties, surface charge, hydrophobicity, molecular weight, production procedure and down stream processing.

Specification

Appearance: White lyophilizate

Solubility: Clear, colorless to light opalescent solution in water (c=20 mg/mL)

Protein (Biuret): 0.15-0.30 mg/mg lyophilizate

Specific activity (10 mg/ml, +37°C, 2-nitrophenyl-β-D-galactopyranoside): ≤0.2 U/mg protein

Aggregated β-Galactose Mutein (TSK 4000): ≤10%

Immunoreactivity (based on ML): 80-120%

Stability: At -15 to -25°C within specification range for 24 months.

Background information

β-Galactose Mutein from *E.coli* is constructed using site-directed mutagenesis of single amino acids in the active site.

Catalog number

11 184 024 103

Pack size

custom fill

Will be supplied as "β-Galactosidase Mutein". Unit of measure is "mg active ingredient".



For further processing only.

Alkaline Phosphatase, EIA Grade

from calf intestine, solution

Application

Marker enzyme for preparing antibody- / antigen-enzyme conjugates incorporated in immunoassay reagents for colorimetric, fluorimetric and luminometric detection.

EC 3.1.3.1

Properties

Nomenclature: Orthophosphoric-monoester phosphohydrolase (alkaline optimum)

Molecular weight: ≥57 kD

Inhibitors: P_i, metal chelating agents, divalent heavy metal ions (e.g., Be²⁺, Zn²⁺), many amino acids (e.g., L-phenylalanine, L-tryptophan, L-cysteine), iodosobenzoate, iodoacet-amide.

Activators: Mg²⁺, Co²⁺, Mn²⁺

pH optimum: 9.8

pH stability: 8.0

Thermal stability: Up to +40°C

Specification

Appearance: Clear, colourless solution in NaCl, 3 mol/L; MgCl₂, 1 mmol/L; ZnCl₂, 0.1 mmol/L; Tea, 30 mmol/L

pH value: 7.0-8.0

Protein (A₂₈₀, 1 mg/ml=1, against water): ≥10 mg/mL

Specific activity (+37°C, 4-NPP): ≥3,000 U/mg

Alkaline Phosphatase (HPLC): ≥95 area% (HPLC profile added to certificate)

Amino groups: 8-16 mol/mol

Carbohydrates, n=2: No limit

Molecular weight: ≥57,000 D

NH₄⁺: ≤30 ppm

pH 5.5 treatment (30 minutes): Corresponds to specification

Origin of bovine intestine: NZL

Stability: At +2 to +8°C within specification range for 15 months.

Background information

Alkaline phosphatase catalyzes the hydrolysis of numerous phosphate esters, such as esters of primary and secondary alcohols, sugar alcohols, cyclic alcohols, phenols and amines. Phosphodiester do not react with Alkaline Phosphatase, EIA Grade. The enzyme hydrolyzes PP_i. The kinetic properties of the enzyme depend on many factors, such as purity of enzyme, concentration of enzyme in the assay, buffer, pH etc.

Catalog number

10 556 602 103

Pack size

custom fill

Will be supplied as "Phosphatase, Alkaline, Calf Intestine". Unit of measure is "g".

For further processing only.

Alkaline Phosphatase, recombinant, highly activefrom calf intestine, expressed in *Pichia pastoris*, solution**Application**

Marker enzyme for preparation of antibody- / antigen-enzyme conjugates incorporated in highly sensitive immunoassay reagents for colorimetric, fluorimetric and luminometric detection. Alkaline Phosphatase is recommended for conjugation via carbohydrate groups (content approximately 30%).

EC 3.1.3.1

Properties**IEP** (IEF, CE): 3.6-4.7**MALDI-TOF MS:**

Total molecular weight: 124±10 kD

Molecular weight protein: 104 kD (=84%)

Molecular weight carbohydrate: 20±10 kD (16±6%)

Accessible N-glycosylation sites: 2/subunit**Branching type:** Higher branched (hybrid) type GlcNAc, Mannose, no NeuAc detected**O-glycosylation sites:** Not detected**Number of isoenzymes:**

Present: 3 (MS)

Based on protein: 1 (MS)

Specification**Appearance:** Clear, colorless solution in NaCl, 3 mol/L; MgCl₂, 5 mmol/L; ZnCl₂, 0.1 mmol/L; Tea, 30 mmol/L, pH approximately 7.6**pH value:** 7.0-8.0**Protein** (A₂₈₀, 1 mg/ml=1, against water): 20±1 mg/mL**Specific activity** (+37°C, 4-NPP): ≥7,000 U/mg**Alkaline Phosphatase** (HPLC): ≥95 area% (HPLC profile added to certificate)**Amino groups:** 5-13 mol/mol**Carbohydrates, n=2:** No limit**Stability:** At +2 to +8°C within specification range for 12 months.**Background information**

Alkaline Phosphatase recombinant, highly active catalyzes the hydrolysis of numerous phosphate esters, such as esters of primary and secondary alcohols, sugar alcohols, cyclic alcohols, phenols and amines. Phosphodiesterases do not react with Alkaline Phosphatase, recombinant. The enzyme hydrolyzes PP_i. The kinetic properties of the enzyme depend on many factors, such as purity of enzyme, concentration of enzyme in the assay, buffer, pH etc.

Catalog number**03 137 031 103****Pack size**

custom fill

Will be supplied as "AP, Yeast, high act., rec., EIA, NaCl". Unit of measure is "g".

For further processing only.

Alkaline Phosphatase, recombinant, highly active, carbohydrate reducedfrom calf intestine, expressed in *Pichia pastoris*, solution**Application**

Marker enzyme for preparation of antibody- / antigen- enzyme conjugates incorporated in highly sensitive immunoassay reagents for colorimetric, fluorimetric and luminometric detection. Recommended for conjugation via amino groups.

EC 3.1.3.1

Properties**IEP** (IEF, CE): 5.2-6.0**MALDI-TOF MS:**

Total molecular weight: 111 kD

Molecular weight protein: 104 kD (= 94%)

Molecular weight carbohydrates: 6.5 kD (= 6%)

Accessible N-glycosylation sites: 2/subunit**Branching type:** Reduced branched type GlcNAc, Mannose no NeuAc detected**O-glycosylation sites:** Not detected**Number of isoenzymes:**

Present: 1 (MS)

Based on protein: 1 (MS)

Specification**Appearance:** Clear, colourless solution in NaCl, 3 mol/L; MgCl₂, 5 mmol/L; ZnCl₂, 0.1 mmol/L; Tea, 30 mmol/L, pH approximately 7.6**pH value:** 7.0-8.0**Protein** (A₂₈₀; 1 mg/ml=1; against water): 20±1 mg/mL**Specific activity** (+37°C, 4-NPP): ≥ 7000 U/mg**Alkaline Phosphatase** (HPLC): ≥95 area% (HPLC profile added to certificate)**Amino groups:** 5-13 mol/mol**Carbohydrates, n=2:** ≤7%**Stability:** At +2 to +8°C within specification range for 12 months.**Background information**

AP recombinant, highly active CR catalyzes the hydrolysis of numerous phosphate esters, such as esters of primary and secondary alcohols, sugar alcohols, cyclic alcohols, phenols and amines. Phosphodiesterases do not react with Alkaline Phosphatase, recombinant. The enzyme hydrolyzes PP_i. The kinetic properties of the enzyme depend on many factors, such as purity of enzyme, concentration of enzyme in the assay, buffer, pH and others. The product contains a significantly reduced carbohydrate moiety.

Catalog number**03 535 452 103****Pack size**

custom fill

Will be supplied as "AP, highly active, recombinant, CR". Unit of measure is "g".

For further processing only.

β-Galactosidasefrom *E. coli*, lyophilizate**Application**

Marker enzyme for the manufacturing of antibody- and antigen-enzyme conjugates incorporated in immunoassays for colorimetric and fluorimetric detection.

EC 3.2.1.23

Properties**Nomenclature:** β-D-galactohydrolase**Molecular weight** (by sequence): 465 kD**Structure:** 4 identical subunits; β-galactosidase contains no carbohydrates**Specification****Appearance:** White lyophilizate**pH value** (c=10 mg/mL, in water): 7.0-8.0**Protein** (Biuret): 0.25-0.5 mg/mg lyophilizate**Activity** (+37°C, 2-NP-β-D-galactoside): ≥120 U/mg lyophilizate**Specific activity** (+37°C, 2-nitrophenyl-β-D-galactopyranoside): ≥300 U/mg protein**Contaminants** (expressed as percentage of β-Galactosidase activity):

β-Fructosidase: ≤0.001

α-Galactosidase: ≤0.001

Glucose-DH: <0.001

α-Glucosidase: ≤0.001

"NADH oxidase": ≤0.001

Na (flame photometric): ≤2,500 ppm**Stability:** At +2 to +8°C: within specification range for 12 months.

Store dry.

Background information

β-Galactosidase hydrolyzes β-D-galactosides.

Catalog number**11 291 963 103****Pack size**

custom fill

Will be supplied as "β-Galactosidase, Lyo.". Unit of measure is "MU".

For further processing only.

β-Galactosidase, recombinant, EIA Gradefrom *E. coli* overproducer, lyophilizate**Application**

Marker enzyme for the manufacturing of antibody- and antigen-enzyme conjugates incorporated in immunoassays for colorimetric and fluorimetric detection.

EC 3.2.1.23

Catalog number**10 570 079 103****Pack size**

custom fill

Will be supplied as "β-Galactosidase, *E. coli*". Unit of measure is "g active ingredient".



For further processing only.

Properties

Nomenclature: β -D-galactohydrolase

Molecular weight (by sequence): 465 kD

Structure: 4 identical subunits, β -galactosidase contains no carbohydrates

Isoelectric point: 4.61

Michaelis constants:

Tris buffer, pH 7.6, +20°C / relation rate:

2-nitrophenyl- β -galactoside: 9.50×10^{-4} mol/L / 1.00

phenyl- β -D-galactoside: 3.23×10^{-3} mol/L / 0.05

lactose: 3.85×10^{-2} mol/L / 0.06

4-nitrophenyl- β -galactoside: 4.45×10^{-4} mol/L / ~ 0.50

Activators: Mg^{2+} and Na^+ (or other monovalent cations) are essential for activity.

pH optimum: 8.0

pH stability: 6.0

Thermal stability: Up to +37°C

Thiol groups: 64 SH groups, approximately 16 of these are accessible for SH-reactive reagents. 4 of these (Cys 76) take part in conjugation.

Specification

Appearance: White lyophilizate, stabilized with phosphate buffer and sucrose

Solubility: Clear, colorless solution in water (c=20 mg/mL)

Protein (Biuret): Approximately 0.25 mg/mg lyophilizate

Specific activity (+37°C, 2-nitrophenyl- β -D-galactopyranoside): ≥ 700 U/mg protein

SH-groups, free (after dialysis): ≥ 12 mol/mol enzyme (corresponds to 465,000 g)

Aggregated β -galactosidase (HPLC): ≤ 3 area% (dimer-part with a molar mass of 0.93×10^6 D)

Stability: At -15 to -25°C within specification range for 24 months. Store under nitrogen.

Background information

β -Galactosidase hydrolyzes β -D-galactosides. Although the enzyme activity with 2-nitrophenyl- β -D-galactoside as substrate is higher than with the 4-isomer, the enzyme reaction with the 4-compound is more sensitive due to a higher absorption coefficient for 4-nitrophenol, ϵ_{405} : 18.5 [$mmol^{-1} \times l \times cm^{-1}$].

Peroxidase (POD), EIA Grade

from horseradish, lyophilizate

Application

Peroxidase (POD), EIA Grade is a marker enzyme enabling peroxidation of reduced dyes in the indicator reaction producing a color, fluorimetric or luminescent derivative of the labeled molecule for further detection and quantification.

EC 1.11.1.7

Properties

Horseradish peroxidase is a 44,173.9 D glycoprotein with 4 lysine residue.

Specification

Appearance: Red-brown lyophilizate

Activity (+25°C, guaiacol, H₂O₂): ≥225 U/mg lyophilizate

Specific Activity (+25°C, ABTS, H₂O₂, pH 5.0): ≥900 U/mg lyophilizate

Purity number (A₄₀₃/A₂₇₅): 3.0-3.5

A₄₀₃ (0.2 mg/ml, against buffer): No limit

Contaminants (expressed as percentage of Peroxidase activity):

ATPase: ≤0.001

Catalase: ≤0.7

Phosphatase, acidic: ≤0.001

Isoenzyme distribution (HPLC): ≥90% (homogeneous with respect to isoenzyme C)

Amino groups: 2-3 mol/mol enzyme

Carbohydrates: 12.0-14.5% (w/w)

Stability: At -15 to -25°C within specification range for 24 months. Keep tightly sealed.

Catalog number

10 815 462 103

Pack size

custom fill

Will be supplied as "Peroxidase (POD) from Horse-radish". Unit of measure is "g".

For further processing only.

Peroxidase (POD), Grade I

from horseradish, lyophilizate

Application

Peroxidase (POD), Grade I is a marker enzyme enabling peroxidation of reduced dyes in the indicator reaction producing a color, fluorimetric or luminescent derivative of the labeled molecule for further detection and quantification.

EC 1.11.1.7

Properties

Horseradish peroxidase is a 44,173.9 D glycoprotein with 4 lysine residue.

Catalog number

10 121 606 103

Pack size

custom fill

Will be supplied as "Peroxidase (POD), Grade I, Horse-radish". Unit of measure is "MU".

For further processing only.

Specification**Appearance:** Red-brown lyophilizate**Solubility:** Clear, red-brown solution in water (c=10 mg/mL)**pH value** (c=10 mg/mL): 6.0-7.0**Activity** (+25°C, guaiacol, H₂O₂): ≥250 U/mg lyophilizate**Purity number** (A₄₀₃/A₂₇₅): 3.0-3.5**Contaminants** (expressed as percentage of Peroxidase activity):

ATPase: ≤0.001

Catalase: ≤0.7

Phosphatase, acidic: ≤0.001

Stability: At +2 to +8 within specification range for 24 months. Keep tightly sealed.**Poly Peroxidase (Poly POD), EIA Grade**

from horseradish, lyophilizate

Application

Poly Peroxidase (Poly POD) is a marker enzyme enabling peroxidation of reduced dyes in the indicator reaction producing a color, fluorimetric or luminescent derivative of the labeled molecule for further detection and quantification.

Product description

Poly Peroxidase is lyophilized in 10 mmol/L potassium phosphate, 50 mmol/l NaCl, 1 mmol/L EDTA, pH 6.1 and saccharose as stabilizer.

EC 1.11.1.7

Properties**Molecular weight:** 0.8 ± 0.2 x 10⁶ D (~ 20 POD-monomers)**Activation:** Is accomplished by MHS (Maleimidohexanoyl-N-hydroxysuccinimide ester) ≥40 MH-groups per Poly POD (MH) are accessible for conjugation with sulfhydryl groups.**Specification****Appearance:** Red-brown lyophilizate**Solubility:** Clear, red-brown solution in water (c= 5 mg/mL)**Specific activity** (+25°C, ABTS): ≥600 U/mg**MH-groups:** ≥2 (mol MH/mol POD)**Stability:** At -60 to -90°C within specification range for 48 months.**Catalog number****11 578 545 103****Pack size**

custom fill

Will be supplied as "Peroxidase, Polymerized (MH)". Unit of measure is "mg active ingredient".



For further processing only.

Chlorophenolred-beta-D-galactopyranoside (CPRG) high pure

sodium salt, powder

ApplicationUse CPRG high pure as a substrate for β -Galactosidase.**Properties****Formula:** $C_{25}H_{21}O_{10}Cl_2SNa$ **Molecular weight:** 607.4 D**Specification****Identity (HPLC):** Corresponds to specification**Appearance:** Orange-red powder**Solubility:** Clear, red colored solution in water (c=20 mg/mL)**A₅₇₈** (c=5 mmol/l water) : ≤ 0.090 **A₆₅₀** (c=5 mmol/l water; turbidity): ≤ 0.030 **Purity (HPLC, 265 nm):** ≥ 97.5 area%**Purity (HPLC, 220 nm):** ≥ 88 area%**Carbinol (HPLC):** ≤ 10 area%**Total (CPRG+Carbinol, HPLC):** ≥ 98 area%**Assay CPRG (A₄₀₅):** 80-110%**Na (flame photometric):** 3-4%**Water (K. Fischer):** $\leq 15\%$ **Thin layer chromatography:** Corresponds to reference**Chlorophenolred, free:** $\leq 0.075\%$ **Galactose, free (enzymatic):** $\leq 5.0\%$ **Reaction rate (β -galactosidase) of sample/CPRG standard:** 80-120%**Reaction rate (β -galactosidase) of sample/2-NP-galactoside:** $\geq 8.5\%$ **Stability:** At -15 to -25°C within specification range for 36 months.

Store dry. Keep tightly closed.

Catalog number**07 930 097 103****Pack size**

custom fill

Will be supplied as "CPRG high pure". Unit of measure is "g".



For further processing only.

Chlorophenolred- β -D-galactopyranoside (CPRG)

sodium salt, powder

ApplicationUse CRPG as a substrate for β -Galactosidase.**CAS:** 99792-79-7**Properties****Formula:** $C_{25}H_{21}O_{10}Cl_2SNa$ **Molecular weight:** 607.4 D**Catalog number****11 379 119 103****Pack size**

custom fill

Will be supplied as "Chlorophenolred-beta;-D-galactopyranoside (CPRG)".Unit of measure is "g".



For further processing only.

Specification**Identity** (HPLC): Corresponds to specification**Appearance:** Orange-red powder**Solubility:** Clear, red colored solution in water (c=20 mg/mL)**A₅₇₈** (c=5 mmol/L water) : ≤0.200**A₆₅₀** (c=5 mmol/L water; turbidity): ≤0.030**Purity** (HPLC, 265 nm): ≥97.5 area%**Purity** (HPLC, 220 nm): ≥88 area%**Carbinol** (HPLC): ≤10 area%**Total** (CPRG+Carbinol, HPLC): ≥98 area%**CPRG** (A₄₀₅): 80-110%**Na** (flame photometric): 3-4%**Water** (K. Fischer): ≤15%**Thin layer chromatography:** Corresponds to reference**Chlorophenolred, free:** ≤0.1%**Galactose, free** (enzymatic): ≤5.0%**Reaction rate** (β-galactosidase) of sample/CPRG standard: 80-120%**Reaction rate** (β-galactosidase) of sample/2-NP-galactoside: ≥8.5%**Stability:** At -15 to -25°C within specification range for 36 months.

Store dry. Keep tightly closed.

3,3',5,5'-Tetramethylbenzidine (TMB)

crystalline powder

Application

Use TMB as substrate solution for horse radish peroxidase (450 nm).

CAS: 54827-17-7**Properties****Formula:** C₁₆H₂₀N₂**Molecular weight:** 240.35 D**Specification****Appearance:** Yellowish to light brown crystalline powder**Melting range:** +168 to +171°C**TMB** (GC): ≥99.5 area%**TMB** (titrimetric, based on dry weight): ≥97.5%**Loss on drying** (for 2 hours at +105°C): ≤1%**Stability:** At +2 to +8°C within specification range for 24 months.

Protect from light.

Catalog number**10 203 700 103****Pack size**

custom fill

Will be supplied as "3,3',5,5'-Tetramethylbenzidine". Unit of measure is "kg".

For further processing only.

4-Aminophenyl Phosphate (pAPP)

disodium salt, powder

Application

Use 4-Aminophenyl Phosphate (pAPP), Disodium Salt as substrate for alkaline phosphatase.

CAS: 75966-16-4**Properties****Formula:** C₆H₇NO₄PNa**Molecular weight:** 211.09 D**Specification****Appearance:** Off white to brownish powder**Solubility:** Clear, fawn to brownish solution in water (c=100 mg/mL)**pH value:** 7.0-8.0**ESI-MS:** 188.011±0.005 D**pAPP (HPLC):** ≥90 area%**Water (K. Fischer):** ≤15%**Na (flame photometric):** 15-21%**p-Nitrophenylphosphate (HPLC):** ≤0.3 area%**Cl:** ≤0.1%**P_i (acid labile):** ≤0.5%**P_i':** ≤1%**Stability:** At +2 to +8°C within specification range for 24 months.

Protect from light.

Catalog number**05 642 965 103****Pack size**

custom fill

Will be supplied as "4-Aminophenyl Phosphate Disodium Salt". Unit of measure is "kg".

For further processing only.

4-Nitrophenyl Phosphate (pNPP)

disodium salt, crystalline powder

Application

Use pNPP as a substrate for alkaline phosphatase.

CAS: 4264-83-9**Properties****Formula:** C₆H₄NO₆PNa₂ × 6 H₂O**Molecular weight:** 371.1 D (pNPP: 219.1 D)**Detection:** at 405 nm**Specification****Appearance:** White to slightly yellow crystalline powder**Solubility:** Clear, colorless to slightly yellow solution in water (c=50 mg/mL)**pH value:** 9±1**Catalog number****10 004 847 103****Pack size**

custom fill

Will be supplied as "4-Nitrophenyl Phosphate, Disodium Salt". Unit of measure is "kg".

Additional formulation: Tablets are available on request.

For further processing only.

4-NPP-Na₂ x 6 H₂O (calculated from value found enzymatically): ≥95%

4-NPP (enzymatic): ≥56%

Na (flame photometric): 13±1%

Water (K. Fischer): 28±3%

4-Nitrophenol (free): ≤0.07%

P_i: ≤ 0.3%

Blank (with TC AP_{opt}, ΔA/30 min): ≤0.015

Reactions rates (AP): 100±5%

Stability: At +2 to +8°C within specification range for 24 months.

Protect from light.

5-Bromo-4-chloro-3-indolyl-phosphate (BCIP)

toluidin, crystalline powder

Application

Use BCIP as precipitating substrate for alkaline phosphatase.

CAS: 6578-06-09

Properties

Formula: C₈H₆NO₄BrClIP x C₇H₉N

Molecular weight: 433.6 D (BCIP: 326.4 D)

Detection: Forms a blue precipitate

Specification

Appearance: White to slightly yellowish microcrystalline powder

BCIP x toluidine (from N): ≥99%

BCIP x toluidine (HPLC): ≥99 area%

N (elementary analysis): 6.39-6.60%

Performance test (incubation with alkaline phosphatase, aerial oxidation): Blue precipitate

Stability: At +2 to +8°C within specification range for 24 months.

Store dry. Protect from light.

Catalog number

10 760 978 103

Pack size

custom fill

Will be supplied as "5-Br-4-Cl-3-indolyl-phosphate". Unit of measure is "g".

For further processing only.

Human Serum

frozen solution

Qualified for **cobas**[®] platforms

Application

Basic matrix for manufacturing calibrators and controls for immunoassays and assays in clinical chemistry.

Product description

Frozen solution from pooled blood donors.

Specification

Appearance of solution: Yellowish slight turbid liquid

Turbidity: ≤0.600 E

pH value (+25°C): 7.0-7.5

Protein (Biuret): ≥65 mg/mL

Cholesterol: ≥135 mg/dL

Triglyceride: 65-206 mg/dL

Ca: ≤2.2 mmol/l

Cholinesterase: ≥4,700 U/L

Creatine kinase: ≤265 U/L

Bioburden: ≤1,000 CFU/mL

Non-reactive in HBsAg: Corresponds to specification

Anti-HIV I+II: Negative

Anti-HCV: Negative

HCV NAT non-reactive: Corresponds to specification

HIV-1 NAT non-reactive: Corresponds to specification

Stability: At -15 to -25°C within specification range for 12 months.

Catalog number

11 758 225 103

Pack size

9 L (samples 100 mL)

Will be supplied as "Serum <-> Human". Unit of measure is "L".



For further processing only.

Trademarks

All mentioned trademarks are owned by or licensed to a Member of the Roche Group

ABTS
COBAS
ELECSYS
LUMIGRADE
STREPTAWELL

Other brands or product names are trademarks of their respective holders

1

1,4-Dithiothreitol (DTT)	26
--------------------------	----

3

3,3',5,5'-Tetramethylbenzidine (TMB)	243
3,5-Dichlorophenolsulfonic Acid	148
3-Hydroxy-1,2,3,4-tetrahydrobenzo[h]quinoline	26

4

4-Aminoantipyrine (4-APP)	148
4-Aminophenyl Phosphate (pAPP)	244
4-Nitrophenyl Phosphate (4-NPP)	150
4-Nitrophenyl Phosphate (pNPP)	244
4-Nitrophenyl- α -D-maltohexaoside	149

5

5-Bromo-4-chloro-3-indolyl-phosphate (BCIP)	245
---	-----

6

6-Phosphogluconolactonase	137
---------------------------	-----

 α

α -Glucosidase	95
α -Ketoglutarate (2-Oxoglutarate) for potassium test	158
α -Ketoglutarate (2-Oxoglutarate)	156
α -Ketoglutarate (2-Oxoglutarate)	157
α -Ketoglutarate (2-Oxoglutarate)	158

 β

β -Galactosidase Mutein	234
β -Galactosidase, recombinant, EIA Grade	238
β -Galactosidase	238
β -Glucuronidase	97

 γ

γ -Glutamyltransferase	103
-------------------------------	-----

A

Acetate-CoA Ligase (Acetyl-CoA Synthetase)	52
Acetyl-Coenzyme A	30
Acid Phosphatase	53
Acyl-CoA Oxidase	53
Acyl-CoA Synthetase (Long-Chain-Fatty-Acid-CoA Ligase)	52
Adenosine-5'-O-(2-thiodiphosphate)	40
Adenosine-5'-O-(3-thiotriphosphate)	40
ADP for potassium test	43
ADP	41
ADP	42
AFP Ag	199
Alanine Aminotransferase (ALT) (GPT)	54

Alanine Aminotransferase (ALT) (GPT)	54
Alcohol Dehydrogenase, chemically modified	56
Alcohol Dehydrogenase	55
Alcohol Dehydrogenase	58
Aldehyde Dehydrogenase	59
Aldose 1-Epimerase (Mutarotase)	60
Alkaline Phosphatase Mutein, recombinant	233
Alkaline Phosphatase, EIA Grade	235
Alkaline Phosphatase, recombinant, highly active, carbohydrate reduced	237
Alkaline Phosphatase, recombinant, highly active	236
Alpha-Biotin-Epsilon-Dig-Lys	195
AMP	43
AMP	44
Aprotinin	24
Ascorbate Oxidase, chemically modified	63
Ascorbate Oxidase	60
Ascorbate Oxidase	62
Aspartate Aminotransferase (AST) (GOT)	64
Aspartate Aminotransferase (AST) (GOT)	65
ATP, Grade II	45
ATP, Grade I	45

B

B2M Ag	201
Benzylidene-4-NP-G7	150
Bovine IgG (PAB\<->R-IgG)	223
Bovine Serum Albumin (BSA), Fraction V	227
Bovine Serum Albumin (BSA), Fraction V	228
Bovine Serum Albumin (BSA), reduced sodium and potassium	229
Bovine Serum Albumin IV	227
Bovine Serum Albumin I	226

C

CA 125 Ag	197
CA 15-3 Ag	195
CA 19-9 Ag	196
CA 72-4 Ag	197
CEA Ag	198
CHAPSO	23
CHAPS	22
Chlorophenolred- β -D-galactopyranoside (CPRG)	242
Chlorophenolred-beta-D-galactopyranoside (CPRG) high pure	242
Cholate	17
Cholesterol Esterase, chemically modified	68
Cholesterol Esterase, Grade II	69
Cholesterol Esterase, Grade I	66
Cholesterol Esterase	70

Cholesterol Esterase	71	Glucose-6-phosphate Dehydrogenase (G6P-DH), chemically modified	88
Cholesterol Oxidase	72	Glucose-6-phosphate Dehydrogenase (G6P-DH)	86
Cholesterol Oxidase	74	Glucose-6-phosphate Dehydrogenase (G6P-DH)	89
Chromogenic Substrate for Lipase	151	Glucose-6-phosphate Dehydrogenase (G6P-DH)	91
Citrate Lyase	75	Glucose-6-phosphate Dehydrogenase (G6P-DH)	92
Citrate Synthase	75	Glucose-6-phosphate Dehydrogenase (G6P-DH)	93
Coenzyme A, Grade II	31	Glucose-6-phosphate Dehydrogenase (G6P-DH)	93
Coenzyme A, Grade I	30	Glucose-6-phosphate Isomerase	94
Coenzyme A, Grade I	31	Glucose-6-phosphate	162
Colipase	76	Glupa-carboxylate	153
Colloidal Gold 20 nm	215	Glutamate Dehydrogenase (NAD(P))	97
Colloidal Gold 40 nm	215	Glutamate Dehydrogenase (NAD(P))	98
cOmplete, EDTA free	24	Glutamate Dehydrogenase (NAD(P))	100
Creatinase	77	Glutamate Dehydrogenase (NAD(P))	101
Creatine Phosphate	159	Glycerol Kinase (GK), concentrated	106
Creatininase	78	Glycerol Kinase (GK)	103
Creatinine Deaminase	79	Glycerol Kinase (GK)	105
CRP Ag	200	Glycerol-3-phosphate Dehydrogenase	107
D			
D(-)-Lactate	160	Glycerol-3-phosphate Oxidase, chemically modified	109
D-Biotin-N-hydroxysuccinimide ester	206	Glycerol-3-phosphate Oxidase	107
D-Biotinoyl-ε-aminocaproic acid-N-hydroxysuccinimide ester	206	Glycylglycine	10
D-Lactate Dehydrogenase (D-LDH), Grade II	117	Guanosine-5'-O-(2-thiodiphosphate)	46
D-Lactate Dehydrogenase (D-LDH), Grade I	115	H	
D-Lactate Dehydrogenase (D-LDH)	115	HAMA Serum, Type II	225
D-Mannitol	27	HAMA Serum, Type I	224
Deoxycholate	18	HBsAg (ad)	202
Di(adenosine-5'-)penta-phosphate	160	HBsAg (ay)	202
Diethanolamine 85%	10	Hepes	11
Dilaurylglycosulfate	17	Hexokinase (HK), chemically modified	112
E			
Ethylidene-4-NP-G7	152	Hexokinase (HK)	110
F			
FAD	32	Hexokinase (HK)	113
Formate Dehydrogenase	80	Human Serum	246
Framework IEP	217	I	
G			
Galactose 1-Dehydrogenase	81	Imidazole	12
Galactose 1-Dehydrogenase	81	K	
Galactose 1-Dehydrogenase	82	Kryptofix 221	28
Glucose Oxidase (GOD), chemically modified	85	L	
Glucose Oxidase (GOD), Grade II	84	L(+)-Alanine	156
Glucose Oxidase (GOD), Grade I	82	L-Lactate Dehydrogenase (L-LDH), chemically modified	121
Glucose-1,6-diphosphate	161	L-Lactate Dehydrogenase (L-LDH)	118
		L-Lactate Dehydrogenase (L-LDH)	119
		L-Lactate Dehydrogenase (L-LDH)	120
		L-Lactate Dehydrogenase (L-LDH)	122
		Lactate 2-Monooxygenase (Lactate oxidase), Grade II	124
		Lactate 2-Monooxygenase (Lactate oxidase), Grade I	122

Lactate 2-Monooxygenase (Lactate oxidase)	125
Lipoprotein Lipase, chemically modified	127
Lipoprotein Lipase	125
Lysozyme	128

M

MAB IgG2b/Fab2a Poly	217
MAB33 IgG1/Fab1 Poly	219
MAB33 IgG1/IgG1 Poly	220
MAB33 IgG1	218
MAB\<AFP\>M-LJ738 IgG	172
MAB\<AFP\>M-TU11 IgG	172
MAB\<C-Peptid\>M-RCPT3F11 IgG	188
MAB\<C-Peptid\>M-RPEP001 IgG	188
MAB\<CEA\>M-TU2 IgG	173
MAB\<CEA\>M-TU3 IgG	174
MAB\<CK-MB\>M-6.12.47 IgG	170
MAB\<CK-MB\>M-7.4.5 IgG	170
MAB\<CK-MM\>Mix	171
MAB\<DD\>M-1.2.57 IgG	174
MAB\<DD\>M-2.1.16 IgG	175
MAB\<Digit\>M-145A41 IgG	190
MAB\<Dig\>M-19-11 IgG	187
MAB\<Ferr\>M-3.170 IgG	176
MAB\<Ferr\>M-4.184 IgG	176
MAB\<FSH\>M-1.303 IgG	177
MAB\<FSH\>M-W3g IgG	178
MAB\<H-S-Amy\>M-Tu66C7 IgG	180
MAB\<H-S-Amy\>M-Tu88E8 IgG	181
MAB\<HCG\>M-INN2 IgG	178
MAB\<HCG\>M-INN22 IgG	179
MAB\<HCT\>M-11.3.59 IgG	189
MAB\<IGE\>M-323 IgG	182
MAB\<IGE\>M-7H8 IgG	182
MAB\<INSULIN\>M-BM1 IgG	183
MAB\<INSULIN\>M-ST3 IgG	184
MAB\<Keratin\>M-AE1 IgG	191
MAB\<Keratin\>M-AE3 IgG	190
MAB\<LH\>M-11412 IgG	184
MAB\<LH\>M-2.406 IgG	185
MAB\<TSH\>M-A8 IgG	186
MAB\<TSH\>M-TU1.20 IgG	186
Malate Dehydrogenase, chemically modified	129
Malate Dehydrogenase, IFCC Quality	130
Malate Dehydrogenase, IFCC Quality	131
Malate Dehydrogenase	129
MES, free acid	15
Multi Analyte Stripe	210

N

N-Acetyl-L-Cysteine	162
N-Carbamoylsarcosine Amidase	66
n-Dodecyl- β -D-maltoside	20
N-Methylhydantoinase (ATP-hydrolyzing)	132
n-Octyl β -D-glucoside	20
NAD(P)H Dehydrogenase (quinone) (Diaphorase)	133
NAD(P)H Dehydrogenase (quinone) (Diaphorase)	134
NAD, Grade II	33
NAD, Grade I	32
NADH, Grade II	35
NADH, Grade II	36
NADH, Grade I	34
NADPH	38
NADP	37
NADP	37
Nitrate Reductase	147

O

Oxalate Oxidase	134
-----------------	-----

P

PAB\<->H-IgG/Fab Poly	221
PAB\<CRP\>S IgG	192
PAB\<E2\>K IgG	192
PAB\<T3\>S IgG	193
PAB\<T4\>S IgG	194
Pefabloc SC (AEBSF)	25
Peroxidase (POD), EIA Grade	240
Peroxidase (POD), Grade II	134
Peroxidase (POD), Grade I	240
Phosphoenolpyruvate (PEP), for potassium test	165
Phosphoenolpyruvate (PEP)	163
Phosphoenolpyruvate (PEP)	164
Phosphogluconate Dehydrogenase (decarboxylating)	136
Pipes	13
Pipes	13
Polidocanol (Thesit)	21
Poly BSA Type II	231
Poly BSA Type I	230
Poly Peroxidase (Poly POD), EIA Grade	241
PSA Ag	200
Pyruvate Kinase	139
Pyruvate Kinase	140
Pyruvate Oxidase	140
Pyruvate	165

R

Rabbit IgG (PAB\<->K-IgG)	222
---------------------------	-----

S

S-Butyrylthiocholine Iodide	166
Sarcosine Oxidase	141
Sheep IgG (PAB\<->S-IgG)	221
Streptavidin Magnetic Particles	210
Streptavidin POD	205
Streptavidin R-Phycoerythrin LumiGrade Reagent	208
Streptavidin R-Phycoerythrin LumiGrade Ultrasensitive Reagent	208
Streptavidin rec. inactive, Poly	233
Streptavidin, recombinant	204
Streptavidin-Mutectin-Sepharose	204
StreptaWell, 384/C1	211
StreptaWell, F8/C1 module	213
StreptaWell, F8/C2 plus, high binding capacity	213
StreptaWell, N-breakapart C8/C1	212
StreptaWell, N-breakapart C8/C2 plus, high binding capacity	212

T

Taurodesoxycholat	19
Thio-NAD	39
Thrombin	143
TOOS	154
Tribromo-hydroxybenzoic acid	154
Triose-phosphate Isomerase	143
Tris	14
Tween 20	21

U

Urease	144
Uricase	145

V

Valinomycin	29
-------------	----